
TENS + Heat Wireless Pain Management Solution

Model 22-033

Intention for Use:

To be used for temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, arm, and leg, due to strain from exercise or normal household and work activities.

It is also intended for symptomatic relief and management of chronic, intractable pain and relief of pain associated with arthritis.

See pg 5 of this manual for further details.


**PLEASE READ THESE INSTRUCTIONS COMPLETELY BEFORE
USING YOUR TENS DEVICE**

INSTRUCTION MANUAL • ENGLISH & ESPAÑOL

STOP! PLEASE ENSURE YOU HAVE ALL OF THE FOLLOWING COMPONENTS BEFORE USING YOUR TENS DEVICE

Electric Stimulator

1 Large, 1 Small
Electrode Pads

USB Cable

Storage Bag

Instruction Manual
Replacement
Electrode Order Form

IF YOU ARE MISSING ANY PARTS, INCLUDING INSERTS OR INSTRUCTION MANUALS, DO NOT RETURN TO PLACE OF PURCHASE. CONTACT CUSTOMER CARE AT 866-326-1313.

INDEX

| | |
|--|-------|
| Introduction/Indications for Use | 4-5 |
| Care and Safety Information | 6-13 |
| Product Features..... | 14-15 |
| Charging the Device | 16-17 |
| Preparing Your Device | 18 |
| Placement & Treatment..... | 19 |
| Pulse Program Descriptions | 20-21 |
| Selecting Your Gel Pad..... | 21 |
| Instructions for Use | 22-25 |
| Caring for Electrode Pads | 26 |
| Cleaning and Storage | 27 |
| Troubleshooting | 28 |
| Device & Label Symbols..... | 29 |
| Electromagnetic Compatibility | 30-33 |
| Product Specifications..... | 34-35 |
| Warranty Information | 36-37 |
| Español..... | 39-76 |

Toll-Free Customer Care Help Line:
1-866-326-1313
Monday – Friday
8:30 a.m. – 4:30 p.m. CST

Distributed by
Veridian Healthcare
1175 Lakeside Drive
Gurnee, IL 60031
www.veridianhealthcare.com

Made in China

#93-1319 08/16
©2016 Veridian Healthcare, LLC

INTRODUCTION

Electronic Stimulator delivers electric pulses and heat generated to the user's skin through the electrodes. The portable and compact device has multiple modes of different pulse frequencies, covering Transcutaneous Electrical Nerve Stimulation (TENS) and Electrical Muscle Stimulation (EMS). It includes operating elements of POWER-ON/OFF button, intensity increase button, intensity decrease button, mode button, timer button, and heat button, and can be attached and detached to the electrode through the four magnetic connectors.

INDICATIONS FOR USE

TENS (Modes 1, 2, 4, 5, 6, 8)

To be used for temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, arm, and leg, due to strain from exercise or normal household and work activities.

It is also intended for symptomatic relief and management of chronic, intractable pain and relief of pain associated with arthritis.

EMS (Modes 1, 3, 7)

To stimulate healthy muscles in order to improve and facilitate muscle performance. To be used for the improvement of muscle tone and firmness, and for strengthening muscles in the arms, abdomen, legs, and buttocks. Not intended for use in any therapy or for the treatment of any medical conditions or diseases.

It is also intended to temporarily increase local blood circulation in the healthy muscles of lower extremities.

IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

NOTE: Read all instructions and inserts included with this device carefully before use. The following basic precautions are needed when using an electrical product.

CAUTION: Failure to read and observe all precautions could result in injury or equipment damage.

Improper care or use of your device may result in injury, damage to the unit or ineffective treatment. Following these instructions will ensure the device's efficacy and long life.

CONTRAINDICATIONS :

- Do not use this device if you have a cardiac pacemaker, implanted defibrillator or other implanted metal or electronic device.
- Do not use this device to treat undiagnosed chronic pain.
- Do not use this device if you have suspected or diagnosed heart problem.
- Do not use this device if you have suspected or diagnosed epilepsy.
- Do not use this device if you have a tendency to bleed internally following an injury.

IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

- Do not use this device during pregnancy.

GENERAL WARNINGS PRIOR TO USE

- If you are under the care of a Physician, consult with your Physician before using this device.
- Long-term effects of this device are not known.
- Do not place the device on or close to your heart.
- Take care when applying device close to the neck; do not place this device on the front or sides of the neck. Severe spasm of the muscles may occur and the contractions may be strong enough to close the airway or cause difficulty in breathing. Stimulation over the neck could also affect hearing or blood pressure.
- Do not apply stimulation across the chest because the introduction of electrical current into the chest may cause rhythm disturbances to the heart.
- Do not place the device on or around your head. The effects of stimulation of the brain are unknown.
- Do not use this device while sleeping.

IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

- Do not use if you feel numbness.
- Do not use this device in or close to water.
- Use only on normal, healthy, clean skin. Do not use on or close to open wounds or rashes, or over swollen, red, infected or inflamed skin.
- Do not apply stimulation over, or in proximity to, cancerous lesions.
- Do not use the device on children or incapacitated persons.
- Consult with your physician before using this device, because the device may cause lethal rhythm disturbances to the heart in susceptible individuals.

GENERAL CAUTIONS PRIOR TO USE

- Read this manual before using this device for the first time.
- Keep this manual available whenever you use this device.
- This device is intended for individual person use only.
- This device is not effective for pain associated with Central Pain Syndromes, such as headaches.

IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

- This device is for pain caused by muscle soreness, and should be placed only around muscles where pain originates.
- The pain may indicate that you have some other health problem. You should know the reason and source of your pain before using these devices.
- The safety of using this device during pregnancy or birth has not been established.
- The effectiveness of this device depends greatly on a person's individual physical condition. It may not always be effective for every user.
- If you have had medical or physical treatment for your muscle pain, consult with your treatment provider before using this device. You should contact your physician prior to using this device following recent surgical procedures. Stimulation may disrupt the healing process.

OPERATING WARNINGS & CAUTIONS

CONSULT WITH YOUR PHYSICIAN IF YOU HAVE ANY OF THE FOLLOWING:

- If you have suspected or diagnosed heart problem.

IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

- If you have suspected or diagnosed epilepsy.
- If you have a tendency to bleed internally following an injury.
- If you recently had surgery, or have ever had surgery on your back.
- If areas of skin lack normal sensations, such as skin that tingles or is numb.
- During menstruation or during pregnancy.
- The unit is intended for the temporary relief of pain caused by muscle soreness. The source or cause of the pain should be addressed by a healthcare professional. Pain could indicate an injury or condition that requires treatment.
- Some people may feel skin irritation or experience a very sensitive feeling in the skin due to electrical stimulation. If this occurs, stop using these devices and consult your physician.
- If skin under the pad feels irritated after using the stimulator for a long period of time, use the stimulator for a shorter period of time.
- Minor redness at stimulation placement is a normal skin reaction. It is not considered a skin irritation, and it will normally disappear within 30 minutes after the electrodes are removed.

IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

If the redness does not disappear after 30 minutes from the removal of electrodes, do not use the stimulator again until after the excessive redness has disappeared.

- Turn off your device if the stimulation feels unpleasant or does not provide pain relief.
- Keep this device out of the reach of children.
- Use this device only with the pads and accessories recommended by the manufacturer.
- Do not use this device when driving, operating machinery or when swimming.
- Before removing the Electrode Pad, be sure to turn it OFF, avoiding unpleasant stimulation.
- TENS is not a substitute for pain medications and other pain management therapies.
- TENS is a symptomatic treatment and, as such, suppresses the sensation of pain that would otherwise serve as a protective mechanism.
- Do not use the unit while sleeping.

IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

- After each use, clean any adhesive residue left on the skin with soap and water.
- The effectiveness of stimulus therapy varies from person to person. It may not be an effective treatment for all users.
- The long-term effects of stimulus therapy are unknown.

STORAGE WARNINGS & CAUTIONS

- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: 14°F – 122°F (-10°C – 50°C); humidity: 30–90% RH.
- Keep the unit out of reach of small children.
- Do not store the device in direct sunlight, dusty or humid environments, or extreme temperatures. Exact figures for appropriate use and storage can be found in the Product Specifications section of this manual.

IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS



CLEANING WARNINGS & CAUTIONS

- Never immerse the unit in water to clean as it may damage the unit.
- Follow the 'Cleaning and Maintenance' portion of this manual for instruction on how to clean and care for your device.

PRODUCT FEATURES


Electrical Stimulator

PRODUCT FEATURES


Large Electrode Pad

Magnetic Electrode
Contact Points

Small Electrode Pad

Magnetic Electrode
Contact Points


CHARGING THE DEVICE

This device includes an internal rechargeable battery within the electrical stimulator unit and USB cable.

A USB wall charger/adaptor is required for charging; not included.

Before first use, the device should be fully charged. Charging time may vary depending on the wall adapter used; 3–6 hours from depleted battery to full charge.

1. Connect the wide end of the USB cable to a wall adapter.
2. Connect the smaller mini USB to the TENS stimulator.
3. During charging the battery symbol on the stimulator will flash. When the device is fully charged, the battery symbol will be solid green.
4. Unplug the charger and stimulator when fully charged; never operate the unit while charging.

CHARGING THE DEVICE

NOTE:

The TENS + Heat device offers 3 user modes, 8 pulse modes and 20 levels of intensity . As such, the power usage and rechargeable battery life will vary depending on the users specific settings during treatment.

For example, a fully charged device utilizing mid-range intensity and only the TENS setting battery life may extend to 20 hours of use. The same setting with the heat function activated may extend to approximately 3 hours of use.

For maximum effectiveness, monitor the battery life symbol and charge the device when the battery level graphic is depleted.




PREPARING YOUR DEVICE

It is necessary to prepare and assemble your device prior to treatment.

NOTE: Save the plastic film for storage after treatment to protect the pads from dirt and debris.

1. While the stimulator device is OFF (the display is NOT illuminated) carefully attach the electrode pads to the control unit via the magnetic connection points.
2. Pull gently to ensure the electrode pad is firmly attached.
3. Carefully remove the plastic film from the electrode pad.
4. You are ready to apply the pad to the affected treatment area.

Do not turn the power ON until the device is in place on the user's skin as indicated in the 'Instructions for Use' section of this manual.


PLACEMENT & TREATMENT

Read all warnings and cautions completely before use.

The electrode pads are applied directly to the skin. Clothing may be worn over the device; clothing should not constrict the device and should allow access to the control unit buttons.

Targeted treatment areas may include:

- shoulder
- elbow
- thigh
- knee
- calf
- ankle
- back of neck
- upper back
- lower back
- feet

Never place device on over the heart, on the front or sides of neck or the face.

PULSE PROGRAM DESCRIPTIONS

| | Pulse Rate | Description | Benefit |
|----|--------------------------|--|---|
| P1 | Combination of the below | Combination of the below | All kinds of pain relief |
| P2 | 69 | Pulse on for 3.4 second and off for 1.6 second | Mild gate pain relief |
| P3 | 12.5–55.5 | Pulse on for 20 second and off for 1 second | Muscle performance improvement with variable pulse rate |
| P4 | 1.2 | Pulse on every 0.85 second | Slow endorphin pain relief |
| P5 | 100 | Pulse on for 10 second and off for 2.5 second | Quick gate pain relief |
| P6 | 100 | Pulse on for 20 second and off for 1 second | Quick gate pain relief and strong pulse amplitude |
| P7 | 20 | Pulse on for 5 second and off for 1 second | Muscle performance improvement |
| P8 | 160 | Pulse on for 10 second and off for 2 second | Quick and deep gate pain relief |

PULSE PROGRAM DESCRIPTIONS


The programs are not specific to treatment at a certain body location or to treat a specific ailment or pain.

The sensations may include—massage, acupuncture, tapping, scraping or a combination of the sensations at variable durations and cycles.

Select the program that provides the best relief for your particular ache and pain.

SELECTING YOUR GEL PAD

This device includes two sizes of gel electrode pads. There is no functional difference in the pad. Select the size for use most appropriate to the size of the treatment area.

INSTRUCTIONS FOR USE


Read all warnings and cautions completely before use. Recommended treatment time is 30-minutes; consult your healthcare professional for a treatment cycle that best suits your injury. Discontinue use immediately if treatment is uncomfortable or painful.

1. Treatment area should be clean before applying the control unit. Note that excessive skin lotion or soap residue may affect the adhesiveness of the pad. Only use on HEALTHY skin areas with no cuts, sores or wounds.
2. Apply the pad to the desired treatment area, taking care to ensure the controller is comfortably within reach of your hand for safety and control during treatment.
3. The electrode pad should be smoothly applied and the gel adhesive should hold the pad firmly in place on the skin.
4. After the device is applied to the desired treatment area, press and hold the **POWER** button for 2-seconds; the display will illuminate.
5. The TENS unit will begin to operate at the lowest intensity level setting for program 1. The HEAT feature will default to off when the unit is first powered on. The default TIMER setting is 30 minutes.

INSTRUCTIONS FOR USE

6. Until you are familiar with use of the unit, begin with the lower intensity levels and gradually progress to more advanced levels as you feel comfortable.

7. **INTENSITY SETTING**—Press and release the '+' button to find your desired therapy sensation level. This device offers 20 intensity levels. Press and release the '-' button to reduce the intensity level.

8. **MODE SETTING**— To change the mode, while the unit is on/operating, press and release the 'M'; the mode will switch to the next in the sequence of 8 pulse types. See page 20 for a detailed description of each pulse type. When changing modes, the intensity will revert to level 0 (lowest sensation).

9. **HEAT SETTING**— To activate the HEAT setting, press the '+/- HEAT'. The HEAT sensation has two levels, press and release the '+/- HEAT' to scroll through the settings: 00—OFF; 01—level 1; 02—level 2.

10. **TIMER SETTING**— Press and release the 'T' to adjust the therapy time in 10-minute increments; the unit will automatically shut off when the timer is done.

INSTRUCTIONS FOR USE

11. The sensation should be strong but comfortable; press the + or – button at any time during the treatment to change the intensity as desired.
12. Press and hold the **POWER** button for 2 seconds to turn the device off at any time during treatment.
13. Be sure the device is OFF (no sensation is felt) before removing the device from your skin. Take care during removal to best protect the adhesive and preserve the life of the electrode pad.
14. Place the plastic cover over the gel side of the pad to prevent them from drying out.
15. To conserve the battery and prevent shocks, the device will automatically power off when not in use.

NOTE:


You may wish to clean the treatment site after use to remove any adhesive residue from the electrode pad.

INSTRUCTIONS FOR USE

DISPLAY LOCK FEATURE:

The display may be locked to prevent the user from accidentally pressing the touch-sensitive buttons.

1. While the device is active, press and release the **POWER** button.
2. A LOCK symbol will appear on the upper left of the display.
3. To release the LOCK, press and release the **POWER** button again. The modes/timer/intensity can now be adjusted.

CARING FOR THE ELECTRODE PADS

Gel Electrode Pads

- The electrode pads can be reused as long as they retain their stickiness. Cover the pads with the included film pieces when not in use to prevent them from drying out
- Storing the pads in a plastic zip bag will help prevent the gel pads from drying out prematurely
- Properly cared for, the electrode pads may last up to 50 uses; the life of the electrode pad is largely dependent upon proper care and the condition of the skin as applied for use
- To help prolong the life of the electrode pad, between uses, moisten the gel side of the pad with a drop of water and allow to dry; replace the film for storage
- If you notice that the pads have become dry, no longer adhere well to the skin, or begin to irritate the skin, they are in need of replacement.

Replacement pads are available through Veridian Healthcare at 866-326-1313, by using the enclosed order form or through many major retailers and ecommerce sites.

CLEANING AND STORAGE

Cleaning

- Unsnap the electrode pad before cleaning the device
- The Electrical Stimulator can be cleaned with a soft, dry cloth
- Never use cleaning agents or excessive water to clean any part of the unit. Do not submerge any part of the unit or let it get wet
- Never disassemble the unit or attempt to repair the control unit

Storage

- Always use the included plastic film to protect the gel electrode from dust and debris
- If the device will not be used for several days or longer, place the film cover over the electrode pads to prevent drying out or getting dirt and debris on them
- Store all components in the included bag while not in use
- Storing the pads in a plastic zip bag will help prevent the gel pads from drying out prematurely
- Protect the unit from mechanical shock or heavy impact
- Avoid exposure to extreme temperatures

TROUBLESHOOTING

| Problem | Solution |
|--|---|
| Display does not illuminate/ is weak / Control unit will not turn on | Recharge the batteries |
| Intensity output is low / no stimulation is occurring | Electrode pad may be dirty; clean and attempt again. Confirm that the electrode snaps are firmly connected to the control unit. Battery power is low, recharge the battery. Device may not be applied properly; ensure the gel pad is making complete contact with the skin. |
| Electrode pads are not adher- ing to the skin | Review the 'Caring for Electrode Pads' section of this manual; pads may need to be replaced |

DEVICE & LABEL SYMBOLS

These symbols may appear on your device, instructions or packaging and may vary by make and model.

| Symbol | Meaning |
|---|---|
|  | Read This Manual/Consult Instructions Before Use |
|  | Caution/consult accompanying documents before use |
|  | Use by Date |
|  | Date of Manufacture |
|  | Batch Code |
|  | Catalogue Number |
|  | Serial Number |

| Symbol | Meaning |
|---|---------------------------------|
|  | Manufacturer |
|  | Temperature Limitation |
|  | Humidity Limitation |
|  | Non-sterile |
|  | Fragile, Handle with Care |
|  | Keep Device Dry |
|  | Product Packaging is Recyclable |

ELECTROMAGNETIC COMPATIBILITY

Table 1

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration - electromagnetic emissions

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

| Emmission Test | Compliance | Electromagnetic Environment Guidance |
|---|------------|---|
| RF emissions CISPR 11 | Group 1 | This monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF emissions CISPR 11 | Class B | This monitor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2 | Class A | |
| Voltage fluctuations/ flicker emissions IEC 61000-3-3 | Complies | |

ELECTROMAGNETIC COMPATIBILITY

Table 2

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

| Immunity Test | IEC 60601 Test Level | Compliance Level | Electromagnetic Environment Guidance |
|--|---|---|--|
| Electrostatic discharge (ESD) IEC 61000-4-2 | ± 6 kV contact ± 8 kV air | ± 6 kV contact ± 8 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %. |
| Electrical fast transient/burst IEC 61000-4-4 | ± 2 kV for power supply lines ± 1 kV for input/output lines | ± 2 kV for power supply lines | Mains power quality should be that of a typical commercial or hospital environment. |
| Surge IEC 61000-4-5 | ± 1 kV line(s) to line(s) ± 2 kV line(s) to earth | ± 1 kV line(s) to line(s) ± 2 kV line(s) to earth | Mains power quality should be that of a typical commercial or hospital environment. |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | <5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s | <5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s | Mains power quality should be that of a typical commercial or hospital environment. If the user of this monitor requires continued operation during power mains interruptions, it is recommended that this monitor be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 3 A/m | 3 A/m | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |

ELECTROMAGNETIC COMPATIBILITY

Table 3

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING Guidance and manufacturer's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

| Immunity Test | IEC 60601 Test Level | Compliance Level | Electromagnetic Environment Guidance |
|-------------------------------|--------------------------|------------------|---|
| Conducted RF IEC 61000-4-6 | 3 Vrms 150 kHz to 80 MHz | 3 Vrms | <p>Portable and mobile RF communications equipment should be used no closer to any part of this monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> $d = 1.2 \sqrt{\frac{P}{f}}$ $d = 1.2 \sqrt{\frac{P}{f}} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3 \sqrt{\frac{P}{f}} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,[*] should be less than the compliance level in each frequency range.^{**}</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p> |
| Radiated RF IEC 61000-4-3 | 3 V/m 80 MHz to 2,5 GHz | 3 V/m | |

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

* Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which this monitor is used exceeds the applicable RF compliance level above, this monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating this monitor.

** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

ELECTROMAGNETIC COMPATIBILITY

Table 4

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING
Recommended separation distances between portable and mobile RF communications equipment and this monitor

This monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this monitor as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter W | Separation distance according to frequency of transmitter m | | |
|--|--|---|--|
| | 150 kHz to 80 MHz $d = 1.2 \sqrt{P}$ | 80 MHz to 800 MHz $d = 1.2 \sqrt{P}$ | 800 MHz to 2,5 GHz $d = 2.3 \sqrt{P}$ |
| 0,01 | 0,12 | 0,12 | 0,23 |
| 0,1 | 0,38 | 0,38 | 0,73 |
| 1 | 1,2 | 1,2 | 2,3 |
| 10 | 3,8 | 3,8 | 7,3 |
| 100 | 12 | 12 | 23 |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

PRODUCT SPECIFICATIONS

| | |
|-------------------------------|---|
| NAME | TENS + Heat Pain Management Solution |
| MODEL NUMBER | 22-033 |
| NUMBER OF MODES | 8 modes |
| TIMER RANGE | 10–60 (minutes) |
| HEAT SETTING LEVELS | Level 1: 98.6°F (37°C); Level 2: 104°F (40°C) |
| WAVEFORM | Biphasic |
| SHAPE | Rectangular |
| CHANNEL | Single Channel |
| MAXIMUM OUTPUT VOLTAGE | 64V@500Ω |
| MAXIMUM OUTPUT CURRENT | 128mA@500Ω |
| PULSE DURATION | 100μSec |
| MAXIMUM FREQUENCY | 160Hz |
| PULSE INTENSITY | 0~84mA +/- 20%; 0~20 stages adjustable |
| AUTOMATIC SHUT-OFF | After completion of set therapy treatment time; 10–60 minutes maximum |
| POWER SOURCE | 3.7V Battery |
| OPERATION ENVIRONMENT | Temperature 50 F - 104 F (10C - 40C); Humidity 30% - 90% |
| STORAGE ENVIRONMENT | Temperature 14F - 122 F (-10 C - 50 C); Humidity 30% - 90% |
| DEVICE DIMENSIONS | 3.5" x 3.03" x .70" |

PRODUCT SPECIFICATIONS

| | |
|-----------------------------|--|
| WEIGHT | Controller weight 1.3 oz |
| INCLUDED | Electrical stimulator, large gel electrode pad, small gel electrode pad, USB cable, instruction manual, storage bag, replacement electrodes order form |
| AVAILABLE SEPARATELY | Replacement electrodes 2-pack, 22-036 |

Specifications are subject to change without notice

1-YEAR LIMITED WARRANTY

Congratulations on your purchase of the TENS + Heat Pain Management Solution. Your TENS + Heat Pain Management Solution is covered by the following limited warranty commencing upon the date of purchase, and subject to the following terms and conditions:

The Warrantor warrants that its TENS + Heat Pain Management Solution will be free from defects in materials and workmanship under normal consumer use for a period of one year for the original purchaser of the product.

Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the TENS + Heat Pain Management System for commercial purposes or subjecting the TENS + Heat Pain Management System to abnormal usage or conditions; or (d) other acts which are not the fault of the Warrantor, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the TENS + Heat Pain Management Solution.

If the TENS + Heat Pain Management Solution fails to conform to this limited warranty, return the TENS + Heat Pain Management Solution postage prepaid to: **Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031.** When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone number. Carefully

1-YEAR LIMITED WARRANTY

package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. At our option, the Warrantor will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will be responsible for parts charges, if any, and repair charges not covered by this limited warranty.

EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL THE WARRANTOR BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

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