## Frequently Asked Questions & Important Information in Understanding Blood Pressure and Your Digital Blood Pressure Monitor

### General Blood Pressure Information

#### What is Blood Pressure?

Blood pressure is the force that is exerted by blood pushing against the walls of the blood vessels throughout the body. The heart is the center of the circulatory system. The heart provides the force that allows the blood to circulate. Blood is forced through blood vessels increasing the pressure during every heartbeat, marking the highest point of pressure in the cycle known as the systolic pressure. The lowest point occurs between beats while the heart rests, this point during the decrease in pressure is referred to as diastolic pressure.

This complete series of events, which occurs in a single heartbeat, is known as the cardiac cycle. Blood pressure is measured in millimeters (mm) of Mercury (Hg). These two measurements are commonly noted with the systolic pressure first, followed by a slash and the diastolic pressure – i.e. 120/80, or ‘120 over 80’.

#### How can home monitoring can help?

Physician office visits can be stressful for a patient. A temporary occurrence of elevated blood pressure at the physician's office is commonly referred to as “White Coat Syndrome”. Home blood pressure monitoring provides an avenue to supplement physician's office measurements.

Home readings can show an accurate indication of blood pressure when taken over a period of time. These records can assist in evaluating your health and in making important decisions in the diagnosis and treatment of a condition when used in conjunction with physician care and evaluation.

#### What is high blood pressure?

Many factors can influence blood pressure measurements – genetics, age, sex, physical activity, anxiety, muscular development, certain medications, altitude or even the time of day. Generally, high blood pressure is blood pressure that is consistently elevated above 140/90 mmHg.

Standards for assessment of high or low blood pressure have been established by the World Health Organization (WHO) as shown on this chart. Veridian Healthcare blood pressure monitors include comparison charts to the WHO standards for your reference.

CATEGORY	SYSTOLIC (mmHg)	DIASTOLIC (mmHg)
<b>Optimal</b>	---<120	<80
<b>Normal</b>	120-129	80-84
<b>High - Normal</b>	130-139	85-89
<b>Stage 1 - Mild</b> Borderline Hypertension	140-159	90-99
<b>Stage 2 - Moderate</b> Hypertension	160-179	100-109
<b>Stage 3 - Severe</b> Hypertension	≥180	≥110
<b>Isolated Systolic Hypertension</b>	≥140	<90
	140-159	<90
	160-179	<90
	≥180	<90

*This chart is only a general guideline. Your physician or healthcare professional will help determine your individual normal blood pressure range.*

### Using a Home Digital Blood Pressure Monitor

#### How accurate is a digital blood pressure monitor?

Veridian Healthcare Digital Blood Pressure Monitors are designed to meet the exacting standards set by the various regulatory bodies that govern all blood pressure monitors including those found in your physician's office.

In order to ensure the accuracy of your measurements, it is strongly suggested to read the complete guidebook and other support materials and follow the recommended operating conditions for the unit including the suitable arm/wrist circumference range for the cuff.

#### Are manual inflate digital blood pressure monitors as accurate as automatic inflate monitors?

Yes. Both styles comply with the same stringent accuracy standards. The only difference is the way the unit is inflated.

#### How often and when should I measure my blood pressure?

It is recommended that you consult with your healthcare professional for the time and frequency that is best suited for you. It is important to take your readings at similar times and conditions when monitoring on a daily basis. This will allow for reliable comparisons of your readings.

#### Why are my home readings different from my doctor's readings?

Many factors affect blood pressure. The anxiousness of a doctor's visit may cause the common condition referred to as “White Coat Syndrome”, which typically produces a higher than “normal” blood pressure reading. Other factors, including the time of day, can influence blood pressure reading results.

#### How long should I wait between measurements?

It is important to wait a minimum of 15 minutes between readings. The wait time may vary depending on individual physiological characteristics.

#### What happens if I do not place the cuff at heart level?

If the cuff is not at heart level, readings will be affected producing inaccurate measurements.

#### Does age affect blood pressure?

Many factors such as age, genetics, sex, altitude, physical activity, muscular development, certain medications or even the time of day can influence blood pressure.

#### How does my blood pressure monitor detect the readings?

The method of measurement that a monitor uses is called the oscillometric method. The monitor detects your blood's movement through the artery in your arm/wrist and converts the movements into a digital reading. The oscillometric method does not require a stethoscope, making the monitor easy to use.

#### Which arm/wrist monitor should I use to take my blood pressure?

Your unit and accompanying detailed instruction guidebook are designed for left arm/wrist use. If you are unable to or should not use your left arm/wrist, modify the instructions for right arm/wrist use. The same arm/wrist should be used for all readings to ensure reliable comparisons.

#### What should I do if my arm circumference is outside of the stated cuff size range?

The prescribed range for your blood pressure monitor is included on the specifications page of your instruction guidebook. Reading accuracy can be affected if an improperly sized cuff is used. Most models offer additional cuff sizes available for order, please consult your instruction manual for further details.

#### What should I do if my wrist circumference is outside of the stated range?

(For wrist cuff models only) The prescribed range for your blood pressure monitor is included on the specifications page of your instruction guidebook. Unit accuracy can be affected if an improperly sized cuff is used. Due to the design of wrist monitors, the cuff is an integral part of the unit and cannot be adapted to other wrist sizes. It is suggested that you consider an arm cuff unit as an alternative for your blood pressure monitoring.

#### What do the error messages mean?

Each monitor has individual error messages that are unit specific. Consult the instruction manual for complete details.

#### Can I use rechargeable batteries?

No, your unit was designed and tested using alkaline batteries.

#### Can I use a blood pressure monitor while exercising?

The oscillometric method of blood pressure monitoring requires quiet, stable conditions. Movement, vibrations or other activity will interfere with the reading and likely cause an error or inaccurate reading. You should always rest a minimum of 15 minutes prior to taking a reading.


#### Can I use a blood pressure monitor in moving vehicles?

No, movement, vibrations or other activity will interfere with the reading and likely cause an error or inaccurate reading.

#### Is a blood pressure monitor waterproof?

Monitors are not waterproof. Refer to detailed guidebook for care and maintenance instructions. It is strongly recommended that the unit be stored and operated within the conditions expressed in the guidebook.

**Toll-Free Customer Care Help Line:**  
**1-866-326-1313**  
**Monday – Friday 8:30 a.m. – 4:30 p.m. CST**

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