## Important Information in Understanding Respiratory Care and Your Nebulizer

# FAQ's

## General Respiratory Care Information

### What is a nebulizer?

A nebulizer is a machine that sprays a fine mist of medicine or moisture into your lungs at normal air pressure. Nebulizers are generally used for short periods several times a day for both a short term illness and chronic illness. The term 'nebulizer' commonly refers to the whole piece of medical equipment but is actually the piece within the medication cup (as illustrated in your Instruction Manual).

### How does nebulizer therapy work?

The physician prescribed medicine is placed in the nebulizer and either a pump or ultrasonic vibrations are used to convert the liquid medication to an aerosol mist, which is then inhaled. This method of delivery is beneficial in the treatment of asthma and other respiratory ailments, as the medicine is administered directly to the bronchial passages of the lungs to provide quick and effective relief.

### What medications can be used with my nebulizer?

Only medications prescribed by a doctor should be used in this machine. Never use non-prescription/over-the-counter medications in a nebulizer. Most common respiratory medicines can be used with your nebulizer, but some cannot. Your healthcare professional or pharmacist can help guide you in the selection of a nebulizer to suit your medical needs.

### How long should a single treatment last?

Treatment time varies due to the amount and type of medication, as well as other factors. Typically a treatment lasts about 7 to 13 minutes.

### What is the difference between a compressor nebulizer and an ultrasonic nebulizer?

A compressor nebulizer use a pump to create air pressure, converting the liquid medication into a mist, which is sprayed into the mouthpiece. An ultrasonic nebulizer creates ultrasonic vibrations to create the mist. The ultrasonic nebulizer creates a smaller and more consistent particle size, allowing the medicine to penetrate deeper into the lungs for greater benefit. Ultrasonic nebulizers are not suitable for all types of medications, your healthcare professional or pharmacist can help guide you in the selection of a nebulizer to suit your medical needs.

### What is the significance of micron sizes and why does it vary between products?

Microns vary as a result of the method of aeration, the compressor's pressure output, and other factors. The ideal particle size is between 1 and 5 microns, which allows the medication to reach the lower airways. Particles from 5 to 10 microns reach the trachea and large bronchi. Particles of 10 microns or more leave the medicine in the mouth and pharynx. Compressors generally have higher micron sizes than ultrasonic units.

## Using a Nebulizer

*Always consult the included Instruction Manual and your physical or respiratory care therapist for detailed instruction on how to best use your nebulizer and best treat your illness.*

### Why does medication remain in the cup after treatment?

Due to surface tension of the liquid medicine, some medication may remain in the nebulizer after treatment is completed. This is normal. If the vapor begins to slow and there is still medication in the cup, gently tap the side of the cup to settle the liquid at the bottom of the cup for nebulization. Once you hear a significant difference in the sound of the nebulizer, and the vapor rate has significantly slowed, you may stop the treatment. Always discard leftover medication.

### Can I use non-prescription cough medicine in my nebulizer?

No. Never use any liquid in the nebulizer other than those prescribed by your doctor. Liquids such as cough medicines or essential oils can harm both the machine and the patient.

### Should any special steps be taken when using the nebulizer with a child?

Children should never use the nebulizer without adult supervision. For use with babies and children the mask should cover the mouth and nose completely to ensure effective inhalation.

### Why did the nebulizer stop producing mist?

This can occur for a number of reasons:

- There may be too much or too little medicine in the nebulizer. The amount of liquid should be between 0.5 ml and 5.0 ml.
- The hose connectors may be blocked and need to be cleaned.
- The nebulizer may not be held upright; always hold the nebulizer vertically during use.
- The medication may not be suitable for nebulization, check with your doctor.
- The filter may be in need of replacement.
- The air hose may be kinked.

### Why does the nebulizer get so noisy?

During normal operation, the motor of the unit produces about 60 decibels of sound. If the unit sounds much noisier than usual, discontinue use and contact Veridian Healthcare.

### What does it mean if the compressor smells like it is burning?

A burning smell may be present if moisture has been allowed into the unit. If this occurs, unplug the unit and wait 30 minutes, then try again. If the burning smell continues, unplug the unit and contact Veridian Healthcare.

### Which provides a better treatment, the mask or the mouthpiece?

Both the mask and the mouthpiece deliver effective treatment. In any case when the nebulizer cannot be held level throughout treatment, such as with children, the mask should be used to prevent the liquid medication from spilling.

### Why does treatment take longer with the mask than with the mouthpiece?

This is normal. Less medication is delivered per breath with the mask because the aerosol is combined with ambient air through the holes in the mask.

### Should the mouthpiece remain in the mouth for the entire treatment?

Yes. Keep the mouthpiece in the mouth during the entire treatment.

### What position should the nebulizer be in during use?

The compressor unit should be placed on a flat, solid surface. The nebulizer should be held vertically throughout treatment to avoid spilling the liquid medicine. Ultrasonic units may be hand-held, consult your individual unit Instruction Manual for further information.

### Can I share the nebulizer with another person?

No. Each user must have their own nebulizer to ensure proper hygiene. The compressor unit may be shared but the nebulizer, mouthpiece, and masks may not.

### Can treatments be done nasally?

Compressor treatments must be done orally. Ultrasonic treatments may be done nasally, consult with your physician for more information regarding your specific medication.

## General Product and Safety Information

### What is the capacity of the medication cup?

The compressor units hold up to 13 ml and the ultrasonic units hold up 8 ml. Larger cups allow more medication to be used without refilling the cup.

### What type of motor does a compressor unit use?

Compressor nebulizers use either a piston compressor or a diaphragm compressor. All Veridian Healthcare nebulizers use piston compressors, which are preferred due to increased durability.

### What is nebulization rate?

Nebulization rate measures the speed with which the medication is converted from a liquid to an aerosol mist. The nebulization rate determines how long your treatment will take.

### What is MMAD?

MMAD stands for Mass Median Aerodynamic Diameter. This is a measurement of the particle size produced by the nebulizer.

## Cleaning, Disinfecting and Maintenance

### How do I clean the compressor unit?

When necessary, wipe the compressor with a damp cloth and allow it to air dry. Never use soap pads or powder cleanser to clean any part of the unit.

### How often should the nebulizer be cleaned and disinfected?

The nebulizer should be cleaned after each use and disinfected each day after the last treatment of the day.

### How do I clean the nebulizer and mouthpiece?

1. Disconnect the air tube, nebulizer and mouthpiece.
2. Twist the top portion of the nebulizer counterclockwise to remove it from the bottom portion.
3. Rinse the nebulizer and the mouthpiece with water.
4. Allow the nebulizer and mouthpiece to air dry, or dry them with a clean, soft towel.
5. Once completely dry, reassemble the nebulizer, mouthpiece and air tube.

### How do I disinfect the nebulizer and mouthpiece?

The nebulizer should be disinfected once daily, after the last use of the day, unless a different disinfection regimen is specified by a healthcare professional. Disinfect the nebulizer as follows:

1. Wash the nebulizer and mouthpiece or mask in warm water, using a mild detergent.
2. Rinse the nebulizer and mouthpiece or mask with hot water.
3. Mix one part white vinegar with three parts distilled water in a quantity sufficient to fully immerse the nebulizer and mouthpiece or mask. Submerge the nebulizer and mouthpiece or mask in the water and vinegar mixture for thirty minutes.
4. Allow the nebulizer and mouthpiece to air dry, or dry them with a clean, soft towel.
5. Once completely dry, reassemble the nebulizer, mouthpiece and air tube.

### Can I wash the air filter?

No, do not attempt to clean the filters or reuse a used filter. Replace the filter every thirty days, or sooner if it turns gray. Use only Veridian Healthcare filters; do not substitute cotton or any other materials for the filters. Never operate the unit without a filter.

### Is it OK if moisture accumulates in the air tube?

Condensation may collect in the hose. If this occurs, the moisture must be completely removed, following these steps:

1. Disconnect the hose from the nebulizer. Keep it plugged into the hose connector on the unit.
2. Turn the power on and allow the compressor to run until all moisture has evaporated.

### How often should the nebulizer be changed?

The nebulizer and mouthpiece parts of the unit should be replaced once a year, whenever they appear dirty or if cracks in the plastic appear. Regular wear on the nebulizer can affect the composition and delivery of the medicine. Regular replacement of the nebulizer is also hygienic.

### How can I obtain replacement parts or accessories?

All units include a standard angled mouthpiece, nebulizer and tubing. Additional styles of mouthpiece, replacement parts and accessories are available for order on the enclosed order form, by calling Customer Care at 866-326-1313 or at your local medical supply dealer.

## Warranties and Support

### What is the warranty for my nebulizer?

Commonly, Veridian Healthcare warranties compressor nebulizers for five years, and ultrasonic nebulizers for three years. Consult the product packaging or instruction manual for detailed information on the warranty for your particular model.

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©2012 Veridian Healthcare, LLC