# SMARTHEART™

---

## Automatic Digital Blood Pressure Wrist Monitor

Model # 01-540

### **INSTRUCTION MANUAL** **ENGLISH AND SPANISH**

Please read this instruction manual  
completely before operating this unit.


**STOP!**  
**PLEASE ENSURE YOU HAVE ALL OF THE**  
**FOLLOWING COMPONENTS BEFORE USING**  
**YOUR DIGITAL BLOOD PRESSURE MONITOR**


Digital Monitor

2 AAA Batteries

Storage Case

Instruction Manual  
Quick Start Guide  
Blood Pressure Reading Tracking Log  
FAQ Information Sheet

**IF YOU ARE MISSING ANY PARTS, INCLUDING INSERTS OR INSTRUCTION MANUALS, DO NOT RETURN TO PLACE OF PURCHASE. CONTACT CUSTOMER CARE AT 866-326-1313.**

---

# INDEX

---

Care & Safety Information .....	4-6
Introduction & Indications for Use .....	7
Blood Pressure Monitor Features .....	8-9
Getting Started .....	10
Battery Installation .....	11
Setting Date & Time .....	12-13
Fitting & Applying Your Cuff .....	14-15
Taking Your Blood Pressure Reading .....	16-17
Interpreting Your Results .....	18-19
Memory Function .....	20-21
Care & Maintenance .....	22-23
Device & Label Symbols .....	24
Display Symbols .....	25
Troubleshooting .....	26
Error Codes .....	27
Technical Alarm Description .....	27
FCC Statement .....	28
Electromagnetic Compatibility .....	29-31
Specifications .....	32
Warranty .....	33
Instrucciones en Español .....	35-67
Reference Page .....	68

---

**Toll-Free Customer Care Help Line: 1-866-326-1313**  
**Monday – Friday 8:30 a.m. – 4:30 p.m. CST**

SmartHeart™  
Manufactured for  
Veridian Healthcare  
1175 Lakeside Drive  
Gurnee, IL 60031

Made in China

#93-1081 08/21  
©2021 Veridian Healthcare, LLC

---

## CARE & SAFETY INFORMATION

---

**NOTE:** Read all instructions carefully before use. The following basic precautions are needed when using an electrical product.

**CAUTION:** Failure to read and observe all precautions could result in personal injury or equipment damage.

Improper care or use of your blood pressure monitor may result in injury, damage to the unit or ineffective treatment. Following these instructions will ensure the blood pressure monitor's efficacy and long life.

### GENERAL CAUTIONS AND WARNINGS

- Measurement results should be discussed with your physician or healthcare professional; never self-diagnose or attempt treatment as this may be dangerous.
- Intended for adult use only; this device is not approved for neonate, infant or child use. Not for use by pregnant women.
- Keep out of reach of children.
- The patient is the intended operator of this device.
- Do not use the unit if it has any damaged parts (including cuff), if it has been submersed in water or dropped.
- If any abnormality occurs, discontinue use until the unit has been examined and repaired.
- Only use included attachments and parts; do not use attachments from other brands or models as these may not be properly calibrated for use with this device and may result in measurement error.
- Do not disassemble the unit or attempt to repair it; substitution of a component different from that supplied may result in measurement error and will void manufacturer warranty.
- Always follow local regulations for proper disposal of the monitor, cuff and batteries.

---

## CARE & SAFETY INFORMATION

---

- Swallowing batteries and battery fluid is extremely dangerous; keep batteries and unit out of reach of children and disabled persons.
- Do not use this device with other diagnostic monitoring equipment on the same limb.
- Do not use near active HF Surgical Equipment and the RF shielded room of an ME system for magneticresonance imaging, where the intensity of the EM Disturbances is high.
- Do not share this device with infected persons where there may be a risk of cross-contamination.
- Consult with your physician before if the cuff may be applied over a wound or inflammation diseases.
- It is inappropriate for people with serious arrhythmia to use this Electronic Sphygmomanometer.
- Consult with your physician before use if any of the following conditions are present: arrhythmia, common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation, arterial sclerosis, poor perfusion, diabetes, pre-eclampsia, renal diseases. inflammatory diseases, intravascular therapy, arteriovenous shunts or if a mastectomy has been performed.
- Motion, trembling and shivering may affect measurement reading.
- For information regarding potential electromagnetic or other interference between digital monitor and other devices and advice on avoidance of such interference, consult with the Electromagnetic Compatibility portion of this manual.
- Monitor should be kept at least 12" from other wireless devices during use.
- If you are allergic to plastic or rubber, consult with your medical professional before using this device.

---

## **OPERATING CAUTIONS AND WARNINGS**

- Prolonged over-inflation may cause congestion, swelling or bruising in some people.
- Too-frequent measurements may cause injury due to blood flow interference.
- Check the circulation of the user before, during and after use.
- If you experience any discomfort or an abnormality occurs, stop using the device immediately; to stop the device during operation, press the **POWER** button and the cuff will automatically deflate.
- Operation outside of stated operating temperature may result in measurement error or device malfunction; operation environment temperature is: 41°F – 104°F (5°C – 40°C); humidity: ≤90% RH.
- Never use this unit while operating a vehicle or in the bathtub or shower.

## **STORAGE CAUTIONS AND WARNINGS**

- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: -4°F – 131°F (-20°C – 55°C); humidity: ≤90% RH.
- Keep the unit out of reach of children.
- Always keep the unit unplugged while not in use; remove the batteries if the unit will not be used for an extended period of time.

## **CLEANING CAUTIONS AND WARNINGS**

- Never immerse the unit in water to clean as it may damage the unit.
- Follow the 'Care & Maintenance' portion of this manual for instruction on how to clean and care for your monitor.

---

## ***INTRODUCTION & INDICATIONS FOR USE***

---

It is recommended that you first seek the advice and recommendation of your physician or healthcare professional when using home diagnostic devices, including blood pressure monitors.

Automatic digital blood pressure monitors use the oscillometric method to electronically measure your blood pressure. The monitor detects your blood's movement through the artery in your arm and converts the movements into a digital reading. The oscillometric method does not require a stethoscope, making the monitor ideal for home use.

Blood pressure readings determined with this device are equivalent to measurements obtained by a trained healthcare professional using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers.

Thank you for purchasing an Automatic Digital Blood Pressure Monitor. With proper care and use, your monitor will provide you with many years of reliable readings.

**PLEASE READ THIS INSTRUCTION MANUAL  
COMPLETELY BEFORE OPERATING THIS UNIT.**


# BLOOD PRESSURE MONITOR FEATURES

**DETAILS OF THE DISPLAY AND INTERPRETATION OF THE SYMBOLS ARE INCLUDED THROUGHOUT THIS MANUAL.**

# BLOOD PRESSURE MONITOR FEATURES

**DETAILS OF THE DIGITAL DISPLAY SCREEN AND INTERPRETATION OF THE SYMBOLS ARE INCLUDED THROUGHOUT THIS MANUAL.**

---

## **GETTING STARTED**

---

### **Before taking a measurement**

1. Your readings should only be interpreted by your physician or healthcare professional with access to your individual medical history. Regular use of a home blood pressure monitor will allow you to track and record your readings for discussion with your physician.
2. Conduct your measurement in a quiet place while seated in a relaxed position. Rest for 15 minutes before taking your reading. If an error occurs or you desire to take a second reading, allow 15 minutes between readings for your blood vessels to return to normal.
3. Avoid smoking, eating, taking medication, alcohol, physical activity or any other stressful activity for 30 minutes prior to taking a reading.
4. Always remove any jewelry or constrictive clothing that may interfere with the cuff placement.
5. Keep yourself and the monitor still during measuring; do not talk during the reading.
6. It is recommended that you take your readings at the same time each day to better monitor any indications in your results.
7. Record your daily measurement on the included chart or some other written document to share with your physician.


---

## BATTERY INSTALLATION

---

### This unit comes complete with 2 AAA batteries.


It is necessary to replace the batteries when the 'Low Battery' symbol appears on the display or when the display does not turn on after the **START/POWER** button is pressed.

The battery compartment is located on the underside of this monitor.

1. Press down and lift the cover.
2. Insert or replace 2 AAA batteries into the battery compartment, ensuring to match the indicated polarity symbols. Always use new batteries.
3. Replace the battery cover.
4. Dispose of batteries according to local disposal and recycling regulations.

It is recommended to remove the batteries if the unit will not be used for an extended period of time.

## SETTING DATE/TIME

IT IS NECESSARY TO SET THE DATE AND TIME FOR THE UNIT EVERY TIME BATTERIES ARE INITIALLY INSTALLED OR REPLACED.

DATE AND TIME WILL DEFAULT TO FACTORY SETTINGS, 12:00/1(M)-1(D), IF YOU DO NOT CHOOSE TO SET THE CORRECT DATA.

The **MEM** and **START/POWER** buttons are located on the right side of the monitor.

The date/time feature is on the bottom left of the display screen, Fig 1. The date and time will alternate every four seconds regardless of the operating mode.

1. When the unit is off, press and hold the **START** and **MEM** button at the same time. The unit will beep once and the month will begin to flash, Fig. 2.
2. Press the **MEM** button to change until the desired number appears. To scroll through the numbers quicker, hold down the **MEM** button as the numbers advance.


Fig. 1

Fig. 2

---

## **SETTING DATE/TIME**

---

3. Press the **START** button to confirm the month entry.
4. Next, the date will blink. Repeat steps 2 and 3 to set the date, then hours, then minutes. This unit does not track the year.
5. After setting the minutes, the unit will automatically exit out of the date/time setting mode.

When the monitor is off, the date and time will still appear on the display.

Time is displayed in 24hr Military Time. This chart explains the time conversion.

### **MILITARY TIME = REGULAR TIME**

1:00 = 1AM	7:00 = 7AM	13:00 = 1PM	19:00 = 7PM
2:00 = 2AM	8:00 = 8AM	14:00 = 2PM	20:00 = 8PM
3:00 = 3AM	9:00 = 9AM	15:00 = 3PM	21:00 = 9PM
4:00 = 4AM	10:00 = 10AM	16:00 = 4PM	22:00 = 10PM
5:00 = 5AM	11:00 = 11AM	17:00 = 5PM	23:00 = 11PM
6:00 = 6AM	12:00 = 12PM	18:00 = 6PM	24:00 = 12AM

---

## **FITTING & APPLYING YOUR CUFF**

---

**PROPER CUFF FIT AND PLACEMENT IS CRITICAL IN OBTAINING ACCURATE BLOOD PRESSURE MEASUREMENTS. PLEASE READ THESE INSTRUCTIONS CAREFULLY AND CONSIDER CONSULTING WITH YOUR PHYSICIAN FOR FURTHER DIRECTION.**

Before applying your blood pressure cuff, be sure your wrist is within the appropriate cuff range:


Wrist Circumference Range Suitable for 5.5" – 7.6" (14 cm – 19.5 cm).

Reminder, it is important to avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading.

If for any reason you are unable to or cannot not use your left wrist, please modify the instructions for cuff application to your right wrist. Your physician can help you identify which wrist is best for you to take measurements from.

### **Applying Your Cuff**

1. Remove any constrictive clothing or jewelry that may interfere with cuff placement.
2. Be seated at a table or desk with your feet flat on the floor.
3. Hold your left arm in front of you with your palm facing upward.

---

## **FITTING & APPLYING YOUR CUFF**

---

4. Apply the cuff to your left wrist so that the digital display face is positioned on the inside area of your wrist facing you, Fig. 1.
5. Adjust the cuff approximately 0.5" from the edge of the head of the ulna bone, Fig. 2.
6. The cuff should fit comfortably, yet snugly around your wrist.

If you are not comfortable with applying your cuff, please seek the assistance of another member of your household or work with your physician to practice the cuff application. Incorrectly applied cuffs may result in inaccurate readings.


Fig. 1


Fig. 2

## TAKING YOUR BLOOD PRESSURE READING

Please read the preceding portions of this manual prior to taking your first reading.

Reminder: It is important to avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading.

1. Rest your elbow on a solid surface with your palm facing upward. Elevate your arm so that the cuff is at the same level as your heart, Fig 1. Be certain that cuff is elevated to the proper level; you may choose to support your hand with a folded towel or book. Relax your left hand. Be seated with your feet flat on the floor and legs uncrossed.
2. Press and release the **START/POWER** button to turn on device. The unit will run a self-test, Fig. 2
3. The display will briefly show a flashing arrow and '0' before the measurement begins, Fig 3.
4. The cuff will automatically inflate to approximately 190mmhg or 40 mmhg higher than the last systolic pressure reading stored in memory.
5. When the inflation has reached optimum level, the display will begin to show the decreasing pressure while you feel the pressure of the cuff decrease.


Fig. 1


Fig. 2

Fig. 3


**YOU CAN STOP THE INFLATION OR DEFLATION PROCESS AT ANY TIME BY PRESSING THE POWER BUTTON.**

## TAKING YOUR BLOOD PRESSURE READING

- When the measurement is complete, your blood pressure measurement and pulse will display separately on the screen alternating repeatedly every 2-3 seconds, Fig 4. The systolic and diastolic readings will display at the same time.
- The hypertension Indicator will indicate your reading range in the right side of the display, Fig 5. For more information see Interpreting Your Results.
- If an irregular heartbeat was detected during the reading, the Irregular heartbeat Detector icon will appear flashing on the display for both the measurement results and pulse results screens, Fig 6. For more information see Interpreting Your Results.
- The reading will automatically be stored in memory, up to 60 readings.
- Select **START** to turn the unit off and conserve energy and battery life. The unit will automatically shut-off and revert to time/date mode after approximately 1 minute.

Systolic/Diastolic Results

Pulse Results

Fig. 4

Fig. 5


Fig. 6

---



## INTERPRETING YOUR RESULTS

---

### Irregular Heartbeat Detector (♥)

An irregular heartbeat is defined as a heartbeat rhythm that is more than 25% slower or 25% faster than the average measurement detected while the monitor is conducting the measurement. When this device detects an irregular rhythm two or more times during the reading, the IHB symbol will appear on the display screen.

This monitor will operate if the icon appears, but the results may be compromised, especially if this symbol appears often. Please consult with your physician or trained healthcare professional for further information regarding an irregular heartbeat and if this symbol appears frequently.

# INTERPRETING YOUR RESULTS

## Hypertension Indicator

This unit features our unique Hypertension Indicator. The color bars on the right side of the monitor display correspond with an icon on the digital display, indicating where the measurement results fall within the World Health Organization standards.

The World Health Organization has established globally accepted standards for the assessment of high or low blood pressure readings. The below chart should be considered only as a guideline, always consult with your physician or health care professional to interpret your individual results. Never adjust your medication(s) dosage, or severely alter your diet or exercise routine without consulting with your physician.

<p><b>Stage 3 Hypertension</b></p> <p>Systolic pressure: <math>\geq 180</math> mmHg or Diastolic pressure: <math>\geq 110</math> mmHg</p>	<p><b>High Normal</b></p> <p>Systolic pressure: <math>\geq 130 \leq 139</math> mmHg or Diastolic pressure: <math>\geq 85 \leq 89</math> mmHg</p>
<p><b>Stage 2 Hypertension</b></p> <p>Systolic pressure: <math>\geq 160 \leq 179</math> mmHg or Diastolic pressure: <math>\geq 100 \leq 109</math> mmHg</p>	<p><b>Normal</b></p> <p>Systolic pressure: <math>\geq 120 \leq 129</math> mmHg or Diastolic pressure: <math>\geq 80 \leq 84</math> mmHg</p>
<p><b>Stage 1 Hypertension</b></p> <p>Systolic pressure: <math>\geq 140 \leq 159</math> mmHg or Diastolic pressure: <math>\geq 90 \leq 99</math> mmHg</p>	<p><b>Optimal</b></p> <p>Systolic pressure: <math>&lt; 120</math> mmHg or Diastolic pressure: <math>&lt; 80</math> mmHg</p>

# MEMORY FUNCTIONS

## Recalling Measurements in Memory

You can recall up to 60 total measurements in memory. These measurements can be shared with your physician or trained healthcare professional.

1. Press and release the **MEM** button. The unit will show how many measurements are stored within the bank, Fig 1.
2. Press the **MEM** button again to scroll through the saved measurements. Measurements will appear on the display from most current to oldest. All results for a given measurement will display, including measurement results, pulse rate, hypertension Indicator, Irregular heartbeat alert (if applicable), and date/time stamp, Fig 2. The measurement results and pulse results will alternate, flashing between measurement and pulse every 2-3 seconds.
3. When the number of readings exceeds 60, the oldest data will be replaced with the new record.
4. Press **START** to turn the monitor off at any time during review of the stored measurements..


Fig. 1


Fig. 2

---

# MEMORY FUNCTIONS

---

## Clearing Measurements From Memory

To delete results in memory bank:

1. Press and release the **MEM** button. The unit will display how many measurements are stored within the bank.
2. Press and hold down the **MEM** button until four horizontal lines flash, Fig. 3. This indicates that all measurements have been erased.
3. The memory bank is now empty.


Fig. 3

---

## CARE & MAINTENANCE

---

Proper care and maintenance of your blood pressure monitor will help ensure long life for the product; improper use of the monitor may void the included warranty.

### Monitor Care & Maintenance

- Do not use any liquids on the monitor; use a soft, dry cloth to clean the monitor as needed.
- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is:  $-4^{\circ}\text{F} - 131^{\circ}\text{F}$  ( $-20^{\circ}\text{C} - 55^{\circ}\text{C}$ ); humidity:  $\leq 90\%$  RH.
- This monitor requires 6 hours to cool from the maximum storage temperature between uses until the monitor is ready for its INTENDED USE when the ambient temperature is  $68^{\circ}\text{F}$  ( $20^{\circ}\text{C}$ ).
- This monitor requires 6 hours to warm from the minimum storage temperature between uses until the monitor is ready for its INTENDED USE when the ambient temperature is  $68^{\circ}\text{F}$  ( $20^{\circ}\text{C}$ ).
- Avoid high temperatures and direct sunlight.
- Always remove the batteries if the unit will not be used for an extended period of time.
- Do not attempt to disassemble your monitor; disassembling your monitor may void the manufacturers warranty.
- Do not attempt to service the monitor while in use or on a patient.
- Do not subject the monitor to strong shocks; take care not to drop the monitor.
- No component can be maintained by user in the monitor. The circuit diagrams, component part lists, descriptions, calibrations instructions, or other information which will assist the user's appropriately qualified technical personnel to repair those parts of equipment which are designated repairable can be supplied.

---

## **CARE & MAINTENANCE**

---

### **Cuff Care & Maintenance**

- The cuff may be spot cleaned with a mild-detergent as needed.
- Never attempt to iron or dry the cuff; if the cuff becomes wet, allow the cuff to air dry completely before use or storage.
- To disinfect the cuff, wipe the inside (skin contact side) of the cuff using a cloth lightly moistened with 70% - 90% ethyl rubbing alcohol and allow to air dry; it is suggested to disinfect the cuff frequently when used in clinical settings or if the cuff is contaminated in any way.
- It is recommended the performance should be checked every 2 years or after a repair. Please contact Customer Care at 866-326-1313.
- It is recommended to clean the cuff every 200 uses.
- The monitor can maintain the safety and performance characteristics for a minimum of 10,000 measurements of three years.
- The cuff integrity is maintained after 1,000 open- close cycles of the closure.

---

## DEVICE & LABEL SYMBOLS

---

These symbols may appear on your device, instructions or packaging and may vary by make and model.

### Symbol | Meaning

	Read This Manual—All included manuals should be read prior to first use. (Background color will be blue; reading symbol will be white)
	Warning—Symbol indicates a warning, prohibition or mandatory action that mitigates a risk that is not necessarily obvious to the device operator.
	Type BF Applied Parts—Indicates that a part of this unit comes in contact with the patient in order to carry out it's intended function; in the case of this device the cuff is the Type BF Applied Part.
	Environment Protection—dispose of this product properly; consult with your local recycling ordinances for proper recycling and disposal.
<b>IP22</b>	IP22—The first characteristic numeral symbol for “Degrees of protection against access to hazardous parts and against solid foreign objects. The second characteristic numeral symbol for “Degrees of protection against ingress of water”.
<b>SN</b>	Serial Number.
	Direct Current.
	Date of Manufacture.

The Electronic Sphygmomanometers corresponds to the below standards: ISO81060-2 : 2013(Non-Invasive Sphygmomanometers- Part 2: Clinical Validation Of Automated Measurement Type).

## DISPLAY SYMBOLS

Symbol	Meaning
<b>SYS</b>	Systolic Blood Pressure—The 'top' blood pressure result
<b>DIA</b>	Diastolic Blood Pressure—The 'bottom' blood pressure result
<b>PUL/min</b>	Pulse Rate—Pulse rate per minute during the measurement
<b>mmHg</b>	Millimeter of Mercury—Unit of blood pressure measurement for this device
	Low Battery—Battery life indicator
	Ready to Inflate Indicator—The cuff is not properly wrapped or has too much air left in it. The arrow will stop flashing when the unit is ready to use
	Irregular Heartbeat Symbol—An irregular heartbeat was detected during the measurement
	Memory—The display shows the measurements stored in memory

# TROUBLESHOOTING

Problem	Possible Cause	Solution
Display shows an abnormal result.	Cuff position was not correct or was not tight enough.	Apply the cuff properly and attempt a new measurement.
	Body posture was not correct during measurement.	Review the 'Taking a Measurement' portion of this manual and attempt a new measurement.
	Speaking, excessive movement, angry, excited or nervous anxiety during measurement.	Wait a period of time and attempt a new measurement after calm has returned; do not speak or move during measurement.
Display is blank when power is on.	Batteries may be expired or installed improperly.	Check the polarity of the batteries and reinstall if necessary; replace batteries with new batteries.
No response when <b>START</b> is pressed or when new batteries have been loaded.	Incorrect operation; possible strong electromagnetic interference.	Remove batteries and wait a minimum of five minutes; reinstall batteries and attempt a new measurement.
Monitor keeps re-inflating.	Circuit is locked.	Remove and reinsert the batteries; attempt new measurement.

## ERROR CODES

Code	Meaning	Corrective Action
E0	Pressure system is unstable before measurement.	Review the 'Taking a Measurement' portion of the manual, wait a minimum of 5 minutes and attempt a new measurement.
E1 / E2	Failed to detect systolic/ diastolic pressure.	
E3/ E4	Pneumatic system leakage or cuff is too loose/ too tight during inflation.	Review the Applying Cuff portion of this manual; re-apply the cuff, wait a minimum of 5 minutes and attempt a new measurement.
E5	The cuff pressure is above 300mmHg.	Check that the cuff is applied properly, wait a minimum of 5 minutes and attempt a new measurement. If the monitor still does not perform or performs abnormally, contact Customer Care at 866-326-1313.
E6	More Than 3 minutes with cuff pressure above 15 mmHg.	
E7	EEPROM accessing error.	
E8	Pressure Sensor parameter error.	
ER A	Device Parameter checking error.	

## TECHNICAL ALARM DESCRIPTION

The monitor will show 'HI' or 'Lo' as technical alarm on LCD with no delay if the determined blood pressure (systolic or diastolic) is outside the rated range specified in part SPECIFICATIONS. In this case, you should consult a physician or check if your operation violated the instructions.

The technical alarm condition (outside the rated range) is preset in the factory and cannot be adjusted or inactivated. This alarm condition is assigned as low priority according to IEC 60601-1-8.

The technical alarm is non-latching and need no reset. The signal displayed on LCD will disappear automatically after about 8 seconds.

---

## **FCC STATEMENT**

---

### **NOTE:**

#### POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for U.S.A. only)

This product has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. The product generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If the product does cause harmful interference to radio or television reception, which can be determined by turning the product on and off, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the product and the receiver.
- Connect the product into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

#### POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for Canada only)

This digital apparatus does not exceed the Class B limits for radio noise emissions from digital apparatus as set out in the interference-causing equipment standard entitled "Digital Apparatus", ICES-003 of the Canadian Department of Communications.

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Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

---

# **ELECTROMAGNETIC COMPATIBILITY**

---

Table 1 EMISSIONS

## **For all ME EQUIPMENT and ME SYSTEMS**

### **Guidance and manufacturer's declaration - electromagnetic emissions**

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

<b>Phenomenon</b>	<b>Compliance</b>	<b>Electromagnetic Environment Guidance</b>
RF emissions	CISPR 11 Group 1, Class B	Home healthcare environment
Harmonic distortion	IEC 61000-3-2 Class A	Home healthcare environment
Voltage fluctuations and flicker	IEC 61000-3-3 Compliance	Home healthcare environment

# ELECTROMAGNETIC COMPATIBILITY

Table 2 — ENCLOSURE PORT  
For all ME EQUIPMENT and ME SYSTEMS

## Guidance and manufacturer's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Phenomenon	Basic EMC Standard	Immunity test levels
		Home healthcare environment
Electrostatic Discharge	IEC 61000-4-2	$\pm 8$ kV contact $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV air
Radiated RF EM field	IEC 61000-4-3	10V/m 80MHz-2.7GHz 80% AM at 1kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	Refer to table 3
Rated power frequency magnetic fields	IEC 61000-4-8	30A/m 50Hz or 60Hz

# ELECTROMAGNETIC COMPATIBILITY

Table 3 — PROXIMITY FIELDS FROM RF WIRELESS COMMUNICATIONS EQUIPMENT

## For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING Guidance and manufacturer's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Test frequency (MHz)	Band (MHz)	Immunity test levels
		Professional healthcare facility environment
385	380-390	Pulse modulation 18Hz, 27V/m
450	430-470	FM, $\pm 5$ kHz deviation, 1kHz sine, 28V/m
710	704-787	Pulse modulation 217Hz, 9V/m
745		
780		
810		
870	800-960	Pulse modulation 18Hz, 28V/m
930		
1720	1700-1990	Pulse modulation 217Hz, 28V/m
1845		
1970		
2450	2400-2570	Pulse modulation 217Hz, 28V/m
5240	5100-5800	Pulse modulation 217Hz, 9V/m
5500		
5785		

---

## PRODUCT SPECIFICATIONS

---

<b>Name</b>	SmartHeart Automatic Blood Pressure Wrist Monitor
<b>Model Number</b>	01-540
<b>Display System</b>	Digital display/LCD
<b>Measuring Method</b>	Oscillometric Method, Automatic Air Inflation and Measurement
<b>Power Source</b>	2 x 1.5V  AAA Batteries
<b>Measuring Range</b>	Pressure: 0-300 mmHg Diastolic: 40-199 mmHg Systolic: 60-260 mmHg Pulse: 40-180 beats/minute
<b>Accuracy</b>	Pressure: $\pm 3$ mmHg; Pulse: $\pm 5\%$
<b>Pressurization</b>	Automatic pressurization by pump
<b>Deflation</b>	Automatic pressure release valve
<b>Memory</b>	60 measurements
<b>Automatic Power-Off</b>	Approximately 1 minute after last button operation
<b>Battery Life</b>	Approximately 270 measurements
<b>Operation Environment</b>	Temperature 41°F – 104°F (5°C – 40°C); Humidity $\leq 90\%$ RH (non-condensing)
<b>Storage &amp; Transport Environment</b>	Temperature -4°F – 131°F (-20°C – 55°C); Humidity $\leq 90\%$ RH (non-condensing)
<b>Monitor Dimensions</b>	3-1/3" x 2-1/2" x 9/10"
<b>Monitor Weight</b>	3.2 oz. (without batteries)
<b>Wrist Circumference Range</b>	Suitable for 5.5" – 7.6" (14 – 19.5 cm)
<b>Accessories</b>	Instruction manual, two AAA batteries, storage case Quick Start Guide, Blood Pressure Log and FAQ

*Specifications are subject to change without notice.*

---

# 1-YEAR LIMITED WARRANTY

---

Congratulations on your purchase of a Digital Blood Pressure Monitor. Your Digital Blood Pressure Monitor is covered by the following limited warranty commencing upon the date of purchase, and subject to the following terms and conditions:

The Warrantor warrants that its Digital Blood Pressure Monitor will be free from defects in materials and workmanship under normal consumer usage for a period of one year for the original purchaser.

Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the Digital Blood Pressure Monitor for commercial purposes or subjecting the Digital Blood Pressure Monitor to abnormal usage or conditions; or (d) other acts which are not the fault of the Warrantor, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the Digital Blood Pressure Monitor.

If the Digital Blood Pressure Monitor fails to conform to this limited warranty, return the Digital Blood Pressure Monitor postage prepaid to: **Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031**. When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone number. Carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. At our option, the Warrantor will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will be responsible for parts charges, if any, and repair charges not covered by this limited warranty.

EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL THE WARRANTOR BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

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