SMARTHEART™

Automatic Blood Pressure Monitor

Model # 01-539


INSTRUCTION MANUAL **ENGLISH AND SPANISH**

Please read this instruction manual completely before operating this unit.


STOP!
**PLEASE ENSURE YOU HAVE ALL OF THE
FOLLOWING COMPONENTS BEFORE USING
YOUR DIGITAL BLOOD PRESSURE MONITOR**


Digital Monitor


Adult Cuff

4 AA
Batteries

Storage Case

Instruction Manual
Quick Start Guide
Blood Pressure Reading Tracking Log
FAQ Information Sheet
Replacement Parts Order Form

**IF YOU ARE MISSING ANY PARTS, INCLUDING INSERTS OR
INSTRUCTION MANUALS, DO NOT RETURN TO PLACE OF
PURCHASE. CONTACT CUSTOMER CARE AT 866-326-1313.**

INDEX

Care & Safety Information	4-6
Introduction & Indications for Use	7
Blood Pressure Monitor Features	8-9
Getting Started	10
Battery Installation	11
Setting the Date/Time	12-13
Fitting & Applying Your Cuff	14-15
Taking Your Blood Pressure Reading	16-17
Interpreting Your Results	18-19
Memory Function	20-21
Care & Maintenance	22-23
Display & Label Symbols	24
Display Symbols	25
Troubleshooting	26
Error Codes	27
Technical Alarm Description	27
FCC Statement	28
Electromagnetic Compatibility	29-31
Product Specifications	32
Warranty	33
Instrucciones en Español	35-67
Reference Page	68

Toll-Free Customer Care Help Line: 1-866-326-1313
Monday – Friday 8:30 a.m. – 4:30 p.m. CST

SmartHeart™
Manufactured for
Veridian Healthcare
1175 Lakeside Drive
Gurnee, IL 60031
www.veridianhealthcare.com

Made in China

#93-1096 02/22
©2022 Veridian Healthcare, LLC

CARE & SAFETY INFORMATION

NOTE: Read all instructions carefully before use. The following basic precautions are needed when using an electrical product.

CAUTION: Failure to read and observe all precautions could result in personal injury or equipment damage.

Improper care or use of your blood pressure monitor may result in injury, damage to the unit or ineffective treatment. Following these instructions will ensure the blood pressure monitor's efficacy and long life.

The blood pressure monitor is not intended to be exposed to the Electromagnetic Interference (EMI) environment, please do not use the blood pressure monitor within the environment of the following device: Magnetic Resonance Imaging (MRI), computerized axial tomography (CT), diathermy, Radio Frequency Identification (RFID), and electromagnetic security systems such as metal detectors.

GENERAL CAUTIONS AND WARNINGS

- Measurement results should be discussed with your physician or healthcare professional; never self-diagnose or attempt treatment as this may be dangerous.
- Intended for adult use only; this device is not approved for neonate, infant or child use. Not for use by pregnant women.
- Keep out of reach of children.
- The patient is the intended operator of this device.
- Do not use the unit if it has any damaged parts (including power cord or plug and cuff), if it has been submerged in water or dropped.
- If any abnormality occurs, discontinue use until the unit has been examined and repaired.
- Only use included attachments and parts; do not use attachments from other brands or models as these may not be properly calibrated for use with this device and may result in measurement error.
- Do not disassemble the unit or attempt to repair it; substitution of a component

CARE & SAFETY INFORMATION

different from that supplied may result in measurement error and will void manufacturer warranty.

- Always follow local regulations for proper disposal of the monitor, cuff and batteries.
- Swallowing batteries and battery fluid is extremely dangerous; keep batteries and unit out of reach of children and disabled persons.
- Do not use this device with other diagnostic monitoring equipment on the same limb.
- Do not use near active HF Surgical Equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of the EM Disturbances is high.
- Do not share this device with infected persons where there may be a risk of cross-contamination.
- Consult with your physician before if the cuff may be applied over a wound or inflammation diseases.
- It is inappropriate for people with serious arrhythmia to use this Electronic Sphygmomanometer.
- Consult with your physician before use if any of the following conditions are present: arrhythmia, common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation, arterial sclerosis, poor perfusion, diabetes, pre-eclampsia, renal diseases, inflammatory diseases, intravascular therapy, arteriovenous shunts or if a mastectomy has been performed.
- Motion, trembling and shivering may affect measurement reading.
- For information regarding potential electromagnetic or other interference between digital monitor and other devices and advice on avoidance of such interference, consult with the Electromagnetic Compatibility portion of this manual.
- Monitor should be kept at least 12" from other wireless devices during use.
- If you are allergic to plastic or rubber, consult with your medical professional before using this device.

CARE & SAFETY INFORMATION

OPERATING CAUTIONS AND WARNINGS

- Prolonged over-inflation may cause congestion, swelling or bruising in some people.
- Too-frequent measurements may cause injury due to blood flow interference.
- Check the circulation of the user before, during and after use.
- If you experience any discomfort or an abnormality occurs, stop using the device immediately; to stop the device during operation, press the **POWER** button and the cuff will automatically deflate.
- Operation outside of stated operating temperature may result in measurement error or device malfunction; operation environment temperature is: 50°F – 104°F (10°C – 40°C); Humidity: ≤90%RH max.
- Never use this unit while operating a vehicle or in the bathtub or shower.

STORAGE CAUTIONS AND WARNINGS

- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: -4°F – 131°F (-20°C – 55°C); Humidity: ≤90%RH max.
- Keep the unit out of reach of children.
- Always keep the unit unplugged while not in use; remove the batteries if the unit will not be used for an extended period of time.

CLEANING CAUTIONS AND WARNINGS

- Never immerse the unit in water to clean as it may damage the unit.
- Follow the 'Care & Maintenance' portion of this manual for instruction on how to clean and care for your monitor.

INTRODUCTION & INDICATIONS FOR USE

It is recommended that you first seek the advice and recommendation of your physician or healthcare professional when using home diagnostic devices, including blood pressure monitors.


Automatic digital blood pressure monitors use the oscillometric method to electronically measure your blood pressure. The monitor detects your blood's movement through the artery in your arm and converts the movements into a digital reading. The oscillometric method does not require a stethoscope, making the monitor ideal for home use.

Blood pressure readings determined with this device are equivalent to measurements obtained by a trained healthcare professional using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers.

Thank you for purchasing an Automatic Blood Pressure Monitor. With proper care and use, your monitor will provide you with many years of reliable readings.


**PLEASE READ THIS INSTRUCTION MANUAL
COMPLETELY BEFORE OPERATING THIS UNIT.**

BLOOD PRESSURE MONITOR FEATURES

DETAILS OF THE DIGITAL DISPLAY SCREEN AND INTERPRETATION OF THE SYMBOLS ARE INCLUDED THROUGHOUT THIS MANUAL.

BLOOD PRESSURE MONITOR FEATURES

GETTING STARTED

Before Taking a Measurement


1. Your readings should only be interpreted by your physician or healthcare professional with access to your individual medical history. Regular use of a home blood pressure monitor will allow you to track and record your readings for discussion with your physician.
2. Conduct your measurement in a quiet place while seated in a relaxed position. Rest for 15 minutes before taking your reading. If an error occurs or you desire to take a second reading, allow 15 minutes between readings for your blood vessels to return to normal.
3. Avoid smoking, eating, taking medication, alcohol, physical activity or any other stressful activity for 30 minutes prior to taking a reading.
4. Always remove any jewelry or constrictive clothing that may interfere with the cuff placement.
5. Keep yourself and the monitor still during measuring; do not talk during the reading.
6. It is recommended that you take your readings at the same time each day to better monitor any indications in your results.
7. Record your daily measurement on the included chart or some other written document to share with your physician.

BATTERY INSTALLATION

This unit requires 4 AA batteries.


It is necessary to replace the batteries when the “Low Battery”  symbol appears on the display or when the display does not turn on after the **START** button is pressed.

The battery compartment is located on the underside of this monitor.

1. Press down and slide off battery cover.
2. Insert or replace 4 AA batteries into the battery compartment, ensuring to match the indicated polarity symbols. Always use new batteries.
3. Replace the battery cover.
4. Dispose of batteries according to local disposal and recycling regulations.

It is recommended to remove the batteries if the unit will not be used for an extended period of time.

SETTING DATE/TIME

IT IS NECESSARY TO SET THE DATE AND TIME FOR THE UNIT EVERY TIME BATTERIES ARE INITIALLY INSTALLED OR REPLACED.

Date and time will default to factory settings, 12:00 a.m. / 1(M)-1(D), if you do not choose to set the correct data.

The date/time feature is on the bottom left of the display screen. This unit stores the date and time, allowing for a date/time stamp on all readings in the memory recall function.

1. When the batteries are installed, the time/date will automatically show on the lower left of the display as 12:00/1(M) -1(D). After set-up time and date will continue to alternate even when the monitor is off and not in use.
2. When the unit is off, press and hold the **START** and **MEM** button at the same time. The month will begin to flash, Fig. 1.
3. Press the **MEM** button to advance the month until the desired number appears. To scroll through the numbers quicker, hold down the **MEM** button as the numbers advance.
4. Press the **START** button to confirm the month entry.
5. Next, the date will blink. Repeat steps 2 and 3 to set the date, then hours, then minutes, Fig 2. This unit does not track the year.
6. After setting the minutes, the unit will automatically exit out of the date/time setting mode and is ready to take a measurement.


Fig. 1




Fig. 2

SETTING DATE/TIME

WHEN THE MONITOR IS OFF, THE DATE/TIME ALTERNATING SCREEN WILL STILL APPEAR ON THE DISPLAY.

Time is displayed in 24hr Military Time. This chart explains the time conversion.

MILITARY TIME = REGULAR TIME


1:00 = 1AM	7:00 = 7AM	13:00 = 1PM	19:00 = 7PM
2:00 = 2AM	8:00 = 8AM	14:00 = 2PM	20:00 = 8PM
3:00 = 3AM	9:00 = 9AM	15:00 = 3PM	21:00 = 9PM
4:00 = 4AM	10:00 = 10AM	16:00 = 4PM	22:00 = 10PM
5:00 = 5AM	11:00 = 11AM	17:00 = 5PM	23:00 = 11PM
6:00 = 6AM	12:00 = 12PM	18:00 = 6PM	24:00 = 12AM

FITTING & APPLYING YOUR CUFF

PROPER CUFF FIT AND PLACEMENT IS CRITICAL IN OBTAINING ACCURATE BLOOD PRESSURE MEASUREMENTS. PLEASE READ THESE INSTRUCTIONS CAREFULLY AND CONSIDER CONSULTING WITH YOUR PHYSICIAN FOR FURTHER DIRECTION.

Before applying your blood pressure cuff, be sure the cuff included with this monitor is the appropriate size.

To determine the proper cuff size, measure the circumference of your left upper arm at the point midway between inside bend of your elbow and your shoulder.

Adult Cuff fits arm circumference 11-3/4" – 16-1/2" (30 – 42 cm)

If for any reason you are unable to or should not use your left arm, please modify the instructions for cuff application to your right arm. Your physician can help you identify which arm is best for you to take measurements from.

Assembling Your Cuff

Your cuff may need to be 'threaded' through the metal bar to form the cuff cylinder for your arm.

1. Lay the open cuff hook and loop side down, size/instruction wording side up on a flat surface.
2. Keeping the words on the inside of the cuff, thread the end of the cuff through the metal bar, Fig 1. The end of the cuff should be able to fold over onto itself for securing via the hook and loop closure, Fig 2.


Fig. 1


Fig. 2

FITTING & APPLYING YOUR CUFF

Applying Your Cuff

Be seated at a table or desk with your feet flat on the floor. To avoid accidentally pulling on or dropping your monitor, the cuff should not be plugged into the monitor until after the cuff is applied to your arm.

The cuff is printed with a useful guide for applying the cuff; please consult the following steps for additional direction.

1. Remove any constrictive clothing or jewelry that may interfere with cuff placement.
2. Slide the cuff cylinder onto your left arm and align the artery mark over the brachial artery on the inside of your arm, Fig 3.

The tube should follow the inside length of the arm, Fig 4.

3. The bottom edge of the cuff should be positioned approximately 0.5" above the elbow joint, Fig 5.
4. When the cuff is closed around the arm, the triangle guide imprint should be within the OK range bar marking; if the triangle is outside of the range, a different sized cuff is needed.
5. The cuff should fit comfortably, yet snugly around your arm, Fig 6. You should be able to insert one finger easily between your arm and the cuff.
6. During use, avoid compression or squeezing of the air hose/tube, as this may cause inflation error or injury due to excessive cuff pressure.


Fig. 3


Fig. 4


Fig. 5


Fig. 6

Fig. 6

TAKING YOUR BLOOD PRESSURE READING

Please read the preceding portions of this manual prior to taking your first reading.

Reminder, it is important to avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading.


Fig. 1

1. Position the monitor on a flat, stable surface with the digital display in view.
2. With the cuff properly positioned on your arm, insert the cuff tubing connector into the port on the left side of your monitor.
3. Rest your elbow on a solid surface with your palm facing upward. Elevate your arm so that the cuff is at the same level as your heart, Fig 1. Relax your left hand.
4. Press the **START** button to turn the unit on; the unit will run a self-test, Fig 2.
5. The display will show the current memory bank in use, U1 or U2, Fig 3. Press the **MEM** button to select a different memory bank if desired. If no selection is made, the unit will default to the memory bank displayed.
6. The display will briefly show a flashing arrow and '0' before the measurement begins, Fig 4.
7. The cuff will automatically inflate to approximately 190mmHg or 40 mmHg higher than the last systolic pressure reading stored in memory.


Fig. 2


Fig. 3

Fig. 4

TAKING YOUR BLOOD PRESSURE READING

YOU CAN STOP THE INFLATION OR DEFLATION PROCESS AT ANY TIME BY PRESSING THE START BUTTON.

- When the inflation has reached optimum level, the display will begin to show the decreasing pressure while you feel the pressure of the cuff decrease.
- When the measurement is complete, your blood pressure measurement and pulse will display separately on the screen alternating repeatedly every 2-3 seconds, Fig 5. The systolic and diastolic readings will display at the same time.
- The Hypertension Indicator will indicate your reading range on the left side of the display, Fig 6. For more information see Interpreting Your Results.
- If an irregular heartbeat was detected during the reading, the Irregular Heartbeat Detector icon will appear flashing on the display for both the measurement results and pulse results screens, Fig 7. For more information see Interpreting Your Results.
- The reading will automatically be stored in the memory bank selected prior to the measurement, up to 60 readings.
- Select **START** to turn the unit off and conserve energy and battery life. The unit will automatically shut-off and revert to time/date mode after approximately 1 minute.
- Disconnect the tubing prior to storage.

Systolic/Diastolic Results

Pulse Results

Fig. 5

Fig. 6




Fig. 7

INTERPRETING YOUR RESULTS

Irregular Heartbeat (♥),

An irregular heartbeat is defined as a heartbeat rhythm that is more than 25% slower or 25% faster than the average measurement detected while the monitor is conducting the measurement. When this device detects an irregular rhythm two or more times during the reading, the IHB symbol will appear on the display.


This monitor will operate if the icon appears, but the results may be compromised, especially if this symbol appears often. Please consult with your physician or trained healthcare professional for further information regarding an irregular heartbeat and if this symbol appears frequently.

INTERPRETING YOUR RESULTS

Hypertension Indicator - WHO Classifications

This unit features our unique Hypertension Indicator. The color bars on the left side of the monitor display correspond with an icon on the digital display, indicating where the measurement results fall within the World Health Organization standards.

The World Health Organization has established globally accepted standards for the assessment of high or low blood pressure readings. The below chart should be considered only as a guideline, always consult with your physician or health care professional to interpret your individual results. Never adjust your medication(s) dosage, or severely alter your diet or exercise routine without consulting with your physician.

The Hypertension Indicator is located on the left side of the display window

<p>Stage 3 Hypertension Systolic pressure: ≥ 180 mmHg or Diastolic pressure: ≥ 110 mmHg</p>	<p>High Normal Systolic pressure: ≥ 130-139 mmHg or Diastolic pressure: ≥ 85-89 mmHg</p>
<p>Stage 2 Hypertension Systolic pressure: ≥ 160-179 mmHg or Diastolic pressure: ≥ 100-109 mmHg</p>	<p>Normal Systolic pressure: ≥ 120-129 mmHg or Diastolic pressure: ≥ 80-84 mmHg</p>
<p>Stage 1 Hypertension Systolic pressure: ≥ 140-159 mmHg or Diastolic pressure: ≥ 90-99 mmHg</p>	<p>Optimal Systolic pressure: < 120 mmHg or Diastolic pressure: < 80 mmHg</p>

MEMORY FUNCTION

Recalling Measurements in Memory

This monitor stores and recalls up to 60 measurements per memory bank, 120 total, plus an average of the last three measurements in each bank. These measurements can be shared with your physician or trained healthcare professional.

1. When the unit is off, press and release the **MEM** button. The display will show the last memory bank used, U1 or U2, and the number of measurements in that bank, Fig 1. Press the **START** button if you wish to select a different memory bank.
2. Press the **MEM** button again to review the saved measurements. The first reading displayed is the average of the last three measurements in that bank, Fig 2. The average blood pressure measurement and pulse will display separately on the screen alternating repeatedly every 2-3 seconds, the A will remain on the display.
3. Press the **MEM** button again to scroll through all of the measurements. Each measurement will indicate the date and time taken, the measurement and pulse reading, Hypertension Indicator results and Irregular Heartbeat if present, Fig 3.
4. After you've scrolled through all the readings, the display will cycle back to the first reading in memory. Press **START** to turn the unit off.

When the number of readings exceeds 60, the oldest data will be replaced with the new record.


Fig. 1


Fig. 2

Fig. 3


Fig. 4

MEMORY FUNCTION

Clearing Measurements From Memory

1. Press and release the **MEM** button. The unit will display how many measurements are stored within the bank.
2. Select the desired user bank to clear by pressing the **START** button to switch between U1 and U2.
3. Press and release the **MEM** button two times to advance to the first measurement in the bank. You cannot clear the memory bank from the 'Average' screen.
4. When the first measurement is on the display (indicated by the 01 at the lower left of the screen), Fig. 1, press and hold the **MEM** button until two horizontal lines flash, Fig. 2. This indicates that all measurements have been erased. Press **START** to turn the unit off.


Fig. 2


Fig. 3

CARE & MAINTENANCE

Proper care and maintenance of your blood pressure monitor will help ensure long life for the product; improper use of the monitor may void the included warranty.

Always disconnect the cuff and tubing from the monitor prior to cleaning and storage.

Monitor Care & Maintenance

- Do not use any liquids on the monitor; use a soft, dry cloth to clean the monitor as needed.
- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: $-4^{\circ}\text{F} - 131^{\circ}\text{F}$ ($-20^{\circ}\text{C} - 55^{\circ}\text{C}$); Humidity: $\leq 90\% \text{RH}$ (non-condensing).
- This monitor requires 6 hours to cool from the maximum storage temperature between uses until the monitor is ready for its INTENDED USE when the ambient temperature is $68^{\circ}\text{F}(20^{\circ}\text{C})$.
- This monitor requires 6 hours to warm from the minimum storage temperature between uses until the monitor is ready for its INTENDED USE when the ambient temperature is $68^{\circ}\text{F}(20^{\circ}\text{C})$.
- Avoid high temperatures and direct sunlight.
- Always remove the batteries if the unit will not be used for an extended period of time.
- Do not attempt to disassemble your monitor; disassembling your monitor may void the manufacturers warranty.
- Do not attempt to service the monitor while in use or on a patient.
- Do not subject the monitor to strong shocks; take care not to drop the monitor.
- No component can be maintained by user in the monitor. The circuit diagrams, component part lists, descriptions, calibrations instructions, or other information which will assist the user's appropriately qualified technical personnel to repair those parts of equipment which are designated repairable can be supplied.

CARE & MAINTENANCE

Cuff Care & Maintenance

- The cuff may be spot cleaned with a mild-detergent as needed.
- To disinfect the cuff, wipe the inside (skin contact side) of the cuff using a cloth lightly moistened with 70% - 90% ethyl rubbing alcohol and allow to air dry; it is suggested to disinfect the cuff frequently when used in clinical settings or if the cuff is contaminated in any way.
- Never attempt to iron or dry the cuff; if the cuff becomes wet, allow the cuff to air dry completely before use or storage.
- Do not attempt to remove the tube and bladder from the nylon cuff as this may damage the cuff.
- Take care not to forcefully bend or kink the tubing during storage.
- It is recommended the performance should be checked every 2 years or after a repair. Please contact Customer Care at 866-326-1313.
- It is recommended to clean the cuff every 200 uses.
- The monitor can maintain the safety and performance characteristics for minimum of 10,000 measurements of three years.
- The cuff integrity is maintained after 1,000 open- close cycles of the closure.

DEVICE & LABEL SYMBOLS

These symbols may appear on your device, instructions or packaging and may vary by make and model.

Symbol | Meaning

	Read This Manual—All included manuals should be read prior to first use. (Background color will be blue; reading symbol will be white)
	Warning—Symbol indicates a warning, prohibition or mandatory action that mitigates a risk that is not necessarily obvious to the device operator.
	Type BF Applied Parts—Indicates that a part of this unit comes in contact with the patient in order to carry out its intended function; in the case of this device the cuff is the Type BF Applied Part.
	Environment Protection—dispose of this product properly; consult with your local recycling ordinances for proper recycling and disposal.
	Keep Dry—This device should be kept dry; never submerge the unit or cuff. Consult with the Care & Maintenance section of this manual for information on cleaning your monitor.
IP20	IP20—The first characteristic numeral symbol for “Degrees of protection against access to hazardous parts and against solid foreign objects. The second characteristic numeral symbol for “Degrees of protection against ingress of water”.
SN	Serial Number.
	Direct Current.
	Date of Manufacture.
	Unsafe.

The Electronic Sphygmomanometers corresponds to the below standards:
ISO81060-2 : 2013(Non-Invasive Sphygmomanometers- Part 2: Clinical Validation Of Automated Measurement Type).

DISPLAY SYMBOLS

Symbol	Meaning
SYS	Systolic Blood Pressure—The 'top' blood pressure result
DIA	Diastolic Blood Pressure—The 'bottom' blood pressure result
PULSE/Min	Pulse Rate—Pulse rate per minute during the measurement
mmHg	Millimeters of Mercury—Unit of blood pressure measurement for this device
	Low Battery—Battery life indicator
	Ready to Inflate Indicator—The cuff is not properly wrapped or has too much air left in it. The arrow will stop flashing when the unit is ready to use
	Irregular Heartbeat Symbol—An irregular heartbeat was detected during the measurement
	Memory—The display shows the measurements stored in memory
	Average—The display shows the average of the last three measurements in memory
	Error—An error has occurred during the measurement
	User 1/User 2—Indicates which user is active

TROUBLESHOOTING

Problem	Possible Cause	Solution
Display shows an abnormal result.	Cuff position was not correct or was not tight enough.	Apply the cuff properly and attempt a new measurement.
	Body posture was not correct during measurement.	Review the 'Taking Your Blood Pressure Reading' portion of this manual and attempt a new measurement.
	Speaking, excessive movement, angry, excited or nervous anxiety during measurement.	Wait a period of time and attempt a new measurement after calm has returned; do not speak or move during measurement.
Display is blank when power is on.	Batteries may be expired or installed improperly.	Check the polarity of the batteries and reinstall if necessary; replace batteries with new batteries.
No response when POWER is pressed or when new batteries have been loaded.	Incorrect operation; possible strong electromagnetic interference.	Remove batteries and wait a minimum of five minutes; reinstall batteries and attempt a new measurement.
	The battery is low.	Replace the batteries.
Monitor keeps re-inflating.	Circuit is locked.	Remove and reinsert the batteries; attempt new measurement.

ERROR CODES

Code	Meaning	Corrective Action
E0	Pressure system is unstable before measurement.	Re-apply the cuff and attempt new measurement.
E1/ E2	Failed to detect systolic/ diastolic pressure.	Sit comfortable and remain still; ensure the cuff is at the same level as your heart.
E3/ E4	Pneumatic system leakage or cuff is too loose/ too tight during inflation.	Re-apply the cuff and attempt new measurement.
E5	The cuff pressure is above 300mmHg.	The measurement range is over 300 mmHg; retake measurement and if error persists, consult your physician.
E6	More than 3 minutes with cuff pressure above 15 mmHg.	Re-apply the cuff and attempt to take a new measurement. Sit comfortable and remain still; attempt new measurement.
E7	EEPROM accessing error.	
E8	Device parameter checking error.	
ER A	Pressure sensor parameter error.	

TECHNICAL ALARM DESCRIPTION

The monitor will show 'HI' or 'Lo' as technical alarm on LCD with no delay if the determined blood pressure (systolic or diastolic) is outside the rated range specified in part SPECIFICATIONS. In this case, you should consult a physician or check if your operation violated the instructions.

The technical alarm condition (outside the rated range) is preset in the factory and cannot be adjusted or inactivated. This alarm condition is assigned as low priority according to IEC 60601-1-8.

The technical alarm is non-latching and need no reset. The signal displayed on LCD will disappear automatically after about 8 seconds.

FCC STATEMENT

NOTE:

POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for U.S.A. only)

This product has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. The product generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If the product does cause harmful interference to radio or television reception, which can be determined by turning the product on and off, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the product and the receiver.
- Connect the product into an outlet on a circuit different from that to which the receiver is connected, if applicable.
- Consult the dealer or an experienced radio/TV technician for help.

POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for Canada only)

This digital apparatus does not exceed the Class B limits for radio noise emissions from digital apparatus as set out in the interference-causing equipment standard entitled "Digital Apparatus", ICES-003 of the Canadian Department of Communications.

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Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

ELECTROMAGNETIC COMPATIBILITY

Table 1 EMISSIONS

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration - electromagnetic emissions

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Phenomenon	Compliance	Electromagnetic Environment Guidance
RF emissions	CISPR 11 Group 1, Class B	Home healthcare environment
Harmonic distortion	IEC 61000-3-2 Class A	Home healthcare environment
Voltage fluctuations and flicker	IEC 61000-3-3 Compliance	Home healthcare environment

ELECTROMAGNETIC COMPATIBILITY

Table 2 — ENCLOSURE PORT
For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Phenomenon	Basic EMC Standard	Immunity test levels
		Home healthcare environment
Electrostatic Discharge	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Radiated RF EM field	IEC 61000-4-3	10V/m 80MHz-2.7GHz 80% AM at 1kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	Refer to table 3
Rated power frequency magnetic fields	IEC 61000-4-8	30A/m 50Hz or 60Hz

ELECTROMAGNETIC COMPATIBILITY

Table 3 — PROXIMITY FIELDS FROM RF WIRELESS COMMUNICATIONS EQUIPMENT

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING Guidance and manufacturer's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Test frequency (MHz)	Band (MHz)	Immunity test levels
		Professional healthcare facility environment
385	380-390	Pulse modulation 18Hz, 27V/m
450	430-470	FM, ± 5 kHz deviation, 1kHz sine, 28V/m
710	704-787	Pulse modulation 217Hz, 9V/m
745		
780		
810		
870	800-960	Pulse modulation 18Hz, 28V/m
930		
1720		
1845		
1970	1700-1990	Pulse modulation 217Hz, 28V/m
2450		
5240		
5500	2400-2570	Pulse modulation 217Hz, 28V/m
5785		
	5100-5800	Pulse modulation 217Hz, 9V/m

PRODUCT SPECIFICATIONS

Name	SmartHeart Automatic Blood Pressure Monitor
Model Number	01-539
Display System	Digital display/LCD
Measuring Method	Oscillometric Method, Automatic Air Inflation and Measurement
Power Source	4 x 1.5V  AA Batteries
Measuring Range	Cuff Pressure: 0-300mmHg, Systolic Pressure: 60-260mmHg, Diastolic Pressure: 40-199mmHg; Pulse: 40-180 beats/minute
Accuracy	Pressure: ± 3 mmHg; Pulse: $\pm 5\%$
Deflation	Automatic pressure release valve
Memory	Built-in memory with 2 banks, 60 measurements per bank
Automatic Power-Off	Approximately 1 minute after last button operation
Battery Life	Approximately 500 measurements
Operation Environment	Temperature 50°F – 104°F (10°C – 40°C); Humidity $\leq 90\%$ RH
Storage Environment	Temperature -4°F – 131°F (-20°C – 55°C); Humidity $\leq 90\%$ RH
Monitor Dimensions	3-7/8" x 5-1/2" x 1-7/8"
Monitor Weight	7.4 oz. (without batteries); 10.6 oz. (with batteries)
Arm Circumference Range	Adult size cuff fits arm circumference: 11-3/4" – 16-1/2" (30 – 42 cm)
Accessories	Instruction Manual, 4 AA Batteries, Quick Start Guide, Storage Case, Replacement Cuff Order Form, Blood Pressure Log, FAQ's
Additional Options	Additional cuffs available on the enclosed order form or by calling Customer Care

Specifications are subject to change without notice

ONE-YEAR LIMITED WARRANTY

Congratulations on your purchase of a Automatic Blood Pressure Monitor. Your Automatic Blood Pressure Monitor is covered by the following limited warranty commencing upon the date of purchase, and subject to the following terms and conditions:

The Warrantor warrants that its Automatic Blood Pressure Monitor will be free from defects in materials and workmanship under normal consumer usage for a period of one year for the original purchaser of the product.

Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the Automatic Blood Pressure Monitor for commercial purposes or subjecting the Automatic Blood Pressure Monitor to abnormal usage or conditions; or (d) other acts which are not the fault of the Warrantor, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the Digital Blood Pressure Monitor.

If the Automatic Blood Pressure Monitor fails to conform to this limited warranty, return the Automatic Blood Pressure Monitor postage prepaid to: **Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031**. When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone number. Carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. At our option, the Warrantor will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will be responsible for parts charges, if any, and repair charges not covered by this limited warranty.

EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL THE WARRANTOR BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

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