
SMARTHEART™

Blood Pressure Monitor

Models # 01-509 series


INSTRUCTION MANUAL **ENGLISH AND SPANISH**

Please read this instruction manual completely before operating this unit.

STOP!
PLEASE ENSURE YOU HAVE ALL OF THE
FOLLOWING COMPONENTS BEFORE USING
YOUR BLOOD PRESSURE MONITOR


Monitor

Storage Case

DC Wall Adapter

Instruction Manual
Quick Start Guide
Blood Pressure Reading Tracking Log
FAQ Information Sheet

IF YOU ARE MISSING ANY PARTS, INCLUDING INSERTS OR INSTRUCTION MANUALS, DO NOT RETURN TO PLACE OF PURCHASE. CONTACT CUSTOMER CARE AT 866-326-1313.

INDEX

Care & Safety Information.....	4-7
Introduction & Indications for Use.....	8
STOP Notification Page.....	9
Blood Pressure Monitor Features.....	10
Getting Started.....	11
Charging the Battery.....	12-13
Setting Date/Time.....	14-15
Fitting & Applying Your Cuff.....	16-17
Taking Your Blood Pressure Reading.....	18-19
Interpreting Your Results.....	20-21
Memory Functions.....	22-23
Care & Maintenance.....	24-25
Device & Label Symbols.....	26
Display Symbols.....	27
Troubleshooting.....	28
Error Codes.....	29
FCC Statement.....	30
Electromagnetic Compatibility.....	31-33
Product Specifications.....	34
Warranty.....	35
Instrucciones en Español.....	37-71
Reference Page.....	72

Toll-Free Customer Care Help Line: 1-866-326-1313
Monday – Friday 8:30 a.m. – 4:30 p.m. CST

SmartHeart™
Manufactured for
Veridian Healthcare
1175 Lakeside Drive
Gurnee, IL 60031

Made in China

#93-1399 03/22
©2022 Veridian Healthcare, LLC

CARE & SAFETY INFORMATION

NOTE: Read all instructions carefully before use. The following basic precautions are needed when using an electrical product.

CAUTION: Failure to read and observe all precautions could result in personal injury or equipment damage.

Improper care or use of your blood pressure monitor may result in injury, damage to the unit or ineffective treatment. Following these instructions will ensure the blood pressure monitor's efficacy and long life.

GENERAL CAUTIONS AND WARNINGS

- Measurement results should be discussed with your physician or healthcare professional; never self-diagnose or attempt treatment as this may be dangerous.
- Intended for adult use only; this device is not approved for infant or child.
- Keep out of reach of children.
- Do not use the unit if it has any damaged parts, if it has been submersed in water or dropped.
- If any abnormality occurs, discontinue use until the unit has been examined and repaired.
- Only use included attachments and parts; do not use attachments from other brands or models as these may not be properly calibrated for use with this device and may result in measurement error.
- Do not share with an infected person.
- Do not disassemble the unit or attempt to repair it; substitution of a component different from that supplied may result in measurement error and will void manufacturer warranty.
- Always follow local regulations for proper disposal of the monitor, cuff and batteries.
- This device is intended for non-invasive measuring and monitoring of arterial blood pressure.
- It is not intended for use on extremities other than the arm or for functions other than obtaining a blood pressure measurement.

CARE & SAFETY INFORMATION

- To verify the calibration of this device, please contact the manufacturer.
- This device may be used only for the purpose described in this booklet. The manufacturer cannot be held liable for damage caused by incorrect application and use.ITROD
- The service life of the cuff may vary by the frequency of washing, skin condition of the user, and adherence to proper storage conditions. The typical service life is 10000 measurements.
- It is recommended that the performance should be checked every 2 years and after maintenance and repair, by retesting at least the requirements in limits of the error of the cuff pressure indication and air leakage (testing at least at 50mmHg and 200mmHg).
- The device doesn't need to be calibrated within two years of reliable service.
- Keep the unit out of reach of infants, young children or pets to avoid inhalation or swallowing of small parts. It is dangerous or even fatal.

OPERATING CAUTIONS AND WARNINGS

- Prolonged over-inflation may cause congestion, swelling or bruising in some people.
- Wait a minimum of 2 minutes between measurements (20-minutes is recommended); excessive frequent measurements may restrict blood circulation.
- Do not use this device while connected to any monitoring equipment on the same limb.
- Consult with your healthcare professional for guidance for use if you have any of the following: intravascular shunt; mastectomy on the arm corresponding to the measurement side of the body.
- Consult with your healthcare professional before using the device in pregnancy including pre-eclampsia, or if diagnosed with arrhythmia or arteriosclerosis.
- If you experience any discomfort or an abnormality occurs, stop using the device immediately; to stop the device during operation, press the **U1 or U2 (POWER)** buttons and the cuff will automatically deflate.

CARE & SAFETY INFORMATION

- Operation outside of stated operating temperature may result in measurement error or device malfunction; operation environment temperature is: 41°F – 104°F (5°C – 40°C); Humidity: 15% to 90% RH (non-condensing).
- Never use this unit while operating a vehicle or in the bathtub or shower.
- To avoid measurement errors, please avoid the condition of strong electromagnetic field radiated interference signal or electrical fast transient/burst signal.
- The blood pressure monitor, its adaptor, and the cuff are suitable for use within the patient environment. If you are allergic to polyester, nylon or plastic, please don't use this device.
- During use, the patient will be in contact with the cuff. The materials of the cuff have been tested and found to comply with requirements of ISO 10993-5:2009 and ISO 10993-10:2010. This cuff will not cause any potential sensation or irritation reaction.
- Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance *d* away from the equipment. The distanced is calculated by the MANUFACTURER from the 80MHz to 5.8 GHz column of Table 4 and Table 9 of IEC 60601-1-2:2014, as appropriate.
- The equipment is not AP/APG equipment and not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- On the rare occasion of a fault causing the cuff to remain fully inflated during measurement, open the cuff immediately. Prolonged high pressure (cuff pressure >300mmHg or constant pressure >15mmHg for more than 3 minutes) applied to the arm may lead to an ecchymosis.
- If the cuff pressure reaches 40 kPa (300 mmHg), the unit will automatically deflate. Should the cuff not deflate when pressures reaches 40 kPa (300 mmHg), detach the cuff from the arm and press any button to stop inflation.
- The operator shall not touch output of batteries /adapter and the patient simultaneously.

CARE & SAFETY INFORMATION

- The plug/adaptor plug pins insulates the device from the main supply. Do not position the device in a position where it is difficult to disconnect from the supply mains to safely terminate operation of ME equipment.

STORAGE CAUTIONS AND WARNINGS

- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: $-4^{\circ}\text{F} - 140^{\circ}\text{F}$ ($-20^{\circ}\text{C} - 60^{\circ}\text{C}$); Humidity: $\leq 93\% \text{RH}$ (non-condensing) at water vapor pressure to 50hPa.
- Keep the unit out of reach of children.
- After charging, unplug the adapter; do not store while plugged in or with adapter attached.
- When not in use, store the device with the adapter in a dry room and protect it against extreme moisture, heat, lint, dust and direct sunlight. Never place any heavy objects on the storage case.

CLEANING CAUTIONS AND WARNINGS

- Never immerse the unit in water to clean as it may damage the unit.
- Follow the 'Care & Maintenance' portion of this manual for instruction on how to clean and care for your monitor.
- Don't use any abrasive or volatile cleaners.

INTRODUCTION & INDICATIONS FOR USE

It is recommended that you first seek the advice and recommendation of your physician or healthcare professional when using home diagnostic devices, including blood pressure monitors.

Blood pressure monitors are for use in measuring blood pressure and heartbeat rate for those with an arm circumference ranging from 8.6" to 16.5" inches (about 22 – 42 cm). It is intended for adult indoor use only.

Blood pressure monitors use the oscillometric method to electronically measure your blood pressure. The monitor detects your blood's movement through the artery in your arm and converts the movements into a digital reading. The oscillometric method does not require a stethoscope, making the monitor ideal for home use.

Blood pressure readings determined with this device are equivalent to measurements obtained by a trained healthcare professional using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers.

Thank you for purchasing a Blood Pressure Monitor. With proper care and use, your monitor will provide you with many years of reliable readings.


**PLEASE READ THIS INSTRUCTION MANUAL
COMPLETELY BEFORE OPERATING THIS UNIT.**

— STOP! —

This device includes an ON/OFF switch located at the top of the monitor.

- **It is necessary to move the switch to the ON position for all functions of this blood pressure monitor.**
- **It is also strongly advised that the switch is left in the ON position if you are a daily or frequent user.** Turning the switch to OFF will result in having to reset the time/date feature and over long periods may erase memory history.
- This device will “power down” after readings with a blank screen and will not activate until the **POWER (U1 or U2)** button(s) are pressed even when the switch is left ON.
- Leaving the device in the ON position may slowly drain the rechargeable battery. Be sure to follow the battery prompts to know when to recharge (pg 12-13 of this manual).

BLOOD PRESSURE MONITOR FEATURES

NOTE: The USER 1 and USER 2 buttons also act as **POWER** buttons for this model

DETAILS OF THE DIGITAL DISPLAY SCREEN AND INTERPRETATION OF THE SYMBOLS ARE INCLUDED THROUGHOUT THIS MANUAL.

GETTING STARTED

Before taking a measurement

1. Your readings should only be interpreted by your physician or healthcare professional with access to your individual medical history. Regular use of a home blood pressure monitor will allow you to track and record your readings for discussion with your physician.
2. Conduct your measurement in a quiet place while seated in a relaxed position. Rest for 15 minutes before taking your reading. If an error occurs or you desire to take a second reading, allow 15 minutes between readings for your blood vessels to return to normal.
3. Avoid smoking, eating, taking medication, alcohol, physical activity or any other stressful activity for 30 minutes prior to taking a reading.
4. Always remove any jewelry or constrictive clothing that may interfere with the cuff placement.
5. Keep yourself and the monitor still during measuring; do not talk during the reading.
6. It is recommended that you take your readings at the same time each day to better monitor any indications in your results.
7. Record your daily measurement on the included chart or some other written document to share with your physician.

Contraindications

1. The device should not be used by any person who may be suspected of, or is pregnant.
2. The device is not suitable for use on patients with implanted, electrical devices, such as cardiac pacemakers, defibrillators.

CHARGING THE BATTERY

This device includes an internal rechargeable battery with DC Wall Adapter.

It is necessary to charge the battery when the unit will not power up, the low battery symbol appears on the display, when the display dims or when the display does not illuminate.



The 'Low Battery' symbol will appear on the lower right of the display.

Before first use, the device should be fully charged. Charging time may take up to 90 minutes from depleted battery.

Never attempt a measurement while the monitor is charging or connected to a wall outlet.

1. The switch at the top of the device must be in the ON position to charge; it will not charge if not in the ON position.
2. Connect the DC cable to the arm monitor and plug into a standard power outlet.
3. When initially connected, the display will illuminate; after several seconds the display will dim but the battery icon will remain. While charging the battery symbol will flash to indicate the current charge level.
4. When the device is fully charged, the battery symbol will stop flashing and remain solid on the display as long as the power cable is plugged in.

CHARGING THE BATTERY

5. Unplug the charger and monitor when fully charged. The 'battery' symbol will not be visible on the display during regular use until it displays as a 'low battery', indicating it is time to charge the monitor.
6. To maintain the settings in the monitor – date, time, stored measurement results – it is recommended that the battery not be allowed to fully deplete whenever possible.
7. A fully charged monitor can take approximately 1 month of daily measurements.

Adapter


Type: BLJ06L060100P-U

Input: 100-240V, 50-60Hz, 0.2Amax


Output: 6V === 1000mA

SETTING DATE/TIME

The date/time feature is on the bottom of the display screen. This unit stores the date, time and year allowing for a date/time stamp on all readings in the memory recall function. The date/time alternate view while in use.

The switch at the top of the device must be on the ON position to set the date and time feature.


HOW TO SET DATE AND TIME

1. When the unit is off (a blank display screen), press and hold the **U1/SET** button for 3-5 seconds to enter clock setting mode.
2. **HOURS/MINUTES**—first the hour will flash on the display. Press the **MEMORY** button to advance the display to the desired hour, press the **U1/SET** button to confirm your selection and advance to the minute setting. Press the **MEMORY** button to advance to the correct minute, press **U1/SET** to confirm.
3. **MONTH**—next the Month will flash on the display. Press the **MEMORY** button to advance the display to the desired month, press the **U1/SET** button to confirm your selection and advance to the date settings.

SETTING DATE/TIME

- DATE**—repeat the same sequence of step 2 to confirm the MONTH.
- YEAR**—repeat the same sequence of step 2 to confirm the YEAR.
- After confirming YEAR, the LCD will display “dONE” and then shut off.

PLEASE SET DATE AND TIME BEFORE INITIAL USE TO ENSURE YOUR READINGS ARE LABELED AND STORED IN MEMORY PROPERLY.


It is strongly advised that the ON/OFF switch is left in the ON position if you are a daily or frequent user. Turning the switch to OFF will result in having to reset the time/date feature and over long periods may erase memory history

FITTING & APPLYING YOUR CUFF

PROPER CUFF FIT AND PLACEMENT IS CRITICAL IN OBTAINING ACCURATE BLOOD PRESSURE MEASUREMENTS. PLEASE READ THESE INSTRUCTIONS CAREFULLY AND CONSIDER CONSULTING WITH YOUR PHYSICIAN FOR FURTHER DIRECTION.

Before applying your blood pressure cuff, be sure your arm is within the appropriate cuff range:

Arm Circumference Range Suitable for
8.6" – 16.5" (22 – 42 cm).

To determine the proper cuff size, measure the circumference of your left upper arm at the point midway between inside bend of your elbow and your shoulder, Fig 1. If your arm is outside of this range, this device may not be suitable for you and your results may be inaccurate.



Reminder, it is important to avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading.

If for any reason you are unable to or should not use your left arm, please modify the instructions for cuff application to your right arm. Your physician can help you identify which arm is best for you to take measurements from.

FITTING & APPLYING YOUR CUFF

Applying Your Cuff

1. Remove any constrictive clothing or jewelry that may interfere with cuff placement.
2. Be seated at a table or desk with your feet flat on the floor.
3. Hold your left arm in front of you with your palm facing upward.
4. Apply the cuff to your left arm so that the digital display face is positioned on the inside area of your arm facing you, and the word 'ELBOW' is centered over the inside bend of your elbow.
5. Adjust the cuff 3/4" - 1" from the bend of your elbow.
6. The cuff should fit comfortably, yet snugly around your arm. You should be able to insert one finger easily between your arm and the cuff.

If you are not comfortable with applying your cuff, please seek the assistance of another member of your household or work with your physician to practice the cuff application. Incorrectly applied cuffs may result in inaccurate readings.

TAKING YOUR BLOOD PRESSURE READING

Please read the preceding portions of this manual prior to taking your first reading.

Reminder, it is important to avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading.

The switch at the top of the device must be in the ON position to take your measurement; it will not take a reading if not in the ON position.

Select the desired user, 1 or 2, before measurement; this selection will determine where the measurement is stored for review.

When the monitor is off, press **U1/SET OR U2** button to POWER on the monitor. Throughout the measurement, the corresponding icon will show on the bottom of the digital display.


Fig. 1

1. Rest your elbow on a solid surface with your palm facing upward. Elevate your arm so that the cuff is at the same level as your heart, Fig 1. Be certain that cuff is elevated to the proper level; you may choose to support your hand with a folded towel or book. Relax your left hand.
2. When the monitor is off, press the **U1/SET** or **U2 (POWER)** button to POWER on the monitor; the unit will run a self test, Fig 2.
3. The display will briefly show '0' before the measurement begins.
4. You will see the numbers increase on the digital display during

Fig. 2

! YOU CAN STOP THE INFLATION OR DEFLATION PROCESS AT ANY TIME BY PRESSING THE POWER (U1/U2) BUTTONS.

TAKING YOUR BLOOD PRESSURE READING

measurement as the cuff inflates.

- When the measurement is complete, the cuff will automatically deflate and your blood pressure measurement and pulse reading results will show on the display screen at the same time, Fig 3.
- The Hypertension Indicator will indicate your reading range on the display, Fig 3. See page 21 for information regarding interpreting the Hypertension Indicator results.
- If an irregular heartbeat was detected during the reading, the Irregular Heartbeat Detection indicator will appear on the display, Fig 4. See page 20 for more information on the Irregular Heartbeat Detection feature.
- The reading will automatically be stored in the USER Memory Bank selected at the start of the reading, up to 250 readings per user.
- Select **U1/SET or U2 (POWER)** button to turn the unit off and conserve energy and battery life.
- The unit will automatically shut-off approximately 1 minute after use.

Fig. 3



Fig. 4

INTERPRETING YOUR RESULTS

Irregular Heartbeat Detection

An irregular heartbeat is defined as a heartbeat rhythm that is excessively slower or faster than the average measurement detected while the monitor is conducting the measurement. When this device detects an irregular rhythm two or more times during the reading, the IHB symbol will appear on the display screen.

This monitor will operate if the icon appears, but the results may be compromised, especially if this symbol appears often. Please consult with your physician or trained healthcare professional for further information regarding an irregular heartbeat and if this symbol appears frequently.

INTERPRETING YOUR RESULTS

Hypertension Indicator

This unit features our unique Hypertension Indicator. The color bars on the left side of the monitor display correspond with an icon on the digital display, indicating where the measurement results fall within the American Heart Association standards.

The American Heart Association has established globally accepted standards for the assessment of high or low blood pressure readings. The below chart should be considered only as a guideline, always consult with your physician or health care professional to interpret your individual results. Never adjust your medication(s) dosage, or severely alter your diet or exercise routine without consulting with your physician.

Bar Color	Category	Systolic (mmHg)		Diastolic (mmHg)
Green	Normal	<120	and	<80
Green	Prehypertension / Elevated	120-139	or	80-89
Yellow	High Blood Pressure, Stage 1 Hypertension	140-159	or	90-99
Orange	High Blood Pressure, Stage 2 Hypertension	160-180	or	100-110
Red	Hypertensive Crisis	>180	or	>110

MEMORY FUNCTIONS

Recalling Measurements in Memory

The switch at the top of the device must be in the ON position to review the memory banks; it will not allow review if not in the ON position.

1. When the monitor is off, press and release the **MEMORY**  button to access the memory bank.
2. Select **U1/SET** or **U2** to review the respective readings in memory.
3. The unit will first display the average of the last three measurements, Fig 1. If there are less than 3 readings in memory, the most recent reading will display.
4. Press the **MEMORY**  button again to begin scrolling through the saved measurements. Measurements will appear on the display from most current to oldest; the memory number alternate in the same position as the date/time stamp.
5. All results for a given measurement will display, including measurement results, pulse rate, Hypertension Indicator and Irregular Heartbeat alert (if applicable). The date/time stamp of the measurement will alternate as well, Fig 2.
6. When the number of readings exceeds 250, the oldest data will be replaced with the new record.



Fig. 1


Fig. 2

MEMORY FUNCTIONS

Clearing Measurements From Memory

To delete all results in a given bank:

1. When the monitor is off, press and release the **MEMORY**  button to access the memory bank feature.
2. Press and hold the **MEMORY**  button for 3 seconds; the display will flash dEL ALL, Fig 1.
3. Press the **MEMORY**  button again to confirm you would like to delete the memory bank, Fig 2; the display will show dONE and turn off.
4. Repeat the same steps for the second memory bank if desired.


Fig. 1

Fig. 2

CARE & MAINTENANCE

Proper care and maintenance of your blood pressure monitor will help ensure long life for the product; improper use of the monitor may void the included warranty.

Before taking a measurement

- At least 30 min required for ME equipment to warm from the minimum storage temperature between uses until it is ready for intended use. At least 30 min required for ME equipment to cool from the maximum storage temperature between uses until it is ready for intended use.
- Before use, make sure the device functions safely and is in proper working condition.

Monitor Care and Maintenance

- Do not use any liquids on the monitor; use a soft, dry cloth to clean the monitor as needed.
- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: $-4^{\circ}\text{F} - 140^{\circ}\text{F}$ ($-20^{\circ}\text{C} - 60^{\circ}\text{C}$); Humidity: 10% to 93%RH (non-condensing).
- Do not attempt to disassemble your monitor; disassembling your monitor may void the manufacturers warranty.
- Do not subject the monitor to strong shocks; take care not to drop the monitor.
- No component of this device can be maintained by the user; the circuit diagrams, component part lists, descriptions, calibration instruction or other information which will assist the user's appropriately qualified technical personnel to repair those parts of equipment which are designated repairable can be supplied.

CARE & MAINTENANCE

Cuff Care and Maintenance

- The cuff should not be submerged or washed in a washing machine, but may be spot cleaned with a mild-detergent as needed.
- Never attempt to iron or dry the cuff; if the cuff becomes wet, allow the cuff to air dry completely before use or storage.
- To disinfect the cuff, wipe the inside (skin contact side) of the cuff using a cloth lightly moistened with 70% - 90% ethyl rubbing alcohol and allow to air dry; it is suggested to disinfect the cuff frequently when used in clinical sittings or if the cuff is contaminated in any way.

DEVICE & LABEL SYMBOLS

These symbols may appear on your device, instructions or packaging and may vary by make and model.

Symbol	Meaning	Symbol	Meaning
	Read This Manual—All included manuals should be read prior to first use.		Keep Dry—This device should be kept dry; never submerge the unit or cuff. Consult with the Care & Maintenance section of this manual for information on cleaning your monitor.
	Warning—Symbol indicates a warning, prohibition or mandatory action that mitigates a risk that is not necessarily obvious to the device operator.		Direct Current.
	Type BF Applied Parts—Indicates that a part of this unit comes in contact with the patient in order to carry out its intended function.		Symbol for "MANUFACTURE DATE".
	Environment Protection—dispose of this product properly; consult with your local recycling ordinances for proper recycling and disposal.		Symbol for indoor use only.
	Serial Number.		Symbol for " Class II Equipment".
		F1	T1A/250V Φ 3.6*10CCC

DISPLAY SYMBOLS

Symbol	Meaning
SYS	Systolic Blood Pressure—The 'top' blood pressure result
DIA	Diastolic Blood Pressure—The 'bottom' blood pressure result
PULSE/MIN	Pulse Rate—Pulse rate per minute during the measurement
mmHg	Milligrams of Mercury—Unit of blood pressure measurement for this device
	Low Battery—Battery life indicator
	Pulse Reading—Appears when pulse reading is displayed
	Irregular Heartbeat Symbol—An irregular heartbeat was detected during the measurement
	User 1/User 2—Indicates which user is active
	Average—The display shows the average of the last three measurements in memory
	Memory—The display shows the measurements stored in memory
	Movement—Excessive movement was detected during reading; retake measurement

TROUBLE SHOOTING

Problem	Possible Cause	Solution
Display shows an abnormal result.	Cuff position was not correct or was not tight enough.	Apply the cuff properly and attempt a new measurement.
	Body posture was not correct during measurement.	Review the 'Taking Your Blood Pressure Reading' portion of this manual and attempt a new measurement.
	Speaking, excessive movement, angry, excited or nervous anxiety during measurement.	Wait a period of time and attempt a new measurement after calm has returned; do not speak or move during measurement.
Display is blank when power is on.	Battery may be dead / without charge.	Recharge the battery; follow the 'charging the battery' instruction portion of this manual.

ERROR CODES

Problem	Symptom	Check This	Remedy
No Power	Display will not light up.	Power is exhausted.	Charge the power.
Low Batteries	Display is dim or shows  +Lo	Power is low.	Charge the power.
Error Messages	E 1	The cuff is not secure.	Refasten the cuff, relax for a moment and remeasure.
	E 2	The cuff is very tight.	Readjust the cuff, not too loose or tight. Remeasure.
	E 3	The pressure of the cuff is excess.	Relax for a moment and remeasure.
	E 10 E 11	The monitor detected motion, talking or the pulse is too poor while measuring.	Relax for a moment and remeasure.
	E 20	The measurement process does not detect the pulse signal.	Loosen the clothing on the arm and then remeasure.
	E 21	The treatment of the measurement failed.	Relax for a moment and remeasure.
	Eexx	A calibration error occurred.	Retake the measurement. If the problem persists, contact the retailer or customer service for further assistance.

**IF AN ERROR CODE PERSISTS OR IS NOT REMEDIED,
PLEASE CONTACT CUSTOMER CARE AT 866-326-1313.**

FCC STATEMENT

NOTE:

POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for U.S.A. only)

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

ELECTROMAGNETIC COMPATIBILITY

Table 1

The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments

Warning: Don't near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation."

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Technical description:

1. All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the excepted service life.
2. Guidance and manufacturer's declaration -electromagnetic emissions and Immunity.

Guidance and manufacture's declaration - electromagnetic emissions

Emissions Test	Compliance
RF emissions CISPR11	Group 1
RF emissions CISPR11	Class B
Harmonic distortion IEC 61000-3-2	Class A
Voltage fluctuations / flicker emissions IEC 61000-3-3	Comply

ELECTROMAGNETIC COMPATIBILITY

Table 2

Guidance and manufacture's declaration - electromagnetic immunity

Immunity Test	IEC 60601-1-2 test level	Compliance level
Electrostatic Discharge IEC 61000-4-2	±8 kV contact ±2 kV, ±4kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV signal input/output 100 kHz repetition frequency	±2 kV for power supply lines ±1 kV signal input/output 100 kHz repetition frequency
Surge IEC 61000-4-5	±0.5 kV, ±1 kV differential mode ±0.5 kV, ±1 kV,±2 kV common mode	±0.5 kV, ±1 kV differential mode ±0.5 kV, ±1 kV,±2 kV common mode
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°.0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°.0 % UT; 250/300 cycle	0 % UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°.0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°.0 % UT; 250/300 cycle
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz
Radiated RF IEC61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz

NOTE UT is the a.c. mains voltage prior to application of the test level.

ELECTROMAGNETIC COMPATIBILITY

Table 3

Guidance and manufacture's declaration - electromagnetic immunity

Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band a (MHz)	Service	Modulation	Maximum Power (w)	Distance (m)	IEC 60601-1-2 Test Level (V/m)	Compliance level (V/m)
	385	380-390	TETRA 400	Pulse modulation 18Hz	1.8	0.3	27	27
	450	380-390	GMRS 460, FRS 460	FM \pm 5kHz deviation 1kHz sine	2	0.3	28	28
	710 745 780	704-787	LTE Band 13, 17	Pulse modulation 217Hz	0.2	0.3	9	9
	810 870 930	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18Hz	2	0.3	28	28
	1720 1845 1970	1700- 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217Hz	2	0.3	28	28
	2450	2400- 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217Hz	2	0.3	28	28
	5240 5240 5785	5100- 5800	WLAN 802.11 a/n	Pulse modulation 217Hz	0.2	0.3	9	9

PRODUCT SPECIFICATIONS

Name	SmartHeart Blood Pressure Monitor
Model Number	01-509 series
Display System	Digital display/LCD
Measuring Method	Oscillometric Method, Automatic Air Inflation & Measurement
Power Source	3.7V 1000mAH Built-in rechargeable li-polymer battery, 6V \equiv 1A AC Adaptor
Rated Cuff Pressure:	Pressure: 0-299 mmHg (0kPa~39.9kPa)
Measurement Pressure:	SYS: 60~230 mmHg (8.0kPa~30.7kPa) DIA: 40~130 mmHg (5.3kPa~17.3kPa) Pulse value: 40-199 beats/minute
Accuracy	Pressure: 41°F – 104°F (5°C – 40°C) within ± 3 mmHg; Pulse: $\pm 5\%$ of reading
Pressurization	Automatic pressurization by pump
Deflation	Automatic pressure release valve
Memory	Built-in memory with 2 banks, 250 measurements per bank
Automatic Power-Off	Approximately 1 minute after last button operation
Battery Life	Approximately 1 month with normal usage
Operation Environment	Temperature 41°F – 104°F (5°C – 40°C); Humidity 15% to 90% RH (non-condensing, but not requiring water vapor) partial pressure greater than 50hPa; Atmospheric pressure range of 700 hPa to 1060 hPa
Storage Environment	Temperature -4°F – 140°F (-20°C – 60°C); Relative humidity range of $\leq 93\%$ RH (non-condensing at water vapor pressure up to 50hPa)
Monitor Dimensions	5.15" x 2.87" x 1.15"
Monitor Weight	10.05 oz.
Arm Circumference Range	Suitable for 8.6" – 16.5" (22 – 42 cm)
Accessories	Instruction manual, DC adapter, Storage Case, Quick Start Guide, FAQ, Blood Pressure Log
Degree of protection	Type BF applied part
Protection against ingress of water	IP22
Software version	A01
Device Classification	Battery Powered Mode: Internally Powered ME Equipment AC Adaptor Powered Mode: Class II ME Equipment

Specifications are subject to change without notice.

5-YEAR LIMITED WARRANTY

Congratulations on your purchase of a Blood Pressure Monitor. Your Blood Pressure Monitor is covered by the following limited warranty commencing upon the date of purchase, and subject to the following terms and conditions:

The Warrantor warrants that its Blood Pressure Monitor will be free from defects in materials and workmanship under normal consumer usage for a period of five years for the original purchaser of the product.

Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the Blood Pressure Monitor for commercial purposes or subjecting the Blood Pressure Monitor to abnormal usage or conditions; or (d) other acts which are not the fault of the Warrantor, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the Blood Pressure Monitor.

If the Blood Pressure Monitor fails to conform to this limited warranty, return the Blood Pressure Monitor postage prepaid to: **Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031**. When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone number. Carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. At our option, the Warrantor will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will be responsible for parts charges, if any, and repair charges not covered by this limited warranty.


EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL THE WARRANTOR BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

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