# ***Premium Pulse Oximeter***

***for Sports and Aviation Enthusiasts***

Model 11-50DP

**PLEASE READ THESE INSTRUCTIONS  
CAREFULLY BEFORE USE.**

**INSTRUCTION MANUAL • ENGLISH & ESPAÑOL**

**Stop!**  
Please ensure you have all of the following  
components before use:

Unit with Quick-  
Release Lanyard

Micro USB  
Cable

Storage Case

Instruction  
Manual

**IF YOU ARE MISSING ANY PARTS, INCLUDING INSERTS OR  
INSTRUCTION MANUALS, DO NOT RETURN TO PLACE OF  
PURCHASE. CONTACT CUSTOMER CARE AT 866-326-1313.**

---

# INDEX

---

Important Safeguards - Care and Safety Information .....	4-5
Introduction .....	6
How Pulse Oximeters Work .....	7
Features of Your Device .....	8
PC Link Software .....	9
Display Symbols .....	10
Display Modes .....	11
Using Your Pulse Oximeter—Measurement .....	12-14
Using Your Pulse Oximeter—Alarms .....	15-17
Using Your Pulse Oximeter—Recording .....	18-19
Quick-Release Lanyard .....	20
Care and Maintenance .....	20
Battery Installation and Replacement .....	21
Troubleshooting .....	22
Device and Label Symbols .....	23
Electromagnetic Compatibility .....	24-27
Specifications .....	28
Warranty .....	29-30
Instrucciones en Español .....	31-60

---

**Toll-Free Customer Care Help Line: 1-866-326-1313**  
**Monday – Friday 8:30 a.m. – 4:30 p.m. CST**

Manufactured for:  
Veridian Healthcare  
1175 Lakeside Drive  
Gurnee, IL 60031  
[www.veridianhealthcare.com](http://www.veridianhealthcare.com)

Made in China

#93-1132 09/22  
©2022 Veridian Healthcare, LLC

---

# IMPORTANT SAFEGUARDS

## CARE & SAFETY INFORMATION

---

**NOTE:** Read all instructions carefully before use. Failure to read and observe all precautions could result in personal injury or equipment damage.

### **GENERAL CAUTIONS AND WARNINGS**

- Keep unit and lanyard out of reach of children as the included lanyard may present an entanglement or choking hazard to small children. Adult supervision required; never leave children unattended with unit or lanyard.
- To avoid electrical shock keep the unit away from water; never immerse the unit in liquid as this will damage the unit.
- Do not handle the unit with wet hands.
- Do not use the unit if it has been submersed in water or dropped.
- If any abnormality in product function occurs, discontinue use until the unit has been examined; discontinue use and return to the manufacturer for review and replacement under the terms of the two-year warranty.
- Do not shake the unit when in operation.
- Do not disassemble the unit or attempt to repair it.

---

# IMPORTANT SAFEGUARDS

## CARE & SAFETY INFORMATION

---

### **OPERATING CAUTIONS AND WARNINGS**

- Never tie or tape Premium Pulse Oximeter to your finger.
- Explosion hazard: Do not use the product in an explosive atmosphere.

### **STORAGE CAUTIONS AND WARNINGS**

- Never put the Premium Pulse Oximeter the unit in direct sunlight, high temperature or humidity.
- Keep the unit out of reach of children.

### **CLEANING CAUTIONS AND WARNINGS**

- Never immerse the Premium Pulse Oximeter in water to clean as it may be damaged.
- Follow the included instructions for cleaning the unit after each use; if the unit has not been used for an extended period, the device should be cleaned prior to use.

---

# INTRODUCTION

---

**Please read all instructions carefully prior to first use.**

Thank you for purchasing a Premium Pulse Oximeter. With proper care and use, your Premium Pulse Oximeter will provide you with many years of reliable service. The Premium Pulse Oximeter measures user blood oxygen saturation (SpO<sub>2</sub>) and pulse rate information for sports and aviation activities in a small, portable, easy-to-use device.

The Premium Pulse Oximeter is intended for use for persons engaging in sporting and aviation activities. Persons who need SpO<sub>2</sub> and pulse rate measurements because of a medical condition should not use the Premium Pulse Oximeter and should consult with their physician.


---

## HOW THE PULSE OXIMETER WORKS

---

The Premium Pulse Oximeter works by shining two light beams into the small vessels or capillaries of the finger to determine the user's blood oxygen saturation percentage and pulse rate that is then displayed on


the Premium Pulse Oximeter's screen. The unit uses the phrase "SpO2" to identify the oxygen saturation level and displays it as a percentage of full capacity (100%). Pulse rate is identified by the phrase "PR bpm" and is displayed as beats per minute.

---

## PRODUCT FEATURES

---


---

## PC LINK SOFTWARE

---

The Premium Pulse Oximeter features the capability of real-time PC monitoring and the ability to download a recorded monitoring session with the included USB cable and software.

Monitoring software can be downloaded at:

[www.veridianhealthcare.com/product-support/](http://www.veridianhealthcare.com/product-support/)

### **SYSTEM REQUIREMENTS:**


- IBM Compatible
- Windows XP, VISTA, 7, 8, 10
- USB Communications Port
- 16MB Memory
- 10MB Free Space on Hard Drive

---

## DISPLAY SYMBOLS

---

When in use, the Pulse Oximeter display reflects a 'live' reading of the wearers current SpO<sub>2</sub> and PR bpm.


---

# DISPLAY MODES

---

## SCREEN DISPLAY ORIENTATION


The Pulse Oximeter display can be adjusted to allow for easy viewing from any angle. While the unit is ON, rotating the finger while the oximeter is on and measuring will automatically change the display orientation, as shown below.


1.

2.

3.

4.

---

## USING YOUR DEVICE— MEASUREMENT

---

The Premium Pulse Oximeter may be used while the user is seated, standing or lying down. The user should not be walking or running during measurements and should take care not to excessively move the arm or hand where the oximeter is attached.

**The Premium Pulse Oximeter features the capability of real-time PC monitoring and the ability to download a recorded monitoring session. The included software instruction guide provides instruction for utilizing these features.**

### PRIOR TO USE

Keeping the sensors clean will help ensure accurate readings. Before and after each use, carefully swab inside the unit with a soft cloth or cotton swab with rubbing alcohol. Do not pour the alcohol directly on or in the unit.

Dirty hands or nail polish (especially dark shades) can also affect reading accuracy. Washing your hands and/or removing any nail polish prior to use to help ensure accurate readings.



---

# USING YOUR DEVICE— MEASUREMENT

---

## USING YOUR PREMIUM PULSE OXIMETER

1. Squeeze the hinged end of the unit to open the measurement side of the unit.
2. Insert your finger nail side up into the unit until resistance is felt. The unit should gently close below the first knuckle of the finger, and the sensors inside the unit should align with the finger nail and pad at the tip of the finger.

**Suitable for use on any finger (excluding the thumb). In our experience, the index or middle finger work best.**

---

## USING YOUR DEVICE— MEASUREMENT

---

3. Press the power button to turn the unit on. The screen will briefly flash a test display, and then begin the measurement.
4. Solid dashes (— —) will appear on the display for a few seconds, then the display will begin to show the oxygen saturation (SpO<sub>2</sub>) and pulse rate (PR bpm) measurements.

Take care to not move excessively during monitoring. These are active readings and may fluctuate slightly as you observe the display, responding to your heart rate (bpm) and other factors.

5. The Premium Pulse Oximeter should only be used for spot-checking blood oxygen saturation and pulse rate and should not remain on the finger for more than a few minutes at a time.
6. Remove the unit from the finger. The unit will automatically turn off after 5 seconds.
7. Clean the unit as directed before storage.

---

## USING YOUR DEVICE— ALARM MODES

---

The Premium Pulse Oximeter features several alarm functions—when the finger is out of position, when the measurement data is outside of the set limits and when the unit has low voltage (low battery). Alarms are set through the device menu.

**NOTE:** When navigating through all features in the following settings, pressing and HOLDING the **POWER** button will enter/change the selected setting; pressing and RELEASING the **POWER** button will toggle through the available options.

1. Press and hold the **POWER** button, the display will show the 'Settings' menu, Fig 1. To toggle within the 'Settings' menu, press and release the **POWER** button until the **Alarm** setting is selected.
2. Press and hold the **POWER** button to enter the Alarm settings menu; the top choice will be **Dir**, Fig 2.


Fig. 1


Fig. 2

**Dir** — Direction; allows user to set maximum and minimum settings for SpO<sub>2</sub> and Pulse Rate wherein the alarm will sound if the results of the ongoing measurement fall above or below the set parameters.

---

## USING YOUR DEVICE— ALARM MODES

---

**Psound** — Pulse Sound/Tone; only set as ON or OFF, when ON the oximeter will sound a tone with every heartbeat as it reads the user's pulse. i.e. beep-beep-beep-beep.

**Alarm** — Indicates whether the alarm function is ON or OFF; all functions—Finger Out, Reading Thresholds, and Low Battery—are ON when this setting is selected as ON, Fig 3.

<b>Dir</b>	<b>set</b>
Psound	off
▶ Alarm	on
Exit	

Fig. 3

3. Toggle the triangle cursor to Dir, press and hold to enter the reading thresholds options.

4. Press and hold the POWER button to enter the SPO2 setting. Within this setting you may select a High and Low point for the SPO2 measurement.

**Dir down**, numbers will go down as you set the digits; **Dir up**, numbers will go up as you set the digits. Press and hold the power button to select the desired setting. Press and release the power button to address each option—High, Low. Toggle to Exit, then press and hold to exit the SPO2 menu.

---

## USING YOUR DEVICE— ALARM MODES

---

5. Repeat the same steps to set the desired alarm threshold for the Pulse Rate (PR) measurement.

<b>Direction</b>	
SpO <sub>2</sub>	Alm
▶ PR	Alm
Exit	

**When the ALARM function is ON, the unit will beep continuously if the measurements are out of the specified range.**

---

## USING YOUR DEVICE— RECORDING

---

The Premium Pulse Oximeter allows for one recording session of up to 24 hours. This recording session can be viewed and saved with the included Pulse Oximeter Tracking Software.

1. Press and hold the **POWER** button, the display will show the 'Settings' menu, Fig 1. To toggle within the 'Settings' menu, press and release the **POWER** button until the Record setting is selected, Fig 1.


Fig. 1

2. Press and hold the **POWER** button to toggle between **OFF** and **ON**, Fig 2; selecting ON will immediately start the measurement recording.


Fig. 2

3. When the device is recording, a small circle will appear in the upper corner of the display screen, Fig 3.


Fig. 3

---

## USING YOUR DEVICE— RECORDING

---

4. When recording continuously for an extended period without touching the **POWER** button or removing the finger from the device, the unit will go into 'sleep' mode, and only the small circle will blink on the display.
5. To turn off the Record option, enter the Settings and set the Record function to off.

**NOTE:** This device will only record a single measurement session; the next recording session will override the previous session.



---

## QUICK-RELEASE LANYARD

---

The lanyard features a quick-release clip for easy release while in use.

**CHOKING HAZARD-** The included lanyard may present an entanglement or choking hazard to small children. Adult supervision required; never leave children unattended with unit or lanyard.


---

## CARE AND MAINTENANCE

---

Before and after each use, carefully swab inside the unit with a soft cloth or cotton swab with rubbing alcohol. Do not pour the alcohol directly on or into the unit.

The unit should be stored in a dry environment, within the ideal temperature range of -40°F–140°F(-40°C–60°C).

---

## BATTERY INSTALLATION AND REPLACEMENT


---

**This unit comes complete with 2 AAA batteries.**

It is necessary to replace the batteries when the Low Battery  symbol appears on the display or when the display does not turn on after the red power button is pressed.

1. Slide off the battery cover.
2. Insert or replace 2 AAA batteries into the battery compartment, ensuring to match the indicated polarity symbols.
3. Replace the battery cover.
4. It is recommended to remove the batteries if the unit will not be used for an extended period of time.

Always use new batteries. Do not mix alkaline, standard or rechargeable batteries. Please dispose of batteries according to local recycling regulations.

---

## **TROUBLESHOOTING**

---

### **SpO2 or PR is not displaying a reading though the power is on**

- Reinsert the finger, remain still

### **SpO2 or PR display readout is unstable and fluctuating highly**

- Finger may not be inserted deep enough
- Finger may be trembling, try to remain still

### **The unit will not turn on**

- Ensure the batteries have been properly installed
- Replace the batteries
- If the unit will not turn on after changing batteries, consult Customer Care

### **Display screen goes blank**

- The unit will automatically shut-off after 5 seconds if no finger is inserted
- Check the batteries and replace if needed

---

## DEVICE AND LABEL SYMBOLS

---

These symbols may appear on your device, instructions or packaging and may vary by make and model:

<b>SpO<sub>2</sub></b>	Blood Oxygen Saturation Value		Low Battery
<b>PRbpm</b>	Pulse Rate (bpm)		Alarm
	Serial Number	--	No Finger Inserted
	Type BP Applied Parts		Warning or Caution
	Environment Protection		Manufacturer
			Manufacturer Date

---

# ELECTROMAGNETIC COMPATIBILITY

---

## GUIDANCE AND MANUFACTURE'S DECLARATION-ELECTROMAGNETIC EMISSION

### FOR ALL EQUIPMENT AND SYSTEMS

Guidance and manufacture's declaration –electromagnetic emission		
The Pulse Oximeter is tended for use in the electromagnetic environment specified below. The user of the Pulse Oximeter should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The 11-50DP Pulse Oximeter uses RF energy only for their internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The 11-50DP Pulse Oximeter is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emission IEC 61000-3-3	Not applicable	

---

# ELECTROMAGNETIC COMPATIBILITY

---

## GUIDANCE AND MANUFACTURE'S DECLARATION-ELECTROMAGNETIC IMMUNITY

### FOR ALL EQUIPMENT AND SYSTEMS

Guidance and manufacture's declaration-electromagnetic immunity			
The Pulse Oximeter is tended for use in the electromagnetic environment specified below. The user of the Pulse Oximeter should assure that it is used in such an environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6KV contact ±8KV air	±6KV contact ±8KV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment

# ELECTROMAGNETIC COMPATIBILITY

## GUIDANCE AND MANUFACTURE'S DECLARATION-ELECTROMAGNETIC IMMUNITY

### FOR EQUIPMENT AND SYSTEMS THAT ARE NOT LIFE-SUPPORTING

Guidance and manufacture's declaration-electromagnetic immunity			
The Pulse Oximeter is tended for use in the electromagnetic environment specified below. The user of the Pulse Oximeter should assure that it is used in such an environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment-guidance
Radiated RF ICE 61000-4-3	3V/m  80MHz to 2.5GHz	3V/m	<p>Portable and mobile RF communication equipment should be used no closer to any part of the 11-50DP Pulse Oximeter, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>recommended separation distance</p> $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P} \quad 80\text{MHz to } 800\text{MHz}$ $d = \left[ \frac{7}{E_1} \right] \sqrt{P} \quad 800\text{MHz to } 2.5\text{GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: </p>
NOTE 1 At 80MHz and 800MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which The 11-50DP Pulse Oximeter is used exceeds the applicable RF compliance level above, the 11-50DP Pulse Oximeter should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the 11-50DP Pulse Oximeter.			
b Over the frequency range 150 KHz to 80 MHz, field strengths should be less than 3V/m.			

# ELECTROMAGNETIC COMPATIBILITY

## RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE EQUIPMENT OR SYSTEM FOR EQUIPMENT OR SYSTEM THAT NOT LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the 11-50DP Pulse Oximeter

The Pulse Oximeter is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Pulse Oximeter can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Pulse Oximeter as recommended below, according to the maximum output power of the communications equipment.

	Separation distance according to frequency of transmitter (m)		
Rated maximum output power of transmitter (W)	150KHz to 80MHz $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	80MHz to 800MHz* $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	800MHz to 2.5GHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

---

## SPECIFICATIONS

---

Name	Premium Pulse Oximeter
Model number	11-50DP
Display system	Digital display/OLED
Power source	2 "AAA" batteries
SpO2 Measuring range	0-100%
SpO2 Accuracy	70-100%, $\pm 2\%$ ; $\leq 69\%$ no definition
Pulse Rate Range	30-250 bpm
Pulse Rate Accuracy	30-250 bpm, $\pm 2\%$
Automatic Shut-off	After 5 seconds of no indication on the sensors
Operation Environment	Temperature 50oF – 104oF (10°C – 40°C); Humidity 30% – 75%  Atmospheric Pressure: 700hPa ~1060hPa
Storage Environment	Temperature -40oF – 140oF (-40°C – 60°C); Humidity 5% – 95%  Atmospheric Pressure: 500hPa ~1060hPa
Unit Dimensions	2-1/4" x 1-1/4" x 1-1/4"
Unit Weight	1.7 oz. (with batteries); 0.9 oz. (without batteries)

*Specifications are subject to change without notice.*

---

## 2-YEAR LIMITED WARRANTY

---

Congratulations on your purchase of a Premium Pulse Oximeter. Your Premium Pulse Oximeter is covered by the following limited warranty commencing upon the date of purchase for a period of two years, and subject to the following terms and conditions:

The warrantor warrants that its Premium Pulse Oximeter will be free from defects in materials and workmanship under normal consumer usage for as long as the original purchaser owns the product.

Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the Premium Pulse Oximeter for commercial purposes or subjecting the Premium Pulse Oximeter to abnormal usage or conditions; or (d) other acts which are not the fault of Veridian, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the Premium Pulse Oximeter.

If the Premium Pulse Oximeter fails to conform to this limited warranty, return the Premium Pulse Oximeter postage prepaid to: **Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031**. When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone number. Carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. At

---

## 2-YEAR LIMITED WARRANTY

---

our option, the warrantor will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will be responsible for parts charges, if any, and repair charges not covered by this limited warranty.

EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL THE WARRANTOR BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

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