


2-Second Digital Thermometer

Model 08-363

THANK YOU FOR PURCHASING A 2-SECOND THERMOMETER. PLEASE READ THESE INSTRUCTIONS CAREFULLY TO ENSURE PROPER USE AND SAFE OPERATION OF YOUR THERMOMETER.



PLEASE READ AND SAVE THESE INSTRUCTIONS

- This product is not intended to diagnose or treat any health problem or disease and should not be used as a substitute for the consultation and advice of a physician or other medical professional.
- Thermometers should only be used under the supervision of an adult.
- Before use, it is recommended that the patient avoid eating, drinking, smoking or exercising for at least twenty minutes prior to taking a measurement. Never walk, run or talk during a measurement.
- If the thermometer has been stored at a very low or high temperature, allow the unit to adjust to room temperature before use.
- Wait a minimum of 1 minute between measurements to allow the unit to readjust to room temperature.

NOTE: Measurement time may vary based on measuring conditions.

NORMAL BODY TEMPERATURE / FEVER INDICATOR

Body temperature can vary from 97°F and 100.0°F and still be considered normal. The medically accepted 'normal' body temperature is 98.6°F. Body temperature is commonly lower in the upon waking than at any point during the rest of an individual's waking hours.

Temperature readings will vary based on the body location point of the reading. Average 'normal' temperature ranges by measurement location site:

Oral 97.0°F–100.0°F (36.1°C–37.8°C)

Rectal 97.3°F–100.4°F (36.3°C–38.0°C)

Axillary 95.4°F–99.0°F (35.2°C–37.2°C)

This thermometer measures the heat flow rate from the user's body. If the heat flow for a given measurement location is within normal range, the thermometer backlight will be green; if the heat flow rate is high, the backlight will be red, indicating a fever MAY be present.


A fever is defined as a body temperature that is elevated above the normal for that person. It is important to determine what is normal for an individual before determining if a fever is present. Tracking an individual's temperature on a consistent basis, at the same body site and at the same time everyday when the person is well or healthy will help determine an individual's normal temperature.

CHANGING MODES - FAHRENHEIT OR CELSIUS OPERATION

This unit defaults to a Fahrenheit measurement when first used and will default to last used mode each time the unit is turned on.


To change modes (the device should be OFF):

- Press and hold the MODE and POWER buttons simultaneously.
- The current scale mode will briefly show on the display and then alternate to the new mode.
- The thermometer is now ready for measurement; either proceed with measurement or press the POWER button to turn the unit off.



MEMORY RECALL

This unit includes a recall of the last measurement taken. When the unit is OFF, press the POWER button to turn the unit on, and release. After the self-test, the unit will briefly show the last reading on the display.



OPERATING INSTRUCTIONS

Remove the LCD protective sticker. Press and release the POWER button to activate the thermometer. The unit will beep and the display will show all of the digital display symbols as it performs a self-test before showing the last reading in memory.

The thermometer will then show the current mode; press the mode button to select the desired measurement location.




The thermometer will default to Oral the first time it is turned; the thermometer will default to the last mode used for all other uses.

NOTE: If using a probe cover, apply before turning device on.

ORAL MEASUREMENTS


- Place the probe under the tongue as near as possible to a heat pocket as noted in the diagram. Readings will be progressively lower the further the probe is from the indicated heat pockets.
- Hold the thermometer in place during the reading; do not bite down on the thermometer. The mouth must remain closed to ensure accurate results. The reading should take approximately two seconds.
- The measurement will start automatically once the probe is in place and will take approximately 2 seconds. A running circle will briefly appear on the display indicating the reading is in progress.
- The unit will sound one long beep when the measurement is complete and the reading will appear on the display.
- Wait a minimum of one minute before attempting a new measurement.
- To turn the unit off, simply press the POWER button. The unit will automatically shut-off after approximately one minute.



CARE AND CLEANING INSTRUCTIONS


Store the unit in the provided protective case when not in use. Do not store the thermometer where it will be exposed to direct sunlight, dust or humidity. Avoid extreme temperatures. Do not attempt to disassemble the unit. If ERR displays and no measurement can be obtained, the unit must be replaced.

Clean the thermometer by wiping with rubbing alcohol (70% isopropyl alcohol) or soap and water. Do not submerge in any liquid or autoclave; this thermometer is NOT waterproof or shockproof. Wipe the thermometer with a dry cloth and allow to air dry.



AXILLARY (UNDERARM) MEASUREMENTS

- Wipe the underarm with a dry towel before application.
- Place the probe tip under the arm so that the tip is touching skin; hold the thermometer parallel, slightly angled upwards as indicated in the diagram. Position the arm across the chest to ensure the tip is completely covered. The reading should take approximately two seconds.
- The measurement will start automatically once the probe is in place and will take approximately 2 seconds. A running circle will briefly appear on the display indicating the reading is in progress.



TROUBLESHOOTING / IDENTIFYING DISPLAY SYMBOLS

Performance may be degraded should one or more of the following occur:

- Operation outside stated temperature and humidity range
- Storage outside state temperature and humidity range
- Mechanical shock

- The unit will sound one long beep when the measurement is complete and the reading will appear on the display.
- Wait a minimum of one minute before attempting a new measurement.
- To turn the unit off, simply press the POWER button. The unit will automatically shut-off after approximately one minute

The thermometer has been designed to display feedback-warnings to ensure proper usage and accurate temperature measurement. In the event of a malfunction or incorrect temperature measurement, the following warnings will be displayed:

Error	Problem/Cause	Solution
H	The measured temperature is above 107.6°F (42°C)	1. Turn the unit OFF and then ON again 2. Retake temperature 3. If still reads H, this indicates that the patient has a high fever; seek medical attention
LO	The measured temperature is below 95°F (35°C)	1. Turn the unit OFF and then ON again 2. Retake temperature 3. If still reads LO, this indicates that the patient has a low body temperature; seek medical attention
EF.1	The thermometer is too cold or too hot to operate	1. Verify that the room temperature within product range 60.8°F–104°F (16°C–40°C) 2. Retake temperature
EF.2	Temperature did not register because of poor contact, patient movement or poor placement	Turn the unit OFF, wait 2 minutes, turn it ON again and retake temperature
EF.3	Unit malfunctioned	Do not use the thermometer – contact Customer Care
LOW	Low battery - Battery icon flashes	Battery is low; enough power remains for approximately 10 measurements
DEAD	Drained battery - Battery icon is the only icon on the display and unit turns off	Battery is completely expired and must be replaced

FCC STATEMENT

IMPORTANT INFORMATION REQUIRED BY THE FCC

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference; and (2) this device must accept any interference received, including interference that may cause undesirable operation.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy, and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna.
Increase the separation between the equipment and receiver.

Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
Consult the dealer or an experienced radio/TV technician for help.

CAUTION: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

ONE-YEAR WARRANTY

Congratulations on your purchase of a 2-Second Digital Thermometer. Your Thermometer is covered by the following limited warranty commencing upon the date of purchase, and subject to the following terms and conditions:

The Warrantor warrants that its Thermometer will be free from defects in materials and workmanship under normal consumer usage for a period of one year for the original purchaser of the Thermometer.

Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the Thermometer for commercial purposes or subjecting the Thermometer to abnormal usage or conditions; or (d) other acts which are not the fault of the Warrantor, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the Thermometer.

If the Digital Thermometer fails to conform to this limited warranty, return the Thermometer to: **Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL 60031.** When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone number. Carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. At our option, the warrantor will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will be responsible for parts charges, if any, and repair charges not covered by this limited warranty.

EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL THE WARRANTOR BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

PROBE COVER ORDER FORM

Disposable Probe Covers 100 Pack (Model 08-322) Qty. ____ @ \$5.00 ea.

Subtotal \$_____
Illinois shipments add 8% Sales Tax \$_____
Shipping, handling & insurance \$_____
ORDER TOTAL \$_____
State _____ Zip _____

SHIP TO ADDRESS:

Name _____
Address _____
City _____
Daytime Telephone No. (_____) _____

BILL TO ADDRESS:

Name _____
Address _____
City _____
Daytime Telephone No. (_____) _____

PAYMENT METHOD:

Check or Money Order VISA MasterCard
 Discover Diners Club American Express

Please allow 2 weeks for delivery; this order form is for shipment within the Continental U.S. and Canada only.

Thank you for your order!

Please send a copy of this form with your payment to:
Veridian Healthcare, LLC
Attn: Customer Service
1175 Lakeside Drive
Gurnee, IL 60031

FAX CREDIT CARD ORDERS TO: 866-480-1717; QUESTIONS? 866-326-1313

Toll-Free Customer Care Help Line: 1-866-326-1313

Monday – Friday 8:30 a.m. – 4:30 p.m. CST

www.veridianhealthcare.com

Distributed by:
Veridian Healthcare
1175 Lakeside Drive
Gurnee, IL 60031
#93-1242 12/20
©2020 Veridian Healthcare, LLC
www.veridianhealthcare.com

Product of Israel
Packed in China
Daytime Telephone No. (_____) _____

