## **TENS + Smart Control Pain Management Solution**

22-035

### ***Intention for Use:***

To be used for temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, arm, and leg, due to strain from exercise or normal household and work activities.

It is also intended for symptomatic relief and management of chronic, intractable pain and relief of pain associated with arthritis.

*See pg 5 of this manual for further details.*


**PLEASE READ THESE INSTRUCTIONS COMPLETELY BEFORE USING YOUR TENS DEVICE**

**INSTRUCTION MANUAL • ENGLISH & ESPAÑOL**

**STOP! PLEASE ENSURE YOU HAVE ALL OF THE FOLLOWING COMPONENTS BEFORE USING YOUR TENS DEVICE**

Electric Stimulator

1 Large, 1 Small  
Electrode Pads

USB Cable

Storage Bag

Instruction Manual  
Replacement  
Electrode Order Form

**IF YOU ARE MISSING ANY PARTS, INCLUDING INSERTS OR INSTRUCTION MANUALS, DO NOT RETURN TO PLACE OF PURCHASE. CONTACT CUSTOMER CARE AT 866-326-1313.**

---

# INDEX

---

|   |       |
|---|-------|
| Introduction and Indications for Use..... | 4-5   |
| Care and Safety Information.....          | 6-13  |
| Product Features.....                     | 14-15 |
| Charging the Device.....                  | 16    |
| Preparing Your Device.....                | 17    |
| Selecting Your Gel Pad.....               | 18    |
| Placement & Treatment.....                | 19    |
| Pulse Program Descriptions.....           | 20-21 |
| TENS+EMS App.....                         | 22    |
| Instructions for Use.....                 | 23-30 |
| Caring for Electrode Pads.....            | 31    |
| Cleaning and Storage.....                 | 32    |
| Troubleshooting.....                      | 33    |
| Device and Label Symbols.....             | 34    |
| Electromagnetic Compatibility.....        | 35-38 |
| Product Specifications.....               | 39-40 |
| Warranty Information.....                 | 41-42 |
| Español.....                              | 43-84 |

---

**Toll-Free Customer Care Help**  
**Line: 1-866-326-1313**  
**Monday – Friday**  
**8:30 a.m. – 4:30 p.m. CST**

Distributed by  
Veridian Healthcare  
1175 Lakeside  
Gurnee, IL 60031  
[www.veridianhealthcare.com](http://www.veridianhealthcare.com)

Made in China

#93-1325 09/17  
©2017 Veridian Healthcare, LLC

---

# INTRODUCTION

---

Electronic Stimulator delivers electric pulses generated to the user's skin through the electrodes. The portable and compact device has multiple modes of different pulse frequencies, covering Transcutaneous Electrical Nerve Stimulation (TENS) and Electrical Muscle Stimulation (EMS). It includes operating elements of POWER ON/OFF button, intensity increase/mode change button, intensity decrease/timer change button, and can be attached and detached to the electrode through the two snap connectors.

---

## INDICATIONS FOR USE

---

### TENS (Modes 1, 2, 4, 5, 6, 8)

To be used for temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, arm, and leg, due to strain from exercise or normal household and work activities.

It is also intended for symptomatic relief and management of chronic, intractable pain and relief of pain associated with arthritis.

### EMS (Modes 1, 3, 7)

To stimulate healthy muscles in order to improve and facilitate muscle performance. To be used for the improvement of muscle tone and firmness, and for strengthening muscles in the arms, abdomen, legs, and buttocks. Not intended for use in any therapy or for the treatment of any medical conditions or diseases.

It is also intended to temporarily increase local blood circulation in the healthy muscles of lower extremities.

---

# IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

---

**NOTE:** Read all instructions and inserts included with this device carefully before use. The following basic precautions are needed when using an electrical product.

**CAUTION:** Failure to read and observe all precautions could result in injury or equipment damage.

Improper care or use of your device may result in injury, damage to the unit or ineffective treatment. Following these instructions will ensure the device's efficacy and long life.

## **CONTRAINDICATIONS :**

Do not use this device if you have a cardiac pacemaker, implanted defibrillator or other implanted metal or electronic device.

Do not use this device to treat undiagnosed chronic pain.

Do not use this device if you have suspected or diagnosed heart problem.

Do not use this device if you have suspected or diagnosed epilepsy.

Do not use this device if you have a tendency to bleed internally following an injury.

Do not use this device if you recently had surgery, or have ever had surgery on your back.

---

# IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

---

Do not use this device if areas of skin lack normal sensations, such as skin that tingles or is numb.

Do not use this device during pregnancy.

## GENERAL WARNINGS PRIOR TO USE

- If you are under the care of a Physician, consult with your Physician before using this device
- Long-term effects of this device are not known
- Do not place the device on or close to your heart
- Take care when applying device close to the neck; do not place this device on the front or sides of the neck. Severe spasm of the muscles may occur and the contractions may be strong enough to close the airway or cause difficulty in breathing. Stimulation over the neck could also affect hearing or blood pressure
- Do not apply stimulation across the chest because the introduction of electrical current into the chest may cause rhythm disturbances to the heart
- Do not place the device on or around your head. The effects of stimulation of the brain are unknown

---

# IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

---

- Do not use this device while sleeping
- Do not use if you feel numbness
- Do not use this device in or close to water
- Use the only on normal, healthy, clean and dry skin. Do not use on or close to open wounds or rashes, or over swollen, red, infected or inflamed skin
- Do not apply stimulation over, or in proximity to, cancerous lesions
- Do not use the device on children or incapacitated persons
- Consult with your physician before using this device, because the device may cause lethal rhythm disturbances to the heart in susceptible individuals

## **GENERAL CAUTIONS PRIOR TO USE**

- Read this manual before using this device for the first time
- Keep this manual available whenever you use this device
- This device is intended for individual person use only
- This device is not effective for pain associated with Central Pain Syndromes, such as headaches

---

## **IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS**

---

- This device is for pain caused by muscle soreness, and should be placed only around muscles where pain originates
- The pain may indicate that you have some other health problem. You should know the reason and source of your pain before using these devices. Do not solely rely on the treatment of this device for pain
- The safety of using this device during pregnancy or birth has not been established
- The effectiveness of this device depends greatly on a person's individual physical condition. It may not always be effective for every user
- If you have had medical or physical treatment for your muscle pain, consult with your treatment provider before using this device. You should contact your physician prior to using this device following recent surgical procedures. Stimulation may disrupt the healing process

---

# IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

---

## OPERATING WARNINGS & CAUTIONS

- CONSULT WITH YOUR PHYSICIAN IF YOU HAVE ANY OF THE FOLLOWING:
  - If you have suspected or diagnosed heart problem
  - If you have suspected or diagnosed epilepsy
  - If you have a tendency to bleed internally following an injury
  - If you recently had surgery, or have ever had surgery on your back
  - If areas of skin lack normal sensations, such as skin that tingles or is numb
  - During menstruation or during pregnancy
  - The unit is intended for the temporary relief of pain caused by muscle soreness. The source or cause of the pain should be addressed by a healthcare professional. Pain could indicate an injury or condition that requires treatment
  - Some people may feel skin irritation or experience a very sensitive feeling in the skin due to electrical stimulation. If this occurs, stop using these devices and consult your physician

---

## **IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS**

---

- If skin under the pad feels irritated after using the stimulator for a long period of time, use the stimulator for a shorter period of time
- Minor redness at stimulation placement is a normal skin reaction. It is not considered a skin irritation, and it will normally disappear within 30 minutes after the electrodes are removed. If the redness does not disappear after 30 minutes from the removal of electrodes, do not use the stimulator again until after the excessive redness has disappeared
- Turn off your device if the stimulation feels unpleasant or does not provide pain relief
- Keep this device out of the reach of children.
- Use this device only with the pads and accessories recommended by the manufacturer
- Do not use this device when driving, operating machinery or when swimming
- Before removing the Electrode Pad, be sure to turn it OFF, avoiding unpleasant stimulation
- TENS is not a substitute for pain medications and other pain management therapies

---

## **IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS**

---

- TENS devices have no curative value
- TENS is a symptomatic treatment and, as such, suppresses the sensation of pain that would otherwise serve as a protective mechanism
- Do not use the unit while sleeping
- After each use, clean any adhesive residue left on the skin with soap and water
- The effectiveness of stimulus therapy varies from person to person. It may not be an effective treatment for all users
- The long-term effects of stimulus therapy are unknown

## **STORAGE WARNINGS & CAUTIONS**

- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: 14°F – 122°F (-10°C – 50°C); humidity: 30-90% RH
- Keep the unit out of reach of small children

---

## **IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS**

---

- Do not store the device in direct sunlight, dusty or humid environments, or extreme temperatures. Exact figures for appropriate use and storage can be found in the Product Specifications section of this manual

### **CLEANING WARNINGS & CAUTIONS**


- Never immerse the unit in water to clean as it may damage the unit
- Follow the 'Cleaning and Maintenance' portion of this manual for instruction on how to clean and care for your device

---


# PRODUCT FEATURES

---

## *Electrical Stimulator*

**Front**


**Back**

---


# PRODUCT FEATURES

---

## *Large Electrode Pad*

## *Small Electrode Pad*

---

## CHARGING THE DEVICE

---

**This device includes a rechargeable lithium battery in the electrical stimulator unit and USB cable.**

Before first use, the device should be charged for approximately 2 hours. Subsequent charging will also take approximately 2 hours.

To charge the Electrical Stimulator Unit:

1. Connect the wide end of the USB cable to a USB wall adapter of 5V and 0.3-2.0A (not included).
2. Connect the smaller mini-USB end to the electrical stimulator.
3. For maximum power, allow to charge fully before use.

### NOTE:

A full charge will power 10 thirty-minute treatment cycles at a medium intensity setting.

During charging the power indicator light will flash green.

When charging is complete the indicator light will be solid green.

Unplug the wall charger when the charging cycle is complete.

---



## PREPARING YOUR DEVICE

---

It is necessary to prepare and assemble your device prior to treatment.

**NOTE: Save the plastic film discs for storage after treatment to protect the pads from dirt and debris.**

1. While the stimulator device is OFF (the green power light is NOT illuminated) carefully attach the electrode pads to the control unit via the snap connection points.
2. Pull gently to ensure the electrode pad is firmly attached.
3. Carefully remove the plastic film from the electrode pad.
4. You are ready to apply the pad to the affected treatment area.


**Do not turn the power ON until the device is in place on the user's skin as indicated in the 'Instructions for Use' section of this manual.**

---

## SELECTING YOUR GEL PAD

---

This device includes two sizes of gel electrode pads. There is no functional difference in the pad. Select the size for use most appropriate to the size of the treatment area.

---

# PLACEMENT & TREATMENT

---


## **Read all warnings and cautions completely before use.**

The electrode pads are applied directly to the skin. Clothing may be worn over the device; clothing should not constrict the device and should allow access to the control unit buttons.

Targeted treatment areas may include:

- shoulder
- elbow
- thigh
- knee
- calf
- ankle
- back of neck
- upper back
- lower back

**Never place device on over the heart, on the front or sides of neck or the face.**

# PULSE PROGRAM DESCRIPTIONS

|    | <b>Pulse Rate</b>        | <b>Description</b>                             | <b>Benefit</b>   |
|----|--------------------------|--|--|
| P1 | Combination of the below | Combination of the below                       | Pain relief and muscle stimulation (shown on pg 5)         |
| P2 | 69                       | Pulse on for 3.4 second and off for 1.6 second | Pain relief with a pulse rate of 69 Hz                     |
| P3 | 12.5–55.5                | Pulse on for 20 second and off for 1 second    | Muscle stimulation with variable pulse rate                |
| P4 | 1.2                      | Pulse on every 0.85 second                     | Pain relief with a pulse rate of 1.2 Hz                    |
| P5 | 100                      | Pulse on for 10 second and off for 2.5 second  | Pain relief with a pulse rate of 100 Hz                    |
| P6 | 100                      | Pulse on for 20 second and off for 1 second    | Pain relief with a pulse rate of 100 Hz and longer on time |
| P7 | 20                       | Pulse on for 5 second and off for 1 second     | Muscle stimulation with a fixed pulse rate of 20 Hz        |
| P8 | 160                      | Pulse on for 10 second and off for 2 second    | Pain relief with a pulse rate of 160 Hz                    |

---


# PULSE PROGRAM DESCRIPTIONS

---

The programs are not specific to treatment at a certain body location or to treat a specific ailment or pain.

The sensations may include—massage, acupuncture, tapping, scraping or a combination of the sensations at variable durations and cycles.

Select the program that provides the best relief for your particular ache and pain.

---

# TENS + EMS APP

---

This device includes a free download of the TENS+EMS app, available in the App store or on Google Play.

This App acts as a remote control to better operate your TENS unit, especially when the device is placed in hard-to-reach areas such as the back.

To download the app, open the appropriate app store via your mobile device and search for:

## TENS+EMS

**NOTE:** there is no 'space' between the words and + symbol.

Your Bluetooth must be ON in order for your mobile device to connect to the TENS unit. Your mobile device will automatically prompt you to accept the connection.

---

# INSTRUCTIONS FOR USE

---

**Read all warnings and cautions completely before use. Recommended treatment time is 30-minutes; consult your healthcare professional for a treatment cycle that best suits your injury. Discontinue use immediately if treatment is uncomfortable or painful.**

1. Treatment area should be clean before applying the control unit. Note that excessive skin lotion or soap residue may affect the adhesiveness of the pad. Only use on HEALTHY skin areas with no cuts, sores or wounds.
2. Apply the pad to the desired treatment area, taking care to ensure the controller is comfortably within reach of your hand for safety and control during treatment.
3. Ensure the electrode pad is smoothly applied and the gel adhesive is holding the pad firmly in place on the skin.

**This device is intended for use with the included TENS+EMS App and can also be used as a stand-alone TENS unit manually without use of the App.**

**Not all features are available unless paired with the App.**

---

# INSTRUCTIONS FOR USE

---

## Instructions for Use WITH the TENS+EMS App

4. After the device is applied to the desired treatment area, turn on the device by pressing the **POWER** button and turn on the App on your phone. You will automatically be prompted to turn on the Bluetooth if not already on. The bottom of the display will indicate 'Bluetooth not connected. Click this BOX to connect.'
5. The TENS unit will begin to operate at the lowest intensity level setting for program 1. Use the slide bar on the App to adjust the intensity to your desired level.
6. The sensation should be strong but comfortable; **Until you are familiar with use of the unit, begin with the lower intensity levels and gradually progress to more advanced levels as you feel comfortable.**
7. To discontinue use prior to the end of treatment time, press **PAUSE (II)** button or slide the power to the 0% position. To resume use, press the **PAUSE (II)** button again.
8. Be sure the device is OFF (no sensation is felt) before removing the device from your skin. Take care during removal to best protect the adhesive and preserve the life of the electrode pad.
9. Place the plastic cover over the gel side of the pad to prevent them from drying out.
10. To conserve the battery and prevent shocks, the device will automatically power off when not in use.

---

# INSTRUCTIONS FOR USE

---

## Instructions for Use WITH the TENS+EMS App —Home Screen

**Treatment Timer**

Press to go to the Timer Setting screen

**Treatment Intensity Slide**

The intensity will default to 0% at the start of treatment and whenever the program mode is changed

**Treatment Program**

Press to go to the Program Mode Screen

**Pause/Stop Treatment**

The timer will pause and remain paused until resumed or the App is closed; treatment will cease if the App is closed.


**NOTE:** If you leave the App, you will have to reconnect the device via Bluetooth before resuming treatment

---

# INSTRUCTIONS FOR USE

---

## Instructions for Use WITH the TENS+EMS App —Timer Setting


---

# INSTRUCTIONS FOR USE

---

## Instructions for Use WITH the TENS+EMS App —Program Mode Screen

### **Return**

Press to return to the main control screen

### **Program Modes**

Select the best program type for your desired treatment; detailed descriptions of each program type can be found on page 20 of this manual


**NOTE:** When changing modes, the intensity as set on the HOME screen will default to 0% and must be adjusted to your desired level.

### **Home/Back/Forward**

Press to return to the home screen or to the other screen options

# INSTRUCTIONS FOR USE

## Instructions for Use WITH the TENS+EMS App —Manual Mode Screen

### Return

Press to return to the main control screen

### Manual Mode

You may manually adjust the following elements to customize the treatment:  
Frequency  
Pulse Width  
Continuous  
Interval

**NOTE:** When changing modes, the intensity as set on the HOME screen will default to 0% and must be adjusted to your desired level.

### Home/Back/Forward

Press to return to the home screen or to the other screen options

---

# INSTRUCTIONS FOR USE

---

## Instructions for Use **WITHOUT** the TENS+EMS App

This device may be used without the TENS+EMS App. For maximum benefit of all features within the App, it is strongly suggested that App be utilized as the “remote control” to operated the device.

1. Follow the previous outlined steps in this manual to prepare your TENS device as shown on page 23 of this manual.
2. Apply the device to the treatment site before turning on”. If I was right, you could not increase the intensity of the device and the device would automatically shut off after the default 30 min if it was not placed on the skin..
3. Press and release the **POWER** symbol button. A GREEN led light will flash briefly and the device will beep to indicate the power is on.
4. The TENS unit will begin to operate at the lowest level setting. The default lowest level mode is a rotation of 7 pulse types at Program 1 as listed on page 20 of this manual.
5. Press and release the ‘+’ button to find your desired therapy sensation level. This device offers 20 intensity levels. Press and release the ‘-’ button to reduce the intensity level. As you press and release the button, the green light will flash and the device will beep 1x for each intensity level adjustment.

---

## INSTRUCTIONS FOR USE

---

6. To change the mode, while the unit is on/operating, press and hold the '+' for 3 seconds; the mode will switch to the next in the sequence of Mode 1 to 8. As you press and hold the button, the green light will illuminate 2x and the device will beep 2x to indicate the change in mode type.
7. The sensation should be strong but comfortable; press the + or - button at any time during the treatment to change the intensity as desired.
8. The device will pulse for several seconds, pause and resume repeatedly until the treatment is finished
9. The device will automatically shut off after 30 minutes of use. To discontinue use prior to the end of treatment time, press and release the **POWER** button to shut the unit off.  
**Do NOT attempt to remove the device during operation.**
10. Be sure the device is OFF (no sensation is felt) before removing the device from your skin. Take care during removal to best protect the adhesive and preserve the life of the electrode pad.
11. Place the plastic cover over the gel pads to prevent them from drying out. Store the device in the included storage case.
12. To conserve the battery and prevent shocks, the device will automatically power off when not in use.

---

# CARING FOR THE ELECTRODE PADS

---

## **Gel electrode pads**

- The electrode pads can be reused as long as they retain their stickiness. Cover the pads with the included film pieces when not in use to prevent them from drying out
- Properly cared for, the electrode pads may last up to 50 uses; the life of the electrode pad is largely dependent upon proper care and the condition of the skin as applied for use
- To help prolong the life of the electrode pad, between uses, moisten the gel side of the pad with a drop of water and allow to dry; replace the film for storage
- If you notice that the pads have become dry, no longer adhere well to the skin, or begin to irritate the skin, they are in need of replacement.

Replacement pads are available through Veridian Healthcare at 866-326-1313, by using the enclosed order form or through many major retailers and ecommerce sites.

---

# CLEANING AND STORAGE

---

## Cleaning

- Unsnap the electrode pad before cleaning the device
- The Electrical Stimulator can be cleaned with a soft, dry cloth
- Never use cleaning agents or excessive water to clean any part of the unit. Do not submerge any part of the unit or let it get wet
- Never disassemble the unit or attempt to repair the control unit

## Storage

- Always use the included plastic film to protect the gel electrode from dust and debris
- If the device will not be used for several days or longer, place the film cover over the electrode pads to prevent drying out or getting dirt and debris on them
- Store all components in the included bag while not in use
- Storing the pads in a plastic zip bag will help prevent the gel pads from drying out prematurely
- Protect the unit from mechanical shock or heavy impact
- Avoid exposure to extreme temperatures

---

# TROUBLESHOOTING

---

| <b>Problem</b>                              | <b>Solution</b>  |
|---|--|
| Stimulation is weak or non-existent         | Ensure the skin is clean and the pad is firmly adhered to the skin<br>Battery is low and must be recharged   |
| Device does not turn on                     | Battery is depleted and must be recharged  |
| Electrode pads are not adhering to the skin | Review the 'Caring for Electrode Pads' section of this manual; pads may need to be replaced  |
| Skin is turning red or irritated            | Stop treatment; wait until skin returns to normal color before attempting again.<br>Test the device on another area; if problem persists, consult your physician.  |
| Bluetooth will not connect                  | Review the 'Caring for Electrode Pads' section of this manual; pads may need to be replaced.<br>Turn on Bluetooth in your phone setting; turn off both the device and App and turn them on again; Switch to operate the device manually. |

# DEVICE & LABEL SYMBOLS

These symbols may appear on your device, instructions or packaging and may vary by make and model.

| Symbol  | Meaning   |
|---|---|
|   | Read This Manual/Consult Instructions Before Use  |
|   | Caution/consult accompanying documents before use |
|  | Use by Date                                       |
|   | Date of Manufacture                               |
|   | Batch Code  |
|   | Catalogue Number                                  |
|  | Serial Number                                     |

| Symbol  | Meaning                         |
|---|---------------------------------|
|  | Manufacturer                    |
|  | Temperature Limitation          |
|  | Humidity Limitation             |
|  | Non-sterile                     |
|  | Fragile, Handle with Care       |
|  | Keep Device Dry                 |
|  | Product Packaging is Recyclable |

---

# ELECTROMAGNETIC COMPATIBILITY

---

Table 1

## For all ME EQUIPMENT and ME SYSTEMS

### Guidance and manufacture's declaration - electromagnetic emissions

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

| <b>Emission Test</b>  | <b>Compliance</b> | <b>Electromagnetic Environment Guidance</b>   |
|---|-------------------|---|
| RF emissions<br>CISPR 11                                    | Group 1           | This monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.   |
| RF emissions<br>CISPR 11                                    | Class B           | This monitor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions<br>IEC 61000-3-2                         | Class A           |   |
| Voltage fluctuations/<br>flicker emissions IEC<br>61000-3-3 | Complies          |   |

# ELECTROMAGNETIC COMPATIBILITY

Table 2

## For all ME EQUIPMENT and ME SYSTEMS

### Guidance and manufacture's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

| Immunity Test  | IEC 60601 Test Level  | Compliance Level  | Electromagnetic Environment Guidance   |
|--|---|---|--|
| Electrostatic discharge (ESD)<br>IEC 61000-4-2   | ± 6 kV contact<br>± 8 kV air  | ± 6 kV contact<br>± 8 kV air  | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.   |
| Electrical fast transient/burst<br>IEC 61000-4-4   | ± 2 kV for power supply lines<br>± 1 kV for input/output lines  | ± 2 kV for power supply lines   | Mains power quality should be that of a typical commercial or hospital environment.  |
| Surge<br>IEC 61000-4-5   | ± 1 kV line(s) to line(s)<br>± 2 kV line(s) to earth  | ± 1 kV line(s) to line(s)<br>± 2 kV line(s) to earth  | Mains power quality should be that of a typical commercial or hospital environment.  |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>IEC 61000-4-11 | <5 % UT<br>(>95 % dip in UT)<br>for 0.5 cycle<br>40 % UT<br>(60 % dip in UT)<br>for 5 cycles<br>70 % UT<br>(30 % dip in UT)<br>for 25 cycles<br><5 % UT<br>(>95 % dip in UT)<br>for 5 s | <5 % UT<br>(>95 % dip in UT)<br>for 0.5 cycle<br>40 % UT<br>(60 % dip in UT)<br>for 5 cycles<br>70 % UT<br>(30 % dip in UT)<br>for 25 cycles<br><5 % UT<br>(>95 % dip in UT)<br>for 5 s | Mains power quality should be that of a typical commercial or hospital environment. If the user of this monitor requires continued operation during power mains interruptions, it is recommended that this monitor be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60 Hz) magnetic field<br>IEC 61000-4-8   | 3 A/m   | 3 A/m   | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.  |

# ELECTROMAGNETIC COMPATIBILITY

Table 3

## For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING Guidance and manufacturer's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

| Immunity Test                 | IEC 60601 Test Level     | Compliance Level | Electromagnetic Environment Guidance   |
|-------------------------------|--------------------------|------------------|--|
| Conducted RF<br>IEC 61000-4-6 | 3 Vrms 150 kHz to 80 MHz | 3 Vrms           | <p>Portable and mobile RF communications equipment should be used no closer to any part of this monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance:</b></p> $d = 1.2 \sqrt{\frac{P}{f}}$ $d = 1.2 \sqrt{\frac{P}{f}} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3 \sqrt{\frac{P}{f}} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,* should be less than the compliance level in each frequency range.** Interference may occur in the vicinity of equipment marked with the following symbol: ((⦿))</p> |
| Radiated RF<br>IEC 61000-4-3  | 3 V/m 80 MHz to 2,5 GHz  | 3 V/m            |  |

**NOTE 1** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

\* Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which this monitor is used exceeds the applicable RF compliance level above, this monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating this monitor.

\*\* Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

# ELECTROMAGNETIC COMPATIBILITY

Table 4

## For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

### Recommended separation distances between portable and mobile RF communications equipment and this monitor

This monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this monitor as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter<br>W | Separation distance according to frequency of transmitter<br>m |   |  |
|--|--|---|--|
|  | 150 kHz to 80 MHz<br>$d = 1.2 \sqrt{P}$                        | 80 MHz to 800 MHz<br>$d = 1.2 \sqrt{P}$ | 800 MHz to 2,5 GHz<br>$d = 2.3 \sqrt{P}$ |
| 0,01   | 0,12   | 0,12                                    | 0,23                                     |
| 0,1  | 0,38   | 0,38                                    | 0,73                                     |
| 1  | 1,2  | 1,2                                     | 2,3                                      |
| 10   | 3,8  | 3,8                                     | 7,3                                      |
| 100  | 12   | 12                                      | 23                                       |

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

---

# PRODUCT SPECIFICATIONS

---

|                               |  |
|-------------------------------|--|
| <b>NAME</b>                   | TENS with APP Pain Management Solution                         |
| <b>MODEL NUMBER</b>           | 22-035   |
| <b>NUMBER OF MODES</b>        | 8 auto pulse modes and 1 manual mode                           |
| <b>TIMER RANGE</b>            | 10-540 (minutes)   |
| <b>WAVEFORM</b>               | Biphasic   |
| <b>SHAPE</b>                  | Rectangular  |
| <b>CHANNEL</b>                | Single Channel   |
| <b>MAXIMUM OUTPUT VOLTAGE</b> | 58V@500Ω   |
| <b>MAXIMUM OUTPUT CURRENT</b> | 116mA@500Ω   |
| <b>PULSE DURATION</b>         | 0-500μSec  |
| <b>MAXIMUM FREQUENCY</b>      | 0-500Hz  |
| <b>PULSE INTENSITY</b>        | 0~116mA +/- 20%  |
| <b>AUTOMATIC SHUT-OFF</b>     | After completion of set therapy treatment time                 |
| <b>POWER SOURCE</b>           | 3.7V Rechargeable Battery                                      |
| <b>OPERATION ENVIRONMENT</b>  | Temperature 50 F - 104 F (10 C - 40C);<br>Humidity 30% - 85%   |
| <b>STORAGE ENVIRONMENT</b>    | Temperature 14 F - 122 F (-10 C - 50 C);<br>Humidity 30% - 90% |
| <b>DEVICE DIMENSIONS</b>      | 2" x 1.75" x 0.375"  |
| <b>WEIGHT</b>                 | Controller weight 0.5 oz                                       |

---

---

# PRODUCT SPECIFICATIONS

---

**INCLUDED** Electrical stimulator, large gel electrode pad, small gel electrode pad, USB cable, instruction manual, storage bag, replacement electrodes order form

---

**AVAILABLE SEPARATELY** Replacement electrodes 2-pack, 22-045

---

*Specifications are subject to change without notice*

---

# 1-YEAR LIMITED WARRANTY

---

Congratulations on your purchase of the TENS+Smart Control. Your TENS+Smart Control is covered by the following limited warranty commencing upon the date of purchase, and subject to the following terms and conditions:

The Warrantor warrants that its TENS+Smart Control will be free from defects in materials and workmanship under normal consumer usage for a period of one year for the original purchaser of the product.

Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the TENS+Smart Control for commercial purposes or subjecting the TENS+Smart Control to abnormal usage or conditions; or (d) other acts which are not the fault of the Warrantor, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the TENS+Smart Control.

If the TENS+Smart Control fails to conform to this limited warranty, return the TENS Pain Management Solution postage prepaid to: **Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031**. When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone number. Carefully package the product to avoid

---

# 1-YEAR LIMITED WARRANTY

---

any damage that may occur while in transit; shipping insurance with returned receipt is recommended. At our option, the Warrantor will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will be responsible for parts charges, if any, and repair charges not covered by this limited warranty.

EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL THE WARRANTOR BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

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