Ear & Forehead Infrared Thermometer

Model 09-342

THESE BRIEF INSTRUCTIONS ALLOW YOU TO QUICKLY BEGIN USING YOUR THERMOMETER. PLEASE REFER TO THE INSTRUCTION MANUAL FOR MORE DETAILED INSTRUCTIONS COVERING THE USE AND FEATURES OF THE THERMOMETER.

MEASURING TEMPERATURE BY EAR


NOTE:

Proper placement and position of the probe is crucial for effective and accurate measurement at the ear.

The ear canal should be 'straight' in order for the sensor probe to work most effectively. Gently pull the ear slightly up and away from the head with the opposite hand that is holding the thermometer

- For children under one year or who cannot sit on their own, lay the child flat and turn the head sideways, then gently pull the ear and place the probe.
- For older children and adults, stand at the shoulder or elbow of the seated patient and pull back the ear and place the probe.

MEASURING TEMPERATURE BY EAR

To measure body temperature at the EAR:

1. Remove the clear cap from the tip of the thermometer, Fig 1.
2. Press and release the **EAR / POWER** button to turn the thermometer on.
3. The display will illuminate and the unit will briefly run a self test, Fig 2.
4. The unit will beep twice, then show the most recent results in memory, Fig 3, and is now ready to take a new measurement.
5. Gently insert the probe into the ear canal until a slight resistance is felt.
6. Press the **EAR** button for a count of one-second and release; you will hear an audible beep. Do not remove the thermometer nor release the button until you hear the beep.
7. The display will show the measurement results, Fig. 4. If the audio is enabled, the results will also be announced aloud.
8. The LED fever alert signal light will illuminate, indicating where the results fall in the fever alert range.
9. After 60-seconds of non-use, the thermometer will automatically beep and shut off.


Fig. 1

Fig. 2


Fig. 3


Fig. 4

FEVER ALERT SIGNAL LIGHT:

GREEN LIGHT— If the results below 99.5°F (37.5°C) indicating no fever is present

AMBER (YELLOW) LIGHT — If the results are between 99.5°F - 100.2°F (37.5°C - 37.9°C) a slight fever is present.

RED LIGHT — If the results are 100.3 °F (38.0°C) and above, a fever is considered present.


MEASURING TEMPERATURE BY FOREHEAD

NOTE:

Proper placement and position of the probe is crucial for effective and accurate measurement.

The forehead should be clean and clear of excessive lotion or oil.

The temple aligns with the eyes on each side.

To measure body temperature at the FOREHEAD:

1. The clear cap should be in place during forehead measurements, Fig 1. The hole in the cap allows the sensor to work effectively, while the rounded surface of the cap helps with the 'swipe' action of taking a measurement at the forehead.
2. Press and release the **EAR / POWER** button to turn the thermometer on.
3. The display will illuminate and the unit will briefly run a self test, Fig 2.
4. The unit will beep twice, then show the most recent results in memory, Fig 3, and is now ready to take a new measurement.
5. While looking at the patient, place the probe tip on one temple. Press and hold the **HEAD** button while slowly sliding the thermometer across the forehead to the opposite temple while keeping contact on the skin.
6. Release the button, the unit will beep, and remove the thermometer from the patient's head. The display will show the measurement results, Fig. 4. If the audio is enabled, the results will also be announced aloud.
7. The LED fever alert signal light will illuminate, indicating where the results fall in the fever alert range.
8. After 60-seconds of non-use, the thermometer will automatically beep and shut off.


Fig. 1

Fig. 2

Fig. 3

MEASURING OBJECT/LIQUID TEMPERATURE

This thermometer can measure objects and liquids from 32°F to 212°F (0°C - 100°C). To measure the temperature of an object or liquid, follow these steps:

1. Press the **EAR / POWER** button to turn the thermometer on. The unit will run a self-test and all symbols on the display will momentarily appear, Fig. 1.
2. The unit will beep twice, then show the most recent results in memory.
3. Press **MEM** and **SET** simultaneously for 3-seconds to enter object measurement mode. The object mode symbol will appear on the display.
4. Aim the infrared lens at the object being measured, holding the thermometer approximately 1.0" away, Fig. 2.
5. Press the **EAR/POWER** button for a count of one-second and release; you will hear an audible beep
6. The display will show the measurement results, Fig. 3. If the audio is enabled, the results will also be announced aloud. The fever light indicator does not activate during object measurements.
7. After 60-seconds of non-use, the thermometer will automatically beep and shut off.


Fig. 1

Fig. 2


Fig. 3

PLEASE NOTE:

- Ear & Forehead temperature readings on this device are equivalent to oral readings.
- For accurate results, wait at least 10 seconds between measurements.
- Illness may be present even without an elevated temperature, particularly in the elderly, in infants under 3 months of age, and those with compromised immune systems. If illness is apparent without a fever, consult with a healthcare professional.
- Certain ear conditions, such as ear canal drainage, ear infection, or prior ear surgery may make ear thermometer readings impossible or inaccurate. If special ear conditions are present, consult with your healthcare professional.
- Never leave a child unattended with a thermometer.

IF YOU HAVE ANY QUESTIONS OR ARE EXPERIENCING DIFFICULTIES, PLEASE CALL THE TOLL-FREE CUSTOMER SERVICE HELP LINE FOR FURTHER ASSISTANCE AT 1-866-326-1313.

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