Instant Ear Thermometer

Frequently
Asked
Questions

FAQ'S

FREQUENTLY ASKED QUESTIONS

TABLE OF CONTENTS

General Product and Safety Information.....	2
Probe Covers.....	3
Accuracy.....	3-4
Cleaning.....	5
Using the Thermometer	5-6
Read-out Variances	7
Battery	8
Warranty.....	8
Instrucciones en Español.....	9-16

Customer Care: **866-326-1313**

GENERAL PRODUCT AND SAFETY INFORMATION

What are the advantages of measuring body temperature at the ear?

The hypothalamus, which is the part of the brain that regulates body temperature, is very close to the eardrum. Ear thermometers, unlike oral, rectal, or underarm thermometers, can detect real-time changes in body temperature. Measuring temperature at the ear is also faster and less intrusive than other methods.

How does the Instant Ear Thermometer work?

The lens captures energy emitted from the eardrum and surrounding tissue and the thermometer converts this energy into a temperature reading.

Is the thermometer dangerous or harmful?

Not at all. The thermometer only receives heat and infrared energy; it does not radiate or emit anything.

Is the thermometer safe for children and infants?

Absolutely. The probe can be fully inserted into children's and infants' ears. The probe is too short to reach the eardrum, and the ear canal is smaller than the thermometer's probe, so it cannot reach any sensitive tissues or structures in the ear. Be sure to follow the instructions, inserting the probe snugly but stopping when resistance is felt.

PROBE COVERS

Should I use probe covers with the Instant Ear Thermometer?

No. The thermometer is calibrated to work without the use of probe covers, for your convenience. Always clean the lens before and after use to maintain accuracy and prevent cross-contamination.

ACCURACY

Is my ear thermometer as accurate as devices used in hospitals?

This thermometer is designed and manufactured to meet the standards used for professional thermometers, and its accuracy has been proven in clinical use.

How do ear thermometer measurements compare to temperature readings taken using other methods?

Ear measurements and oral measurements are comparable; the normal range of temperatures is 97.0° F to 99.0° F (36.1° C – 37.2° C). Rectal temperatures are typically 1.0° F (0.5° C) higher, while underarm temperatures are 1.0° F (0.5° C) lower.

ACCURACY

Why do multiple measurements produce different readings?

Many factors may contribute to differing results. To prevent different readings, clean the lens before and after each use, use the same ear for each measurement, position the probe in the same way each time, and wait at least 10 seconds between readings.

Is the thermometer still accurate after being dropped or subjected to another heavy impact?

Your thermometer performs a self-test each time it is turned on. As long as the self-test display (all indicators shown on the display at once) is shown when the thermometer is activated, an accurate reading should occur if proper technique is used.

CLEANING

How often should I clean the thermometer?

The lens should be cleaned before and after each use, to eliminate earwax and fingerprints and ensure the accuracy of your readings. Use an alcohol swab, or a soft cotton swab moistened with 70% isopropyl alcohol, to gently cleanse the probe lens and thermometer casing. Never use a paper towel to clean your thermometer, as it may scratch the lens. Do not use cleaning agents other than alcohol to clean the thermometer; abrasives, thinners or gasoline will damage the unit. Do not immerse the thermometer or allow liquid to enter the casing.

Is the thermometer waterproof?

The probe tip is waterproof, but other parts of the unit (such as the display, battery cover, or buttons) are vulnerable to water damage. Never immerse the thermometer or allow liquid to enter the casing.

USING THE THERMOMETER

Can the thermometer be used on infants? What about the elderly?

The thermometer is safe and accurate for all ages, and is especially beneficial for use in infants and children due to its non-invasive nature and the speed with which a reading can be taken. Keep in mind, illness may be present even without an elevated temperature, particularly in the elderly, in infants under 3 months of age, and those with compromised immune systems. If illness is apparent without a fever, consult with a healthcare professional. Never leave children unsupervised with a thermometer.

USING THE THERMOMETER

Which ear should I use when taking measurements?

There is generally no significant difference between ear temperatures, although special conditions such as infection, ear canal drainage, or previous ear surgery may make one ear more suitable than the other. For best results, always use the same ear to take measurements for the same patient.

Why do different ears in the same patient produce different readings?

Individual variations or a differing amount of wax or dirt in the ears may cause different readings. The position of the probe may also produce differing results. For best results, always use the same ear to take measurements for the same patient.

Can the Instant Ear Thermometer be used to measure temperature at sites on the body other than the ear?

No. This thermometer is designed specifically for ear temperature measurements.

USING THE THERMOMETER

Can I use the thermometer to take the temperature of an animal?

No. The thermometer is meant for human use only.

Can I use the thermometer to determine ovulation by measuring basal body temperature?

No. This unit is not appropriate for basal readings.

READ-OUT VARIANCES

Can earwax affect the reading?

While a small amount of earwax in the ear does not have an effect, a larger amount may affect results. Earwax on the lens may cause an inaccurate reading. Therefore, clean the lens before and after each use.

Can an ear thermometer be used when the patient has an ear infection?

To prevent unnecessary pain or discomfort, we recommend that the measurement be taken in an unaffected ear or with another type of thermometer at another measurement site.

Will ambient air temperatures affect results?

An accurate reading can be produced in an air temperature anywhere between 60.8° F and 104.0° F (16.0° C – 40.0° C).

BATTERY

What type of battery is used in the unit?

One CR2032 lithium battery is used in the thermometer.

How long should the battery last?

The battery will last for approximately 5,000 measurements.

WARRANTY

What is the warranty on the Instant Ear Thermometer?

The warrantor guarantees that its Instant Ear Thermometer will be free from manufacturing defects under normal use for as long as the original purchaser owns the product. This warranty covers only normal use, and does not apply to use in any clinical or commercial applications. This warranty does not cover batteries or other power sources that may be provided with or used with the Digital Thermometer. This warranty is voided if the Digital Thermometer product is misused or abused in any manner.

If the Instant Ear Thermometer fails to operate during the time the original purchaser owns it, return it postage prepaid to: Repair Department, 1175 Lakeside Drive, Gurnee, IL 60031. The warrantor will repair or replace the defective unit at the option of the warrantor. Repair or replacement of the defective unit, at warrantor's option, is the sole remedy under this warranty. For complete detailed warranty information, refer to the unit instruction manual.

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©2017 Veridian Healthcare