
# SMARTHEART™

## Pedometer with Body Fat Function

Model # 19-006

### Function

- |                 |                         |
|-----------------|-------------------------|
| 1. LCD Display  | 6. Set/Reset Button     |
| 2. Up           | 7. Sensor Plates        |
| 3. Down         | 8. Battery Cover Screws |
| 4. Mode         | 9. Belt Clip            |
| 5. Start Button |                         |

### Features

- Step counter up to 99999 steps
- Distance traveled up to 999.99 KM/ 999.99 Miles
- Calories burnt up to 9999.9 KCal
- Determine body fat level fat by human body electric resistance
- Five groups user body fat data memory for five users
- Real-time clock
- Walking Timer function
- Metric and Imperial system selection

### Mode & Up/Down Key Functions Table

Time	Step Timer	Pedometer	Body Fat
▲ 12/24 HRS	▲ Step Distance Calorie Stride	▲ Memory ▼ Personal Data ▼ Weight	

### SETTING UP YOUR PEDIOMETER

Prior to use, set the date and time and user specifications to help ensure accurate measurements and results.

#### Time Setting

1. Press any button to wake unit from sleep (blank screen) mode.
2. Press and release **MODE** to cycle through available screens to the real time clock screen.
3. Press and hold **SET** until the digits on display flash.
4. Press **UP / DOWN** to scroll to current HOUR; press **SET** to confirm.
5. Press **UP / DOWN** to scroll to current MINUTES; press **SET** to confirm.

#### 12/24 Hour Format Setting

- When the display shows the current time, press **UP** to change to 12-Hour or 24-Hour format

#### Setting the Stride

To calculate distance traveled and calories burned accurately, you must first determine your stride length. The default stride is 40" / 40 cm.

1. Press and release **MODE** to cycle through available screens until the display shows Step.
2. Press **UP/DOWN** until the word "Stride" shows up on the display.
3. Press and hold **SET** until the digits on display start to flash to adjust the Stride length in inches (cm)
4. Press **UP / DOWN** to scroll to scroll to the correct length; press **SET** to confirm.

To determine your average stride length: select a starting point, walk 10 steps. Measure the distance traveled from starting point to tip of forward foot and divide by 10. This is your stride length.

**NOTE:** Stride Range: 10" - 72" (30cm - 180cm)

#### Setting User Body Weight

To ensure the calorie calculations are correct you will need to input the correct weight.

1. Press **MODE** until the top of the display shows "Step".
2. Press **UP** or **DOWN** until the display shows "KG" or "LB".
3. Press the **SET** for inputting your body weight. The digits on the display will start to flash ##.
4. Press **UP** or **DOWN** to scroll to your correct weight.
5. Press **SET** to confirm setting.

**NOTE:** Range: 22lb - 500lb (10kg - 250kg); Default weight: 110lb/50kg

#### Setting User Body Details for Body Fat

To ensure the Fat Scan is correct you will need to input the correct weight.

At normal time display, press **MODE** twice to enter Body Fat Mode.

1. Press **UP**, you will see 5 different user groups memory; M1 - M5.
2. Choose one of the 5 data groups to represent your personal data.
3. Press and hold **SET** to enter personal data setting mode.
4. Press **SET** to shift settings in sequence: Weight/Height/ Age/Gender/Confirm Setting.
5. Press **UP** or **DOWN** to set the desired values.

#### Disclaimer:

These calculations are not intended to take the place of the advice and instruction from your personal health provider. Please consult with your healthcare professional before starting an exercise program.

### USING YOUR PEDIOMETER

#### Step Counter Function

Steps are detected via the movement of the waist. For accurate readings, mount the unit closest to the center of the body.

- Press **MODE** until the top of the display shows the word "Step". The step counter will automatically start counting once motion is detected.
- Your counter will continue even if unit appears 'asleep' / no display is shown. Press any button to turn on the display.

**NOTE:** 5-Step error correction- Pedometer must detect 5 consecutive steps before steps begin being counted. This is to eliminate step counts from accidental movement.

#### Step Timer Function

The pedometer can track the actual time you spend walking with the step timer function.

- Under step timer mode, the timer will start counting once walking motion is detected.
- Press **MODE** and scroll to the walking figure to see the amount of time spent walking.
- The timer will show your cumulative walking time for a day - ie taking a morning walk for fitness, your daily activity, and other home routine will give you the total walking minutes for the entire day.

#### Distance Traveled Function

The pedometer can track the actual distance walked with the distance traveled function.

- Press **MODE** until the top of the display shows the word "ML" (KM).
- The distance traveled counter will automatically start counting once motion is detected.
- The distance traveled will show your cumulative walking distance for an entire day.

#### Calorie Calculation

The pedometer can track the calories burned throughout your active day based on your tracked steps.

- Press **MODE** until the top of the display shows "STEP".
- Press **UP/DOWN** until the display shows "KCAL" on the bottom right hand side of the display
- The pedometer will automatically start calculating once motion is detected.

To reset any value - ie step counter, timer, distance or calorie counter - press and hold the **SET** button until the digits revert back to zero.

**NOTE:** Resetting the timer back to zero also resets the step count, calorie calculation and distance traveled back to zero.

#### Metric and Imperial

The pedometer can be set to either Imperial or Metric system.

- Under any mode, press and hold the **MODE** for 3-4 seconds.
- Press **MODE** to "Body Fat" mode and check the top line indicator to see if the unit is in imperial or metric format.

#### Body Fat Scanning Feature

- Go to body fat mode.
- Press **UP** button to select the data group memory (M1-M5) which contains your personal data.
- Press **START** button to start scanning process.
- Put 2 thumbs on the two sensor plates for around 4 seconds and the result will show on the display. It will show your body fat percentage.

	Lean	Fit	Acceptable	Obese	Overweight
Men	5 - 9.9%	10 - 14.4%	14.5 - 24.9%	25 - 34.9%	35 - 50%
Woman	5 - 9.9%	10 - 19.9%	20 - 29.9%	30 - 34.9%	35 - 50%

Error: Men and Women >50%; Men & Women <5%.

If "ERR" message shows up on the display, it means poor skin contact on the sensor plates. Please try again.

**NOTE:** "ERR" message persists, the skin on your thumbs may be too dry. You will need to moisten your thumbs with little water to enable contact on the sensor plates.

#### Battery Replacement

Battery included: One (1) CR-2032

Change the battery when display is dim or won't activate.

1. Carefully unscrew the 2 screws on the top left and right of the LCD Display.
2. Take the button cell battery "CR-2032" out and replace with fresh button cell battery.
3. Close the cover and secure the unit with the 2 screws.

#### Auto Shut Off / Sleep

If no motion is detected for 1 minute, the unit will automatically go to sleep. The value of the pedometer will still be retained in memory. It will resume once motion is detected.

#### Care and Maintenance

- Do not disassemble. Doing so could cause damage and affect the accuracy of the pedometer measurement.
- Replace battery with care. Do not touch other pedometer parts. Doing so could damage pedometer function.
- Do not get wet; avoid excessive perspiration.

#### Troubleshooting

Problems	Cause & Remedy
Not all your steps were detected.	Mount the pedometer properly according to the manual and ensure your walking method is correct.
LCD is black.	The pedometer has been put under direct sunlight or operated in high temperature too long. It will recover normal status under the shade.

Distributed by:  
Veridian Healthcare  
1175 Lakeside Drive  
Gurnee, IL 60031  
www.veridianhealthcare.com

Made in China  
#93-1554 06/22  
©2022 Veridian Healthcare, LLC

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