# SMARTHEART™

---

## Automatic Digital Blood Pressure Arm Monitor

Model # 01-563


### **INSTRUCTION MANUAL** **ENGLISH AND SPANISH**

Please read this instruction manual completely before operating this unit.


**STOP!**  
**PLEASE ENSURE YOU HAVE ALL OF THE**  
**FOLLOWING COMPONENTS BEFORE USING**  
**YOUR DIGITAL BLOOD PRESSURE MONITOR**


Digital Monitor

Wide-Range Cuff

4 AA  
Batteries

Instruction Manual  
Quick Start Guide  
Blood Pressure Reading Tracking Log  
FAQ Information Sheet  
Replacement Parts Order Form

IF YOU ARE MISSING ANY PARTS, INCLUDING INSERTS OR INSTRUCTION MANUALS, DO NOT RETURN TO PLACE OF PURCHASE. CONTACT CUSTOMER CARE AT 866-326-1313.

---

# INDEX

---

Care & Safety Information .....	4-8
Introduction & Indications for Use .....	9
Blood Pressure Monitor Features .....	10-11
Getting Started .....	12
Battery Installation .....	13
Setting Date & Time .....	14-15
Fitting & Applying Your Cuff .....	16-17
Taking Your Blood Pressure Reading .....	18-19
Interpreting Your Results .....	20-21
Memory Function .....	22-23
Care & Maintenance .....	24-25
Device & Label Symbols .....	26
Display Symbols .....	27
Troubleshooting .....	28
Error Codes .....	29
FCC Statement .....	30
Electromagnetic Compatibility .....	31-34
Specifications .....	35
Warranty .....	36
Instrucciones en Español .....	37-72

---

**Toll-Free Customer Care Help Line: 1-866-326-1313**  
**Monday – Friday 8:30 a.m. – 4:30 p.m. CST**

SmartHeart™  
Manufactured for  
Veridian Healthcare  
1175 Lakeside Drive  
Gurnee, IL 60031

Made in China

#93-1562 10/22  
©2022 Veridian Healthcare, LLC

---

## CARE & SAFETY INFORMATION

---

**NOTE:** Read all instructions carefully before use. The following basic precautions are needed when using an electrical product.

**CAUTION:** Failure to read and observe all precautions could result in personal injury or equipment damage.

Improper care or use of your blood pressure monitor may result in injury, damage to the unit or ineffective treatment. Following these instructions will ensure the blood pressure monitor's efficacy and long life.

### **WARNINGS**

- In order to ensure comparable values, always measure your blood pressure at the same time of day.
- Before every measurement, relax for about five minutes.
- If you want to perform several measurements on the same person, wait five minutes between each measurement.
- Do not take a measurement within 30 minutes of eating, drinking, smoking or exercising.
- Repeat the measurement if you are unsure of the measured value.
- The measured values taken by you are for your information only – they are not substitute for a medical examination. Discuss the measured values with your doctor and never base any medical decisions on them (e.g. medicines and their dosages).
- Using the blood pressure monitor outside your home environment or whilst on the move (e.g. whilst travelling in a car, ambulance or helicopter, or whilst undertaking physical activity such as playing sport) can influence the measurement accuracy and cause incorrect measurements.
- Do not use the blood pressure monitor on newborns or patients with pre-eclampsia. We recommend consulting a doctor before using the blood pressure monitor during pregnancy.
- Cardiovascular diseases may lead to incorrect measurements or have a detrimental effect on measurement accuracy. The same also applies to very low

---

## CARE & SAFETY INFORMATION

---

blood pressure, diabetes, circulatory disorders and arrhythmias as well as chills or shaking.

- This device is not intended for use by people (including children) with restricted physical, sensory or mental skills or a lack of experience and/or a lack of knowledge, unless they are supervised by a person who is responsible for their safety or are instructed by such a person in how to use the device. Supervise children around the device to ensure they do not play with it.
- The blood pressure monitor must not be used in connection with a high-frequency surgical unit.
- Only use the device on people who have the specified upper arm measurement for the device.
- Please note that when inflating, the functions of the limb in question may be impaired.
- During the blood pressure measurement, the blood circulation must not be stopped for an unnecessarily long time. If the device malfunctions, remove the cuff from the arm.
- Avoid any mechanical restriction, compression or bending of the cuff line.
- Do not allow sustained pressure in the cuff or frequent measurements. The resulting restriction of the blood flow may cause injury.
- Make sure that the cuff is not placed on an arm in which the arteries or veins are undergoing medical treatment, e.g. intravascular access or therapy, or an arteriovenous (AV) shunt.
- Do not use the cuff on people who have undergone a mastectomy.
- Do not place the cuff over wounds as this may cause further injury.
- Place the cuff on your upper arm only. Do not place the cuff on other parts of the body.

### CAUTIONS

- You can either use the blood pressure monitor with batteries or with a mains

---

## CARE & SAFETY INFORMATION

---

part. Please note that data transfer and data storage is only possible when your blood pressure monitor is supplied with power. As soon as the batteries are empty or the mains part is disconnected from the power supply, the blood pressure monitor loses the date and time.

- To conserve the batteries, the blood pressure monitor switches off automatically if no buttons are pressed for one minute.
- The device is only intended for the purpose described in these instructions for use. The manufacturer is not liable for damage resulting from improper or careless use.
- The blood pressure monitor is made from precision and electronic components. The accuracy of the measured values and service life of the device depend on its careful handling:
  - Protect the device from impacts, humidity, dirt, marked temperature fluctuations and direct sunlight.
  - Do not drop the device.
  - Do not use the device in the vicinity of strong electromagnetic fields and keep it away from radio systems or mobile telephones.
  - Only use the cuff included with the delivery or original replacement cuffs. Otherwise incorrect measured values will be recorded.
- We recommend that the batteries be removed if the device will not be used for a prolonged period of time.

### **WARNING ON HANDLING BATTERIES**

- If your skin or eyes come into contact with battery fluid, rinse the affected areas with water and seek medical assistance.
-  **Choking hazard!** Small children may swallow and choke on batteries. Store batteries out of the reach of small children.
- Observe the plus (+) and minus (-) polarity signs while inserting the batteries in the battery housing.

---

## CARE & SAFETY INFORMATION

---

- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Protect batteries from excessive heat.
- **⚠ Risk of explosion!** Do not throw batteries into a fire.
- Do not charge or short-circuit batteries.
- If the device will not be used for some time, take the batteries out of the battery compartment.
- Use identical or equivalent battery types only.
- Always replace all batteries at the same time.
- Do not use rechargeable batteries!
- Do not disassemble, open or crush the batteries.

### **i** INSTRUCTIONS FOR REPAIRS AND DISPOSAL

- Batteries do not belong in household waste. Please dispose of empty batteries at the collection points intended for this purpose.
- Do not open the device. Failure to comply will invalidate the warranty.
- Do not repair or adjust the device yourself. Proper operation can no longer be guaranteed in this case.
- Repairs must only be carried out by Customer Services or authorized retailers. Before making a claim, please check the batteries first and replace them if necessary.
- For environmental reasons, do not dispose of the device in the household waste at the end of its useful life. Dispose of the device at a suitable local collection or recycling point in your country. Dispose of the device in accordance with local Legislation on Waste Electrical and Electronic Equipment. If you have any questions, please contact the local authorities responsible for waste disposal.

---

## CARE & SAFETY INFORMATION

---

### **NOTES ON ELECTROMAGNETIC COMPATIBILITY**

- The device is suitable for use in all environments listed in these instructions for use, including domestic environments.
- The use of the device may be limited in the presence of electromagnetic disturbances. This could result in issues such as error messages or the failure of the display/device.
- Avoid using this device directly next to other devices or stacked on top of other devices, as this could lead to faulty operation. If, however, it is necessary to use the device in the manner stated, this device as well as the other devices must be monitored to ensure they are working properly.
- The use of accessories other than those specified or provided by the manufacturer of this device can lead to an increase in electromagnetic emissions or a decrease in the device's electromagnetic immunity; this can result in faulty operation.
- Failure to comply with the above can impair the performance of the device.

---

## ***INTRODUCTION & INDICATIONS FOR USE***

---

It is recommended that you first seek the advice and recommendation of your physician or healthcare professional when using home diagnostic devices, including blood pressure monitors.


Automatic digital blood pressure monitors use the oscillometric method to electronically measure your blood pressure. The monitor detects your blood's movement through the artery in your arm and converts the movements into a digital reading. The oscillometric method does not require a stethoscope, making the monitor ideal for home use.

Blood pressure readings determined with this device are equivalent to measurements obtained by a trained healthcare professional using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers.


Thank you for purchasing an Automatic Digital Blood Pressure Monitor. With proper care and use, your monitor will provide you with many years of reliable readings.

**PLEASE READ THIS INSTRUCTION MANUAL  
COMPLETELY BEFORE OPERATING THIS UNIT.**

# BLOOD PRESSURE MONITOR FEATURES

# BLOOD PRESSURE MONITOR FEATURES

**DETAILS OF THE DIGITAL DISPLAY SCREEN AND INTERPRETATION OF THE SYMBOLS ARE INCLUDED THROUGHOUT THIS MANUAL.**

---

## **GETTING STARTED**

---

### **Before Taking a Measurement**

1. Your readings should only be interpreted by your physician or healthcare professional with access to your individual medical history. Regular use of a home blood pressure monitor will allow you to track and record your readings for discussion with your physician.
2. Conduct your measurement in a quiet place while seated in a relaxed position. Rest for 15 minutes before taking your reading. If an error occurs or you desire to take a second reading, allow 15 minutes between readings for your blood vessels to return to normal.
3. Avoid smoking, eating, taking medication, alcohol, physical activity or any other stressful activity for 30 minutes prior to taking a reading.
4. Always remove any jewelry or constrictive clothing that may interfere with the cuff placement.
5. Keep yourself and the monitor still during measuring; do not talk during the reading.
6. It is recommended that you take your readings at the same time each day to better monitor any indications in your results.
7. Record your daily measurement on the included chart or some other written document to share with your physician.

---

## BATTERY INSTALLATION

---

**This unit requires 4 AA batteries.**


It is necessary to replace the batteries when the “Low Battery” symbol appears on the display or when the display does not turn on after the **POWER** button is pressed.

The battery compartment is located on the underside of this monitor.

The ‘Low Battery’ symbol will appear on the upper left of the display.

1. Depress the clip on the battery compartment and remove the cover.
2. Insert or replace 4 AA batteries into the battery compartment, ensuring to match the indicated polarity symbols. Always use new batteries.
3. Replace the battery cover.
4. Dispose of batteries according to local disposal and recycling regulations.

It is recommended to remove the batteries if the unit will not be used for an extended period of time.

## SETTING DATE/TIME

**IT IS NECESSARY TO SET THE DATE AND TIME FOR THE UNIT EVERY TIME BATTERIES ARE INITIALLY INSTALLED OR REPLACED**

It is essential to set the date and time correctly. Otherwise, you will not be able to save your measured values correctly with a date and time and access them again later

The **M1**, **M2** and **POWER** buttons are located below the display.


1. When the batteries are first installed, the device will automatically prompt you to set the Date/Time feature.\*  
*\*NOTE: If the batteries have already been installed—when the unit is off (a blank display screen), press and hold the **POWER** button for 5-7 seconds to make the unit go into clock setting mode.*
2. **Military/ Standard Clock**- When you first enter into the clock setting mode you have the option to choose between a 24 or 12 hour clock. Press the **M1** button to toggle between the two settings. Once a selection is made press the **POWER** button to advance to the next step.
3. **YEAR**—next the year will flash on the display. Press the **M** button to advance the display to the desired year, press the **POWER** button to confirm your selection and advance to the date settings.

---

## SETTING DATE/TIME

---

- MONTH/DATE**—repeat the same sequence of step 2 to confirm the MONTH and DATE setting.
- HOURS/MINUTES**—repeat the same sequence of step 2 to confirm the HOUR and MINUTE setting.
- After the Date/Time setting function is complete, the monitor is ready for use.


## FITTING & APPLYING YOUR CUFF

**PROPER CUFF FIT AND PLACEMENT IS CRITICAL IN OBTAINING ACCURATE BLOOD PRESSURE MEASUREMENTS. PLEASE READ THESE INSTRUCTIONS CAREFULLY AND CONSIDER CONSULTING WITH YOUR PHYSICIAN FOR FURTHER DIRECTION.**

Before applying your blood pressure cuff, be sure you have selected the appropriate size cuff.

To determine the proper cuff size, measure the circumference of your left upper arm at the point midway between inside bend of your elbow and your shoulder.

**Wide-Range Adult Cuff** fits arm circumference 8.6" – 16.5" (22–42 cm)

If for any reason you are unable to or should not use your left arm, please modify the instructions for cuff application to your right arm. Your physician can help you identify which arm is best for you to take measurements from.

### Assembling Your Cuff

Your cuff may need to be 'threaded' through the metal bar to form the cuff cylinder for your arm.

1. Lay the open cuff hook and loop side down, size/instruction wording side up on a flat surface.
2. Keeping the words on the inside of the cuff, thread the end of the cuff through the metal bar, Fig 1. The end of the cuff should be able to fold over onto itself for securing via the hook and loop closure, Fig 2.


Fig. 1


Fig. 2

# FITTING & APPLYING YOUR CUFF

## Applying Your Cuff

Be seated at a table or desk with your feet flat on the floor. To avoid accidentally pulling on or dropping your monitor, the cuff should not be plugged into the monitor until after the cuff is applied to your arm.

The cuff is printed with a useful guide for applying the cuff; please consult the following steps for additional direction.

1. Remove any constrictive clothing or jewelry that may interfere with cuff placement.
2. Slide the cuff cylinder onto your left arm and align the artery mark over the brachial artery on the inside of your arm, Fig 3.

The tube should follow the inside length of the arm, Fig 4.

3. The bottom edge of the cuff should be positioned approximately 0.5" above the elbow joint, Fig 5.
4. When the cuff is closed around the arm, the triangle guide imprint should be within the OK range bar marking; if the triangle is outside of the range, a different sized cuff is needed.
5. The cuff should fit comfortably, yet snugly around your arm, Fig 6. You should be able to insert one finger easily between your arm and the cuff.
6. During use, avoid compression or squeezing of the air hose/tube, as this may cause inflation error or injury due to excessive cuff pressure.


Fig. 3


Fig. 4


Fig. 5


Fig. 6

## TAKING YOUR BLOOD PRESSURE READING

Please read the preceding portions of this manual prior to taking your first reading.

Reminder, it is important to avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading.

1. Position the monitor on a flat, stable surface with the digital display in view. Be seated with your feet flat on the floor and legs uncrossed.
2. After applying the cuff insert the cuff tubing connector into the port on the left side of your monitor.
3. Rest your elbow on a solid surface with your palm facing upward. Elevate your arm so that the cuff is at the same level as your heart, Fig 1. Relax your left hand.
4. Press and release the **POWER** button to turn the unit on; the unit will briefly flash a self-test, Fig 2.
5. The cuff will automatically inflate; you will see the numbers increase on the display screen and feel the cuff inflate.


Fig. 1


Fig. 2

### PROPER CUFF PLACEMENT ICON

After a few moments, the screen will show the OK Cuff symbol, indicating that the cuff is applied properly. If this icon does not show, the measurement will halt and the error E3 will show on the display. If an error occurs, wait a few minutes, reapply the cuff and attempt a new measurement.

## TAKING YOUR BLOOD PRESSURE READING

**YOU CAN STOP THE INFLATION OR DEFLATION PROCESS AT ANY TIME BY PRESSING THE POWER BUTTON.**

- When the measurement is complete, the cuff will automatically deflate, then your blood pressure measurement and pulse reading results will show on the display screen at the same time, Fig 3.
- The measurement results will default to the memory bank in use. To select a memory bank, while the results show on the display, press **M1** or **M2** to store the results.
- The Hypertension Indicator will indicate your reading range on the left display in comparison to the color bars on the face of the monitor, Fig 4. For more information see Interpreting Your Results.
- If an irregular heartbeat was detected during the reading, the Irregular Heartbeat Detection indicator will appear on the display, Fig 5. For more information see Interpreting Your Results.
- The reading will automatically be stored in memory, up to 60 readings per user bank.
- Select **POWER** to turn the unit off and conserve energy and battery life.
- The unit will automatically shut-off approximately 1 minute after use.
- Disconnect the tubing prior to storage.


Fig. 3


Fig. 5




Fig. 6

# INTERPRETING YOUR RESULTS

## Irregular Heartbeat Indicator

An irregular heartbeat is defined as a heartbeat rhythm that is more than 25% slower or 25% faster than the average measurement detected while the monitor is conducting the measurement. When this device detects an irregular rhythm two or more times during the reading, the IHB symbol will appear on the display.

This monitor will operate if the icon appears, but the results may be compromised, especially if this symbol appears often. Please consult with your physician or trained healthcare professional for further information regarding an irregular heartbeat and if this symbol appears frequently.

# INTERPRETING YOUR RESULTS

## Hypertension Indicator - WHO Classifications

This unit features our unique Hypertension Indicator. The color bars on the left side of the monitor display correspond with an icon on the digital display, indicating where the measurement results fall within the World Health Organization standards.

The World Health Organization has established globally accepted standards for the assessment of high or low blood pressure readings. The below chart should be considered only as a guideline, always consult with your physician or health care professional to interpret your individual results. Never adjust your medication(s) dosage, or severely alter your diet or exercise routine without consulting with your physician.

### Stage 3 Hypertension

Systolic pressure:  $\geq 180$  mmHg

or

Diastolic pressure:  $\geq 110$  mmHg

### Stage 2 Hypertension

Systolic pressure:  $\geq 160 \leq 179$  mmHg

or

Diastolic pressure:  $\geq 100 \leq 109$  mmHg

### Stage 1 Hypertension

Systolic pressure:  $\geq 140 \leq 159$  mmHg

or

Diastolic pressure:  $\geq 90 \leq 99$  mmHg

### High Normal

Systolic pressure:  $\geq 130 \leq 139$  mmHg

or

Diastolic pressure:  $\geq 85 \leq 89$  mmHg

### Normal

Systolic pressure:  $\geq 120 \leq 129$  mmHg

or

Diastolic pressure:  $\geq 80 \leq 84$  mmHg

### Optimal

Systolic pressure:  $< 120$  mmHg

or

Diastolic pressure:  $< 80$  mmHg

# MEMORY FUNCTION

## Recalling Measurements in Memory

This monitor stores and recalls up to 60 readings per user. These measurements can be shared with your physician or trained healthcare professional.

This unit can recall:

- the average of the last 3 measurement
- all morning measurements (5:00 a.m. thru 9 a.m.)
- all evening measurements (6 p.m. thru 8 p.m.)
- all measurements in a given memory bank sequentially

1. When the unit is OFF, press and release the **M1** or **M2** button to select the memory bank for review.
2. The unit will first display the average of the last three measurements, Fig 1.
3. Press the **M1** or **M2** button again to display the average of all AM measurements, then PM measurements.
4. Press **M1** or **M2** again to begin scrolling through all saved measurements. Measurements will appear on the display from most current to oldest; the memory number will appear flashing at the bottom of the screen next to the measurement that was taken.
5. All results for a given measurement will display, including measurement results, pulse rate, Hypertension Indicator and Irregular Heartbeat alert (if applicable). The date/time stamp is at the bottom of the display, Fig 2.


Fig. 1


Fig. 2

---

## MEMORY FUNCTION


---

- When the number of readings exceeds 60, the oldest data will be replaced with the new record.
- Press **POWER** to turn the monitor off at any time during review of the stored measurements.

### Clearing Measurements From Memory

To delete all results in a given bank:

- When the monitor is off, press and hold the **M1** or **M2** button to access the memory bank feature.
- Continue to hold the **M** button; after several seconds, the display will show CL and 00.
- The memory bank is now empty.

---

## CARE & MAINTENANCE

---

Proper care and maintenance of your blood pressure monitor will help ensure long life for the product; improper use of the monitor may void the included warranty.

Always disconnect the cuff and tubing from the monitor prior to cleaning and storage.

### Monitor Care & Maintenance

- Do not use any liquids on the monitor; use a soft, dry cloth to clean the monitor as needed.
- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: -4°F – 122°F (-20°C – 50°C); Humidity: <90%RH (non- condensing).
- This monitor requires 6 hours to cool from the maximum storage temperature between uses until the monitor is ready for its INTENDED USE when the ambient temperature is 68°F(20 °C).
- This monitor requires 6 hours to warm from the minimum storage temperature between uses until the monitor is ready for its INTENDED USE when the ambient temperature is 68°F(20 °C).
- Avoid high temperatures and direct sunlight.
- Always remove the batteries if the unit will not be used for an extended period of time.
- Do not attempt to disassemble your monitor; disassembling your monitor may void the manufacturers warranty.
- Do not attempt to service the monitor while in use or on a patient.
- Do not subject the monitor to strong shocks; take care not to drop the monitor.
- No component can be maintained by user in the monitor. The circuit diagrams, component part lists, descriptions, calibrations instructions, or other information which will assist the user's appropriately qualified technical personnel to repair those parts of equipment which are designated repairable can be supplied.

---

# CARE & MAINTENANCE

---

## Cuff Care & Maintenance

- The cuff may be spot cleaned with a mild-detergent as needed.
- To disinfect the cuff, wipe the inside (skin contact side) of the cuff using a cloth lightly moistened with 70% - 90% ethyl rubbing alcohol and allow to air dry; it is suggested to disinfect the cuff frequently when used in clinical settings or if the cuff is contaminated in any way.
- Never attempt to iron or dry the cuff; if the cuff becomes wet, allow the cuff to air dry completely before use or storage.
- Do not attempt to remove the tube and bladder from the nylon cuff as this may damage the cuff.
- Take care not to forcefully bend or kink the tubing during storage.
- It is recommended the performance should be checked every 2 years or after a repair. Please contact Customer Care at 866-326-1313.
- It is recommended to clean the cuff every 200 uses.
- The monitor can maintain the safety and performance characteristics for minimum of 10,000 measurements of three years.
- The cuff integrity is maintained after 1,000 open- close cycles of the closure.

## DEVICE & LABEL SYMBOLS

These symbols may appear on your device, instructions or packaging and may vary by make and model.

### Symbol | Meaning

	Read This Manual—All included manuals should be read prior to first use. (Background color will be blue; reading symbol will be white)
	Warning—Symbol indicates a warning, prohibition or mandatory action that mitigates a risk that is not necessarily obvious to the device operator.
	Type BF Applied Parts—Indicates that a part of this unit comes in contact with the patient in order to carry out its intended function; in the case of this device the cuff is the Type BF Applied Part.
	Environment Protection—dispose of this product properly; consult with your local recycling ordinances for proper recycling and disposal.
	Keep Dry—This device should be kept dry; never submerge the unit or cuff. Consult with the Care & Maintenance section of this manual for information on cleaning your monitor.
<b>IP21</b>	IP21 —Protection against solid foreign objects and against harmful effects due to ingress of water.
<b>SN</b>	Serial Number.
	Direct Current.
	Manufacture; Date of Manufacture.

The Electronic Sphygmomanometers corresponds to the below standards:  
ISO81060-2 : 2013(Non-Invasive Sphygmomanometers- Part 2: Clinical Validation Of Automated Measurement Type).

## DISPLAY SYMBOLS

Symbol	Meaning
<b>SYS</b>	Systolic Blood Pressure
<b>DIA</b>	Diastolic Blood Pressure
<b>PUL/min</b>	Pulse Rate per Minute
<b>mmHg</b>	Millimeters of Mercury
	Low Battery
	Pulse Rate Results
	Irregular Heartbeat Symbol
	User 1 / User 2
<b>(M)</b>	Memory Storage
<b>(M) A</b>	Last 3-Reading Average
<b>(M) AM</b>	Average of all A.M. (morning) Measurements
<b>(M) PM</b>	Average of all P.M. (evening) Measurements
	Correct Cuff Wrap Icon

# TROUBLESHOOTING

Problem	Possible Cause	Solution
Display shows an abnormal result.	Cuff position was not correct or was not tight enough.	Apply the cuff properly and attempt a new measurement.
	Body posture was not correct during measurement.	Review the 'Taking Your Blood Pressure Reading' portion of this manual and attempt a new measurement.
	Speaking, excessive movement, angry, excited or nervous anxiety during measurement.	Wait a period of time and attempt a new measurement after calm has returned; do not speak or move during measurement.
Display is blank when power is on.	Batteries may be expired or installed improperly.	Check the polarity of the batteries and reinstall if necessary; replace batteries with new batteries.
No response when <b>POWER</b> is pressed or when new batteries have been loaded.	Incorrect operation; possible strong electromagnetic interference.	Remove batteries and wait a minimum of five minutes; reinstall batteries and attempt a new measurement.
	The battery is low.	Replace the batteries.

# ERROR CODES

Code	Possible Causes	Corrective Action
Monitor cannot be switched on	Batteries are missing, inserted incorrectly or flat.	Check batteries and insert four identical new batteries if necessary.
Cuff will not inflate	Cuff connector plug is incorrectly positioned in the connector socket of the monitor.	Check the connection between the cuff connector plug and socket.
E1	Pulse could not be detected correctly.	Check whether the cuff has been applied correctly. Do not talk or move during the measuring procedure.
E2	The pulse could not be detected.	Check whether the cuff has been positioned correctly. Do not talk or move during the measuring procedure.
E3	The cuff was applied too firmly or too loosely.	Apply the cuff so that two fingers fit between the cuff and the upper arm.  Air tube is not correctly connected to the monitor. Check that the connector plug is correctly positioned. If this error occurs frequently, you should use a new cuff.
E4	Errors occur during the measurement; excessive movement	Please repeat the measurement. Make sure that the air tube is properly inserted.
E5	Inflation pressure is higher than 300 mmHg.	Please rest for at least 1 minute and take the measurement again.
E6	There is a system error.	Contact customer service if this error message appears.
	The batteries are almost empty.	Replace batteries.
Implausible measure values	Implausible measured values often occur due to inappropriate handling of the device or mistakes during the measuring procedure.	Please repeat the measurement according to the instructions.

Switch the device off if an error symbol appears. Check for all possible causes and note the instructions for self-measurement. Relax for a minute and then take the measurement again.

In such cases, repeat the measurement. Ensure that you do not move or speak. If necessary, reinsert or replace the batteries.

---

## **FCC STATEMENT**

---

### **NOTE:**

#### POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for U.S.A. only)

This product has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. The product generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If the product does cause harmful interference to radio or television reception, which can be determined by turning the product on and off, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the product and the receiver.
- Connect the product into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

#### POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for Canada only)

This digital apparatus does not exceed the Class B limits for radio noise emissions from digital apparatus as set out in the interference-causing equipment standard entitled "Digital Apparatus", ICES-003 of the Canadian Department of Communications.

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Table 1 EMISSIONS

## For all ME EQUIPMENT and ME SYSTEMS

### Guidance and manufacturer's declaration - electromagnetic emissions

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

<b>Emission test</b>	<b>Compliance</b>	<b>Electromagnetic Environment Guidance</b>
RF emissions CISPR 11	Group 1	The Blood pressure monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Blood pressure monitor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

# ELECTROMAGNETIC COMPATIBILITY

Table 2 — IMMUNITY

## For all ME EQUIPMENT and ME SYSTEMS

### Guidance and manufacturer's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance levels	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	$\pm 8$ kV contact $\pm 15$ kV air	$\pm 8$ kV contact $\pm 15$ kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50Hz) magnetic Field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

# ELECTROMAGNETIC COMPATIBILITY

Table 3

## For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING Guidance and manufacturer's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	NA	<p>Portable and mobile RF communications equipment should be used no closer to any part of this monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance:</b></p> $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3 \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,* should be less than the compliance level in each frequency range.**</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>
Radiated RF IEC 61000-4-3	10 V/m	10 V/m	

**NOTE 1** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

\* Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which this monitor is used exceeds the applicable RF compliance level above, this monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating this monitor.

\*\* Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

# ELECTROMAGNETIC COMPATIBILITY

Table 4

## For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING Recommended separation distances between portable and mobile RF communications equipment and this monitor

This monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this monitor as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m			
	150 kHz to 80 MHz (out ISM and amateur radio bands) $d = 1.2 \sqrt{P}$	150 kHz to 80 MHz (out ISM and amateur radio bands) $d = 2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3 \sqrt{P}$
0,01	0.12	0.2	0.12	0.23
0,1	0.38	0.632	0.38	0.73
1	1.2	2	1.2	2.3
10	3.8	6.32	3.8	7.3
100	12	20	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

---

## PRODUCT SPECIFICATIONS

---

<b>Name</b>	Automatic Digital Blood Pressure Arm Monitor
<b>Model Number</b>	01-563
<b>Display System</b>	Digital display/LCD
<b>Measuring Method</b>	Oscillometric Method, Automatic Air Inflation and Measurement
<b>Power Source</b>	4 x 1.5V  AA Batteries
<b>Measuring Range</b>	Cuff Pressure: 0-300 mmHg, Systolic pressure:50-280mmHg, Diastolic Pressure: 30-200mmHg; Pulse: 40-199 beats/minute
<b>Accuracy</b>	Pressure: +/-5 mmHg at 120-80mmHg; Pulse: ±5%
<b>Deflation</b>	Automatic pressure release valve
<b>Memory</b>	Built-in memory with 2 banks, 60 measurements per bank
<b>Automatic Power-Off</b>	Approximately 1 minute after last button operation
<b>Battery Life</b>	>300 cycles
<b>Operation Environment</b>	Temperature 50°F – 104°F (10°C – 40°C); Humidity <90% RH max
<b>Storage and Transport Environment</b>	Temperature -4°F – 131°F (-20°C – 55°C); Humidity, <90% RH max
<b>Monitor Dimensions</b>	5.28" x 3.94" x 2.36" (134mm x 100mm x 60mm)
<b>Monitor Weight</b>	7.8 oz. (without batteries/cuff); 4.5 oz. (cuff only)
<b>Arm Circumference Range</b>	Wide Range size cuff fits arm circumference: 8.6" – 16.5" (22 – 42 cm)
<b>Accessories</b>	Instruction Manual, Quick Start Guide, Storage Case, Replacement Cuff Order Form, Blood Pressure Log and FAQ
<b>Additional Options</b>	Additional cuffs available on the enclosed order form or by calling Customer Care

*Specifications are subject to change without notice.*

---

## 1-YEAR LIMITED WARRANTY

---

Congratulations on your purchase of a Digital Blood Pressure Monitor. Your Digital Blood Pressure Monitor is covered by the following limited warranty commencing upon the date of purchase, and subject to the following terms and conditions:

The Warrantor warrants that its Digital Blood Pressure Monitor will be free from defects in materials and workmanship under normal consumer usage for a period of one year for the original purchaser.

Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the Digital Blood Pressure Monitor for commercial purposes or subjecting the Digital Blood Pressure Monitor to abnormal usage or conditions; or (d) other acts which are not the fault of the Warrantor, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the Digital Blood Pressure Monitor.

If the Digital Blood Pressure Monitor fails to conform to this limited warranty, return the Digital Blood Pressure Monitor postage prepaid to: **Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031**. When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone number. Carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. At our option, the Warrantor will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will be responsible for parts charges, if any, and repair charges not covered by this limited warranty.

EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL THE WARRANTOR BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

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