
# **TENS Wired Pain Management Solution**

Model 22-040

***Intention for Use:***

***The TENS Wired Pain Management Solution is intended for temporary relief of pain associated with sore and aching muscles due to strain from exercise or normal household and work activities.***

**PLEASE READ THESE INSTRUCTIONS COMPLETELY BEFORE USING YOUR TENS DEVICE**

**INSTRUCTION MANUAL • ENGLISH & ESPAÑOL**

**STOP! PLEASE ENSURE YOU HAVE ALL OF THE FOLLOWING COMPONENTS BEFORE USING YOUR TENS DEVICE**

Control Unit

2 Sets Electrode Pads (4 total)

2 Electrode Wires

Plastic Storage Board

USB Cable

Storage Bag

Instruction Manual  
Replacement  
Electrode Order Form

**IF YOU ARE MISSING ANY PARTS, INCLUDING INSERTS OR INSTRUCTION MANUALS, DO NOT RETURN TO PLACE OF PURCHASE. CONTACT CUSTOMER CARE AT 866-326-1313.**

---

# INDEX

---

Care & Safety Information .....	4-11
Product Features.....	12-13
Charging the Device .....	14
Preparing Your Device .....	15
Placement & Treatment.....	16-17
Pulse Program Descriptions.....	18-19
Instructions for Use .....	20-22
Caring for Electrode Pads .....	23
Cleaning & Storage .....	24
Troubleshooting .....	25
Device & Label Symbols.....	26
FCC Statement.....	27
Safety Test Standards.....	28
Electromagnetic Compatibility .....	29-32
Product Specifications.....	33
Warranty Information .....	34-35
Español.....	37-72

---

**Toll-Free Customer Care Help Line: 1-866-326-1313**  
**Monday – Friday 8:30 a.m. – 4:30 p.m. CST**

Distributed by  
Veridian Healthcare  
1175 Lakeside  
Gurnee, IL 60031  
[www.veridianhealthcare.com](http://www.veridianhealthcare.com)

Made in China

#93-1278 12/20  
©2020 Veridian Healthcare, LLC

---

# IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

---

**NOTE:** Read all instructions and inserts included with this device carefully before use. The following basic precautions are needed when using an electrical product.

**CAUTION:** Failure to read and observe all precautions could result in injury or equipment damage.

Improper care or use of your device may result in injury, damage to the unit or ineffective treatment. Following these instructions will ensure the device's efficacy and long life.

## **CONTRAINDICATIONS :**

Do not use this device if you have a cardiac pacemaker, implanted defibrillator or other implanted metal or electronic device.

Do not use this device to treat undiagnosed chronic pain.

Do not use this device if you have suspected or diagnosed heart problem.

Do not use this device if you have suspected or diagnosed epilepsy.

Do not use this device if you have a tendency to bleed internally following an injury.

Do not use this device if you recently had surgery, or have ever had surgery on your back.

---

# IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

---

Do not use this device if areas of skin lack normal sensations, such as skin that tingles or is numb.

Do not use this device during pregnancy.

## GENERAL WARNINGS PRIOR TO USE

- If you are under the care of a Physician, consult with your Physician before using this device.
- Long-term effects of this device are not known.
- Do not place the device on or close to your heart.
- Take care when applying device close to the neck; do not place this device on the front or sides of the neck. Severe spasm of the muscles may occur and the contractions may be strong enough to close the airway or cause difficulty in breathing. Stimulation over the neck could also affect hearing or blood pressure.
- Do not apply stimulation across the chest because the introduction of electrical current into the chest may cause rhythm disturbances to the heart.
- Do not place the device on or around your head. The effects of stimulation of the brain are unknown.

---

# IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

---

- Do not use this device while sleeping.
- Do not use if you feel numbness.
- Do not use this device in or close to water.
- Use the only on normal, healthy, clean and dry skin. Do not use on or close to open wounds or rashes, or over swollen, red, infected or inflamed skin.
- Do not apply stimulation over, or in proximity to, cancerous lesions.
- Do not use the device on children or incapacitated persons.
- Consult with your physician before using this device, because the device may cause lethal rhythm disturbances to the heart in susceptible individuals.

## **GENERAL CAUTIONS PRIOR TO USE**

- Read this manual before using this device for the first time.
- Keep this manual available whenever you use this device.
- This device is intended for individual person use only.
- This device is not effective for pain associated with Central Pain Syndromes, such as headaches.

---

## **IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS**

---

- This device is for pain caused by muscle soreness, and should be placed only around muscles where pain originates.
- The pain may indicate that you have some other health problem. You should know the reason and source of your pain before using these devices. Do not solely rely on the treatment of this device for pain.
- The safety of using this device during pregnancy or birth has not been established.
- The effectiveness of this device depends greatly on a person's individual physical condition. It may not always be effective for every user.
- If you have had medical or physical treatment for your muscle pain, consult with your treatment provider before using this device. You should contact your physician prior to using this device following recent surgical procedures. Stimulation may disrupt the healing process.

---

# IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

---

## OPERATING WARNINGS & CAUTIONS

CONSULT WITH YOUR PHYSICIAN IF YOU HAVE ANY OF THE FOLLOWING:

- If you have suspected or diagnosed heart problem.
- If you have suspected or diagnosed epilepsy.
- If you have a tendency to bleed internally following an injury.
- If you recently had surgery, or have ever had surgery on your back.
- If areas of skin lack normal sensations, such as skin that tingles or is numb.
- During menstruation or during pregnancy.
- The unit is intended for the temporary relief of pain caused by muscle soreness. The source or cause of the pain should be addressed by a healthcare professional. Pain could indicate an injury or condition that requires treatment.
- Some people may feel skin irritation or experience a very sensitive feeling in the skin due to electrical stimulation. If this occurs, stop using these devices and consult your physician.

---

## **IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS**

---

- If skin under the pad feels irritated after using the stimulator for a long period of time, use the stimulator for a shorter period of time.
- Minor redness at stimulation placement is a normal skin reaction. It is not considered a skin irritation, and it will normally disappear within 30 minutes after the electrodes are removed. If the redness does not disappear after 30 minutes from the removal of electrodes, do not use the stimulator again until after the excessive redness has disappeared.
- Turn off your device if the stimulation feels unpleasant or does not provide pain relief.
- Keep this device out of the reach of children.
- Use this device only with the pads and accessories recommended by the manufacturer.
- Do not use this device when driving, operating machinery or when swimming.
- Before removing the Electrode Pad, be sure to turn it OFF, avoiding unpleasant stimulation.
- TENS is not a substitute for pain medications and other pain management therapies.

---

## **IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS**

---

- TENS devices have no curative value.
- TENS is a symptomatic treatment and, as such, suppresses the sensation of pain that would otherwise serve as a protective mechanism.
- Do not use the unit while sleeping.
- After each use, clean any adhesive residue left on the skin with soap and water.
- The effectiveness of stimulus therapy varies from person to person. It may not be an effective treatment for all users.
- The long-term effects of stimulus therapy are unknown.

## **STORAGE WARNINGS & CAUTIONS**

- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: 14°F – 140°F (-10°C – 60°C); humidity: 30-95% RH.
- Keep the unit out of reach of small children.

---

## **IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS**

---


- Do not store the device in direct sunlight, dusty or humid environments, or extreme temperatures. Exact figures for appropriate use and storage can be found in the Product Specifications section of this manual.

### **CLEANING WARNINGS & CAUTIONS**

- Never immerse the unit in water to clean as it may damage the unit.
- Follow the 'Cleaning and Maintenance' portion of this manual for instruction on how to clean and care for your device.

# PRODUCT FEATURES

## Control Unit


---

# PRODUCT FEATURES

---

## *Electrode Pads and Connection Wires*

---

## CHARGING THE DEVICE

---


**This device includes a rechargeable battery and USB cord.**

Before first use, the battery should be charged for approximately 1 hour. Subsequent charging will also take approximately 1 hour.

To charge the device:

1. Connect the wide end of the USB cable to a wall adapter (not included).
2. Connect the smaller end to the control device.
3. While charging the battery icon  on the upper right of the display will flash.
4. When the device is fully charged the battery icon will remain solid and cease flashing.

**NOTE:** A full charge will power six 60-minute treatment cycles before needing to be recharged.

---



## PREPARING YOUR DEVICE

---

It is necessary to prepare and assemble your device prior to treatment.

**NOTE: Save the plastic film discs for storage after treatment to protect the pads from dirt and debris.**

1. While the device is OFF (the display is NOT illuminated) carefully attach the electrode pads to the control unit via the snap connection points.
2. Pull gently to ensure the electrode pad is firmly attached.
3. Carefully remove the plastic film from the electrode pad.
4. You are ready to apply the pad to the affected treatment area.

**Do not turn the power ON until the device is in place on the user's skin as indicated in the 'Instructions for Use' section of this manual.**

---

# PLACEMENT & TREATMENT

---


## **Read all warnings and cautions completely before use.**

The electrode pads are applied directly to the skin. Clothing may be worn over the device; clothing should not constrict the device or wires and should allow access to the control unit buttons.

Targeted treatment areas may include:

- shoulder
- elbow
- thigh
- knee
- calf
- ankle
- back of neck
- upper back
- lower back

**Never place device on over the heart, on the front or sides of neck or the face.**



---

# PLACEMENT & TREATMENT

---

The Wired TENS should be applied on alternating sides of the area where you are experiencing pain or may be applied directly on the location dependent upon the size of the treatment area.

Two or four pads may be used for treatment. It is necessary to use BOTH pads on a given lead wire (white and red caps) to ensure a circuit is created when the device is activated or the pulse sensation will not be felt.

# PULSE PROGRAM DESCRIPTIONS

	<b>Pulse Rate</b>	<b>Description</b>	<b>Benefit</b>
P1	65	Pulse on for 3.4 second and off for 1.6 second	Mild gate pain relief
P2	13-56	Pulse on for 20 second and off for 1 second	Muscle performance improvement with variable pulse rate
P3	1.2	Pulse on every 850 ms	Slow endorphin pain relief
P4	100	Pulse on for 3.6 second and off for 1 second	Quick gate pain relief and strong pulse amplitude
P5	100	Pulse on for 10 second and off for 2.5 second	Quick gate pain relief
P6	Combina- tion P1-P6	Combination of the above	All kinds of pain relief
P7	20	Pulse on for 8.6 second and off for 2 second	Muscle performance improvement
P8	165	Pulse on for 12 second and off for 2 second	Quick and deep gate pain relief

---

## PULSE PROGRAM DESCRIPTIONS

---

The programs are not specific to treatment at a certain body location or to treat a specific ailment or pain.

The sensations may include—massage, acupuncture, tapping, scraping or a combination of the sensations at variable durations and cycles.

Select the program that provides the best relief for your particular ache and pain.

---

# INSTRUCTIONS FOR USE

---

**Read all warnings and cautions completely before use. Recommended treatment time is 20-minutes; consult your healthcare professional for a treatment cycle that best suits your injury. Discontinue use immediately if treatment is uncomfortable or painful.**

1. Treatment area should be clean before applying the control unit. Note that excessive skin lotion or soap residue may affect the adhesiveness of the pad. Only use on HEALTHY skin areas with no cuts, sores or wounds.
2. Apply the pads to the desired treatment area, taking care to ensure the controller is comfortably within reach of your hand for safety and control during treatment.
3. Ensure the electrode pads are smoothly applied and the gel adhesive is holding the pad firmly in place on the skin.
4. The power switch is located at the top of the control unit; slide the switch to the on position and the digital display will illuminate.

On/Off Switch

---

## INSTRUCTIONS FOR USE

---

5. The TENS unit will begin to operate at the lowest intensity level setting for program 1.
6. Press and release the '+' button to find your desired therapy sensation level. This device offers 18 intensity levels. Press and release the '-' button to reduce the intensity level.
7. To change the mode, while the unit is on/operating, press and release the 'M'; the mode will switch to the next in the sequence of 8 pulse types. The current mode in use will flash on the display. See page 18 for a detailed description of each pulse type. When changing modes, the intensity will revert to level 1.
8. The sensation should be strong but comfortable; press the + or - button at any time during the treatment to change the intensity as desired.

**NOTE: Until you are familiar with use of the unit, begin with the lower intensity levels and gradually progress to more advanced levels as you feel comfortable.**

---

## INSTRUCTIONS FOR USE

---

9. The timer default treatment time is 20-minutes. To adjust the time in 10 minute increments, press and release the 'T' button to scroll to the desired time up to 60 minutes.  
To discontinue use prior to the end of treatment time, press Play/Pause (>||) button or slide the power to the OFF position.
10. Be sure the device is OFF (no sensation is felt) before removing the device from your skin. Take care during removal to best protect the adhesive and preserve the life of the electrode pad.
11. Place the plastic covers over the gel side of the pads to prevent them from drying out.
12. To conserve the battery and prevent shocks, the device will automatically power off when not in use.

**Note:**

You may wish to clean the treatment site after use to remove any adhesive residue from the electrode pad.

---

# CARING FOR THE ELECTRODE PADS

---

## **Gel electrode pads**

- The electrode pads can be reused as long as they retain their stickiness. Cover the pads with the included film pieces when not in use to prevent them from drying out.
- Properly cared for, the electrode pads may last up to 50 uses; the life of the electrode pad is largely dependent upon proper care and the condition of the skin as applied for use.
- To help prolong the life of the electrode pad, between uses, moisten the gel side of the pad with a drop of water and allow to dry; replace the film for storage.
- If you notice that the pads have become dry, no longer adhere well to the skin, or begin to irritate the skin, they are in need of replacement.

Replacement pads are available through Veridian Healthcare at 866-326-1313, by using the enclosed order form or through many major retailers and ecommerce sites.

---

# CLEANING & STORAGE


---

## Cleaning

- Unplug the electrode wires before cleaning the device.
- The Control Unit can be cleaned with a soft, dry cloth.
- Never use cleaning agents or excessive water to clean any part of the unit. Do not submerge any part of the unit or let it get wet.
- Never disassemble the unit or attempt to repair the control unit.

## Storage

- During a period of frequent therapy treatments, use the included plastic storage board piece to place the electrode pads and wrap the wires to avoid tangling.
- If the device will not be used for several days or longer, place the film cover over the electrode pads to prevent drying out or getting dirt and debris on them.
- Store all components in the included bag while not in use.
- Storing the pads in a plastic zip bag will help prevent the gel pads from drying out prematurely.
- Protect the unit from mechanical shock or heavy impact.
- Avoid exposure to extreme temperatures.

---

# TROUBLESHOOTING

---

<b>Problem</b>	<b>Solution</b>
Display does not illuminate/ is weak / Control unit will not turn on.	Recharge the battery.
Intensity output is low / no stimulation is occurring.	Electrode pad may be dirty; clean and attempt again. Confirm that the electrode wires are firmly connected to the control unit. Battery power is low, recharge the battery. Device may not be applied properly; ensure BOTH pads from the same wire are adhered to the body.
Electrode pads are not adher- ing to the skin.	Review the 'Caring for Electrode Pads' section of this manual; pads may need to be replaced.

# DEVICE & LABEL SYMBOLS

These symbols may appear on your device, instructions or packaging and may vary by make and model.

Symbol	Meaning	Symbol	Meaning
	Read This Manual—All included manuals should be read prior to first use.		Keep Device Dry
	Date of Manufacture		Product Packaging is Recyclable
	Batch Code		Type BF Applied Parts—Indicates that a part of this unit comes in contact with the patient in order to carry out its intended function; in the case of this device the cuff is the Type BF Applied Part.
	Serial Number		Environment Protection—dispose of this product properly; consult with your local recycling ordinances for proper recycling and disposal.
	Manufacturer	<b>IP22</b>	IP code of the device
	Fragile, Handle with Care		
	This Way Up		

---

# FCC STATEMENT

---

## NOTE:

### POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for U.S.A. only)

The subject device has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

The product generates, uses, and can radiate radio frequency energy and, if not installed and used accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that the interference will not occur in a particular installation. If the product does cause harmful interference to radio or television reception, which can be determined by turning the product on or off, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the product and the receiver.
- Consult the dealer or an experienced radio/TV technician for help.
- Connect the product into an outlet on a circuit different from that to which the receiver is connected.

Changes or modifications to this product not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

---

# SAFETY TEST STANDARDS

---

Medical Devices Directive 93/42/EEC

IEC 60601-1:2005+A1:2012/EN 60601-1:2006 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance

IEC60601-1-2 2014/EN60601-1-2 2014 Medical electrical equipment - Part 1-2: General requirements for safety - Collateral standard: Electromagnetic compatibility - Requirements and tests

IEC 60601-2-10:2012/EN 60601-2-10:2000+A1:2001 Medical electrical equipment - Part 2-10: Particular requirements for the safety of nerve and muscle stimulators

IEC60601-1-11:2015 Medical electrical equipment -- Part 1-11: General requirements for basic safety and essential performance -- Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home health care environment.

EN ISO 15223-1:2016 Medical devices -- Symbols to be used with medical device labels, labelling and information to be supplied -- Part 1: General requirements

EN 1041 Information supplied by the manufacturer with medical devices

IEC/60601-1-6/ EN 60601-1-6 Medical electrical equipment – Part1-6: General requirements for basic safety and essential performance – Collateral standard: Usability

IEC 60601-1-11/ EN 60601-1-11 Medical electrical equipment – Part 1-11: General requirements for basic safety and essential performance – Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in home health care environment

IEC 62304/ EN 62304 Medical device software - Software life-cycle processes

IEC 62366/ EN 62366 Medical devices - Application of usability engineering to medical devices

ISO 10993-1 Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process

---

# ELECTROMAGNETIC COMPATIBILITY

---

- 1) This product needs special precautions regarding electromagnetic compatibility (EMC) and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile radio frequency (RF) communications equipment.
- 2) Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.
- 3) Caution: This unit has been thoroughly tested and inspected to assure proper performance and operation!
- 4) Caution: This machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

Table 1

## For all ME EQUIPMENT and ME SYSTEMS

### Guidance and manufacturer's declaration - electromagnetic emissions

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Emmission Test	Compliance	Electromagnetic Environment Guidance
RF emissions CISPR 11	Group 1	This device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

# ELECTROMAGNETIC COMPATIBILITY

Table 2

## For all ME EQUIPMENT and ME SYSTEMS

### Guidance and manufacture's declaration - electromagnetic immunity

This device is intended for use in the electromagnetic environment specified below.

The customer or the user of this device should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of this monitor requires continued operation during power mains interruptions, it is recommended that this monitor be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

**NOTE** UT is the a.c. mains voltage prior to application of the test level.

# ELECTROMAGNETIC COMPATIBILITY

Table 3

## For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING Guidance and manufacturer's declaration - electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of this monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. <b>Recommended separation distance:</b> $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2,5 GHz  Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,* should be less than the compliance level in each frequency range.** Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	

**NOTE 1** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

\* Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which this monitor is used exceeds the applicable RF compliance level above, this monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating this device.

\*\* Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

# ELECTROMAGNETIC COMPATIBILITY

Table 4

**For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING**  
**Recommended separation distances between portable and mobile RF communications equipment and this monitor**

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this monitor as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

---

# PRODUCT SPECIFICATIONS

---

<b>NAME</b>	TENS Wired Pain Management Solution
<b>MODEL NUMBER</b>	22-040
<b>CHANNEL</b>	Dual Channel
<b>OUTPUT VOLTAGE</b>	0~42V +/- 20% (at 500Ohm load)
<b>PULSE INTENSITY</b>	0~84mA +/- 20%; 0~18 stages adjustable
<b>PULSE WIDTH</b>	0~100 uS
<b>PULSE RATE</b>	0 ~ 200Hz
<b>AUTOMATIC SHUT-OFF</b>	After completion of set therapy treatment time; 10–60 minutes maximum
<b>POWER SOURCE</b>	Rechargeable Lithium Polymer Battery
<b>OPERATION ENVIRONMENT</b>	Temperature 50 F - 104 F (-10 C - 40C); Humidity 40% - 90%
<b>STORAGE ENVIRONMENT</b>	Temperature 14 F - 140 F (-10 C - 60 C); Humidity 30% - 95%
<b>DEVICE DIMENSIONS</b>	3-1/2" x 1-5/8" x 3/8"
<b>WEIGHT</b>	Controller weight 1.3 oz
<b>INCLUDED</b>	Control unit, electrode wire cable with 2 connectors, 2x electrode pads set, plastic storage board, USB cable, instruction manual, storage bag, replacement electrodes order form
<b>AVAILABLE SEPARATELY</b>	Replacement electrodes 4-pack, 22-046
<b>EDITION NUMBER</b>	Edition V1.0

*Specifications are subject to change without notice*

---

# 1-YEAR LIMITED WARRANTY

---

Congratulations on your purchase of the TENS Wired Pain Management Solution. Your TENS Wired Pain Management Solution is covered by the following limited warranty commencing upon the date of purchase, and subject to the following terms and conditions:

The Warrantor warrants that its TENS Wired Pain Management Solution will be free from defects in materials and workmanship under normal consumer use for a period of one year for the original purchaser of the product.

Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the TENS Wired Pain Management System for commercial purposes or subjecting the TENS Wired Pain Management System to abnormal usage or conditions; or (d) other acts which are not the fault of the Warrantor, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the TENS Wired Pain Management Solution.

If the TENS Wired Pain Management Solution fails to conform to this limited warranty, return the TENS Wired Pain Management Solution postage prepaid to: **Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031.** When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone number. Carefully

---

# 1-YEAR LIMITED WARRANTY

---

package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. At our option, the Warrantor will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will be responsible for parts charges, if any, and repair charges not covered by this limited warranty.

EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL THE WARRANTOR BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

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