SMARTHEART™

Talking Blood Pressure Monitor

Model # 01-521


INSTRUCTION MANUAL ENGLISH, ESPAÑOL AND LE FRANÇAIS

Please read this instruction manual completely before operating this unit.

STOP!
PLEASE ENSURE YOU HAVE ALL OF THE
FOLLOWING COMPONENTS BEFORE USING
YOUR DIGITAL BLOOD PRESSURE MONITOR


Digital Monitor


Wide-Range
Arm Cuff

Storage Case

4 AA Batteries

Instruction Manual
Quick Start Guide
Blood Pressure Reading Tracking Log
FAQ Information Sheet

IF YOU ARE MISSING ANY PARTS, INCLUDING INSERTS OR INSTRUCTION MANUALS, DO NOT RETURN TO PLACE OF PURCHASE. CONTACT CUSTOMER CARE AT 866-326-1313.

INDEX

Care & Safety Information.....	4-5
Introduction & Indications for Use.....	6
Digital Blood Pressure Monitor Features.....	7
Getting Started.....	8
Battery Installation.....	9
Setting Memory Bank / User.....	10
Setting Date and Time.....	10
Setting Language and Volume.....	11
Applying Your Blood Pressure Cuff.....	12-13
Taking Your Blood Pressure Reading.....	14-15
Interpreting your Results.....	16-17
Memory Functions.....	18
Care and Maintenance.....	19
Device and Label Symbols.....	20
Description of Display Symbols.....	21
Troubleshooting.....	22
FCC Statement.....	23
Electromagnetic Compatibility.....	24-27
Specifications.....	28
Warranty.....	29
Instrucciones en Español.....	31-60
Instructions en Français.....	61-92

Toll-Free Customer Care Help Line:
1-866-326-1313
Monday – Friday
8:30 a.m. – 4:30 p.m. CST

SmartHeart™
Manufactured for
Veridian Healthcare
1175 Lakeside Drive
Gurnee, IL 60031

Made in China

#93-1398 08/18
©2018 Veridian Healthcare, LLC

CARE & SAFETY INFORMATION

NOTE: Read all instructions carefully before use. The following basic precautions are needed when using an electrical product.

CAUTION: Failure to read and observe all precautions could result in personal injury or equipment damage.

Improper care or use of your blood pressure monitor may result in injury, damage to the unit or ineffective treatment. Following these instructions will ensure the blood pressure monitor's efficacy and long life.

GENERAL CAUTIONS AND WARNINGS

- Measurement results should be discussed with your physician or healthcare professional; never self-diagnose or attempt treatment as this may be dangerous
- Intended for adult use only; this device is not approved for infant or child use
- Keep out of reach of children
- Do not use the unit if it has any damaged parts (including power cord or plug and cuff), if it has been submersed in water or dropped
- If any abnormality occurs, discontinue use until the unit has been examined and repaired
- Only use included attachments and parts; do not use attachments from other brands or models as these may not be properly calibrated for use with this device and may result in measurement error
- Do not disassemble the unit or attempt to repair it; substitution of a component different from that supplied may result in measurement error and will void manufacturer warranty
- Always follow local regulations for proper disposal of the monitor, cuff and batteries

CARE & SAFETY INFORMATION

OPERATING CAUTIONS AND WARNINGS

- Prolonged overinflation may cause congestion, swelling or bruising in some people
- If you experience any discomfort or an abnormality occurs, stop using the device immediately; to stop the device during operation, press the **START/STOP** button and the cuff will automatically deflate
- Operation outside of stated operating temperature may result in measurement error or device malfunction; operation environment temperature is: 50°F – 104°F (10°C – 40°C); Humidity: 15% – 93% RH max
- Never use this unit while operating a vehicle or in the bathtub or shower

STORAGE CAUTIONS AND WARNINGS

- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: -13°F – 158°F (-25°C – 70°C); Humidity: ≤93%
- Keep the unit out of reach of small children
- Always keep the unit unplugged while not in use; remove the batteries if the unit will not be used for an extended period of time

CLEANING CAUTIONS AND WARNINGS

- Never immerse the unit in water to clean as it may damage the unit
- Follow the 'Cleaning and Maintenance' portion of this manual for instruction on how to clean and care for your monitor

INTRODUCTION & INDICATIONS FOR USE

It is recommended that you first seek the advice and recommendation of your physician or healthcare professional when using home diagnostic devices, including blood pressure monitors.


Automatic digital blood pressure monitors use the oscillometric method to electronically measure your blood pressure. The monitor detects your blood's movement through the artery in your arm and converts the movements into a digital reading. The oscillometric method does not require a stethoscope, making the monitor ideal for home use.

Blood pressure readings determined with this device are equivalent to measurements obtained by a trained healthcare professional using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers.

Thank you for purchasing a SmartHeart Automatic Digital Blood Pressure Monitor. With proper care and use, your monitor will provide you with many years of reliable readings.

**PLEASE READ THIS INSTRUCTION MANUAL
COMPLETELY BEFORE OPERATING THIS UNIT.**

BLOOD PRESSURE MONITOR FEATURES

DETAILS OF THE DIGITAL DISPLAY SCREEN AND INTERPRETATION OF THE SYMBOLS ARE INCLUDED THROUGHOUT THIS MANUAL.

GETTING STARTED

Before Taking a Measurement

1. Your readings should only be interpreted by your physician or healthcare professional with access to your individual medical history. Regular use of a home blood pressure monitor will allow you to track and record your readings for discussion with your physician.
2. Conduct your measurement in a quiet place while seated in a relaxed position. Rest for 15 minutes before taking your reading. If an error occurs or you desire to take a second reading, allow 15 minutes between readings for your blood vessels to return to normal.
3. Avoid smoking, eating, taking medication, alcohol, physical activity or any other stressful activity for 30 minutes prior to taking a reading.
4. Always remove any jewelry or constrictive clothing that may interfere with the cuff placement.
5. Keep yourself and the monitor still during measuring; do not talk during the reading.
6. It is recommended that you take your readings at the same time each day to better monitor any indications in your results.
7. Record your daily measurement on the included chart or some other written document to share with your physician.

BATTERY INSTALLATION

This unit comes complete with 4 AA batteries.

It is necessary to replace the batteries when the Low Battery symbol appears on the display, Fig 1, or when the display does not turn on after the **START/STOP** button is pressed.


Fig. 1

The battery compartment is located on the underside of this monitor.

1. Depress the clip in the battery compartment cover and remove the cover.
2. Insert or replace 4 AA batteries into the battery compartment, ensuring to match the indicated polarity symbols, Fig 2. Always use new batteries.
3. Replace the battery cover.
4. Dispose of batteries according to local disposal and recycling regulations.



Fig. 2

It is recommended to remove the batteries if the unit will not be used for an extended period of time.

SETTING MEMORY BANK / USER

This unit offers 4 memory banks, allowing for tracking of multiple user results in the household with the same monitor.

1. With power off, press **SET** button to enter the system set-up functions.
2. The User / Memory Bank in use will show on the upper right of the display screen.
3. Press the **M** button to scroll through the users until the desired user is shown, 1, 2, 3, or 4. The User selected will default to the selected user every time the device is powered on until a new user is selected.
4. After selection, press **START/STOP** to exit set-up, or press **SET** to confirm and move to the Date and Time setting feature.


SETTING DATE AND TIME

IT IS NECESSARY TO SET THE DATE AND TIME FOR THE UNIT EVERY TIME BATTERIES ARE INITIALLY INSTALLED OR REPLACED.

IF YOU CHOOSE NOT TO SET THE DATE AND TIME, YOU MAY PRESS THE START/STOP BUTTON TO FORWARD THROUGH THE TIME AND DATE FEATURE WITHOUT SETTING A PREFERENCE.

After setting the Memory Bank, the screen will advance to setting the Date and Time.

1. The month will blink on the lower left of the display; press **M** to advance to the desired Month. Press **SET** to confirm.
2. Repeat Step 1 to confirm the Day, Hour and then Minutes.
3. After selection, press **START/STOP** to exit set-up, or press **SET** to confirm and move to the Language Selection feature.

SETTING LANGUAGE AND VOLUME

This unit offers 3 language settings for audio prompts and reading results - English, Spanish, or French - with volume control, or the audio option can be turned completely OFF.

After setting the Date and Time, the screen will advance to setting the Language and Volume.


LANGUAGE:

1. The display will show L1 and announce 'English'.
2. Press the **M** button to scroll through the options until the desired language is shown, **L1** (English), **L2** (Spanish), **L3** (French), or **OFF**.
NOTE: The language selected will default to the selected language every time the device is powered on until a new language is selected.
3. Press **SET** to confirm the selected language. Proceed to setting the volume.

VOLUME:

4. The volume level bars will show on the left side of the display.
5. Press the **M** button to increase the volume; you will hear an audio prompt for each level. Continue to press **M** until the desired volume is reached.
6. After selection, press **SET** to confirm and then **START/STOP** to exit set-up and prepare for your measurement.


FITTING AND APPLYING YOUR CUFF

PROPER CUFF FIT AND PLACEMENT IS CRITICAL IN OBTAINING ACCURATE BLOOD PRESSURE MEASUREMENTS. PLEASE READ THESE INSTRUCTIONS CAREFULLY AND CONSIDER CONSULTING WITH YOUR PHYSICIAN FOR FURTHER DIRECTION.

Before applying your blood pressure cuff, be sure you have selected the appropriate size cuff.

To determine the proper cuff size, measure the circumference of your left upper arm at the point midway between inside bend of your elbow and your shoulder.

Wide-Range Cuff fits arm circumference 8.6" – 16.5" (22–42 cm)

If for any reason you are unable to or should not use your left arm, please modify the instructions for cuff application to your right arm. Your physician can help you identify which arm is best for you to take measurements from.

Assembling Your Cuff

Your cuff may need to be 'threaded' through the metal bar to form the cuff cylinder for your arm.

1. Lay the open cuff hook and loop side down, size/instruction wording side up on a flat surface.
2. Keeping the words on the inside of the cuff, thread the end of the cuff through the metal bar, Fig 1. The end of the cuff should be able to fold over onto itself for securing via the hook and loop closure, Fig 2.


Fig. 1


Fig. 2

FITTING AND APPLYING YOUR CUFF

Applying Your Cuff

Be seated at a table or desk with your feet flat on the floor. To avoid accidentally pulling on or dropping your monitor, the cuff should not be plugged into the monitor until after the cuff is applied to your arm.

1. Remove any constrictive clothing or jewelry that may interfere with cuff placement.
2. Slide the cuff cylinder onto your left arm and align the air hose over the main artery on the inside of your arm, Fig 3. For proper positioning, the tube should follow the inside length of the arm, Fig 4.
3. The bottom edge of the cuff should be positioned approximately one inch above the elbow joint, Fig 5.
4. The INDEX mark arrow should be within the range bar marking as shown, Fig 6; if the INDEX arrow is outside of the range, a different sized cuff is needed.

The cuff should fit comfortably, yet snugly around your arm. You should be able to insert one finger easily between your arm and the cuff.


Fig. 3


Fig. 4


Fig. 5


Fig. 6

TAKING YOUR BLOOD PRESSURE READING

Please read the preceding portions of this manual prior to taking your first reading.

Reminder, it is important to avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading.


Fig. 1

Getting Started

1. Position the monitor on a flat, stable surface with the digital display in view.
2. Insert the cuff tubing connector into the port on the left side of your monitor.
3. Rest your elbow on a solid surface with your palm facing upward. Elevate your arm so that the cuff is at the same level as your heart, Fig 1. Relax your left hand.
4. Press the **START/STOP** button to turn the unit on; the unit will briefly flash a self-test, Fig 2. A voice tone will indicate when the unit is ready for testing.


Fig. 2

NOTE: Unit will not function if residual air from previous testing is present in cuff. The LCD may flash  until pressure is stabilized.

5. The cuff will automatically inflate to approximately 180mmHg; you will see the numbers go up on the display screen and the heart icon may flash, Fig 3.


Fig. 3

TAKING YOUR BLOOD PRESSURE READING

YOU CAN STOP THE INFLATION OR DEFLATION PROCESS AT ANY TIME BY PRESSING THE START/STOP BUTTON.

- When the inflation has reached optimum level, the display will begin to show the decreasing pressure while you feel the pressure of the cuff decrease.
- When the measurement is complete, your blood pressure measurement and pulse will display simultaneously on the screen with a voice broadcast.
- The Hypertension Indicator will indicate your reading range in the left side of the display, Fig 4. See page 16 for a chart interpreting the Hypertension Indicator results.
- If an irregular heartbeat was detected during the reading, the Irregular Heartbeat Detector indicator will appear on the display, Fig 5. See page 15 for more information on the Irregular Heartbeat Detector feature.
- The reading will automatically be stored in memory bank selected prior to the reading, up to 99 readings.
- Select **START/STOP** to turn the unit off and conserve energy and battery life. The unit will automatically shut-off after 3-minutes.
- Disconnect the cuff and tubing prior to storage.


Fig. 4




Fig. 5

INTERPRETING YOUR BLOOD PRESSURE READING

Irregular Heartbeat Detector (♥)

An irregular heartbeat is defined as a heartbeat rhythm that is more than 25% slower or 25% faster than the average measurement detected while the monitor is conducting the measurement. When this device detects an irregular rhythm two or more times during the reading, the IHB symbol will appear on the display screen.

This monitor will operate if the icon appears, but the results may be compromised, especially if this symbol appears often. Please consult with your physician or trained healthcare professional for further information regarding an irregular heartbeat and if this symbol appears frequently.

INTERPRETING YOUR BLOOD PRESSURE READING

Hypertension Indicator

This unit features our unique Hypertension Indicator. The color bars on the left side of the monitor display correspond with an icon on the digital display, indicating where the measurement results fall within the World Health Organization standards.

The World Health Organization has established globally accepted standards for the assessment of high or low blood pressure readings. The below chart should be considered only as a guideline, always consult with your physician or health care professional to interpret your individual results. Never adjust your medication(s) dosage, or severely alter your diet or exercise routine without consulting with your physician.

 **Stage 3 Hypertension**
Systolic pressure: ≥ 180 mmHg
or
Diastolic pressure: ≥ 110 mmHg

 **High Normal**
Systolic pressure: ≥ 130 -139 mmHg
or
Diastolic pressure: ≥ 85 -89 mmHg

 **Stage 2 Hypertension**
Systolic pressure: ≥ 160 -179 mmHg
or
Diastolic pressure: ≥ 100 -109 mmHg

 **Normal**
Systolic pressure: ≥ 120 -129 mmHg
or
Diastolic pressure: ≥ 80 -84 mmHg

 **Stage 1 Hypertension**
Systolic pressure: ≥ 140 -159 mmHg
or
Diastolic pressure: ≥ 90 -99 mmHg

 **Optimal**
Systolic pressure: < 120 mmHg
or
Diastolic pressure: < 80 mmHg

MEMORY FUNCTIONS

Recalling Measurements in Memory

This monitor stores and recalls up to 99 measurements, plus an average of the last three measurements in each memory bank / per user.

Select the user to review; see page 10 for instructions on how to change the User.

All results will be announced aloud in the selected memory bank if the voice feature is On.

1. Press and release the **M** button. The unit will first display the average of the last three measurements, Fig 1.
2. Press the **M** button again to begin scrolling through the saved measurements. Measurements will appear on the display from most current to oldest; the memory number will appear on the lower right corner of the display, Fig 2. Press **M** again to proceed to the result, Fig 3.
3. Press the **START/STOP** button to turn the unit off.

Clearing Measurements from Memory

Memory for a specific bank may be entirely deleted while in Memory Recall mode.

1. Press and hold the **SET** button for approximately 3 seconds to delete all memory records from the selected group with voice broadcast "Memory Clear" and then transfer into testing mode, Fig. 3.
2. Press the **START/STOP** button to turn the unit off.

NOTE: Memory can not be recovered once it has been deleted.


Fig. 1


Fig. 2


Fig. 3

CARE AND MAINTENANCE

Proper care and maintenance of your blood pressure monitor will help ensure long life for the product; improper use of the monitor may void the included warranty.

Always disconnect the cuff and tubing from the monitor prior to cleaning and storage.

Monitor Care and Maintenance

- Do not use any liquids on the monitor; use a soft, dry cloth to clean the monitor as needed
- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: $-13^{\circ}\text{F} - 158^{\circ}\text{F}$ ($-25^{\circ}\text{C} - 70^{\circ}\text{C}$); Humidity: $\leq 93\%$
- Always keep the unit unplugged while not in use; remove the batteries if the unit will not be used for an extended period of time
- Do not attempt to disassemble your monitor; disassembling your monitor may void the manufacturers warranty
- Do not subject the monitor to strong shocks; take care not to drop the monitor

Cuff Care and Maintenance

- The cuff should not be submerged or washed in a washing machine, but may be spot cleaned with a mild-detergent as needed
- Never attempt to iron or dry the cuff; if the cuff becomes wet, allow the cuff to air dry completely before use or storage
- Do not attempt to remove the tube and bladder from the nylon cuff as this may damage the cuff
- Take care not to forcefully bend or kink the tubing during storage

DEVICE & LABEL SYMBOLS

These symbols may appear on your device, instructions or packaging and may vary by make and model.

Symbol	Meaning
	Read This Manual—All included manuals should be read prior to first use
	Warning—Symbol indicates a warning, prohibition or mandatory action that mitigates a risk that is not necessarily obvious to the device operator
	Type BF Applied Parts—Indicates that a part of this unit comes in contact with the patient in order to carry out its intended function
	Environment Protection—dispose of this product properly; consult with your local recycling ordinances for proper recycling and disposal
	Serial Number

Symbol	Meaning
	Keep Dry—This device should be kept dry; never submerge the unit or cuff. Consult with the Care and Maintenance section of this manual for information on cleaning your monitor
	Symbol for “MANUFACTURE DATE”
	Symbol for “MANUFACTURER”
	Keep off Sunlight
	Prohibited
	The product conforms to the requirements of the EC Directive MDD (93/42/EEC) on medical devices.

DESCRIPTION OF DISPLAY SYMBOLS

Symbol	Meaning
SYS	Systolic Blood Pressure—The 'top' blood pressure result
DIA	Diastolic Blood Pressure—The 'bottom' blood pressure result
PULSE	Pulse Rate—Pulse rate per minute during the measurement
mmHg	Millimeters of Mercury—Unit of blood pressure measurement for this device
	Low Battery—Battery life indicator
	Pulse Reading—Appears when pulse reading is displayed.
	Irregular Heartbeat Symbol—An irregular heartbeat was detected during the measurement
AVG	Average—The display shows the average of the last three measurements in memory
	User / Memory Bank —Indicates which memory bank is in use for reading storage and memory recall review

TROUBLESHOOTING

Problem	Possible Cause	Solution
Blood pressure results are not within typical range	Cuff is too tight or not properly positioned on the arm.	Firmly reposition cuff approximately 1–2cm (1/2") above the elbow joint.
	Inaccurate test results due to body movement or monitor movement.	Sit in a relaxed position with arm placed near heart. Avoid speaking or moving body parts while testing. Make sure the monitor unit is placed in a stationary position throughout the testing period.
"Err" displayed	Cuff fails to inflate properly.	Make sure hose is properly fastened to cuff and monitor unit.
	Improper operation.	Read user manual carefully and re-test properly.
	Pressurization is over cuff rated pressure 300mmHg	Read user manual carefully and re-test properly.

FCC STATEMENT

NOTE:

POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for U.S.A. only)

This product has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. The product generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If the product does cause harmful interference to radio or television reception, which can be determined by turning the product on and off, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the product and the receiver.
- Connect the product into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for Canada only)

This digital apparatus does not exceed the Class B limits for radio noise emissions from digital apparatus as set out in the interference-causing equipment standard entitled "Digital Apparatus", ICES-003 of the Canadian Department of Communications.

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Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

ELECTROMAGNETIC COMPATIBILITY

Table 1 EMISSIONS

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration - electromagnetic emissions

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Phenomenon	Compliance	Electromagnetic Environment Guidance
RF emissions	CISPR 11 Group 1, Class B	Home healthcare environment
Harmonic distortion	IEC 61000-3-2 Class A	Home healthcare environment
Voltage fluctuations and flicker	IEC 61000-3-3 Compliance	Home healthcare environment

ELECTROMAGNETIC COMPATIBILITY

Table 2 — ENCLOSURE PORT
For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Phenomenon	Basic EMC Standard	Immunity test levels
		Home healthcare environment
Electrostatic Discharge	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Radiated RF EM field	IEC 61000-4-3	10V/m 80MHz-2.7GHz 80% AM at 1kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	Refer to table 3
Rated power frequency magnetic fields	IEC 61000-4-8	30A/m 50Hz or 60Hz

ELECTROMAGNETIC COMPATIBILITY

Table 3 — PROXIMITY FIELDS FROM RF WIRELESS COMMUNICATIONS EQUIPMENT

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING Guidance and manufacturer's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Test frequency (MHz)	Band (MHz)	Immunity test levels
		Professional healthcare facility environment
385	380-390	Pulse modulation 18Hz, 27V/m
450	430-470	FM, ± 5 kHz deviation, 1kHz sine, 28V/m
710	704-787	Pulse modulation 217Hz, 9V/m
745		
780		
810		
870	800-960	Pulse modulation 18Hz, 28V/m
930		
1720		
1845		
1970	1700-1990	Pulse modulation 217Hz, 28V/m
2450		
5240		
5500	2400-2570	Pulse modulation 217Hz, 28V/m
5785		
	5100-5800	Pulse modulation 217Hz, 9V/m

ELECTROMAGNETIC COMPATIBILITY

Table 4

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and this monitor

This monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this monitor as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \sqrt[3]{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt[3]{P}$	800 MHz to 2,5 GHz $d = 2.3 \sqrt[3]{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

PRODUCT SPECIFICATIONS

Name	SmartHeart Premium Talking Digital Blood Pressure Arm Monitor
Model Number	01-521
Display System	Digital display/LCD
Measuring Method	Oscillometric Method, Automatic Air Inflation and Measurement
Power Source	4  AA Batteries
Measuring Range	Pressure: 0-300 mmHg; Pulse: 30-180 beats/minute
Accuracy	Pressure: ± 3 mmHg (0.4kPa); Pulse: $\pm 5\%$ of reading
Deflation	Automatic pressure release valve
Memory	Built-in memory enabling display of up to 396 measurements 4 Groups with 99/user
Automatic Power-Off	Approximately 3 minute after last button operation
Battery Life	Approximately 2 months with 3 tests per day
Operation Environment	Temperature 50°F – 104°F (10°C – 40°C); Humidity 15 – 93% RH max
Storage Environment	Temperature -13°F – 158°F (-25°C – 70°C); Humidity $\leq 93\%$ RH max
Monitor Dimensions	5.83"L x 3.94"W x 2.21"H
Monitor Weight	13.47 oz. (without batteries)
Arm Circumference Range	Wide-Range size cuff fits arm circumference: 8.6" – 16.5" (22 – 42 cm)
Accessories	Storage pouch, Instruction Manual, four AA alkaline batteries, Quick Start Guide, Blood Pressure Log

Specifications are subject to change without notice

FIVE-YEAR LIMITED WARRANTY

Congratulations on your purchase of a Digital Blood Pressure Monitor. Your Digital Blood Pressure Monitor is covered by the following limited warranty commencing upon the date of purchase, and subject to the following terms and conditions:

The Warrantor warrants that its Digital Blood Pressure Monitor will be free from defects in materials and workmanship under normal consumer usage for a period of five years for the original purchaser of the product.

Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the Digital Blood Pressure Monitor for commercial purposes or subjecting the Digital Blood Pressure Monitor to abnormal usage or conditions; or (d) other acts which are not the fault of the Warrantor, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the Digital Blood Pressure Monitor.

If the Digital Blood Pressure Monitor fails to conform to this limited warranty, return the Digital Blood Pressure Monitor postage prepaid to: **Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031**. When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone number. Carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. At our option, the Warrantor will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will be responsible for parts charges, if any, and repair charges not covered by this limited warranty.

EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL THE WARRANTOR BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

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