# SMARTHEART™

---

## Premium Talking Blood Pressure Wrist Monitor

Model # 01-527


### **INSTRUCTION MANUAL** **ENGLISH AND SPANISH**

Please read this instruction manual  
completely before operating this unit.

**STOP!**  
**PLEASE ENSURE YOU HAVE ALL OF THE**  
**FOLLOWING COMPONENTS BEFORE USING**  
**YOUR DIGITAL BLOOD PRESSURE MONITOR**


Digital Monitor

2 AAA Batteries

Storage Case

Instruction Manual  
Quick Start Guide  
Blood Pressure Reading Tracking Log  
FAQ Information Sheet

**IF YOU ARE MISSING ANY PARTS, INCLUDING  
INSERTS OR INSTRUCTION MANUALS, DO NOT  
RETURN TO PLACE OF PURCHASE. CONTACT  
CUSTOMER CARE AT 866-326-1313.**

---

# INDEX

---

Care & Safety Information .....	4-5
Introduction & Indications for Use .....	6
Digital Blood Pressure Monitor Features .....	7
Getting Started .....	8
Battery Installation .....	9
Setting Memory Bank / User .....	10
Setting Date/Time .....	10
Setting Language & Volume .....	11
Applying Your Blood Pressure Cuff .....	12-13
Taking Your Blood Pressure Reading .....	14-15
Interpreting Your Results .....	16-17
Memory Functions .....	18
Care & Maintenance .....	19
Device & Label Symbols .....	20
Description of Display Symbols .....	21
Troubleshooting .....	22
FCC Statement .....	23
Electromagnetic Compatibility .....	24-27
Specifications .....	28
Warranty .....	29
Instrucciones en Español .....	31-60
Instructions en Français .....	61-92

---

**Toll-Free Customer Care Help Line:**  
**1-866-326-1313**  
**Monday – Friday**  
**8:30 a.m. – 4:30 p.m. CST**

SmartHeart™  
Manufactured for  
Veridian Healthcare  
1175 Lakeside Drive  
Gurnee, IL 60031

Made in China

#93-1402 11/19  
©2019 Veridian Healthcare, LLC

---

## CARE & SAFETY INFORMATION

---

**NOTE:** Read all instructions carefully before use. The following basic precautions are needed when using an electrical product.

**CAUTION:** Failure to read and observe all precautions could result in personal injury or equipment damage.

Improper care or use of your blood pressure monitor may result in injury, damage to the unit or ineffective treatment. Following these instructions will ensure the blood pressure monitor's efficacy and long life.

### GENERAL CAUTIONS AND WARNINGS

- Measurement results should be discussed with your physician or healthcare professional; never self-diagnose or attempt treatment as this may be dangerous.
- Intended for adult use only; this device is not approved for infant or child.
- Keep out of reach of children.
- Do not use the unit if it has any damaged parts (including power cord or plug and cuff), if it has been submersed in water or dropped.
- If any abnormality occurs, discontinue use until the unit has been examined and repaired.
- Only use included attachments and parts; do not use attachments from other brands or models as these may not be properly calibrated for use with this device and may result in measurement error.
- Do not share the cuff with an infected person.
- Do not disassemble the unit or attempt to repair it; substitution of a component different from that supplied may result in measurement error and will void manufacturer warranty.
- Always follow local regulations for proper disposal of the monitor, cuff and batteries.

### OPERATING CAUTIONS AND WARNINGS

- Prolonged over-inflation may cause congestion, swelling or bruising in some people.

---

## CARE & SAFETY INFORMATION

---

- Wait a minimum of 2 minutes between measurements (20-minutes is recommended); excessive frequent measurements may restrict blood circulation.
- Do not use this device while connected to any monitoring equipment on the same limb.
- Consult with your healthcare professional for guidance for use if you have any of the following: intravascular shunt; mastectomy on the arm corresponding to the measurement side of the body.
- Consult with your healthcare professional before using the device in pregnancy including pre-eclampsia, or if diagnosed with arrhythmia or arteriosclerosis.
- If you experience any discomfort or an abnormality occurs, stop using the device immediately; to stop the device during operation, press the **START STOP** button and the cuff will automatically deflate.
- Operation outside of stated operating temperature may result in measurement error or device malfunction; operation environment temperature is: 50°F – 104°F (10°C – 40°C); Humidity: 15% – 93%RH (non-condensing).
- Never use this unit while operating a vehicle or in the bathtub or shower.

### STORAGE CAUTIONS AND WARNINGS

- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: -13°F – 158°F (-25°C – 70°C); Humidity: ≤93%RH (non-condensing).
- Keep the unit out of reach of small children.
- Always keep the unit unplugged while not in use; remove the batteries if the unit will not be used for an extended period of time.

### CLEANING CAUTIONS AND WARNINGS

- Never immerse the unit in water to clean as it may damage the unit.
- Follow the 'Cleaning and Maintenance' portion of this manual for instruction on how to clean and care for your monitor.

---

## **INTRODUCTION & INDICATIONS FOR USE**

---

It is recommended that you first seek the advice and recommendation of your physician or healthcare professional when using home diagnostic devices, including blood pressure monitors.

The Blood Pressure Monitor is digital monitors intended for use in measuring blood pressure and heartbeat rate with wrist circumference ranging from 5.3" to 8.5" inches (about 13.5 – 21.5 cm). It is intended for adult indoor use only.

Automatic digital blood pressure monitors use the oscillometric method to electronically measure your blood pressure. The monitor detects your blood's movement through the artery in your arm and converts the movements into a digital reading. The oscillometric method does not require a stethoscope, making the monitor ideal for home use.

Blood pressure readings determined with this device are equivalent to measurements obtained by a trained healthcare professional using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers.

Thank you for purchasing an Automatic Digital Blood Pressure Monitor. With proper care and use, your monitor will provide you with many years of reliable readings.

**PLEASE READ THIS INSTRUCTION MANUAL  
COMPLETELY BEFORE OPERATING THIS UNIT.**

# BLOOD PRESSURE MONITOR FEATURES

Digital Display

**DETAILS OF THE DIGITAL DISPLAY SCREEN AND INTERPRETATION OF THE SYMBOLS ARE INCLUDED THROUGHOUT THIS MANUAL.**

---

## GETTING STARTED

---

### Before taking a measurement

1. Your readings should only be interpreted by your physician or healthcare professional with access to your individual medical history. Regular use of a home blood pressure monitor will allow you to track and record your readings for discussion with your physician.
2. Conduct your measurement in a quiet place while seated in a relaxed position. Rest for 15 minutes before taking your reading. If an error occurs or you desire to take a second reading, allow 15 minutes between readings for your blood vessels to return to normal.
3. Avoid smoking, eating, taking medication, alcohol, physical activity or any other stressful activity for 30 minutes prior to taking a reading.
4. Always remove any jewelry or constrictive clothing that may interfere with the cuff placement.
5. Keep yourself and the monitor still during measuring; do not talk during the reading.
6. It is recommended that you take your readings at the same time each day to better monitor any indications in your results.
7. Record your daily measurement on the included chart or some other written document to share with your physician.

---

## BATTERY INSTALLATION

---


### This unit comes complete with 2 AAA batteries.

It is necessary to replace the batteries when the Low Battery symbol appears on the display or when the display does not turn on after the **START STOP** button is pressed.

The unit will broadcast "Low Battery" when battery life is depleting and unable to inflate cuff.

The battery compartment is located on the bottom of this monitor.

1. The cover is located on the bottom of the monitor, under the START/STOP button. Carefully release the cover to open the battery compartment.
2. Insert or replace 2 AAA batteries into the battery compartment, ensuring to match the indicated polarity symbols, Fig 2. Always use new batteries.
3. Replace the battery cover.
4. Dispose of batteries according to local disposal and recycling regulations.

The 'Low Battery' symbol will appear in the upper left of the display.

It is recommended to remove the batteries if the unit will not be used for an extended period of time.

---

## SETTING MEMORY BANK / USER

---

This unit offers 4 memory banks, allowing for tracking of multiple user results in the household with the same monitor.

1. With power off, press **SET** button to enter the system set-up functions.
2. The User / Memory Bank in use will show on the upper right of the display screen.
3. Press the **M** button to scroll through the users until the desired user is shown, 1, 2, 3, or 4. The User selected will default to the selected user every time the device is powered on until a new user is selected.
4. After selection, press **START/STOP** to exit set-up, or press **SET** to confirm and move to the Date and Time setting feature.

---

## SETTING DATE AND TIME

---

**IT IS NECESSARY TO SET THE DATE AND TIME FOR THE UNIT EVERY TIME BATTERIES ARE INITIALLY INSTALLED OR REPLACED.**

**IF YOU CHOOSE NOT TO SET THE DATE AND TIME, YOU MAY PRESS THE START/STOP BUTTON TO FORWARD THROUGH THE TIME AND DATE FEATURE WITHOUT SETTING A PREFERENCE.**

After setting the Memory Bank, the screen will advance to setting the Date and Time.

1. The month will blink on the lower left of the display; press **M** to advance to the desired Month. Press **SET** to confirm.
2. Repeat Step 1 to confirm the Day, Hour and then Minutes.
3. After selection, press **START/STOP** to exit set-up, or press **SET** to confirm and move to the Language Selection feature.

---

## SETTING LANGUAGE & VOLUME

---

This unit offers 3 language settings for audio prompts and reading results - English, Spanish, or French - with volume control, or the audio option can be turned completely OFF.

After setting the Date and Time, the screen will advance to setting the Language and Volume.


### LANGUAGE:

1. The display will show L1 and announce 'English'.
2. Press the **M** button to scroll through the options until the desired language is shown, **L1** (English), **L2** (Spanish), **L3** (French), or **OFF**.  
NOTE: The language selected will default to the selected language every time the device is powered on until a new language is selected.
3. Press **SET** to confirm the selected language. Proceed to setting the volume.

### VOLUME:

4. The volume level bars will show on the left side of the display.
5. Press the **M** button to increase the volume; you will hear an audio prompt for each level. Continue to press **M** until the desired volume is reached.
6. After selection, press **SET** to confirm and then **START/STOP** to exit set-up and prepare for your measurement.

---

## FITTING & APPLYING YOUR CUFF

---

**PROPER CUFF FIT AND PLACEMENT IS CRITICAL IN OBTAINING ACCURATE BLOOD PRESSURE MEASUREMENTS. PLEASE READ THESE INSTRUCTIONS CAREFULLY AND CONSIDER CONSULTING WITH YOUR PHYSICIAN FOR FURTHER DIRECTION.**


Before applying your blood pressure cuff, be sure your wrist is within the appropriate cuff range: Wrist Circumference Range Suitable for 5.3" – 8.5" (13.5 – 21.5 cm).

Reminder, it is important to avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading.

If for any reason you are unable to or should not use your left wrist, please modify the instructions for cuff application to your right wrist. Your physician can help you identify which wrist is best for you to take measurements from.

### Applying Your Cuff

1. Remove any constrictive clothing or jewelry that may interfere with cuff placement.
2. Be seated at a table or desk with your feet flat on the floor.

## FITTING & APPLYING YOUR CUFF

3. Hold your left arm in front of you with your palm facing upward.
4. Apply the preformed cuff to your left wrist so that the digital display face is positioned on the inside area of your wrist facing you, Fig. 1.
5. Adjust the cuff 1/4" - 1/2" from the edge of the head of the ulna bone, Fig. 2.
6. The cuff should fit comfortably, yet snugly around your wrist.

If you are not comfortable with applying your cuff, please seek the assistance of another member of your household or work with your physician to practice the cuff application. Incorrectly applied cuffs may result in inaccurate readings.


Fig. 1



Fig. 2

# TAKING YOUR BLOOD PRESSURE READING

Please read the preceding portions of this manual prior to taking your first reading.

Reminder, it is important to avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading.

## Getting Started

1. Position the monitor on a flat, stable surface with the digital display in view.
2. Rest your elbow on a solid surface with your palm facing upward. Elevate your arm so that the cuff is at the same level as your heart, Fig 1. Relax your left hand.
3. Press the **START/STOP** button to turn the unit on; the unit will briefly flash a self-test, Fig 2. A voice tone will indicate when the unit is ready for testing.


Fig. 2

**NOTE:** Unit will not function if residual air from previous testing is present in cuff. The LCD may flash  until pressure is stabilized.

5. The cuff will automatically inflate to approximately 180mmHg; you will see the numbers go up on the display screen and the heart icon may flash, Fig 3.

Fig. 3

# TAKING YOUR BLOOD PRESSURE READING

**YOU CAN STOP THE INFLATION OR DEFLATION PROCESS AT ANY TIME BY PRESSING THE START/STOP BUTTON.**

- When the inflation has reached optimum level, the display will begin to show the decreasing pressure while you feel the pressure of the cuff decrease.
- When the measurement is complete, your blood pressure measurement and pulse will display simultaneously on the screen with a voice broadcast.
- The Hypertension Indicator will indicate your reading range in the left side of the display, Fig 4. See page 16 for a chart interpreting the Hypertension Indicator results.
- If an irregular heartbeat was detected during the reading, the Irregular Heartbeat Detector indicator will appear on the display, Fig 5. See page 15 for more information on the Irregular Heartbeat Detector feature.
- The reading will automatically be stored in memory bank selected prior to the reading, up to 99 readings.
- Select **START/STOP** to turn the unit off and conserve energy and battery life. The unit will automatically shut-off after 3-minutes.
- Disconnect the cuff and tubing prior to storage.


Fig. 4




Fig. 5

# INTERPRETING YOUR BLOOD PRESSURE READING

## Irregular Heartbeat Detector **IHB**

An irregular heartbeat is defined as a heartbeat rhythm that is more than 25% slower or 25% faster than the average measurement detected while the monitor is conducting the measurement. When this device detects an irregular rhythm two or more times during the reading, the IHB symbol will appear on the display screen.


This monitor will operate if the icon appears, but the results may be compromised, especially if this symbol appears often. Please consult with your physician or trained healthcare professional for further information regarding an irregular heartbeat and if this symbol appears frequently.

# INTERPRETING YOUR BLOOD PRESSURE READING

## Hypertension Indicator

This unit features our unique Hypertension Indicator. The color bars on the left side of the monitor display correspond with an icon on the digital display, indicating where the measurement results fall within the World Health Organization standards.

The World Health Organization has established globally accepted standards for the assessment of high or low blood pressure readings. The below chart should be considered only as a guideline, always consult with your physician or health care professional to interpret your individual results. Never adjust your medication(s) dosage, or severely alter your diet or exercise routine without consulting with your physician.

### Stage 3 Hypertension

Systolic pressure:  $\geq 180$  mmHg  
or  
Diastolic pressure:  $\geq 110$  mmHg

### Stage 2 Hypertension

Systolic pressure:  $\geq 160 \leq 179$  mmHg  
or  
Diastolic pressure:  $\geq 100 \leq 109$  mmHg

### Stage 1 Hypertension

Systolic pressure:  $\geq 140 \leq 159$  mmHg  
or  
Diastolic pressure:  $\geq 90 \leq 99$  mmHg

### High Normal

Systolic pressure:  $\geq 130 \leq 139$  mmHg  
or  
Diastolic pressure:  $\geq 85 \leq 89$  mmHg

### Normal

Systolic pressure:  $\geq 120 \leq 129$  mmHg  
or  
Diastolic pressure:  $\geq 80 \leq 84$  mmHg

### Optimal

Systolic pressure:  $< 120$  mmHg  
or  
Diastolic pressure:  $< 80$  mmHg

# MEMORY FUNCTIONS

## Recalling Measurements in Memory

This monitor stores and recalls up to 99 measurements, plus an average of the last three measurements in each memory bank / per user.

Select the user to review; see page 10 for instructions on how to change the User.

All results will be announced aloud in the selected memory bank if the voice feature is On.

1. Press and release the **M** button. The unit will first display the average of the last three measurements, Fig 1.
2. Press the **M** button again to begin scrolling through the saved measurements. Measurements will appear on the display from most current to oldest; the memory number will appear on the lower right corner of the display, Fig 2. Press M again to proceed to the result, Fig 3.
3. Press the **START/STOP** button to turn the unit off.

Fig. 1


Fig. 2

## Clearing Measurements from Memory

Memory for a specific bank may be entirely deleted while in Memory Recall mode.

1. Press and hold the **SET** button for approximately 3 seconds to delete all memory records from the selected group with voice broadcast "Memory Clear" and then transfer into testing mode, Fig. 3.
2. Press the **START/STOP** button to turn the unit off.


Fig. 3

**NOTE:** Memory can not be recovered once it has been deleted.

---

## CARE & MAINTENANCE

---

Proper care and maintenance of your blood pressure monitor will help ensure long life for the product; improper use of the monitor may void the included warranty.

### Monitor Care & Maintenance

- Do not use any liquids on the monitor; use a soft, dry cloth to clean the monitor as needed.
- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: -13°F – 158°F (-25°C – 70°C); Humidity: ≤93%.
- Always keep the unit unplugged while not in use; remove the batteries if the unit will not be used for an extended period of time.
- Do not attempt to disassemble your monitor; disassembling your monitor may void the manufacturers warranty.
- Do not subject the monitor to strong shocks; take care not to drop the monitor.

### Cuff Care & Maintenance

- The cuff should not be submerged or washed in a washing machine, but may be spot cleaned with a mild-detergent as needed.
- Never attempt to iron or dry the cuff; if the cuff becomes wet, allow the cuff to air dry completely before use or storage.
- Take care not to forcefully bend or kink the tubing during storage.

## DEVICE & LABEL SYMBOLS

These symbols may appear on your device, instructions or packaging and may vary by make and model.

Symbol	Meaning	Symbol	Meaning
	Read This Manual—All included manuals should be read prior to first use		Keep Dry—This device should be kept dry; never submerge the unit or cuff. Consult with the Care and Maintenance section of this manual for information on cleaning your monitor
	Warning—Symbol indicates a warning, prohibition or mandatory action that mitigates a risk that is not necessarily obvious to the device operator		Symbol for "MANUFACTURE DATE"
	Type BF Applied Parts—Indicates that a part of this unit comes in contact with the patient in order to carry out it's intended function		Symbol for "MANUFACTURER"
	Environment Protection—dispose of this product properly; consult with your local recycling ordinances for proper recycling and disposal		Keep off Sunlight
			Prohibited
	Serial Number		The product conforms to the requirements of the EC Directive MDD (93/42/EEC) on medical devices.

## DESCRIPTION OF DISPLAY SYMBOLS

Symbol	Meaning
<b>SYS</b>	Systolic Blood Pressure—The 'top' blood pressure result
<b>DIA</b>	Diastolic Blood Pressure—The 'bottom' blood pressure result
<b>PULSE</b>	Pulse Rate—Pulse rate per minute during the measurement
<b>mmHg</b>	Millimeters of Mercury—Unit of blood pressure measurement for this device
	Low Battery—Battery life indicator
	Pulse Reading—Appears when pulse reading is displayed.
<b>IHB</b>	Irregular Heartbeat Symbol—An irregular heartbeat was detected during the measurement
<b>Ave</b>	Average—The display shows the average of the last three measurements in memory
	User / Memory Bank —Indicates which memory bank is in use for reading storage and memory recall review

## TROUBLE SHOOTING

<b>Problem</b>	<b>Possible Cause</b>	<b>Solution</b>
Blood pressure results are not within typical range.	Cuff is too tight or not properly positioned on the wrist.	Adjust the cuff 1/4" - 1/2" from the edge of the head of the ulna bone.
	Inaccurate test results due to body movement or monitor movement.	Sit in a relaxed position with arm placed near heart. Avoid speaking or moving body parts while testing. Make sure the monitor unit is placed in a stationary position throughout the testing period.
"Err" displayed.	Cuff fails to inflate properly.	Make sure hose is properly fastened to cuff and monitor unit.
	Improper operation.	Read user manual carefully and re-test properly.
	Pressurization is over cuff rated pressure 300mmHg.	Read user manual carefully and re-test properly.

---

## **FCC STATEMENT**

---

### **NOTE:**

#### POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for U.S.A. only)

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

#### POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for Canada only)

This digital apparatus does not exceed the Class B limits for radio noise emissions from digital apparatus as set out in the interference-causing equipment standard entitled "Digital Apparatus", ICES-003 of the Canadian Department of Communications.

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Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

# ELECTROMAGNETIC COMPATIBILITY

Table 1

## For all ME EQUIPMENT and ME SYSTEMS

### Guidance and manufacture's declaration - electromagnetic emissions

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

<b>Emission Test</b>	<b>Compliance</b>	<b>Electromagnetic Environment Guidance</b>
RF emissions CISPR 11	Group 1	This monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	This monitor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

# ELECTROMAGNETIC COMPATIBILITY

Table 2

## For all ME EQUIPMENT and ME SYSTEMS

### Guidance and manufacture's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of this monitor requires continued operation during power mains interruptions, it is recommended that this monitor be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

# ELECTROMAGNETIC COMPATIBILITY

Table 3

## For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING Guidance and manufacturer's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of this monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance:</b></p> $d = 1.2 \sqrt[3]{\frac{P}{f}}$ <p><math>d = 1.2 \sqrt[3]{\frac{P}{f}}</math> 80 MHz to 800 MHz</p> <p><math>d = 2.3 \sqrt[3]{\frac{P}{f}}</math> 800 MHz to 2,5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,* should be less than the compliance level in each frequency range.**</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	

**NOTE 1** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

\* Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which this monitor is used exceeds the applicable RF compliance level above, this monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating this monitor.

\*\* Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

# ELECTROMAGNETIC COMPATIBILITY

Table 4

## For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

### Recommended separation distances between portable and mobile RF communications equipment and this monitor

This monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this monitor as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

---

# PRODUCT SPECIFICATIONS

---

<b>Name</b>	Premium Talking Blood Pressure Wrist Monitor
<b>Model Number</b>	01-527
<b>Display System</b>	Digital display/LCD
<b>Measuring Method</b>	Oscillometric Method, Automatic Air Inflation and Measurement
<b>Power Source</b>	2 x 1.5V  AAA Batteries
<b>Measuring Range</b>	Pressure: 0-300 mmHg; Pulse: 30-180 beats/minute
<b>Accuracy</b>	Pressure: 41°F – 104°F (5°C – 40°C) within $\pm 3$ mmHg; Pulse: $\pm 5\%$ of reading
<b>Pressurization</b>	Automatic pressurization by pump
<b>Deflation</b>	Automatic pressure release valve
<b>Memory</b>	Built-in memory with 4 banks, 396 readings 99 measurements per bank
<b>Automatic Power-Off</b>	Approximately 3 minutes after last button operation
<b>Battery Life</b>	Approximately 2 months at 3 tests per day
<b>Operation Environment</b>	Temperature 50°F – 104°F (10°C – 40°C); Humidity 15% – 93% RH (non-condensing)
<b>Storage Environment</b>	Temperature -13°F – 158°F (-25°C – 70°C); Humidity $\leq 93\%$ RH (non-condensing)
<b>Monitor Dimensions</b>	3.31" L x 2.52" W x 1.14" H
<b>Monitor Weight</b>	4.06 oz. (without batteries)
<b>Wrist Circumference Range</b>	Suitable for 5.3" – 8.5" (13.5 – 21.5 cm)
<b>Accessories</b>	Instruction manual, two AAA batteries, storage case Quick Start Guide, Blood Pressure Log
<b>Degree of protection</b>	Type BF applied part
<b>Protection against ingress of water</b>	IP22

*Specifications are subject to change without notice*

---

## 5-YEAR LIMITED WARRANTY

---

Congratulations on your purchase of a Digital Blood Pressure Monitor. Your Digital Blood Pressure Monitor is covered by the following limited warranty commencing upon the date of purchase, and subject to the following terms and conditions:

The Warrantor warrants that its Digital Blood Pressure Monitor will be free from defects in materials and workmanship under normal consumer usage for a period of five years for the original purchaser of the product.

Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the Digital Blood Pressure Monitor for commercial purposes or subjecting the Digital Blood Pressure Monitor to abnormal usage or conditions; or (d) other acts which are not the fault of the Warrantor, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the Digital Blood Pressure Monitor.

If the Digital Blood Pressure Monitor fails to conform to this limited warranty, return the Digital Blood Pressure Monitor postage prepaid to: **Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031**. When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone number. Carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. At our option, the Warrantor will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will be responsible for parts charges, if any, and repair charges not covered by this limited warranty.

EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL THE WARRANTOR BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

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