# instruction manual


## PULSE OXIMETER

MODEL 302-3110B

Please read this instruction manual completely before operating this unit.

**STOP!**  
**PLEASE ENSURE YOU HAVE ALL OF  
THE FOLLOWING COMPONENTS BEFORE  
USING YOUR PULSE OXIMETER.**

Pulse Oximeter  
with Lanyard

2 AAA  
Batteries  
(Preinstalled)

Instruction  
Manual

**IF YOU ARE MISSING ANY PARTS,  
INCLUDING INSERTS OR INSTRUCTION  
MANUALS, DO NOT RETURN TO PLACE  
OF PURCHASE. CONTACT CUSTOMER  
CARE AT 866-326-1313.**

---

## INDEX

---

|                                     |       |
|-------------------------------------|-------|
| Care & Safety Information .....     | 4-5   |
| Introduction .....                  | 6     |
| Pulse Oximeter Features .....       | 7     |
| How Pulse Oximeters Work .....      | 8     |
| Trackstar Wellness App .....        | 9     |
| How to Use Pulse Oximeter .....     | 10-12 |
| Pulse Oximeter Menu Functions ..... | 13-15 |
| Battery Installation .....          | 15    |
| Lanyard Information.....            | 16    |
| Care & Maintenance .....            | 17-19 |
| Troubleshooting.....                | 20    |
| FCC.....                            | 21    |
| Electromagnetic Compatibility.....  | 22-27 |
| Specifications .....                | 28    |
| Warranty.....                       | 29    |
| Instrucciones en Español .....      | 30-59 |

---

Toll-Free Customer Care Help Line:  
**1-866-326-1313**

Monday – Friday **8:30 a.m. – 4:30 p.m. CST**  
[www.veridianhealthcare.com](http://www.veridianhealthcare.com)

Trackstar™

Manufactured for  
Veridian Healthcare  
1175 Lakeside Drive  
Gurnee, IL 60031

Made in China

#93-1682 05/25

©2025 Veridian Healthcare, LLC

---

## CARE & SAFETY INFORMATION

---

### **WARNINGS**

- For general wellness use only
- It is not a medical device and should not be used for clinical diagnosis or monitoring of medical conditions
- Do not use the Oximeter during sports activities
- Do not use the Oximeter for continuous monitoring of patients.
- Keep unit and lanyard out of reach of children as the included lanyard may present an entanglement or choking hazard to small children. Adult supervision required; never leave children unattended with unit or lanyard
- To avoid electrical shock keep the unit away from water; never immerse the unit in liquid as this will damage the unit
- Do not handle the unit with wet hands
- Do not use the unit if it has been submersed in water or dropped
- If any abnormality in product function occurs, discontinue use until the unit has been examined; discontinue use and return to the manufacturer for review and replacement under the terms of the two-year warranty
- Do not shake the unit when in operation

---

## CARE & SAFETY INFORMATION

---

- Do not disassemble the unit or attempt to repair it

### **OPERATING CAUTIONS AND WARNINGS**

- Never tie or tape Pulse Oximeter to your finger
- Explosion hazard: do not use the product in an explosive atmosphere

### **STORAGE CAUTIONS AND WARNINGS**

- Never store the Pulse Oximeter the unit in direct sunlight, high temperature or humidity
- Keep the unit out of reach of small children

### **CLEANING CAUTIONS AND WARNINGS**

- Never immerse the device in water to clean as it may be damaged
- Follow the included instructions for cleaning the unit after each use; if the unit has not been used for an extended period, the device should be cleaned prior to use

---

## INTRODUCTION

---

**Please read all instructions carefully prior to first use.**

Thank you for purchasing a Pulse Oximeter. With proper care and use, your Pulse Oximeter will provide you with many years of reliable service.

The pulse Oximeter is a reusable device and intended for spot-checking of oxygen saturation and pulse rate for use with the finger of adult in home environments.

It is for general wellness use only.

It is not a medical device and should not be used for clinical diagnosis or monitoring of medical conditions.

It is not intended to be used under motion or low perfusion scenarios.

Some studies have indicated that certain physical traits, including, but not limited to darker skin pigmentation, finger thickness and peripheral perfusion may create variations in readings that make the oxygen level appear higher than it actually is. Users should establish a baseline value when healthy and at rest. It is recommended that users consult a licensed healthcare provider on how to interpret user readings based on physical traits or characteristics.


**PLEASE READ THIS INSTRUCTION MANUAL  
COMPLETELY BEFORE OPERATING THIS UNIT.**

---

---

## PULSE OXIMETER FEATURES

---


---

## HOW PULSE OXIMETERS WORK

---

The Pulse Oximeter works by shining two light beams into the small vessels or capillaries of the finger to determine the user's the blood oxygen saturation percentage and pulse rate that is then displayed on the Pulse Oximeter's screen. The unit uses the phrase "SpO2" to identify the oxygen saturation level and displays it as a percentage of full capacity (100%). Pulse rate is identified by the phrase "PR bpm" and is displayed as beats per minute.

Light Beam Sensors

---

## TRACKSTAR WELLNESS APP

---

Installing the Trackstar Wellness app is fast and easy.

Simply download the app, prepare your device, and start your health journey. Consult your healthcare professional before starting a new health regime, changing medication and for aid in interpreting your personal results.

1. Download the **Trackstar Wellness** app from the App Store (iOS) or Google Play (Android).
2. Launch the app and follow the instructions for account profile set up.
3. Pair your device following the on-screen directions.

**THIS PULSE OXIMETER OPERATES WITH OR WITHOUT THE TRACKSTAR APP. RESULTS HISTORY CAN ONLY BE ACCESSED THROUGH THE APP.**

---

## INSTRUCTIONS FOR USE

---

### APP INSTALLATION AND PAIRING:

If you received this oximeter as part of a healthcare maintenance program, directly from your insurance or healthcare provider, or other health management program, please consult with the program paperwork prior to use and set up. For all other uses, see the cover of this manual for information on the app set up and installation prior to use.

- Ensure your phone bluetooth is enabled, device is on and the app is open to establish connection.
- When device is first powered on, the bluetooth symbol will flash on the LED display indicating pairing is in process.
- When the device is paired, the bluetooth symbol will be solid (not blinking).
- It is not necessary to sync after every measurement; this device will store up to 50 measurements. Sync regularly to avoid losing data.

---

## **HOW TO USE YOUR PULSE OXIMETER**

---

### **USING YOUR PULSE OXIMETER**

The Pulse Oximeter may be used while the user is seated, standing or lying down. The user should not be walking or running during measurements and should take care not to excessively move the arm or hand where the oximeter is attached.

### **PRIOR TO USE**

Keep the sensors clean to help ensure accurate readings. Before and after each use, carefully swab inside the unit with a soft cloth or cotton swab with rubbing alcohol. Do not pour the alcohol directly on or into the unit.

Dirty hands or nail polish (especially dark shades) can also affect reading accuracy. Washing your hands and/or removing any nail polish prior to use to help ensure accurate readings.

**Suitable for use on any finger (excluding the thumb). In most instances, the index or middle finger work best.**

---

## HOW TO USE YOUR PULSE OXIMETER

---

1. Squeeze the hinged end of the unit to open the measurement side of the unit.
2. Insert finger nail side up into the unit until resistance is felt. The unit should gently close below the first knuckle of the finger, and the sensors inside the unit should align with the finger nail and pad at the tip of the finger.
3. Press the power button to turn the unit on.
4. After about 8-seconds, the display will begin to show the oxygen saturation (SpO<sub>2</sub>) and pulse rate (BPM) measurements. Take care to not move excessively during monitoring. These are active readings and may fluctuate slightly as you observe the display, responding to your heart rate (bps) and other factors.
6. Remove the unit from the finger. The unit will automatically turn off after 8 seconds.
7. Clean the unit as directed before storage.

---

## **PULSE OXIMETER MENU FUNCTIONS**

---

### **MENU SETTINGS OPERATION**

- After turning on the Oximeter, press and hold the power button for about 2 seconds. The Oximeter will display the settings menu. Short press to move “\*” to a corresponding option, and hold the button to confirm your selection.
- Quickly press the power button to rotate the display.

### **GENERAL SET UP**

- Short press to move “\*” to the corresponding setting, and hold the button to set alarm or sound to ON or OFF.
- When the alarm is set to ON and the measured Spo2 or PR Values go beyond the upper or lower limits, the Oximeter will sound an alert.
- When alarm is set to OFF and the measured values go beyond the limit, the Oximeter will not sound an alert.
- When sound is set to ON, the Oximeter plays a sound indicating the pulse rate during the measurement.

---

## PULSE OXIMETER MENU FUNCTIONS

---

- When sound is set to OFF, no pulse rate sound is played.
- Press the button to select a brightness level ranging from 1 to 5. The greater the value, the greater the brightness of the screen.
- While the “\*” is on the default option, hold the button to restore factory settings.

### ALARM SETUP

- Press the button to switch between options. On this screen, you can set the upper limit and lower limit of the spo2 and PR alarms.
- Select the corresponding option and hold the button to change the upper or lower limit.
- Move “\*” to the Exit option and hold the button to return to Measurement Mode.

|                                                                               |
|-------------------------------------------------------------------------------|
| MENU<br>GENERAL SETUP ENTER<br>ALARM SETUP ENTER<br>DATA REVIEW ENTER<br>EXIT |
|-------------------------------------------------------------------------------|

|                                                                                               |
|-----------------------------------------------------------------------------------------------|
| GENERAL SETUP<br>ALARM ON<br>BEEP OFF<br>DEMO PLAY<br>BRIGHTNESS 3<br>DEFAULT RESTORE<br>EXIT |
|-----------------------------------------------------------------------------------------------|

|                                                                                           |
|-------------------------------------------------------------------------------------------|
| ALARM SETUP<br>SPO2 ALM HI 100<br>SPO2 ALM LO 94<br>PR ALM HI 120<br>PR ALM LO 50<br>EXIT |
|-------------------------------------------------------------------------------------------|

---


## BATTERY INSTALLATION

---

**This unit comes complete with 2 AAA batteries preinstalled.**

It is necessary to replace the batteries when the Low Battery symbol appears on the display or when the display does not turn on after the power button is pressed.

1. Slide off the battery cover.
2. Insert or replace 2 AAA batteries into the battery compartment, ensuring to match the indicated polarity symbols.
3. Replace the battery cover.
4. It is recommended to remove the batteries if the unit will not be used for an extended period of time.

Always use new batteries. Do not mix alkaline, standard or rechargeable batteries. Please dispose of batteries according to local recycling regulations.


---

## LANYARD

---

To thread the lanyard onto the device, insert the thin end of the lanyard through the hang hole; pass the thicker end through the thin end and tighten.

**CHOKING HAZARD**-The included lanyard may present an entanglement or choking hazard to small children. Adult supervision required; never leave children unattended with unit or lanyard.

---

## **CARE & MAINTENANCE**

---

### **CLEANING**

Your equipment should be cleaned on a regular basis. If there is heavy pollution or lots of dust and sand in your place, the equipment should be cleaned more frequently. Before cleaning the equipment, consult your hospital's regulations for cleaning the equipment. Recommended cleaning agents are:

- a) Mild soap (diluted).
- b) Ethanol (70%).

To clean your equipment, follow these rules:

- a) Shut down the pulse Oximeter;
- b) Clean the display screen using a soft, clean cloth dampened with a glass cleaner;
- c) Clean the exterior surface of the equipment and probe using a soft cloth dampened with the cleaner;
- d) wipe off all the cleaning solution with a dry cloth after cleaning if necessary;
- e) Dry your equipment in a ventilated, cool place.

---

## **CARE & MAINTENANCE**

---

### **CLEANING CAUTIONS**

- Always dilute according the manufacturer's instructions or use lowest possible concentration.
- Do not immerse part of the equipment in the liquid.
- Do not pour liquid onto the equipment or accessories.
- Never use abrasive materials (such as steel wool or silver polish), or erosive cleaners (such as acetone or acetone-based cleaners).
- If you spill liquid onto the equipment, contact us or your service personnel.

### **MAINTENANCE**

The equipment's design life expectancy is about 2 years, keep your equipment and accessories free of dust and dirt, and follow these rules:

- A. Please clean the equipment before use. Remove the batteries inside the battery cassette if the equipment will not be operated for a long time;
- B. Replace the batteries in time when

---

## CARE & MAINTENANCE

---

the battery voltage indicate lamps were empty;

- C. It is recommended that the equipment should be kept in a dry environment with no corrosive gases and good ventilation any time. The moisture and high-light environments will affect its lifetime and even might damage the equipment.
- D. It is best to preserve the product in a place where the temperature is between -20°C to 60°C and the relative humidity is less than 95%.
- E. The packed equipment can be transported by ordinary conveyance. The equipment not be transported mixed with toxic, harmful, corrosive materials.

### DISINFECTION

Clean the pulse oximeter before disinfecting it. The recommended disinfectant is ethanol 70%. Disinfection steps are the same as cleaning.

### DISPOSAL

Dispose of the pulse oximeter in accordance with local environment and waste disposal laws and regulations.

---

## **TROUBLESHOOTING**

---

### **SpO2 or PR is not displaying a reading though the power is on**

- The user's oxygen level is too low to be measured; consult your physician if the measurement cannot be read
- Finger size is too big or small; select different finger for measurement
- Finger may not inserted deep enough
- Finger may be trembling; select different finger for measurement
- Excessive ambient light may be affecting reading; change location and attempt new measurement

### **The unit will not turn on**

- Ensure the batteries have been properly installed and replace if needed
- If the unit will not turn on after checking batteries, consult Customer Care

### **Display screen goes blank**

- The unit will automatically shut-off after 8 seconds if no finger is inserted
- Check the batteries and replace if needed

---

## FCC STATEMENT

---

FCC ID: 2A949-YM2022

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Caution: The user is cautioned that changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC Radiation Exposure Statement:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating.

---

## ELECTROMAGNETIC COMPATIBILITY

---

The equipment complies with the requirement of standard EN 60601-1-2:2014  
“Electromagnetic Compatibility - Medical Electrical Equipment”.

### GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSION

This model is intended for use in the electromagnetic environment specified below. The customer or the user of this model should assure that it is used in such an environment.

| Emissions test                                             | Compliance     | Electromagnetic environment-guidance                                                                                                                                                                                            |
|------------------------------------------------------------|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| RF emissions<br>CISPR 11                                   | Group 1        | The oximeter uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.                                               |
| RF emissions<br>CISPR 11                                   | Class B        | The oximeter model is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions<br>IEC 61000 -3-2                       | Not applicable |                                                                                                                                                                                                                                 |
| Voltage fluctuations / flicker emissions<br>IEC 61000 -3-3 | Not applicable |                                                                                                                                                                                                                                 |

---

# ELECTROMAGNETIC COMPATIBILITY

---

The equipment complies with the requirement of standard EN 60601-1-2:2014 "Electromagnetic Compatibility - Medical Electrical Equipment".

## GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSION

This model is intended for use in the electromagnetic environment specified below. The customer or the user of this model should assure that it is used in such an environment.

| Emissions test                                             | Compliance     | Electromagnetic environment-guidance                                                                                                                                                                                            |
|------------------------------------------------------------|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| RF emissions<br>CISPR 11                                   | Group 1        | The oximeter uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.                                               |
| RF emissions<br>CISPR 11                                   | Class B        | The oximeter model is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions<br>IEC 61000 -3-2                       | Not applicable |                                                                                                                                                                                                                                 |
| Voltage fluctuations / flicker emissions<br>IEC 61000 -3-3 | Not applicable |                                                                                                                                                                                                                                 |

# ELECTROMAGNETIC COMPATIBILITY

Guidance and manufacture's declaration –electromagnetic emission

The Pulse Oximeter is tended for use in the electromagnetic environment specified below. The user of the Pulse Oximeter should assure that it is used in such an environment.

| Immunity Test                                    | IEC60601 test level                                                                                      | Compliance Level                                                           | Electromagnetic environment- guidance                                                                                                          |
|--------------------------------------------------|----------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic discharge (ESD)<br>IEC 61000-4-2   | $\pm 8$ kVcontact<br><br>$\pm 2$ kV, $+4$ kV,<br>$\pm 8$ kV, $\pm 15$ kV air                             | $\pm 8$ kV contact<br><br>$\pm 2$ kV, $+4$ kV, $\pm 8$ kV, $\pm 15$ kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %. |
| Electrostatic transient / burst<br>IEC 61000-4-4 | $\pm 2$ kV for power supply lines<br>100 kHz repetition frequency<br>$\pm 1$ kV for input / output lines | N/A                                                                        | N/A                                                                                                                                            |
| Surge<br>IEC 61000-4-5                           | $+0.5$ kV, $\pm 1$ kV differential mode line-line                                                        | N/A                                                                        | N/A                                                                                                                                            |

# ELECTROMAGNETIC COMPATIBILITY

|                                                                                                                   |                                                                                                                                                                                                                                                             |                        |                                                                                                                                                |
|-------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>Voltage dips, short interruptions and voltage variations on power supply input lines</p> <p>IEC 61000-4-11</p> | <p>0 %UT (100 % dip in UT) for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°</p> <p>0 %UT (100 % dip in UT) for 1 cycle at 0°</p> <p>70 %UT (30 % dip in UT) for 25/30 cycle at 0°</p> <p>0 %UT (100 % dip in UT) for 250/300 cycle at 0°</p> | <p>N/A</p>             | <p>N/A</p>                                                                                                                                     |
| <p>Power frequency (50/60Hz) magnetic field</p> <p>IEC 61000-4-8</p>                                              | <p>30 A/m, 50/60Hz</p>                                                                                                                                                                                                                                      | <p>30 A/m, 50/60Hz</p> | <p>Power frequency magnetic fields should be at levels characteristic of atypical location in atypical commercial or hospital environment.</p> |

NOTE: UT is the a. c. mains voltage prior to application of the test level

# ELECTROMAGNETIC COMPATIBILITY

Guidance and manufacture's declaration –electromagnetic emission

The Pulse Oximeter is tended for use in the electromagnetic environment specified below. The user of the Pulse Oximeter should assure that it is used in such an environment.

| Immunity Test                 | IEC60601 test level                                                   | Compliance Level | Electromagnetic environment- guidance                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|-------------------------------|-----------------------------------------------------------------------|------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Conducted RF<br>IEC 61000-4-6 | 3 Vrms 150 kHz to 80 MHz<br>6 Vrms 150 kHz to 80MHz outside ISM bands | N/A              | Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to                                                                                                                                                                                                                                                                                                                                                                              |
| Radiated RF<br>IEC 61000-4-3  | 10 V/m<br>80 MHz to 2.7 GHz                                           | 10V/m            | <p>Recommended separation distance</p> $d \left[ \frac{3.5}{P_1} \right] \sqrt{P}$ $d \left[ \frac{3.5}{E_1} \right] \sqrt{P}$ <p>80MHz to 800MHz</p> $d \left[ \frac{7}{E_1} \right] \sqrt{P}$ <p>80MHz to 2.7GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts(W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range b</p> |

---

# ELECTROMAGNETIC COMPATIBILITY

---

|               |                   |       |                                                                                       |
|---------------|-------------------|-------|---------------------------------------------------------------------------------------|
| Radiated RF   | 10 V/m            | 10V/m | Interference may occur in the vicinity of equipment marked with the following symbol: |
| IEC 61000-4-3 | 80 MHz to 2.7 GHz |       |                                                                                       |

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

---

a) The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz; 3,5 MHz to 4,0 MHz; 5,3 MHz to 5,4 MHz; 7 MHz to 7,3 MHz; 10,1 MHz to 10,15 MHz; 14 MHz to 14,2 MHz; 18,07 MHz to 18,17 MHz; 21,0 MHz to 21,4 MHz; 24,89 MHz to 24,99 MHz; 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

b) The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.

c) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model is used exceeds the applicable RF compliance level above, the model should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the YM402.

d) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

---

## PRODUCT SPECIFICATIONS

---

|                            |                                                                                                                   |
|----------------------------|-------------------------------------------------------------------------------------------------------------------|
| NAME                       | Pulse Oximeter                                                                                                    |
| MODEL NUMBER               | 302-3110B                                                                                                         |
| DISPLAY SYSTEM             | OLED/TFT                                                                                                          |
| POWER SOURCE               | 2 “AAA” batteries                                                                                                 |
| SPO2 MEASURING RANGE       | 35-100%                                                                                                           |
| SPO2 ACCURACY              | 70-100%, $\pm 2\%$ ; 0%-69% unspecified                                                                           |
| PULSE RATE MEASURING RANGE | 30-250 bpm                                                                                                        |
| PULSE RATE ACCURACY        | 30-250 bpm, $\pm 3$ bpm                                                                                           |
| PR RESOLUTION              | 1 BPM                                                                                                             |
| AUTOMATIC SHUT-OFF         | After 8 seconds of no indication<br>on the sensors                                                                |
| OPERATION ENVIRONMENT      | Temperature 50°F – 104°F (10°C – 40°C)<br>Humidity 15-95% noncondensing<br>Atmospheric Pressure: 70kPa -106kPa    |
| STORAGE ENVIRONMENT        | Temperature -4°F – 140°F (-20°C – 60°C)<br>Humidity 10-95% noncondensing<br>Atmospheric Pressure: 50kPa -107.4kPa |
| UNIT DIMENSIONS            | 62 x 35 x 34 mm                                                                                                   |
| UNIT WEIGHT                | 30g (without batteries)                                                                                           |

*Specifications are subject to change without notice*

---

## ONE-YEAR LIMITED WARRANTY

---

Congratulations on your purchase of a Pulse Oximeter. Your Pulse Oximeter is covered by the following limited warranty commencing upon the date of purchase for a period of one year, and subject to the following terms and conditions:

The warrantor warrants that its Pulse Oximeter will be free from defects in materials and workmanship under normal consumer usage for one year as the original purchaser owns the product.

Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the Pulse Oximeter for commercial purposes or subjecting the Pulse Oximeter to abnormal usage or conditions; or (d) other acts which are not the fault of the warrantor, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the Pulse Oximeter.

If the Pulse Oximeter fails to conform to this limited warranty, return the Pulse Oximeter to: **Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031**. When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone number. Carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. At our option, The warrantor will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will be responsible for parts charges, if any, and repair charges not covered by this limited warranty.

EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL THE WARRANTOR BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

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## measure • track • share

Setting up your health monitor is fast and easy. Simply download the app, prepare your device, and start your health journey. Consult your healthcare professional before starting a new health regime, changing medication and for aid in interpreting your personal results.

1. Download the **Trackstar Wellness** app from the App Store (iOS) or Google Play (Android).
2. Launch the app and follow the instructions for account profile set up.
3. Consult the included instruction manual for detailed instructions, tips, cautions and cleaning direction.