instruction manual

INFRARED THERMOMETER

MODEL 302-6110B


Please read this instruction manual completely before operating this unit.


STOP!
PLEASE ENSURE YOU HAVE ALL OF THE
FOLLOWING COMPONENTS BEFORE USING
YOUR DIGITAL THERMOMETER

Digital
Thermometer

Two (2) AAA
Batteries
Preinstalled

Instruction Manual
APP Quick Start
Guide

INTENTION FOR USE

This infrared forehead thermometer is intended for the intermittent measurement of human body temperature from the skin surface of forehead. The device can be re-used by people of all ages for home use and clinical use.

Intended use: The infrared forehead thermometer is used to measure body temperature by measuring forehead.

INDEX

Warnings & Cautions	4-5
Body Temperature & Fever	6
Tips for Measuring	6
Product Features.....	7
Display Symbols	8
Trackstar Wellness App	9
Battery Installation	10
Date/Time Setting	11
Measurement Mode Options	12
Switching Between F/C	12
Night/Sleeping Silent Option	12
Measuring Body Temperature	13-14
Surface / Object Temperature	14
Memory Recall.....	15
Care & Safety	16
Cleaning & Storage.....	17
Troubleshooting.....	18
Device & Label Symbols.....	19
FCC Information	20
Electromagnetic Chart.....	21-26
Specifications	27-28
Warranty	29
Instructions in Spanish	31-59
Trackstar Wellness App	60

Toll-Free
Customer Care Help Line:
1-866-326-1313
Monday - Friday
8:30 a.m. - 4:30 p.m. CST
www.veridianhealthcare.com

trackstar™
Manufactured for:
Veridian Healthcare
1175 Lakeside Drive
Gurnee, IL 60031
www.veridianhealthcare.com
Made in China
#93-1730 12/25
©2025 Veridian Healthcare, LLC

WARNINGS & CAUTIONS

NOTE: Read all instructions carefully before use. The following basic precautions are needed when using an electrical product.

CAUTION: Failure to read and observe all precautions could result in personal injury or equipment damage.

GENERAL CAUTIONS AND WARNINGS

- Use as directed in this manual; misuse may void the warranty.
- This thermometer may be used for professional purposes or for personal home use.
- Measurement results should be discussed with your physician or healthcare professional; never self-diagnose, attempt treatment or alter medications this may be dangerous.
- The patient is the intended operator of this device.
- Wait approximately 30 minutes before measurement after strenuous physical activity, bathing, spending time outdoors or other activity which may artificially affect temperatures.
- Operation outside of stated operating temperature may result in measurement error or device malfunction; operation environment temperature is 41°F -104°F (5°C - 40°C); Humidity: 15-85%RH
- If the thermometer is stored outside of operating temperature range, allow at least 30 minutes in correct range before taking a measurement.
- Keep in a clean, dry area; do not store the thermometer in direct sunlight, dusty or humid environments, or extreme temperatures.
- Clean the lens before use; do not touch the lens directly with your fingers.
- Do not expose this thermometer to electric shocks.
- Do not drop this thermometer.

WARNINGS & CAUTIONS

- Follow the included cleaning and storage instructions; do not submerge the thermometer or wet.
- Do not use this thermometer if it is damaged.
- Do not disassemble the unit or attempt to repair it; substitution of a component different from that supplied may result in measurement error and will void manufacturer warranty.
- Remove the batteries if you do not wish to use the unit for an extended period of time.
- Do not attempt to service the ME equipment while in use with a patient.
- Portable and mobile RF communications can affect the device. The device needs special precaution regarding EMC according to the EMC information provided in the accompany documents ME equipment should not be cleaned and disinfected while in use.
- Keep the battery away from children. Remove battery from the device when not in operation for a long time.

BODY TEMPERATURE & FEVER


Body temperature can vary from 97°F(36.1°C) and 99°F(37.2°C) and still be considered normal. The medically accepted 'normal' body temperature is 98.6°F(37°C). Body temperature is commonly lower upon waking than at any point during the rest of an individual's active hours.

A fever is defined as a body temperature that is elevated above the normal for that person. It is important to determine what is normal for an individual before determining if a fever is present. Tracking an individual's temperature on a consistent basis, at the same body site and at the same time every day will help determine an individual's normal temperature.


TIPS FOR MEASURING BODY TEMPERATURE

- Wait at least one minute between consecutive readings, to ensure accuracy.
- Attempting to take temperature readings from sites on the body other than the forehead may produce inaccurate results.
- The patient should remain still while the reading is being taken.
- If an error message appears, see the Troubleshooting section of this manual for instructions on how to resolve the problem.
- Readings taken while asleep should not be compared directly to readings taken while awake, as body temperature while asleep is typically lower.
- Do not take body temperature readings within 30 minutes of being outdoors, exercising or bathing.
- Level of activity, time of day, and other factors may affect body temperature. Keep records of your individual body temperature to use as a reference when judging a fever.

PRODUCT FEATURES

DISPLAY SYMBOLS

OBJECT
MODE

BLUETOOTH

FOREHEAD
MODE

MEMORY
MODE

POSITION
SENSOR ICON

DATE/TIME
INDICATOR

SPEAKER ON /
SPEAKER OFF

BATTERY
STRENGTH

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TRACKSTAR WELLNESS APP

Installing the Trackstar Wellness app is fast and easy.

Simply download the app, prepare your device, and start your health journey. Consult your healthcare professional before starting a new health regime, changing medication and for aid in interpreting your personal results.

1. Download the **Trackstar Wellness** app from the App Store (iOS) or Google Play (Android).
2. Launch the app and follow the instructions for account profile set up.
3. Pair your thermometer following the on-screen directions.

THIS THERMOMETER WORKS WITH OR WITHOUT THE TRACKSTAR APP. SOME FUNCTIONS, SUCH AS DATA STORAGE AND RETRIEVAL, ARE ONLY ACCESSIBLE THROUGH THE APP.


BATTERY INSTALLATION

THIS UNIT INCLUDES 2 AAA BATTERIES, PREINSTALLED.

When using for the first time, remove the pull strip to activate the batteries.

It is necessary to replace the batteries when the low battery symbol  appears on the display or the thermometer will not power on.

1. Slide open the battery compartment on the backside of the device.
2. Insert or replace 2 AAA batteries into the battery compartment, ensuring to match the indicated polarity symbols. Always use new batteries.
3. Replace the battery cover.
4. Dispose of batteries according to local disposal and recycling regulations.

It is recommended to remove the batteries if the unit will not be used for an extended period of time.

DATE / TIME SETTING

The Date / Time setting option is automatically set when this thermometer is paired with an App or assistive bluetooth device. Users can set the date / time manually at any time. The sequence for setting this function requires Fahrenheit and Celsius selection as the first step.

1. To switch between Fahrenheit and Celsius, while the unit is OFF press and hold the side **SETTING** button. The icon on the display will confirm what setting is in use. Press and release **POWER** to change modes. Press **SETTING** to confirm
2. The current time format will flash on the display, AM/PM (12:00 hour) or Military Time (24:00 hour). Press **POWER** to advance to the select format, press **SETTING** to confirm.
3. The current hour will flash on the display. Press **POWER** to advance to the hour. Press **SETTING** to confirm.
4. Repeat the same steps to set the Minute, Year, Month, Date:
POWER to advance to desired minute; **SETTING** to confirm
POWER to advance to desired year; **SETTING** to confirm
POWER to advance to desired month; **SETTING** to confirm
POWER to advance to desired date; **SETTING** to confirm

After the date is set, the display will show OFF.

The thermometer is now ready for measurement.

MEASUREMENT MODE OPTIONS

The Infrared Thermometer offers a selection of settings for use on people or objects.

To toggle between modes, when the device is ON, press and hold the **SETTING** button to select the desired measurement type. The symbol indicating the **PERSON** or **OBJECT** will appear on the display.

- **PERSON** — used for body temperature measurements taken at the forehead.
- **OBJECT/SURFACE** — used for checking baby bottles, bathwater, food or similar items. Always use caution when checking the temperature of food or beverages.

Follow instructions for use in this manual for proper use.

SWITCHING BETWEEN °F / °C

Your thermometer can display results in either degrees Fahrenheit (°F) or degrees Celsius (°C).

To switch between Fahrenheit and Celsius, while the unit is OFF press and hold the side **SETTING** button. The icon on the display will confirm what setting is in use. Press and release **POWER** to change modes.

NIGHT / SLEEPING SILENT OPTION

Night Mode is a convenient way to conduct measurements without an audible beep indicating the thermometer has powered on or a measurement is complete.

1. While the thermometer is powered ON, press and release the **SETTING** button on the left side of the thermometer.
2. The speaker icon on the upper right of the display screen will show a speaker symbol when ON, and a speaker with line when the speaker is OFF.

MEASURING BODY TEMPERATURE

To measure body temperature using your Infrared Thermometer, follow the below steps.

This thermometer offers three sensor positions; swivel the head to the best position for your use. When the probe is face down into the thermometer, the display is blank and the unit is OFF.

1. Rotate the probe to the desired position; the display will activate and run a self-test; all symbols on the display will momentarily appear, Fig. 1. If the sensor is already in the desired position, press the **POWER** button to turn the thermometer on.


Fig. 1

2. The display will show dashed lines, the selected mode, and the temperature mode ($^{\circ}\text{F}/^{\circ}\text{C}$). The thermometer is ready for measurement.
3. Press and release the **POWER** button while aiming the infrared lens at the individual's forehead, holding the thermometer approximately 1.0" away.
4. This thermometer includes a positioning sensor; if the sensor is too far, the positioning sensor icon will appear; move the unit closer to the center of the forehead.
5. Within one second, the thermometer will beep and the results will appear on the display, Fig. 2. The display will illuminate, indicating if a fever may be present:
GREEN — $< 99.1^{\circ}\text{F}$ (37.3°C); no fever is present
YELLOW — 99.1°F - 99.9°F (37.3°C - 37.7°C); mild fever is present
RED — $\geq 100.0^{\circ}\text{F}$ ($\geq 37.8^{\circ}\text{C}$); high fever is present
6. The Bluetooth symbol flashing indicates that the Bluetooth is not connected, while the Bluetooth symbol often appears to indicate that the Bluetooth is connected.
7. Rotate the probe to the original position, the thermometer will turn off. The thermometer will automatically turn off after approximately 1 minute.


Fig. 2

MEASURING BODY TEMPERATURE

NOTE: Proper placement and position of the sensor is crucial for effective and accurate measurement.

The forehead should be clean and clear of excessive lotion or oil and hair should be moved to the side.

MEASURING SURFACE / OBJECT TEMPERATURE

The Infrared Thermometer can measure objects and liquids from 32°F to 212°F (0°C - 100°C). To measure the temperature of an object or liquid using your Infrared Thermometer, follow these steps:



CAUTION: Surface mode measures the temperature at the surface of the object measured. Use common sense caution when testing food, drink or bathwater to protect against possible burns.

1. Rotate the probe to the desired position; the display will activate and run a self-test; all symbols on the display will momentarily appear. If the sensor is already in the desired position, press the **POWER** button to turn the thermometer on.
2. Press and hold the **SETTING** button for approximately 3-seconds; the mode on the display will alternate from PERSON to OBJECT icon.
3. Aim the infrared lens at the object, holding the thermometer approximately 1.0" away.
4. The thermometer will emit a single beep and the temperature results will be displayed. In object mode, the display will default to green backlight only.
5. Rotate the probe to the original position and the thermometer will turn off. The thermometer will automatically turn off after approximately 1 minute.

MEMORY RECALL

This Infrared Thermometer automatically stores the last 30 readings in memory, both person and object modes. To view the readings stored in memory, follow these steps:

1. When the thermometer is ON, press and hold the **SETTING:** and **POWER** buttons simultaneously.
2. The memory icon will appear on the lower right of the display; the display will show a 1 and alternate to the most recent reading in memory including the date/time of measurement.

3. Press the **SETTING:** button to scroll through the readings from most recent to oldest.
4. Press **POWER** at any time to exit memory recall and take a new measurement, or rotate the sensor head to the down position.

CARE & SAFETY

Improper care or use of your thermometer may result in injury, damage to the unit or inaccurate readings, which can lead to a misdiagnosed fever. Following these instructions will ensure the thermometer's accuracy and long life.

- Never allow the thermometer to be used without adult supervision.
- The infrared technology used in this thermometer is safe and is not harmful in any way. However, we recommend avoiding looking into the lens during use.
- Clean the lens before and after each use.
- Do not touch the lens with fingers.
- Dropping the thermometer or causing other heavy impact may damage it.
- Do not store the thermometer in direct sunlight, dusty or humid environments, or extreme temperatures.
- Do not disassemble the thermometer or attempt to repair it.
- If a reading cannot be taken, and the troubleshooting steps do not fix the problem, the thermometer must be replaced.

CLEANING & STORAGE

The lens should be cleaned before and after each use, to eliminate dirt and fingerprints and ensure the accuracy of your readings. Proper storage will prevent degradation and damage to the unit.

CLEANING:

- An alcohol swab, or a soft cotton cloth moistened with 70% alcohol, should be used to gently cleanse the lens. Do not use water to wash the thermometer.
- The body of the thermometer can be cleaned with a soft, dry cloth.
- Never use a paper towel to clean the probe tip, as it may scratch the lens.
- Do not immerse the thermometer or allow liquid to enter the casing.
- Do not use cleaning agents other than alcohol to clean the thermometer; abrasives, thinners or gasoline will damage the unit.
- Do not use the thermometer again until it is completely dry and allowed to remain untouched for at least 30 minutes.

STORAGE:

- Store the unit in an area away from sunlight, dust and humidity.
- Protect the thermometer from mechanical shock or heavy impact.
- Avoid extended exposure to extreme temperatures. The thermometer should be stored in an environment between -13°F and 131°F (-25°C—55°C), and with relative humidity of 15% - 95% or lower.
- Remove the batteries during extended storage.

TROUBLESHOOTING

Performance may be degraded should one or more of the following occur:

- Operation outside stated temperature and humidity range
- Storage outside state temperature and humidity range
- Mechanical shock
- Patient temperature is below ambient temperature

The following chart indicates error messages that may appear on your display.

Error Message	Problem	Solution
	Thermometer is not functioning	Remove battery for one minute and replace. Reinstall and retry. If error persists, contact Customer Care.
	Ambient (room) temperature is not within the range of 41°F-104°F (5°C-40°C)	Allow thermometer to rest in a room at ambient temperature for at least 30-minutes and retry.
	Temperature results are too high: Forehead over 109.4°F (43.0°C); Object mode over 212°F (100°C)	Read 'Tips for Taking Temperature' closely and attempt new measurement
	Temperature results are too low: Forehead lower 93.2°F (34.0°C); Object mode over 32°F (0°C)	Allow thermometer to rest in a room at ambient temperature for at least 30-minutes and retry.
  	Solid battery Flashing battery Empty battery	Use as normal Replace batteries soon Replace batteries immediately

DEVICE & LABEL SYMBOLS

These symbols may appear on your device, instructions or packaging and may vary by make and model.

Symbol	Explanation	Symbol	Explanation
	Caution		Atmospheric Pressure Limitation
	Direct Current		Storage / Transport Humidity Limit; 15% - 95% RH
	Batch Code		Date of Manufacture
	Storage / Transport Temp Limit; -13°F-131°F (-25°C-55°C)		Manufacturer
	Type BF Applied Parts	IP22	Degree of protection against harmful ingress of water and particulate matter; 2 Protected and against solid foreign objects of 12.5mm and greater; 2 Protection against vertically falling water drops when ENCLOSURE tilted up to 15°
	Read This Manual (Symbol is in Blue)		
	Recycling Symbol		
	DISPOSAL: Do not dispose of this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.		Unique Device Identifier

FCC STATEMENT

IMPORTANT INFORMATION REQUIRED BY THE FCC

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference; and (2) this device must accept any interference received, including interference that may cause undesirable operation.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy, and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna.

Increase the separation between the equipment and receiver.

Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

Consult the dealer or an experienced radio/TV technician for help.

CAUTION: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

INTRODUCTION

This appliance conforms to the following standards: ASTM E1965-98 Standard Specification for Infrared Thermometers for Intermittent Determination of Patient Temperature, ISO 80601-2-56 Medical electrical equipment —Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement, IEC 60601-1-11 Medical electrical equipment —Part 1-11: General requirements for basic safety and essential performance –Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of IEC 60601-1-2(EMC) , AAMI ANSI ES60601-1(Safety) standards.

ELECTROMAGNETIC COMPATIBILITY

The device satisfies the EMC requirements of the international standard IEC 60601-1-2. The requirements are satisfied under the conditions described in the table below. The device is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the instructions for use. Portable and mobile HF communications equipment can affect the device. Use of the unit in conjunction with non-approved accessories can affect the device negatively and alter the electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

Table 1

Guidance and manufacturer's declaration -electromagnetic emissions and Immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Emissions Test	Compliance
Conducted emissions CISPR 11	Not Applicable
Radiated emissions CISPR 11	Group 1 Group B
Harmonic emissions IEC 61000-3-2	Not Applicable
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not Applicable

ELECTROMAGNETIC COMPATIBILITY

Table 2

Guidance and manufacture's declaration - electromagnetic emissions

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical transient/burst IEC 61000-4-4	±2 kV for power supply lines 100kHz repetition frequency ±1 kV for input/output lines	Not Applicable
Surge IEC 61000-4-5	±0.5 kV, ±1 kV differential mode line-line	Not Applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; (100% dip in UT) for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT; (100% dip in UT) for 1 cycle at 0° 70 % UT ; (30% dip in UT) for 25/30 cycles at 0° 0 % UT (100% dip in UT) for 250/300 cycle at 0°	Not Applicable
Power frequency (50/60 Hz) Magnetic field IEC 61000-4-8	30 A/m, 50/60Hz	30 A/m, 50/60Hz
Proximity magnetic fields	IEC 61000-4-39	See Table 3
NOTE: UT is the a.c. mains voltage prior to application of the test level.		

ELECTROMAGNETIC COMPATIBILITY

Table 3

Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields

Test Frequency	Modulation	IMMUNITY TEST LEVEL(A/m)
30 kHz a)	CW	8
134,2 kHz	Pulse Modulation b) 2,1 kHz b)	65 c)
13,56 MHz	Pulse Modulation b) 50 kHz b)	7,5 c)
a) This test applicable only to ME EQUIPMENT and ME SYSTEMS intended for use in the HOME HEALTHCARE ENVIRONMENT. b) The carrier shall be modulated using a 50% duty cycle square wave signal. c) r.m.s., before modulation is applied.		

Table 4

Guidance and manufacture's declaration - electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity Test Level	IEC 60601 Test Level	Compliance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms 150 kHz to 80 MHz outside ISM bands a	Not Applicable
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7GHz	10 V/m

ELECTROMAGNETIC COMPATIBILITY

Table 4 -- continued

WARNINGS!

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

b The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.

c Field strengths from fixed transmitters, such as base stations for radio(-cellular/ cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

d Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

ELECTROMAGNETIC COMPATIBILITY

Table 5

Recommended separation distances between RF wireless communications equipment

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between RF wireless communications equipment and the device as recommended below, according to the maximum output power of the communications equipment.

Frequency Compliance MHz	Maximum Power W	Distance	IEC 60601 Test Level	Level
385	1,8	0,3	27	27
450	2	0,3	28	28
710	0,2	0,3	9	9
745				
780				
810	2	0,3	28	28
870				
930				
1720	2	0,3	28	28
1845				
1970				
2450	2	0,3	28	28
5240	0,2	0,3	9	9
5500				
5785				

NOTE 1: The guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people

ELECTROMAGNETIC COMPATIBILITY

WARNINGS!

- This device should not be used in the vicinity or on the top of other electronic equipment such as cell phone, transceiver or radio control products. If you have to do so, the device should be observed to verify normal operation.
- The use of accessories and power cord other than those specified, with the exception of cables sold by the manufacturer of the equipment or system as replacement parts for internal components, may result in increased emissions or decreased immunity of the equipment or system.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Portable and mobile RF communications can affect the device. The device needs special pre-cautions regarding EMC according to the EMC information provided in the accompany documents.
- Do not use the devices in the MR environment.
- The Operator should not use the system and should inform the customer service, if the ESSENTIAL PERFORMANCE is lost or degraded due to EM DISTURBANCES.

PRECAUTION:

The performance of the device may be degraded should one or more of the following occur:

- Operation outside the manufacturer's stated temperature and humidity range.
- Storage outside the manufacturer's stated temperature and humidity range.
- Mechanical shock (for example, drop test) or degraded sensor.
- Patient temperature is below ambient temperature.

SPECIFICATIONS

Name	Infrared Thermometer
Model Number	302-6110B
Display System	Display with 0.1° F (0.1° C) resolution
Measuring Site	Forehead
Reference Body Site	Oral (This thermometer converts the forehead temperature to display its “oral equivalent.”)
Mode Operation	Adjusted Mode
Lab	Forehead mode:
Accuracy	± 0.4°F (0.2°C) during 95.9°F-107.6°F (35.5°C-42.0°C) at 59.0°F-95.0°F (15°C-35°C) operating temperature range ±0.5°F (0.3°C) for other measuring and operating temperature range Object mode: ±4% or ±4°F (2°C) whichever is greater
Measure time	1 Second
Operating Temperature	41°F-104°F (5°C-40°C) 15%-85%RH, non-condensing Atmospheric Pressure: 70kPa ~ 106kPa
Storage/transport	-13°F-131°F (-25°C~ 55°F), 15%-95%RH, non-condensing
Temperature	Atmospheric Pressure: 70kPa ~ 106kPa
Clinical	0-1 year:

SPECIFICATIONS

Accuracy	Clinical bias: - 0.31°F -0.17°C; Clinical repeatability: 0.11°F (0.06°C); Limits of agreement: 0.43°F (0.24°C) 1- 5 years: Clinical bias: - 0.31°F (-0.17°C) ; Clinical repeatability: 0.13°F (0.07°C); Limits of agreement: 0.40°F (0.22°C) over 5 years: Clinical bias: - 0.31°F (-0.17°C) ; Clinical repeatability: 0.36F (0.20°C); Limits of agreement: 0.47°F (0.26°C)
Shock	Withstands drop of 3 feet
Unit Dimensions	1.62" x 1.01" x 6.05" (41.4 x 25.9 x 153.9mm)
Weight	Approx.102 grams (with batteries)
Battery	2 x AAA size; DC3V
Battery Life	Approx. 6000 readings
Applied Part	Probe
Expected service life/Shelf life	Three years
Ingress protecting rating	IP22
Accessories	Instruction manual, two AAA size batteries
Contraindications	None
Bluetooth technical parameters	Modulation mode: GFSK Frequency range: 2400-2483.5MHz Bandwidth occupied: ≤2MHz Transmit power: ≤20db

Specifications are subject to change without notice.

ONE-YEAR LIMITED WARRANTY

The warrantor guarantees that its Infrared Thermometer will be free from manufacturing defects for a period of one year from the date of purchase. This warranty does not cover batteries or other power sources that may be provided with or used with the Infrared Thermometer. This warranty is voided if the Infrared Thermometer is misused or abused in any manner.

If the Infrared Thermometer fails to operate during the time the original purchaser owns it, return it postage prepaid to: **Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031**. When returning a product, please also include your name, address, phone number and brief description of the specific problem with your unit. Take care to carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. The warrantor will repair or replace the defective unit at our option. Repair or replacement of the defective unit, at the warrantor's option, is the sole remedy under this warranty.

ANY IMPLIED WARRANTIES WHICH THE PURCHASER MAY HAVE ARE LIMITED IN DURATION TO THE TIME THAT THE ORIGINAL CONSUMER PURCHASER OWNS THE PRODUCT. Some states do not allow limitations or how long an implied warranty lasts, so the above limitation may not apply to you.

This warranty constitutes the warrantor's only reasonability and obligation to repair and/or replace materials or components, or refund the purchase price. The warrantor will not be responsible for any indirect, incidental, special consequential or punitive damages or other loss, including, but not limited to damage to or loss of other property or liable to the purchaser for any amount in excess of the cost of repair and/or replacement of the unit, or the purchaser for any amount in excess of the cost of repair and/or replacement of the unit, or the purchase price of the unit. Some states do not allow the exclusive or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

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

De lunes a viernes
8:30 a.m. - 4:30 p.m. CST
www.veridianhealthcare.com

trackstar™
Manufactured for
Veridian Healthcare
1175 Lakeside Drive
Gurnee, IL 60031

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measure • track • share

Setting up your health monitor is fast and easy. Simply download the app, prepare your device, and start your health journey. Consult your healthcare professional before starting a new health regime, changing medication and for aid in interpreting your personal results.

1. Download the **Trackstar Wellness** app from the App Store (iOS) or Google Play (Android).
2. Launch the app and follow the instructions for account profile set up.
3. Consult the included instruction manual for detailed instructions, tips, cautions and cleaning direction.