

GE Blood Glucose Monitoring System 333

Instructions for Use

1. Preface

Thank you for selecting the GE Blood Glucose Monitoring System 333. This instruction provides all the information you need to operate this product for accurate test results. Please read this entire instruction before you start any testing. If you have any questions or concerns, please contact the GE Diabetes Customer Support Center toll free at 1-866-613-7085 or ask your questions by the mobile app. We will make every effort to assist you.

2. Intended Use

The GE Blood Glucose Monitoring System 333 is intended to be used for the quantitative measurement of glucose (sugar) in fresh capillary whole blood samples drawn from the fingertips, forearm or palm. The GE Blood Glucose Monitoring System 333 is intended to be used by a single person and should not be shared. The GE Blood Glucose Monitoring System 333 is intended for self testing outside the body (in vitro diagnostic use) by people with diabetes at home as an aid to monitor the effectiveness of diabetes control. The GE Blood Glucose Monitoring System 333 should not be used for the diagnosis of, or screening for diabetes or for neonatal use. Alternative site testing should be done only during steady - state times (when glucose is not changing rapidly).

The GE Blood Glucose Test Strips 333 are for use with the GE Blood Glucose Meter 333 to quantitatively measure glucose (sugar) in fresh capillary whole blood samples drawn from the fingertips, forearm or palm.

GE Control Solution 333 are for use with the GE Blood Glucose Meter 333 and the GE Blood Glucose Test Strip 333 to check that the meter and test strips are working together properly and that the test is performing correctly.

3. Precaution

- Before using the GE Blood Glucose Monitoring System 333 to test your blood glucose, please read all of the instructions and conduct all of the tests including the quality control test.
- Please perform the quality control test regularly to make sure the test results are accurate.
- The GE Blood Glucose Meter 333 can only be used with GE Blood Glucose Test Strips 333. Other brands' test strips should not be used under any circumstances. The use of other brands' strips may give inaccurate results.
- If the GE Blood Glucose Meter 333 and Test Strips are exposed to a substantial change in temperature, please wait 45 minutes before measurement.
- Follow all environmental protection regulations when disposing of batteries.
- The GE Blood Glucose Meter 333 will be set as mg/dL by default when sold in United States. If meter display shows "mmol/L" during setting or test, please contact GE Diabetes Customer Support Center. Use of the wrong unit of measurement may cause incorrect treatment.
- Please connect the Blood Glucose Meter 333 with the mobile app regularly to manage and review your daily glucose test results.
- Avoid contact with dripping or splashing liquids.
- The minimum blood sample size to test using the GE Blood Glucose Monitoring System 333 is 0.75 µL.
- Blood sample size above 3.0 µL might contaminate the test strip port and the meter.
- Blood sample size below 0.75 µL may cause an inaccurate result or may prevent a meter reading. An Er4 reading will be displayed if the sample size is too small. In this case, repeat the test with a new test strip.


Sample Size Example



Important Safety Notes:

- The GE Blood Glucose Meter 333 and Lancing Device are for single patient use. Do not share them with anyone including other family members.
- All parts of the kit are considered biohazardous and can potentially transmit infectious diseases, even after following the cleaning and disinfecting procedures.
- Users should wash their hands thoroughly with soap and water before and after handling the meter, lancing device, or test strips. Please refer to the section "15. Cleaning and Disinfecting Procedures".

4. The GE 333 Blood Glucose Meter




mem	Indicates a test result stored in memory		Indicate when to apply the blood sample
AVG	Indicates the average result		Appears after you insert test strip into meter
CS	Indicates a control solution test result	Error	Appears when an error occurs
mmol/L mg/dL	Unit of test result	88:88	Current date under time mode or testing date under memory mode
888	Test result		Warns when the operational temperature limit is exceeded during testing
	Warns when the batteries are low or must be replaced	AM	Indicates the time in 12H format
Day	Indicates the time period of test results	88:88	Displays current time under time mode or testing time under memory mode
	Indicates that the Bluetooth is turned on		
	Indicates meal markers or other conditions with test		
EN OK ES NO	Manufacturing use only		

5. The GE Blood Glucose Test Strip 333 and Battery Replacement


Inserting the GE Blood Glucose Test Strip 333:

1. Hold the GE Blood Glucose Test Strip 333 between your thumb and middle finger with the view window "V" facing up.
2. Put your forefinger on the side of the strip as shown.
3. Insert the GE Blood Glucose Test Strip 333 into test strip port until it clicks and firmly stops.



Removing the GE Blood Glucose Test Strip 333:

1. Hold the GE Blood Glucose Test Strip 333 as shown.
2. Rotate the GE Blood Glucose Test Strip 333 counterclockwise and pull up simultaneously.
3. Take the GE Blood Glucose Test Strip 333 out of the test strip port. Please follow your healthcare professional's direction and discard used strips properly.



PRECAUTION

- Close the GE Blood Glucose Test Strip 333 vial immediately after removing a test strip.
- Do not reuse GE Blood Glucose Test Strips 333.
- Do not use expired GE Blood Glucose Test Strips 333.
- When you open a new vial of GE Blood Glucose Test Strips 333, record the date on the vial. Discard the vial of test strips after 3 months from opening.
- Store the GE Blood Glucose Test Strip 333 between 39-86°F (4-30°C) and in a location (10-90% relative humidity). Do not expose to direct sunlight or heat.
- For detailed information, please refer to the GE Blood Glucose Test Strip 333 Package Insert.
- No suggest to test blood sugar when the flight is strictly restrict the Bluetooth connection. Please turn off you meter to avoid the auto-on of Bluetooth.

Your GE Blood Glucose Meter 333 comes with two AAA, batteries installed. Two new batteries will provide power to perform approximately 800 tests under normal use with one data transfer. Press the main button or insert a strip to activate your meter.

1. Turn the GE Blood Glucose Meter 333 over. Press and slide the battery cover to open.
2. Install the batteries. Be sure to put the batteries in the correct direction.
3. Slide the battery cover back until it snaps into place.
4. The GE Blood Glucose Meter 333 will enter Self-Testing Mode automatically when replacing the battery (all symbols will appear on the screen).



5. Press any button to exit the self-test and enter Setting Mode.
6. You must set the time and date when replacing the batteries. See the following "Setting Up Your Meter - Setting the Date, Time, and Volume". Test results are still stored in the memory.

CAUTION

- Please follow the local regulation to properly discard a used battery.

6. Setting Up Your Meter - Setting the Date, Time and Volume

You can enter the Setting Mode two ways.

1. Replace the Battery
After removing the battery, press the main button several times until there is no signal on screen, then follow the battery installation steps to replace batteries. The GE Blood Glucose Meter 333 will perform a self-test. Press the main button to exit the self-test and enter the Setting Mode.
2. With Battery Inserted
Press the main button to turn on the GE Blood Glucose Meter 333. Hold down the main button for 7 seconds. During this time the screen will go blank until you hear a beep. After the beep, the meter will turn on into the Setting Mode. The display screen will show setting data.

NOTE

- Use the Left and Right buttons to select your setting data and press the main button to confirm each selection. After confirming all of the settings, you will return to the Time Mode.

1. Year Setting

With the year format blinking, press the Left or Right button until you see the current year. Press the Main button to confirm. Once the year is confirmed, the Month Setting will appear.

2. Month setting

With the month blinking, press the Left or Right button until you see the current month. Press the Main button to confirm. Once the month is confirmed, the Day Setting will appear.

3. Day setting

With the day blinking, press the Left or Right button until you see the current day. Press the Main button to confirm. Once the day is confirmed, the Time Setting will appear.

4. Time format 12/24H selection

With the time format blinking, press the Left or Right button to switch between 12H and 24H. Press the Main button to confirm. Once confirmed, the Hour Setting will appear.

5. Hour setting

With the hour blinking, press the Left or Right button until you see the current hour. Press the Main button to confirm. Once the hour is confirmed, the Minute Setting will appear.

6. Minute setting

With the minutes blinking, press the Left or Right button until you see the current minute. Press the Main button to confirm. Once the minute is confirmed, the Volume Setting will appear.

7. Volume Setting

With the volume blinking, press the Left or Right button to turn the volume on or off. Press the Main button to confirm and finish the settings.

8. Ending Setting

After volume setting, for confirmation you will hear a beep (if volume is turned on). All the settings are saved and completed and will return to time screen.

NOTE

- If you do not change any settings during Meter Set-up for over 2 minutes, the GE Blood Glucose Meter 333 will leave setting mode and power off automatically.

Turning On / Off the Meter

1. How to turn on the GE Blood Glucose Meter 333: Press the Main button or Insert a test strip.
2. Manual Power off: To power off the GE Blood Glucose Meter 333, press and hold the Main button for 4 seconds.
3. Auto Power off: The GE Blood Glucose Meter 333 will power off automatically after 2 minutes if no buttons are pressed or no strip is inserted.

7. Performing a Blood Glucose Test

Please refer to the GE Blood Glucose Test Strip 333 Insert or GE Blood Glucose Monitoring System 333 User's Manual for detail procedure.

CAUTION

- Do not apply your blood drop to the GE Blood Glucose Meter 333 is performing an internal test and will display "Lo" and "Error" if you apply blood too soon. If this occurs, please repeat the test with a new test strip.
- Record the opening date of a new test strip vial. Discard the vial of test strips 3 months after opening.
- Always keep the metal contacts of the test strip port clean. If any dust or impurities are present, please clean with a small, soft brush.
- The GE Blood Glucose Monitoring System 333 and GE Lancing Device are intended for a single user only. Do not share these devices with anyone, including other family members.
- All parts of this kit are considered biohazards and can potentially transmit infectious diseases, even after you have performed the cleaning and disinfecting procedure.
- Users should wash hands thoroughly with soap and water after handling the meter, lancing device, and test strips.



8. Setting Markers

According to actual life status, you can set meal markers with test result.

You have to perform meal marker setting immediately before the meter displays the reading of your test result.

1. After test, the test result is displayed on the screen with time, date and flash "Before".
2. Press the Right or Left button to choose among "None" (no marker), "Before" (before meal), "After" (after meal), the marker will be recorded when continuing to blood test with needed one.


	Before meal: You execute this blood glucose test before a meal.
	After meal: You execute this blood glucose test after a meal.
	(Blank) no marker: regular test or definition by user.

NOTE

- If markers has disappeared or not properly displayed, please contact the GE Diabetes Customer Support Center.

9. Auto code

The GE Blood Glucose Meter 333 will automatically detect the Code number on the strip. You do not have to manually change the Code key or adjust the code number.



10. "Lo" and "Hi" readings


The GE Blood Glucose Meter 333 displays results between 20 and 600 mg/dL.

If the test result is below 20 mg/dL, "Lo" will appear on the screen. Please repeat your test with a new test strip.

If you still get a "Lo" result, contact your healthcare professional.

If the test result is above 600 mg/dL, "Hi" will appear on the screen. Please repeat your test with a new test strip.

If you still get a "Hi" result, contact your healthcare professional.



11. About Quality Control Test

What is a Quality Control Test?

To ensure proper monitoring function, it is necessary to regularly perform a quality control test. Use one of GE Control Solution 333 when testing your GE Blood Glucose Monitoring System 333 in the Control Solution Mode. If the test result is within the Control Solution Range printed on the strip vial label, the GE Blood Glucose Monitoring System 333 passes the Quality Control Test and your GE Blood Glucose Monitoring System 333 is working properly.

Control Solution Range:

Example of Control Solution Range printed on your test strip vial label.

	Control Solution Range mg/dL
	Level 2: 107-145
	Level 4: 315-426

NOTE

- Level 2 and Level 4 of GE Control Solution 333 are available for GE Blood Glucose Monitoring System 333 but must be purchased separately. If you want to purchase Level 2 or Level 4, please contact the GE Diabetes Customer Support Center.

When should a Quality Control Test Be Performed ?

- To ensure that your GE Blood Glucose Meter 333 and GE Blood Glucose Test Strip 333 are working properly.
- To confirm that you are following the correct testing procedures.
- To prepare for your initial blood glucose test.
- To check the GE Blood Glucose Test Strip 333 when you open a new vial of strips.
- To check your GE Blood Glucose Meter 333 after it has been dropped, damaged or exposed to liquids.
- If you suspect that your test results are inaccurate, or if your test results are not consistent with the way you feel.
- To practice testing.

Required Items for Quality Control Tests

To perform a quality control test, prepare the items below :

- GE Blood Glucose Meter 333
- GE Blood Glucose Test Strips 333
- GE Control Solution 333

12. Performing a Quality Control Test

1. Take one GE Blood Glucose Test Strip 333 from the vial and close the vial cap immediately.
2. Insert the GE Blood Glucose Test Strip 333 horizontally (not vertically) with the colored view window facing up into the test strip port.
3. While the blood drop icon is flashing on the display window, press and hold the Main button for at least 3 seconds until the "CS" symbol appears.
4. You will see the blinking "Lo" icon and "CS" icon on the screen, prompting you to apply GE Control Solution 333.
5. Shake the bottle of GE Control Solution 333 well before opening the cap. Place the cap on a flat surface.
6. Place a drop of control solution onto the top of the cap.
7. Gently touch the sample port of the strip with the control solution from the top of the cap.
8. The screen will display the count time starting from 5 (you will hear a beep if the volume is turned on).
9. Tightly replace the cap on the GE Control Solution 333 bottle.
10. The control solution result will appear. Compare your Quality Control Test result to the Control Solution Range printed on the GE Blood Glucose Test Strip 333 vial label.



CAUTION


- Your Control Solution Test results will not be included in the average calculations, however, they can still be recalled and viewed. The Control Solution Test result will be shown with the "CS" icon on the screen.
- The Control Solution Test should be conducted between 59-104°F (15-40°C) and in a location with 10-90% RH.
- Before "Lo" and "CS" appear, do not touch the control solution to the sample port on the strip. The GE Blood Glucose Meter 333 is performing an internal check. Touching the control solution to the sample port before prompted will result in

an error message: "Error" and " " and be accompanied by beeps (if volume is turned on).
 - Do not drip the control solution on to the sample port of the test strip directly. The reagent on the strip might leak into the bottle of control solution and may cause the degeneration of the control solution. This could contaminate the meter via the test strip port.
 - Keep the test strip port clean and dry. Clean immediately if the test strip port is stained or is overly exposed to moisture.
 - Do not touch the tip of the control solution bottle. If the tip is touched, clean with paper and pour out 2 drops of control solution.



13. Recalling Test Results and Average

The GE Blood Glucose Meter 333 is able to automatically store a maximum of 500 test results with time and date. If your meter has stored 500 results, the newest test result will replace the oldest one.



To recall your test memory, start the meter without a test strip inserted.

1. Press the Main button to switch from the Time Mode to Memory Mode. You will see the "mem" icon in the upper left corner of the display. The latest test result will appear with the number "1".
2. Use either the Left or Right button on the side to review all previous results with date and time. Results display from the most recent (Sequence no. "1") to the oldest (Sequence no. "500") in the lower right hand corner of the screen.
3. The Quality Control Test result can be recalled from the saved data. When " " appears under the data, the test was performed using the GE Control Solution 333. The control solution test result will not be used when recalling test average results.
4. Quick Searching: To view test results in sequential order automatically, enter the Memory Mode. Press the Right or Left button for 2 seconds. (The Right button is for viewing from the most current to the oldest test results. The Left button is for viewing from the oldest to the most current test results). Release the button to stop at any particular test result.

The GE Blood Glucose meter 333 provides several average test results. View the 1-day, 7-days, 14-days, 30-days, 60-days and 90-days average test results for better blood glucose monitoring.

1. Press the Main button twice to switch the screen to Average Mode.
2. On the average screen, use either the Left or Right button to view 1-day, 7-days, 14-days, 30-days, 60-days or 90-days averages. (Figure a)



- NOTE**
- The Average function requires that the correct time and date is set. Test results must exist during the desired time interval. For example: To get a 14-day average on 1/30/15, you must have test results dated between 1/17 and 1/30. If no test results exist during that time frame you will not receive an average.
 - The "Lo", "Hi" results, control solution test results will be excluded from average calculations.

14. Maintenance

Keep your meter and test strip free of dust, water or any other liquid. Store the meter in the carrying case when not in use. If your meter is dropped or damaged, perform a quality control test with the control solution before performing a blood glucose test.

Cleaning Meter

Clean the outside of the meter with disinfecting wipes listed on section 15. Do not get the meter test strip port wet.

Cleaning Test Strip Port

If your test strip port is stained with blood, control solution or any liquid, please use disinfecting wipes listed on section 15 to clean it immediately. Do not immerse in water. Perform a quality control test to ensure the GE Blood Glucose Meter 333 is working properly.

15. Cleaning and Disinfecting Procedures

Indirect transmission of Human Immunodeficiency Virus (HIV), Hepatitis B Virus (HBV) and Hepatitis C Virus (HCV) during the delivery of healthcare services has been increasingly reported. Persons using blood glucose monitoring systems have been identified as one risk group due to the shared use of fingerstick (lancing) devices and point of care blood testing devices. The cleaning procedure is to remove dust, blood and body fluid from the surface and should be performed whenever the meter or lancing device is visibly dirty. Performing the cleaning procedure once per week is recommended. The disinfecting procedure is necessary to kill pathogens such as HBV on the housing materials of the device. The cleaning and disinfecting procedure should be performed at least once per week. The meter and lancing device should be cleaned before being disinfected.

If the meter is being operated by a second person who is providing testing assistance to the user, the meter and lancing device should be cleaned and disinfected prior to use by the second person. NOTE: the cleaning procedure can only remove visible contaminants from surfaces. Only the disinfecting procedure can eliminate non-visible pathogens. CavWipes DISINFECTING TOWELETTES may be used to clean and disinfect the meter and lancing device. Users may purchase the disinfecting wipes from the manufacturer (Metrex), distributors or major on line retail sites such as www.amazon.com and www.ebay.com. You can contact Metrex customer service at 1-800-841-1428 and info@metrex.com (Monday through Friday, between 6 AM and 4 PM PST)

- NOTE**
- Clean and disinfect the outside of the device only. Do not remove battery cover when cleaning and disinfecting.

To clean the meter:

1. Thoroughly wipe the entire surface of the meter with disinfecting wipes listed above to clean any possible dirt, dust, blood and other body fluids.

To disinfect the meter:

2. Take another disinfecting wipe and wipe the meter thoroughly. (Note: All blood and body fluids should be cleaned from surface before performing the disinfecting procedure)
3. Allow the surface to remain wet for 2 minute.
4. Allow to air dry.


- NOTE**
- Your GE Blood Glucose Meter 333 has been tested to ensure that there is no change in the performance or external materials of the device after 550 cleaning cycles and 550 disinfecting cycles. The testing simulates 2 cleaning and disinfecting cycles per week over the typical life of the meter (5 years).

To clean the lancing device:

1. Thoroughly wipe the entire surface of the lancing device with disinfecting wipes listed above to clean any possible dirt, dust, blood and other body fluids.

To disinfect the lancing device:

2. Take another disinfecting wipe and wipe the Lancing device thoroughly. (Note: All blood and body fluids should be cleaned from surface before performing disinfecting procedure)
3. Allow the surface to remain wet for 2 minute.



- NOTE**
- Your GE Lancing Device has been tested to ensure that there is no change in the performance of the device after 550 cleaning cycles and 550 disinfecting cycles. The testing simulates 2 cleaning and disinfecting cycles per week over the typical life of the meter (5 years).

- CAUTION**
- Users should wash their hands thoroughly with soap and water after handling the meter, lancing device or test strips.
 - Please examine your LCD screen, test strip port, buttons and surface of your meter and lancing device after cleaning and disinfecting cycles. Stop using the meter and/or lancing device if any of the following occur:
 - Thin, silver streaks appear on the screen,
 - The screen becomes cracked, soft, dissolved, brittle or swollen.
 - You are unable to turn on/off your meter, operate the left/ right button, the lancing device release button or depth adjustable cap.
 - You are unable to enter meter settings, function modes or recall your testing results.
- If you have any questions or concerns, please contact your GE Blood Glucose Monitoring System 333 authorized representative or call the GE Diabetes Support Center.

16. Error Messages and Troubleshooting


Temperature Error

1. In order to get accurate testing results, please test between 50-104°F (10-40°C).
2. When the temperature is below 50°F (10°C) or over 104°F (40°C), the meter will not function and the "Error" symbol will blink.
3. If the GE Blood Glucose Meter 333 and Test Strips are exposed to a substantial change in temperature, please wait 45 minutes before measurement.




Battery Error

1. The " " symbol blinks when the battery power is low. Please change the battery as soon as possible. You can still test.
2. The " " and "Error" symbols blink when the battery is too low. Meter cannot test. Please change the battery immediately.



Strip Error


1. When the test strip is inserted incorrectly, the meter will beep 4 times while the " " symbol flashes. Do not apply a sample on the test strip. Please re-insert the unused test strip correctly (refer to section "5. The GE Blood Glucose Test Strip 333 and Battery Replacement" for more information).
2. If, after performing the required steps, the "Error" and flashing " " symbols appear on the screen, it means the GE Blood Glucose Test Strip 333 has been inserted incorrectly more than twice. Please, re-insert the unused test strip again after reviewing the steps on how to insert a test strip correctly or call Customer Service for support. This error message can also indicate that you may have used the wrong test strips. Please check the test strip vial to ensure you are using the GE Blood Glucose Test strips 333.



Sampling Error

Please do not apply the blood to the sample port of the test strip before the meter displays " ". Discard the test strip if the meter shows "Error".

- Er1-** The inserted test strip has been used or damaged. Please use a new test strip.
Er2- Th meter has malfunctioned. Do a Quality Control Test or reinstall the batteries to check if the meter works properly.
Er3- Th signal transmission is disrupted. Repeat the test using a new test strip.
Er4- The blood sample volume is insufficient. Repeat the test using a new test strip.
Er5- An issue calibrating the Meter has occurred. Please follow the steps below:
- (1) Remove the test strip from the meter.
 - (2) Turn off the meter (press and hold the Main button for 3 seconds).
 - (3) Press the main button to turn ON the meter.



If after turning on the Meter, you do not see Er5, your meter is functioning properly and able to perform a test. If one of the above error messages still appears, please contact GE Diabetes Customer Service at 1-866-613-7085.

Meter Malfunction

If the GE Blood Glucose Meter 333 will not turn on, please follow the steps below:

1. Open the battery cover and remove the batteries.
2. Wait for 5 minutes and reinsert the batteries as described as instructed on section "5. The GE Blood Glucose Test Strip 333 and Battery Replacement".

The meter should work normally after finishing the above steps. If the meter still does not work, please contact GE Diabetes Customer Support Center at 1-866-613-7085.

CAUTION

A blood sample should ONLY be applied to the test strip after the test strip has been inserted correctly and the Meter's screen is showing the image of the test strip and a blood drop flashing symbol. If the blood drop flashing symbol is NOT showing on the screen, do not apply a sample to the test strip. Please re-insert the unused test strip correctly. The screen MUST show the image of the test strip accompanied by the blood drop symbol flashing (this symbol take 3 seconds to show after correctly inserting the test strip) in order for a sample to be applied. Please consult your User Manual and/or contact Customer Service for support on how to correctly insert a GE Blood Glucose Test strip 333.

17. Limitations

- Blood glucose meter readings may be significantly lower than "true glucose levels" in a hyperglycemic-hyposmolar state, with or without ketosis.
- Caution is advised when glucose values are below 50 mg/dL or above 250 mg/dL. Consult a Physician as soon as possible, if values in this range are obtained.
- Healthcare professionals should evaluate their technique and their patients' technique regarding the use of the GE Blood Glucose Monitoring System 333 regularly. To accomplish this, it is recommended that blood glucose monitoring results be compared with a concurrently obtained laboratory measurement on the same blood sample. A proven clinical laboratory method employing hexokinase or glucose oxidase should be used as the comparative method.
- Hands and fingers contaminated with sugar from foods or beverages may cause falsely elevated results.
- The results of blood glucose measurements are different for measurements with whole blood and plasma.
- Storage of strips near bleach as well as bleach containing products will affect the results of the GE Blood Glucose Test Strips 333.
- Inaccurate test results may be obtained at altitudes greater than 10,000 feet (3,048 meters) above sea level.
- Hematocrits below 20% may cause higher results. Hematocrits above 60% may cause lower results.
- High concentrations of Uric acid ≥ 10 mg/dL, Cholesterol ≥ 600 mg/dL, and Ascorbic acid (Vitamin C) ≥ 5 mg/dL may interfere with the glucose test causing inaccurate test results.
- Do not perform the blood glucose test at temperatures below 10°C (50°F) or above 40°C (104°F), nor below 10% or above 90% relative humidity.
- GE Blood Glucose Test Strips 333 are designed for use with capillary whole blood samples. Do not use serum or plasma samples.
- Not for screening or diagnosis of diabetes mellitus.
- Not for use on critically ill patients, dehydrated patients or hyper-osmolar patients.
- Not for neonatal use.
- Alternate site testing should not be used to calibrate continuous glucose monitoring systems (CGMs).
- Alternate site testing should only be done during steady-state times (when glucose is not changing rapidly).
- Results from alternate site testing should not be used in insulin dose calculation.

18. Bluetooth function and data transmission

1. Bluetooth default is with " " symbol when mobile App turn on, it will automatically connect with your mobile device after user finished the pairing. Once connected with the App, test results will start syncing with the connected mobile device. When the sync is completed, the meter Bluetooth will disconnect automatically.



- 1) Bluetooth pairing. The symbol of Bluetooth " " appears.
- 2) Data uploading. The GE meter 333 Plus meter is now uploading test results to your mobile device.

NOTE

1. The meter will automatically turn off the Bluetooth after all the blood glucose readings have synced. In this case, the " " icon will disappear on the screen.
2. Animated " " symbol on the screen means the meter is transferring data.
3. No suggest to test blood sugar when the flight is strictly restrict the Bluetooth connection. Please turn off you meter to avoid the auto-on of Bluetooth.
4. When the meter and App cannot connect to each other successfully, please reboot the Bluetooth on both sides and try again. If the issues still exist, please contact your local distributors for further assistance.

3. Rightest CARE App (May not be available in your country)

"RIGHTEST Care" App is designed for diabetes patients which assistants to provide digital and mobility system for your blood glucose analysis. Quick understanding and managed your blood glucose efficiently. Search "RIGHTEST Care" with iOS or Android system. You may start with the function from pairing the glucose meter to App. More instruction of App, please click "Learn and Assist" for "Operation Guide" on the App. For more information, please visit the official website: <http://www.bionime.com.tw>

19. FEDERAL COMMUNICATIONS COMMISSION (FCC) STATEMENT

15.21
 You are cautioned that changes or modifications not expressly approved by the part responsible for compliance could void the user's authority to operate the equipment.
 15.105(b)
 This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with Part 15 of the FCC Rules.

Operation is subject to the following two conditions:

- 1) this device may not cause harmful interference and
- 2) this device must accept any interference received, including interference that may cause undesired operation of the device.

20. Specification

Measurement Technology	Oxidase Electrochemical Sensor
Sample	Capillary whole blood
Minimum Sample Volume	0.75 microliters
Measuring Range	20-600 mg/dL
Test Time	5 seconds
Test Strip Code Setting	Autocode
Memory Capacity	500 blood glucose test results with date and time
Power Saving	Turns off automatically 2 minutes after last user action. To turn off manually, press themain button " " for 3 seconds.
Operating Temperature	50 - 104°F (10 - 40°C)
Operating Relative Humidity	10 - 90%
Hematocrit	20 - 60%
Power Supply	Two AAA Batteries
Battery Life	Approximately 800 tests, each test with one data transfer
Meter Dimension	85mm x 57mm x 22.5 mm
Meter Weight	81.0 ± 5g with batteries
Monitor	LCD
Display Area	39 mm x 39.5 mm
Meter Storage/ Transportation Conditions	14 - 140°F (-10 - 60°C)
Test Strip Storage/ Transportation Conditions	39 - 86°F (4 - 30°C), 10 - 90% relative humidity

21. Warranty

Bionime Corporation warrants that this product will be free from defects in materials and workmanship for five years from the date of purchase. This warranty does not apply to the performance of a GE Blood Glucose Meter 333 that has been altered, misused, tampered with or abused in any way. This warranty applies only to the original purchaser of the GE Blood Glucose Monitoring System 333. Please complete and return the enclosed warranty card. Different models have different specifications. This warranty applies only to the GE Blood Glucose Monitoring System 333; other models are not covered with this warranty card. GE is a trademark of General Electric Company and is used under license by Bionime Corporation, No. 100, Sec. 2, Daqing St., South Dist., Taichung City 40242, Taiwan

NOTE

- During blood glucose measurement, the GE Blood Glucose Meter 333 itself may come into contact with blood. All parts of the GE Blood Glucose Monitoring System 333 are considered biohazardous and can potentially transmit infectious diseases. Please follow your healthcare professional's direction to properly dispose of the used GE Blood Glucose Monitoring System 333 after removing the batteries.

22. Customer Service

Please review all of the instructions to make sure you are performing the steps correctly. If you have any questions or concerns, please contact the GE Diabetes Support Center toll free at 1-866-613-7085 or visit our web site: www.gediabetes.com. (Monday through Friday 8:00 AM to 5:00 PM PST. If you have questions or need assistance outside the operational days and times, please contact your healthcare provider.)

23. Component Manufacturer Information

Blood Glucose Meter, Test Strip, Control Solution
 Manufacturer: Bionime Corporation
 No. 100, Sec. 2, Daqing St., South Dist., Taichung City 40242, Taiwan



Rev. Date:2019-01

Manufacturer:
 BIONIME CORPORATION
 No. 100, Sec. 2, Daqing St., South Dist., Taichung City 40242, Taiwan
 Made in Taiwan

GE Diabetes Customer Support Center:
 1450 E. Spruce Street, Bldg. #B,
 Ontario, CA 91761, USA
 Tel: 1-866-613-7085
www.gediabetes.com

GE is a trademark of General Electric Company and is used under license by Bionime Corporation, No. 100, Sec. 2, Daqing St., South Dist., Taichung City 40242, Taiwan