
# 30-Second Thermometer

Model: 59-309

**THANK YOU FOR PURCHASING A 30-SECOND THERMOMETER. PLEASE READ THESE INSTRUCTIONS CAREFULLY TO ENSURE PROPER USE AND SAFE OPERATION OF YOUR THERMOMETER.**

**PLEASE READ AND SAVE THESE INSTRUCTIONS**

- This product is not intended to diagnose or treat any health problem or disease and should not be used as a substitute for the consultation and advice of a physician or other medical professional.
- Thermometers should only be used under the supervision of an adult.
- Before use, it is recommended that the patient avoid eating, drinking, smoking or exercising for at least twenty minutes prior to taking a measurement. Never walk, run or talk during a measurement.

**NORMAL BODY TEMPERATURE / FEVER INDICATOR**

Body temperature can vary from 97°F and 100.0°F and still be considered normal. The medically accepted 'normal' body temperature is 98.6°F. Body temperature is commonly lower in the upon waking than at any point during the rest of an individual's waking hours.

Temperature readings will vary based on the body location point of the reading. Average 'normal' temperature ranges by measurement location site:

Oral	97.0°F–100.0°F (36.1°C–37.8°C)
Rectal	97.3°F–100.4°F (36.3°C–38.0°C)
Axillary	95.4°F–99.0°F (35.2°C–37.2°C)

A fever is defined as a body temperature that is elevated above the normal for that person. It is important to determine what is normal for an individual before determining if a fever is present. Tracking an individual's temperature on a consistent basis, at the same body site and at the same time everyday when the person is well or healthy will help determine an individual's normal temperature.

**CHANGING MODES - FAHRENHEIT OR CELSIUS OPERATION**

This unit takes measurements in both Fahrenheit and Celsius. While the unit is OFF, press and hold the POWER button. After the brief self-test, the unit will show the last reading in either °F or °C. Continue to hold the POWER button, and display will change the mode and display the new temperature reading mode. Release the POWER button and begin your reading.

**MEMORY RECALL**

This unit includes a recall of the last measurement taken. When the unit is OFF, press the **POWER** button to turn the unit on, and release. After the self-test, the unit will briefly show the last reading on the display.


**OPERATING INSTRUCTIONS**

Press and release the **POWER** button to activate the thermometer. The unit will beep and the display will show all of the digital display symbols as it performs a self-test before showing the last reading in memory.

NOTE: If using a probe cover, apply before turning device on.

**ORAL MEASUREMENTS — Approximate measurement time 30-40 seconds**

1. Place the probe under the tongue as near as possible to a heat pocket as noted in the diagram. Readings will be progressively lower the further the probe is from the indicated heat pockets.
2. Hold the thermometer in place during the reading; do not bite down on the thermometer. The mouth must remain closed to ensure accurate results.
3. During the measurement the °F(°C) will flash while the numbers go up on the display. When the reading is complete and the peak temperature has been reached, the °F will cease flashing and the results will appear on the display.
4. If the reading is below 100°F(37.7°C), the unit will sound ten long beeps; if the reading is at or above 100°F(37.7°C), the unit will sound ten sets of three short beeps.
5. Wait a minimum of one minute before attempting a new measurement.
6. To turn the unit off, simply press the **POWER** button. The unit will automatically shut-off after approximately ten minutes.

**AXILLARY (UNDERARM) MEASUREMENTS — Approximate measurement time 40-50 seconds**

1. Wipe the underarm with a dry towel before application.
2. Place the probe tip under the arm so that the tip is touching skin; hold the thermometer parallel, slightly angled upwards as indicated in the diagram. Position the arm across the chest to ensure the tip is completely covered.
3. During the measurement the °F(°C) will flash while the numbers go up on the display. When the reading is complete and the peak temperature has been reached, the °F will cease flashing and the results will appear on the display.
4. If the reading is below 100°F(37.7°C), the unit will sound ten long beeps; if the reading is at or above 100°F(37.7°C), the unit will sound ten sets of three short beeps.
5. Wait a minimum of one minute before attempting a new measurement.
6. To turn the unit off, simply press the **POWER** button. The unit will automatically shut-off after approximately ten minutes.

**RECTAL MEASUREMENTS — Approximate measurement time 20-30 seconds**

1. It is recommended that a probe cover is used for rectal measurements; apply a water-soluble lubricant for easier insertion.

2. Lay the patient on their side or stomach with their legs hanging down.
3. Gently insert the tip of the probe NO MORE THAN .5 INCH into the rectum. Stop immediately if you detect resistance.
4. Hold the thermometer in place during the reading.
5. During the measurement the °F(°C) will flash while the numbers go up on the display. When the reading is complete and the peak temperature has been reached, the °F will cease flashing and the results will appear on the display.
6. If the reading is below 100°F(37.7°C), the unit will sound ten long beeps; if the reading is at or above 100°F(37.7°C), the unit will sound ten sets of three short beeps.
7. Wait a minimum of one minute before attempting a new measurement.
8. To turn the unit off, simply press the **POWER** button. The unit will automatically shut-off after approximately ten minutes.

**FEVERLINE INDICATOR**


This thermometer includes our Feverline Indicator, the colored bars next to the digital display. After the measurement, the small triangle will appear with the results, offering a quick visual cue where the results fall—no fever, slight fever or fever. If the result is below 96.4°F (35.8°C), then there will be no indicator.

Temperature	Color	Result
96.4°F≤99.0°F (35.8°C≤37.2°C)	Green	No Fever Present
99.0°F≤100.0°F (37.2°C≤37.8°C)	Yellow	Slight Fever
100.0°F≤109.9°F (37.8°C≤42.9°C)	Red	Fever

**REPLACING THE BATTERY**

This thermometer includes 1 LR41 alkaline or SR41 silver oxide battery. If the battery is low or expired, the battery symbol will appear on the display; replace the battery as soon as possible.

1. Remove the battery cap by pulling gently.
2. Looking at the unit from the back side (ratings label face up) carefully slide the battery tray out so the entire battery is exposed.
3. Remove the battery with a non-metal tool. Do not attempt to remove the thin plastic battery cover bar.
4. Replace the battery with a LR41 alkaline or SR41 silver oxide type battery with the positive (+) side facing up.
5. Carefully replace the battery cap.
6. Do not dispose of used batteries in the trash. Consult your local recycling center for the proper way to dispose of batteries in your area. Keep used batteries away from children, pets and excessive heat.

**PROBE COVER USE**

Disposable, hygienic probe covers are recommended for use with stick thermometers to prevent contamination. Probe covers are intended for one-time use only; a new probe cover should be used with each application. Probe covers can be ordered using the attached form or at your local retailer.

**CARE AND CLEANING INSTRUCTIONS**

Store the unit in the provided protective case when not in use. Do not store the thermometer where it will be exposed to direct sunlight, dust or humidity. Avoid extreme temperatures. Do not attempt to disassemble the unit.

Before and after use—clean the thermometer by wiping with rubbing alcohol (70% isopropyl alcohol) diluted with water, or wash with cold, soapy water. Use only a mild, non-abrasive soap. Do not submerge in any liquid or autoclave. Wipe the thermometer with a dry cloth and allow to air dry.

Never sterilize by boiling, steam gas or autoclave.

**TROUBLESHOOTING / ERROR CODES**

Performance may be degraded should one or more of the following occur:

- Operation outside stated temperature and humidity range
- Storage outside state temperature and humidity range
- Mechanical shock
- Patient temperature is below ambient temperature

Error	Problem/Cause	Solution
Hi °F(°C)	The measured temperature is above 109.9°F (42.9°C)	Results are above range; wait several minutes and attempt new measurement.
Lo °F(°C)	The measured temperature is below 90.0°F (32°C)	Results are above range; wait several minutes and attempt new measurement.
Err	Unit malfunctioned	Attempt new measurement; if error persists contact Customer Care.
Blank Display	Unit will not power on	Battery is depleted; battery polarity may be incorrect; check battery and replace with a new battery if necessary.

**PRODUCT SPECIFICATIONS**

- BATTERY:** 1 LR41 or SR41 battery
- BATTERY LIFE:** Approximately 24 months; minimum of 1000 measurements
- AUTOMATIC SHUT-OFF:** Approximately after 10 minutes of non-use
- MEASUREMENT RANGE:** 90.0°F – 109.9°F (32° – 42.9°C)
- ACCURACY:** ±0.2°F; 95.9°F - 107.6°F (±0.1°C, 35.5°C - 42.0°C); ±0.4°F (0.2°C) outside of that range
- OPERATING TEMPERATURE:** 41°F – 104.0°F (5°C – 40°C); Humidity 15-95% (non condensing)
- STORAGE TEMPERATURE:** -4.0°F – 131.0°F (-20°C – 55°C); Humidity 15-95% (non condensing)
- MEASUREMENT SCALE:** Fahrenheit and Celsius
- MEMORY:** Recall of last measurement taken
- DISPLAY:** Liquid crystal display
- DIMENSIONS:** 4.9" X 1.7" X .35"
- WEIGHT:** .02 ounces with batteries

Specifications are subject to change without notice.

**FCC STATEMENT**

IMPORTANT INFORMATION REQUIRED BY THE FCC

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference; and (2) this device must accept any interference received, including interference that may cause undesirable operation.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy, and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

CAUTION: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**LIFETIME LIMITED WARRANTY**

Congratulations on your purchase of a 30-Second Digital Thermometer. Your Thermometer is covered by the following limited warranty commencing upon the date of purchase, and subject to the following terms and conditions:

The Warrantor warrants that its Thermometer will be free from defects in materials and workmanship under normal consumer usage for the life of the product for the original purchaser of the Thermometer. Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the Thermometer for commercial purposes or subjecting the Thermometer to abnormal usage or conditions; or (d) other acts which are not the fault of the Warrantor, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the Thermometer.

If the Digital Thermometer fails to conform to this limited warranty, return the Thermometer to: **Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031**. When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone number. Carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. At our option, the warrantor will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will be responsible for parts charges, if any, and repair charges not covered by this limited warranty.

EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL THE WARRANTOR BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

**INSTRUCTIONS • ENGLISH**

**PROBE COVER ORDER FORM**


Disposable Probe Covers 100 Pack (Model 08-322) Qty. \_\_\_\_\_ @ \$5.00 ea. \$ \_\_\_\_\_

SUBTOTAL \$ \_\_\_\_\_

Illinois shipments add 8% Sales Tax \$ \_\_\_\_\_

Shipping, handling & insurance \$ 5.00

ORDER TOTAL \$ \_\_\_\_\_

**PAYMENT METHOD:**

Check or Money Order  VISA  MasterCard

Discover  Diners Club  American Express

Card No. \_\_\_\_\_ Security Code \_\_\_\_\_

Expiration Date: \_\_\_\_\_

Signature \_\_\_\_\_

Print Name \_\_\_\_\_

Please allow 2 weeks for delivery; this order form is for shipment within the Continental U.S. and Canada only.

**Thank you for your order!**  
Please send a copy of this form with your payment to:  
Veridian Healthcare, LLC  
Attn: Customer Service  
1175 Lakeside Drive  
Gurnee, IL 60031

**SHIP TO ADDRESS:**

Name \_\_\_\_\_

Address \_\_\_\_\_

City \_\_\_\_\_ State \_\_\_\_\_ Zip \_\_\_\_\_

Daytime Telephone No. (\_\_\_\_) \_\_\_\_\_

**BILL TO ADDRESS: If different from SHIP TO address**

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## PROBE COVER ORDER FORM

Disposable Probe Covers 100 Pack (Model 08-322) Qty. \_\_\_\_\_ @ \$5.00 ea. \$ \_\_\_\_\_  
 SUBTOTAL \$ \_\_\_\_\_  
 Illinois shipments add 8% Sales Tax \$ \_\_\_\_\_  
 Shipping, handling & insurance \$ 5.00  
 ORDER TOTAL \$ \_\_\_\_\_

**PAYMENT METHOD:**  
 Check or Money Order  VISA  MasterCard  
 Discover  Diners Club  American Express  
 Card No. \_\_\_\_\_ Security Code \_\_\_\_\_  
 Expiration Date: \_\_\_\_\_  
 Signature \_\_\_\_\_  
 Print Name \_\_\_\_\_

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### Thank you for your order!

Please send a copy of this form with your payment to:

Veridian Healthcare, LLC  
 Attn: Customer Service  
 1175 Lakeside Drive  
 Gurnee, IL 60031

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FAX CREDIT CARD ORDERS TO: 866-480-1717; QUESTIONS? 866-326-1313