SMARTHEART®

Talking Blood Pressure Arm Monitor

Model # 01-746T


INSTRUCTION MANUAL **ENGLISH AND SPANISH**

Please read this instruction manual
completely before operating this unit.


STOP!
PLEASE ENSURE YOU HAVE ALL OF
THE FOLLOWING COMPONENTS
BEFORE USING YOUR DIGITAL
BLOOD PRESSURE MONITOR.


Digital Monitor

Wide-Range Cuff

4 AAA
Batteries

Instruction Manual
Quick Start Guide
Blood Pressure Reading Tracking Log
FAQ Information Sheet

**IF YOU ARE MISSING ANY PARTS, INCLUDING
INSERTS OR INSTRUCTION MANUALS, DO NOT
RETURN TO PLACE OF PURCHASE. CONTACT
CUSTOMER CARE AT 866-326-1313.**

INDEX

Care & Safety Information	4-6
Introduction & Indications for Use	7
Blood Pressure Monitor Features	8-9
Getting Started	10
Setting Language/Volume	11
Setting Date/Time	12
Battery Installation	13
Fitting & Applying Your Cuff	14-15
Taking Your Blood Pressure Reading	16-17
Interpreting Your Results	18-19
Memory Function	20-21
Care & Maintenance	22-23
Display Symbols	24
Troubleshooting	25
Error Codes	26
Device & Label Symbols	27
FCC Statement	28
Electromagnetic Compatibility	29-33
Specifications	34
Warranty	35
Instrucciones en Español	37-71
Reference Page	72

Toll-Free Customer Care Help
Line: 1-866-326-1313
Monday – Friday
8:30 a.m. – 4:30 p.m. CST

SmartHeart®
Manufactured for
Veridian Healthcare
1175 Lakeside Drive
Gurnee, IL 60031

Made in China

#93-1629 12/2024
©2024 Veridian Healthcare, LLC

CARE & SAFETY INFORMATION

NOTE: Read all instructions carefully before use. The following basic precautions are needed when using an electrical product.

CAUTION: Failure to read and observe all precautions could result in personal injury or equipment damage.

Improper care or use of your blood pressure monitor may result in injury, damage to the unit or ineffective treatment. Following these instructions will ensure the blood pressure monitor's efficacy and long life.

GENERAL CAUTIONS AND WARNINGS

- Measurement results should be discussed with your physician or healthcare professional; never self-diagnose or attempt treatment as this may be dangerous.
- The Electronic Sphygmomanometer is intended to measure the systolic, diastolic and mean blood pressure as well as pulse rate via non-invasive oscillometric technique in which an inflatable cuff is wrapped around the upper arm. It can be used on adult individuals.
- Keep out of reach of children.
- The patient is the intended operator of this device.
- Do not use the unit if it has any damaged parts (including power cord or plug and cuff), if it has been submersed in water or dropped.
- If any abnormality occurs, discontinue use until the unit has been examined and repaired.
- Only use included attachments and parts; do not use attachments from other brands or models as these may not be properly calibrated for use with this device and may result in measurement error.
- Do not disassemble the unit or attempt to repair it; substitution of a component different from that supplied may result in measurement error and will void manufacturer warranty.
- Always follow local regulations for proper disposal of the monitor, cuff and batteries.
- Swallowing batteries and battery fluid is extremely dangerous; keep

CARE & SAFETY INFORMATION

batteries and unit out of reach of children and disabled persons.

- Do not use this device with other diagnostic monitoring equipment on the same limb.
- Do not use near active HF Surgical Equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of the EM Disturbances is high.
- Do not share this device with infected persons where there may be a risk of cross-contamination.
- Consult with your physician before if the cuff may be applied over a wound or inflammation diseases.
- It is inappropriate for people with serious arrhythmia to use this Electronic Sphygmomanometer.
- Consult with your physician before use if any of the following conditions are present: arrhythmia, common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation, arterial sclerosis, poor perfusion, diabetes, pre-eclampsia, renal diseases. Inflammatory diseases, intravascular therapy, arteriovenous shunts or if a mastectomy has been performed.
- Motion, trembling and shivering may affect measurement reading.
- For information regarding potential electromagnetic or other interference between digital monitor and other devices and advice on avoidance of such interference, consult with the Electromagnetic Compatibility portion of this manual.
- The blood pressure monitor is not intended to be exposed to the Electromagnetic Interference (EMI) environment.  Please do not use the blood pressure monitor within the environment of the following device: Magnetic Resonance Imaging (MRI), computerized axial tomography (CT), diathermy, Radio Frequency Identification (RFID), and electromagnetic security systems such as metal detectors.
- Monitor should be kept at least 12" from other wireless devices during use.
- If you are allergic to plastic or rubber, consult with your medical professional before using this device.

CARE & SAFETY INFORMATION

OPERATING CAUTIONS AND WARNINGS

- Prolonged over-inflation may cause congestion, swelling or bruising in some people.
- Too-frequent measurements may cause injury due to blood flow interference.
- Check the circulation of the user before, during and after use.
- If you experience any discomfort or an abnormality occurs, stop using the device immediately; to stop the device during operation, press the **POWER** button and the cuff will automatically deflate.
- Operation outside of stated operating temperature may result in measurement error or device malfunction; operation environment temperature is: 41°F – 104°F (5°C – 40°C); Humidity: 15~90% RH max (no condensation).
- Never use this unit while operating a vehicle or in the bathtub or shower.

STORAGE CAUTIONS AND WARNINGS

- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: -4°F – 131°F (-20°C – 55°C); Humidity, 10~93% RH max.
- Keep the unit out of reach of children.
- Always keep the unit unplugged while not in use; remove the batteries if the unit will not be used for an extended period of time.

CLEANING CAUTIONS AND WARNINGS

- Never immerse the unit in water to clean as it may damage the unit.
- Follow the 'Care & Maintenance' portion of this manual for instruction on how to clean and care for your monitor.

INTRODUCTION & INDICATIONS FOR USE

It is recommended that you first seek the advice and recommendation of your physician or healthcare professional when using home diagnostic devices, including blood pressure monitors.

Automatic digital blood pressure monitors use the oscillometric method to electronically measure your blood pressure. The monitor detects your blood's movement through the artery in your arm and converts the movements into a digital reading. The oscillometric method does not require a stethoscope, making the monitor ideal for home use.


Thank you for purchasing an Automatic Digital Blood Pressure Monitor. With proper care and use, your monitor will provide you with many years of reliable readings.

**PLEASE READ THIS INSTRUCTION
MANUAL COMPLETELY BEFORE
OPERATING THIS UNIT.**

BLOOD PRESSURE MONITOR FEATURES

BLOOD PRESSURE MONITOR FEATURES

DETAILS OF THE DIGITAL DISPLAY SCREEN AND INTERPRETATION OF THE SYMBOLS ARE INCLUDED THROUGHOUT THIS MANUAL.

GETTING STARTED

Before Taking a Measurement

1. Your readings should only be interpreted by your physician or healthcare professional with access to your individual medical history. Regular use of a home blood pressure monitor will allow you to track and record your readings for discussion with your physician.
2. Conduct your measurement in a quiet place while seated in a relaxed position. Rest for 15 minutes before taking your reading. If an error occurs or you desire to take a second reading, allow 15 minutes between readings for your blood vessels to return to normal.
3. Avoid smoking, eating, taking medication, alcohol, physical activity or any other stressful activity for 30 minutes prior to taking a reading.
4. Always remove any jewelry or constrictive clothing that may interfere with the cuff placement.
5. Keep yourself and the monitor still during measuring; do not talk during the reading.
6. It is recommended that you take your readings at the same time each day to better monitor any indications in your results.
7. Record your daily measurement on the included chart or some other written document to share with your physician.

BATTERY INSTALLATION


This unit requires 4 AAA batteries.

It is necessary to replace the batteries when the battery life indicator shows a "Low Battery" or when the display does not turn on after the **POWER** button is pressed.

The battery compartment is located on the underside of this monitor.

1. Remove the cover to the battery compartment.
2. Insert or replace 4 AAA batteries into the battery compartment, ensuring to match the indicated polarity symbols. Always use new batteries.
3. Replace the battery cover.
4. Dispose of batteries according to local disposal and recycling regulations.

It is recommended to remove the batteries if the unit will not be used for an extended period of time.

The 'Low Battery' symbol will appear on the lower right of the display.

SETTING LANGUAGE/VOLUME

IT IS NECESSARY TO SET PREFERENCES - LANGUAGE, VOLUME, DATE, TIME - WHEN FIRST SET UP AND WHEN BATTERIES HAVE BEEN REPLACED.


USER BUTTON

MEMORY BUTTON


LANGUAGE SELECTION

1. When the unit is **OFF**, press and hold the **USER** button to enter the Audio Volume selection.
2. **LANGUAGE**—The current language in use will flash on the display EN (English) or ES (Spanish). Press the **MEMORY** button to toggle between selections. Press the **USER** button to confirm.

AUDIO VOLUME ON / OFF


3. **VOLUME**—VL will flash VL (volume) 1, 2 or 0; pres **MEMORY** to toggle between selections; press the **USER** button to confirm.
0 – No volume / audio **OFF**.

SETTING DATE/TIME

DATE TIME SETTING

4. **YEAR**—first the year will flash on the display. Press the **MEMORY** to advance the display to the desired year, press the **USER** button to confirm your selection and advance to the date settings.
5. **MONTH/DATE**—repeat the same sequence of step 2 to confirm the MONTH and DATE setting.
6. **HOURS/MINUTES**—repeat the same sequence of step 2 to confirm the HOUR and MINUTE setting.
7. After the Date/Time setting function is complete, the monitor is ready for use. The date/time feature is on the top of the display screen. The date and time will alternate while the device is in use.
8. This unit stores the date and time, including year allowing for a date/time stamp on all readings in the memory recall function.


FITTING & APPLYING YOUR CUFF

PROPER CUFF FIT AND PLACEMENT IS CRITICAL IN OBTAINING ACCURATE BLOOD PRESSURE MEASUREMENTS. PLEASE READ THESE INSTRUCTIONS CAREFULLY AND CONSIDER CONSULTING WITH YOUR PHYSICIAN FOR FURTHER DIRECTION.

Before applying your blood pressure cuff, be sure you have selected the appropriate size cuff.

To determine the proper cuff size, measure the circumference of your left upper arm at the point midway between inside bend of your elbow and your shoulder.

Wide-Range Adult Cuff fits arm circumference 8.6" – 16.5" (22–42 cm)

If for any reason you are unable to or should not use your left arm, please modify the instructions for cuff application to your right arm. Your physician can help you identify which arm is best for you to take measurements from.

Assembling Your Cuff

Your cuff may need to be 'threaded' through the metal bar to form the cuff cylinder for your arm.

1. Lay the open cuff hook and loop side down, size/instruction wording side up on a flat surface.
2. Keeping the words on the inside of the cuff, thread the end of the cuff through the metal bar, Fig 1. The end of the cuff should be able to fold over onto itself for securing via the hook and loop closure, Fig 2.


Fig. 1


Fig. 2

FITTING & APPLYING YOUR CUFF

Applying Your Cuff

Be seated at a table or desk with your feet flat on the floor. To avoid accidentally pulling on or dropping your monitor, the cuff should not be plugged into the monitor until after the cuff is applied to your arm.

The cuff is printed with a useful guide for applying the cuff; please consult the following steps for additional direction.

1. Remove any constrictive clothing or jewelry that may interfere with cuff placement.
2. Slide the cuff cylinder onto your left arm and align the artery mark over the brachial artery on the inside of your arm, Fig 3.


Fig. 3

The tube should follow the inside length of the arm, Fig 4.

3. The bottom edge of the cuff should be positioned approximately 0.5" above the elbow joint, Fig 5.
4. When the cuff is closed around the arm, the triangle guide imprint should be within the OK range bar marking; if the triangle is outside of the range, a different sized cuff is needed.


Fig. 4

5. The cuff should fit comfortably, yet snugly around your arm, Fig 6. You should be able to insert one finger easily between your arm and the cuff.
6. During use, avoid compression or squeezing of the air hose/tube, as this may cause inflation error or injury due to excessive cuff pressure.


Fig. 5


Fig. 6

TAKING YOUR BLOOD PRESSURE READING

Please read the preceding portions of this manual prior to taking your first reading.

Reminder, it is important to avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading.

Select the desired user, 1 or 2, before measurement.


Fig. 1

When the unit is off, press and release the **USER** button to toggle between users. Press **POWER** to confirm selection. All readings will be stored in the selected user memory bank.

1. Position the monitor on a flat, stable surface with the digital display in view. Be seated with your feet flat on the floor and legs uncrossed.
2. After applying the cuff insert the cuff tubing connector into the port on the left side of your monitor.
3. Rest your elbow on a solid surface with your palm facing upward. Elevate your arm so that the cuff is at the same level as your heart, Fig 1. Relax your left hand.
4. Press and release the **POWER** button to turn the unit on; the unit will briefly flash a self-test, Fig 2.


Fig. 2

YOU CAN STOP THE INFLATION OR DEFLATION PROCESS AT ANY TIME BY PRESSING THE POWER BUTTON.

EXCESSIVE MOVEMENT symbol – Icon will flash if excessive movement is detected during reading. Remain still and start a new measurement.

TAKING YOUR BLOOD PRESSURE READING

5. The cuff will automatically inflate; you will see the numbers increase on the display screen and feel the cuff inflate.
6. When the measurement is complete, the cuff will automatically deflate, then your blood pressure measurement and pulse reading results will appear on the display screen at the same time, Fig 3.
7. The Hypertension Indicator will indicate where your reading falls within hypertension standards, Fig 4. For more information see Interpreting Your Results.
8. If an irregular heartbeat was detected during the reading, the Irregular Heartbeat Detection indicator will appear on the display, Fig 5. For more information see Interpreting Your Results.
9. The reading will automatically be stored in memory, up to 199 readings.
10. Select **POWER** to turn the unit off and conserve energy and battery life. The unit will automatically shut-off approximately 1 minute after use.
11. Disconnect the tubing prior to storage.


Fig. 3

Fig. 4


Fig. 5

OK symbol – indicates the cuff placement was correct during measurement.

UP ARROW / INFLATION symbol – Will appear during cuff inflation.



DOWN ARROW / INFLATION symbol – Will appear during cuff deflation; may appear before measurement while air is released from the cuff.

INTERPRETING YOUR RESULTS

Irregular Heartbeat Indicator ♥)))

An irregular heartbeat is defined as a heartbeat rhythm that is more than 25% slower or 25% faster than the average measurement detected while the monitor is conducting the measurement. When this device detects an irregular rhythm two or more times during the reading, the IHB symbol will appear on the display.

This monitor will operate if the icon appears, but the results may be compromised, especially if this symbol appears often. Please consult with your physician or trained healthcare professional for further information regarding an irregular heartbeat and if this symbol appears frequently.

INTERPRETING YOUR RESULTS

Hypertension Indicator

This unit features our unique Hypertension Indicator. The color bars on the right side of the monitor display correspond with an icon on the digital display, indicating where the measurement results fall within the World Health Organization standards.

The World Health Organization has established globally accepted standards for the assessment of high or low blood pressure readings. The below chart should be considered only as a guideline, always consult with your physician or health care professional to interpret your individual results. Never adjust your medication(s) dosage, or severely alter your diet or exercise routine without consulting with your physician.

GREEN BAR

SYS: 90~139 mmHg or DIA: 60~89 mmHg

YELLOW BAR

SYS 140~159 mmHg or DIA 90~99 mmHg,

SYS <90 mmHg and DIA < 60 mmHg

RED BAR

SYS $\geq$ 160 mmHg or $\geq$ 100 mmHg

MEMORY FUNCTION

Recalling Measurements in Memory

Select the desired user, 1 or 2, before memory review.

When the unit is off, press and release the **USER** button to toggle between users. Press **POWER** to confirm selection.

This monitor stores and recalls up to 199 readings per user. These results can be shared with your physician or trained healthcare professional.

1. When the unit is OFF, press and release the **MEMORY** button to enter the memory recall function.
2. The display will first show the average of the last three measurements.
3. Press and release **MEMORY** to begin scrolling through the saved measurements.
4. Measurements will appear on the display from most current to oldest; the memory number will appear at the bottom left of the display.
5. All results for a given measurement will display, including measurement results, pulse rate, Hypertension Indicator and Irregular Heartbeat alert (if applicable), Fig 2.
6. When the number of readings exceeds 199, the oldest data will be replaced with the new record
7. Press **POWER** to turn the monitor off at any time during review of the stored measurements.


Fig. 1


Fig. 2

MEMORY FUNCTION

Clearing Measurements From Memory

To delete all results in memory:

1. Select the desired memory bank to delete.
2. Press **MEMORY** to enter Memory Recall mode.
3. Press and hold **MEMORY** for several seconds; the display will show NO on the lower right of the display.
4. The memory bank is now empty.

CARE & MAINTENANCE

Proper care and maintenance of your blood pressure monitor will help ensure long life for the product; improper use of the monitor may void the included warranty.

Always disconnect the cuff and tubing from the monitor prior to cleaning and storage.

Monitor Care & Maintenance

- Do not use any liquids on the monitor; use a soft, dry cloth to clean the monitor as needed.
- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: -4°F – 131°F (-20°C – 55°C); Humidity: 10%~93%RH (non- condensing).
- This monitor requires 6 hours to cool from the maximum storage temperature between uses until the monitor is ready for its INTENDED USE when the ambient temperature is 68°F(20 °C).
- This monitor requires 6 hours to warm from the minimum storage temperature between uses until the monitor is ready for its INTENDED USE when the ambient temperature is 68°F(20 °C).
- Avoid high temperatures and direct sunlight.
- Always remove the batteries if the unit will not be used for an extended period of time.
- Do not attempt to disassemble your monitor; disassembling your monitor may void the manufacturers warranty.
- Do not attempt to service the monitor while in use or on a patient.
- Do not subject the monitor to strong shocks; take care not to drop the monitor.

CARE & MAINTENANCE

- No component can be maintained by user in the monitor. The circuit diagrams, component part lists, descriptions, calibrations instructions, or other information which will assist the user's appropriately qualified technical personnel to repair those parts of equipment which are designated repairable can be supplied.

Cuff Care & Maintenance

- The cuff may be spot cleaned with a mild-detergent as needed.
- To disinfect the cuff, wipe the inside (skin contact side) of the cuff using a cloth lightly moistened with 70% - 90% ethyl rubbing alcohol and allow to air dry; it is suggested to disinfect the cuff frequently when used in clinical settings or if the cuff is contaminated in any way.
- Never attempt to iron or dry the cuff; if the cuff becomes wet, allow the cuff to air dry completely before use or storage.
- Do not attempt to remove the tube and bladder from the nylon cuff as this may damage the cuff.
- Take care not to forcefully bend or kink the tubing during storage.
- It is recommended the performance should be checked every 2 years or after a repair. Please contact Customer Care at 866-326-1313.
- It is recommended to clean the cuff every 200 uses.
- The monitor can maintain the safety and performance characteristics for minimum of 10,000 measurements of three years.
- The cuff integrity is maintained after 1,000 open- close cycles of the closure.

DISPLAY SYMBOLS

Symbol	Meaning
SYS	Systolic Blood Pressure—The 'top' blood pressure result
DIA	Diastolic Blood Pressure—The 'bottom' blood pressure result
PUL/min	Pulse Rate—Pulse rate per minute during the measurement
mmHg	Millimeters of Mercury—Unit of blood pressure measurement for this device
	Battery—Battery life indicator
	Irregular Heartbeat Symbol—An irregular heartbeat was detected during the measurement
	Memory—The display shows the measurements stored in memory
	Excessive Movement—User has moved excessively during measurement and results may be compromised
	Cuff Placement Indicator—OK shows when cuff is wrapped correctly during measurement
	User Icon—Indicates which user is in use and which memory bank reading results will be stored.
	UP/DOWN Arrows—Will show on display during inflation or deflation of the cuff

TROUBLESHOOTING

Problem	Possible Cause	Solution
Unit does not power on	Batteries may be expired or installed improperly.	Check the polarity of the batteries and reinstall if necessary; replace batteries with new batteries.
Cuff fails to inflate or inflate fully	Air tube may not be fully inserted or may be twisted	Ensure tube is inserted firmly into monitor; ensure tube is not creased or compromised for smooth air flow
	Air tube may be broken or leaking	Contact Customer Care to order a replacement cuff
Error codes persist; won't allow measurement	Excessive movement or talking during reading	Review the instructions for use and attempt new measurement; do not move or talk during readings
Air leakage of the cuff	Cuff may be too loose	Reapply the cuff and attempt new measurement
	The airbag may be damaged	Contact Customer Care to order a replacement cuff

ERROR CODES

Code	Meaning
Er U	Pressure failed to reach 30 mmHg in 12 seconds
Er H	Inflation reaches 295 mmHg and deflates automatically after 20 ms
Er 1	Pulse rate was not detected properly
Er 2	Excessive movement resulting in measurement failing to complete
Er 3	Measurement result is abnormal
Er 23	SYS value is lower than 57 mmHg
Er 24	SYS value is higher than 255 mmHg
Er 25	DIA value is lower than 25 mmHg
Er 26	DIA value is higher than 195 mmHg

DEVICE & LABEL SYMBOLS

These symbols may appear on your device, instructions or packaging and may vary by make and model.

Symbol | Meaning

	Read This Manual—All included manuals should be read prior to first use. (Background color will be blue; reading symbol will be white)
	Warning—Symbol indicates a warning, prohibition or mandatory action that mitigates a risk that is not necessarily obvious to the device operator.
	Type BF Applied Parts—Indicates that a part of this unit comes in contact with the patient in order to carry out it's intended function; in the case of this device the cuff is the Type BF Applied Part.
	Environment Protection—dispose of this product properly; consult with your local recycling ordinances for proper recycling and disposal.
IP21	Enclosure protection grade
SN	Serial Number.
	Direct Current.
	Date of Manufacture.
	Unsafe

The Electronic Sphygmomanometers corresponds to the below standards:
ISO81060-2 : 2013(Non-Invasive Sphygmomanometers- Part 2: Clinical Validation Of Automated Measurement Type).

FCC STATEMENT

NOTE:

POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for U.S.A. only)

This product has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. The product generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If the product does cause harmful interference to radio or television reception, which can be determined by turning the product on and off, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the product and the receiver.
- Connect the product into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

ELECTROMAGNETIC COMPATIBILITY

Appendix 1 EMC Information

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration - electromagnetic emissions

The Arm Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Arm Blood Pressure Monitor should assure that it is used in such an environment.

Phenomenon	Compliance	Electromagnetic Environment Guidance
RF emissions CISPR 11	Group 1	The Arm Blood Pressure Monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Arm Blood Pressure Monitor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC61000-3-2	N.A.	
Voltage fluctuations/flicker emissions IEC61000-3-3	N.A.	

ELECTROMAGNETIC COMPATIBILITY

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration - electromagnetic emissions

The Arm Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Arm Blood Pressure Monitor should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrical fast transient/burst IEC 61000-4-4	Not Applicable	Not Applicable
Surge IEC 61000-4-5	Not Applicable	Not Applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Not Applicable	Not Applicable
Power frequency Magnetic field IEC 61000-4-8	30A/m, 50/60Hz	30A/m, 50/60Hz
Conducted RF IEC 61000-4-6	Not Applicable	Not Applicable
Radiated RF IEC 61000-4-3	10 V/m 80 MHz - 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz - 2,7 GHz 80 % AM at 1 kHz

NOTE: UT is the a.c. mains voltage prior to application of the test level

ELECTROMAGNETIC COMPATIBILITY

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration - electromagnetic emissions

The Arm Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Arm Blood Pressure Monitor should assure that it is used in such an environment.

Radiated RF0-4.3 (Test specifications for ENCLOSURE/PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Max Power (W)	Distance (m)	IEC 60601-1-2 Test Level (V/m)	Compliance level (V/m)
	385	380-390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27	27
	450	460-470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28	28
	710 745 780	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9	9
	810 870 930	800-960	GSM 800/900, TETRA 800, DEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28	28

ELECTROMAGNETIC COMPATIBILITY

Radiated RF0-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Max Power (W)	Distance (m)	IEC 60601-1-2 Test Level (V/m)	Compliance level (V/m)
	1720 1845 1970	1700- 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28	28
	2450	2400- 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28	28
	5240 5500 5785	5100- 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9	9

ELECTROMAGNETIC COMPATIBILITY

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration - electromagnetic emissions

The Arm Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Arm Blood Pressure Monitor should assure that it is used in such an environment.

Radiated RF IEC61000-4-39 (Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields)	Test Frequency	Moduation	IEC 60601-1-2 Test Level (A/m)	Compliance level (A/m)
	30 kHz	CW	8	8
	134.2 kHz	Pulse modulation 2.1 kHz	65	65
	13.56 kHz	Pulse modulation 50 kHz	7.5	7.5

NOTE:

- Optional AC adaptor should comply with the requirement of IEC 60601-1 standard.
- Use only the exclusive AC adapter specified by authorized dealers. Other AC adapter may vary in output voltage and polarities and may represent a risk on your life and damaging the device.

PRODUCT SPECIFICATIONS

Name	Talking Blood Pressure Arm Monitor
Model Number	01-746T
Display System	Digital display/LCD
Measuring Method	Oscillometric Method, Automatic Air Inflation and Measurement
Power Source	4 x 1.5V  AAA Batteries;
Measuring Range	Systolic pressure: 57~255mmHg, Diastolic Pressure: 25~195mmHg; Pulse: 40-199 beats/minute
Accuracy	Pressure: ± 3 mmHg; Pulse Rate: $\pm 5\%$
Deflation	Automatic pressure release valve
Memory	2 memory banks; 199 measurements each
Automatic Power-Off	Approximately 1 minutes after last button operation
Battery Life	Approximate measurements: 460 in non-talking mode; 290 English; 260 Spanish
Operation Environment	Temperature 41°F – 104°F (5°C – 40°C); Humidity 15~90% RH max (no condensation)
Storage and Transport Environment	Temperature -4°F – 131°F (-20°C – 55°C); Humidity 10~93% RH max (no condensation)
Monitor Dimensions	3.85" x 3.75" x 2.75" (118 x 98mm x 62.5mm)
Monitor Weight	7.8 oz. (220g) (without batteries)
Arm Circumference Range	Wide Range size cuff fits arm circumference: 8.6" – 16.5" (22–42 cm)
Accessories	Instruction Manual, Quick Start Guide, Blood Pressure Log, FAQ

Specifications are subject to change without notice.

1-YEAR LIMITED WARRANTY

Congratulations on your purchase of a Digital Blood Pressure Monitor. Your Digital Blood Pressure Monitor is covered by the following limited warranty commencing upon the date of purchase, and subject to the following terms and conditions:

The Warrantor warrants that its Digital Blood Pressure Monitor will be free from defects in materials and workmanship under normal consumer usage for a period of one year for the original purchaser.

Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the Digital Blood Pressure Monitor for commercial purposes or subjecting the Digital Blood Pressure Monitor to abnormal usage or conditions; or (d) other acts which are not the fault of the Warrantor, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the Digital Blood Pressure Monitor.

If the Digital Blood Pressure Monitor fails to conform to this limited warranty, return the Digital Blood Pressure Monitor postage prepaid to: **Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031**. When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone number. Carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. At our option, the Warrantor will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will be responsible for parts charges, if any, and repair charges not covered by this limited warranty.

EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL THE WARRANTOR BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

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