Upper Arm Talking Blood Pressure Monitor

Model UAM-901 Wide Range Cuff

INSTRUCTION MANUAL

ENGLISH & SPANISH


Please read this instruction manual
completely before operating this unit.


STOP!
**PLEASE ENSURE YOU HAVE ALL OF THE
FOLLOWING COMPONENTS BEFORE USING
YOUR BLOOD PRESSURE MONITOR.**

Monitor

Wide-Range Cuff

Four (4) AAA
Batteries

Instruction Manual

Toll-Free Customer Care Help Line:

1-866-326-1313

Monday – Friday

8:30 a.m. – 4:30 p.m. CST

www.veridianhealthcare.com

INDEX

Introduction	4
Care and Safety	5-10
Blood Pressure Monitor Features	11-12
Battery Installation	13
Setting Date & Time	14-15
Setting Language	16
Fitting & Applying Cuff	17-18
Taking Your BP Reading	19-20
Memory Function	21-22
Care & Maintenance	23-24
Hypertension Indicator	25
Irregular Pulse Rate Indicator	26
Troubleshooting	27
Device & Symbols	28
FCC Statement	29
Electromagnetic Compatibility	30-32
Specifications	33
Warranty	34
Español	35-70

Zewa®
Manufactured for
Veridian Healthcare
1175 Lakeside Drive
Gurnee, IL 60031

Made in China

#93-1646 11/25
©2025 Veridian Healthcare, LLC

INTRODUCTION

GENERAL DESCRIPTION

Thank you for selecting the Blood Pressure Monitor. The monitor features blood pressure measurement, pulse rate measurement and the result storage.

Readings taken by the monitor are equivalent to those obtained by a trained observer using the cuff and stethoscope auscultation method.

This manual contains important safety and care information, and provides step by step instructions for using the product.

Read the manual thoroughly before using the product.

INDICATIONS FOR USE

This Blood Pressure Monitor is digital monitor intended for use in measuring blood pressure and pulse rate with arm circumference ranging from:

Wide Range 8.6" – 16.5" (22 cm – 42 cm).

It is intended for adult indoor use only.

INDICATIONS FOR USE

- The device should not be used by any person who may be suspected of, or is pregnant.
- The device is not suitable for use on patients with implanted, electrical devices, such as cardiac pacemakers, defibrillators.

MEASUREMENT PRINCIPLE

This product uses the Oscillometric Measuring method to detect blood pressure. Before every measurement, the unit establishes a "zero pressure" equivalent to the air pressure. Then it starts inflating the arm cuff, meanwhile, the unit detects pressure oscillations generated by beat-to-beat pulsatile, which is used to determine the systolic and diastolic pressure, and also pulse rate.

CARE & SAFETY INFORMATION

NOTE: Read all instructions carefully before use. The following basic precautions are needed when using an electrical product.

CAUTION: Failure to read and observe all precautions could result in personal injury or equipment damage.

PRECAUTIONS

- This device is intended for adult use in homes only.
- The device is not suitable for use on neonatal patients, pregnant women, patients with implanted, electronic devices, patients with pre-eclampsia, premature ventricular beats, atrial fibrillation, peripheral, arterial disease and patients undergoing intravascular therapy or arterio-venous shunt or people who received a mastectomy. Please consult your doctor prior to using the unit if you suffer from illnesses.
- The device is not suitable for measuring the blood pressure of children. Ask your doctor before using it on older children.
- The device is not intended for patient transport outside a healthcare facility.
- The device is not intended for public use.
- This device is intended for no-invasive measuring and monitoring of arterial blood pressure. It is not intended for use on extremities other than the arm or for functions other than obtaining a blood pressure measurement.
- Do not confuse self-monitoring with self-diagnosis. This unit allows you to monitor your blood pressure. Do not begin or end medical treatment without asking a physician for treatment advice.
- If you are taking medication, consult your physician to determine the most appropriate time to measure your blood pressure. Never change a prescribed medication without consulting your physician.
- Do not take any therapeutic measures on the basis of a self measurement. Never alter the dose of a medicine prescribed

CARE & SAFETY INFORMATION

by a doctor. Consult your doctor if you have any question about your blood pressure.

- When the device was used to measure patients who have common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation, the best result may occur with deviation. Please consult your physician about the result.
- Don't kink the connection tube during use, otherwise, the cuff pressure may continuously increase which can prevent blood flow and result in harmful injury to the PATIENT.
- When using this device, please pay attention to the following situation which may interrupt blood flow and influence blood circulation of the patient, thus cause harmful injury to the patient: connection tubing kinking too frequent and consecutive multiple measurements; the application of the cuff and its pressurization on any arm where intravascular access or therapy, or an arterio-venous (A-V) shunt, is present; inflating the cuff on the side of a mastectomy.
- **Warning:** Do not apply the cuff over a wound; otherwise it can cause further injury.
- Do not inflate the cuff on the same limb which other monitoring ME equipment is applied around simultaneously, because this could cause temporary loss of function of those simultaneously-used monitoring ME equipment.
- On the rare occasion of a fault causing the cuff to remain fully inflated during measurement, open the cuff immediately. Prolonged high pressure (cuff pressure > 300mmHg or constant pressure > 15mmHg for more than 3 minutes) applied to the arm may lead to an ecchymosis.
- Please check that operation of the device does not result in prolonged impairment of patient blood circulation.
- When measurement, please avoid compression or restriction of the connection tubing.
- The device cannot be used with HF surgical equipment at the same time.

CARE & SAFETY INFORMATION

- The ACCOMPANYING DOCUMENT shall disclose that the SPHYGMOMANOMETER was clinically investigated according to the requirements of ISO 81060-2:2013.
- To verify the calibration of the AUTOMATED SPHYGMOMANOMETER, please contact the manufacturer.
- This device is contraindicated for any female who may be suspected of, or is pregnant. Besides providing inaccurate readings, the effects of this device on the fetus are unknown.
- Too frequent and consecutive measurements could cause disturbances in blood circulation and injuries.
- This unit is not suitable for continuous monitoring during medical emergencies or operations. Otherwise, the patient's arm and fingers will become anaesthetic, swollen and even purple due to a lack of blood.
- When not in use, store the device with the adapter in a dry room and protect it against extreme moisture, heat, lint, dust and direct sunlight. Never place any heavy objects on the storage case.
- This device may be used only for the purpose described in this booklet. The manufacturer cannot be held liable for damage caused by incorrect application.
- This device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in this booklet.
- The equipment is not AP/APG equipment and not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- **Warning:** No servicing/maintenance while the ME equipment is in use.
- The patient is an intended operator.
- The patient can measure data and change batteries under normal circumstances and maintain the device and its accessories according to the user manual.

CARE & SAFETY INFORMATION

- To avoid measurement errors, please avoid the condition of strong electromagnetic field radiated interference signal or electrical fast transient/burst signal.
- The blood pressure monitor, and the cuff are suitable for use within the patient environment. If you are allergic to polyester, nylon or plastic, please don't use this device.
- During use, the patient will be in contact with the cuff. The materials of the cuff have been tested and found to comply with requirements of ISO 10993-5:2009 and ISO 10993-10:2010. It will not cause any potential sensitization or irritation reaction.
- Adaptor is specified as a part of ME EQUIPMENT.
- If you experience discomfort during a measurement, such as pain in the arm or other complaints, press the START/STOP button to release the air immediately from the cuff. Loosen the cuff and remove it from your arm.
- If the cuff pressure reaches 40 kPa (300 mmHg), the unit will automatically deflate. Should the cuff not deflate when pressures reaches 40 kPa (300 mmHg), detach the cuff from the arm and press the START/STOP button to stop inflation.
- Before use, make sure the device functions safely and is in proper working condition. Check the device, do not use the device if it is damaged in any way. The continuous use of a damaged unit may cause injury, improper results, or serious danger.
- Do not wash the cuff in a washing machine or dishwasher!
- The service life of the cuff may vary by the frequency of washing, skin condition, and storage state. The typical service life is 10000 times.
- It is recommended that the performance should be checked every 2 years and after maintenance and repair, by retesting at least the requirements in limits of the error of the cuff pressure indication and air leakage (testing at least at 50mmHg and 200mmHg).

CARE & SAFETY INFORMATION

- Please dispose of ACCESSORIES, detachable parts, and the ME EQUIPMENT according to the local guidelines.
- Manufacturer will make available on request circuit diagrams, component part lists, descriptions, calibration instructions, etc., to assist to service personnel in parts repair.
- The plug/adaptor plug pins insulates the device from the main supply. Do not position the device in a position where it is difficult to disconnect from the supply mains to safely terminate operation of ME equipment.
- The operator shall not touch output of batteries /adaptor and the patient simultaneously.
- **Cleaning:** Dust environment may affect the performance of the unit. Please use the soft cloth to clean the whole unit before and after use. Don't use any abrasive or volatile cleaners.
- The device doesn't need to be calibrated within two years of reliable service.
- If you have any problems with this device, such as setting up, maintaining or using, please contact Customer Care. Don't open or repair the device by yourself in the event of malfunctions. The device must only be serviced, repaired and opened by individuals at authorized sales/service centers.
- Keep the unit out of reach of infants, young children or pets to avoid inhalation or swallowing of small parts. It is dangerous or even fatal.
- Be careful to strangulation due to cables and hoses, particularly due to excessive length.
- At least 30 min required for ME equipment to warm from the minimum storage temperature between uses until it is ready for intended use. At least 30 min required for ME equipment to cool from the maximum storage temperature between uses until it is ready for intended use.
- This equipment needs to be installed and put into service in accordance with the information provided in the ACCOMPA-


CARE & SAFETY INFORMATION

NYING DOCUMENTS;

- Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance d away from the equipment. The distance d is calculated by the MANUFACTURER from the 80MHz to 5.8 GHz column of Table 4 and Table 9 of IEC 60601-1-2:2014, as appropriate.
- Please use ACCESSORIES and detachable parts specified/authorized by MANUFACTURE. Otherwise, it may cause damage to the unit or danger to the user/patients.
- There is no luer lock connectors are used in the construction of tubing, there is a possibility that they might be inadvertently connected to intravascular fluid systems, allowing air to be pumped into a blood vessel.
- Please use the device under the environment which was provided in the user manual. Otherwise, the performance and lifetime of the device will be impacted and reduced.

BLOOD PRESSURE MONITOR FEATURES


MONITOR COMPONENTS

**PLEASE READ THIS
INSTRUCTION MANUAL COMPLETELY
BEFORE OPERATING THIS UNIT.**

BLOOD PRESSURE MONITOR FEATURES

LCD DISPLAY SIGNAL


Symbol	Description	Explanation
SYS	Systolic blood pressure	High blood pressure
DIA	Diastolic blood pressure	Low blood pressure
PUL /min	Pulse display	Pulse in beats per minute
AVG	Average value	The average value of blood pressure
1/3	Memory	Indicate it is in the memory mode and which group of memory it is.
	Motion indicator	Motion may result in an inaccurate measurement
mmHg	mmHg	Measurement Unit of the blood pressure (1mmHg=0.133kPa)
	Low battery	Batteries are low and need to be replaced
	Irregular pulse rate	Blood pressure monitor is detecting an irregular pulse rate during measurement
I	Blood pressure level indicator	Indicate the blood pressure level
M D Y 00/00/00 PM	Current Time	Month/Day/Year, Hour/Minute
	Pulse rate	Blood pressure monitor is detecting a pulse rate during measurement
	User 1/ User 2	Start measurement for User 1 or User 2
	The voice switch	Indicating voice switch

BATTERY INSTALLATION

INSTALLING AND REPLACING THE BATTERIES

1. Open the battery cover.
2. Install the batteries by matching the correct polarity, as shown.
3. Replace the cover.

Replace the batteries whenever the below happens.


- The shows
- The display dims
- The display does not light up

CAUTION

- Do not use new and used batteries together.
- Do not use different types of batteries together.
- Do not dispose the batteries in fire. Batteries may explode or leak.
- Remove batteries if the device is not likely to be used for some time.
- Worn batteries are harmful to the environment. Do not dispose with daily garbage.
- Remove the old batteries from the device following your local recycling guidelines.

CHOICE OF POWER SUPPLY

1. Battery powered mode:
6VDC 4xAAA batteries
2. AC adaptor powered mode:
5V 1A (Not included)
(Please use the AC adapter and USB cable which is authorized by the manufacturer!)
Please unplug the adaptor to depart from using utility power.

CAUTION


In order to get the best effect and protect your monitor, please use the right batteries and special power adapter which complies with local safety standard.

SETTING DATE & TIME


SETTING DATE, TIME AND MEASUREMENT UNIT

It is important to set the clock before using your blood pressure monitor, so that a time stamp can be assigned to each record that is stored in the memory. (The setting range of the year: 2018—2058 time: 12 H)


1. When the monitor is off, press **SET** button shortly it will display the time. Then press **SET** button for 3 seconds to enter the mode for year setting.


2. Press **MEM** button to change the [YEAR]. Each press will increase the numeral by one in a cycling manner.

3. When you get the right year, press **SET** button to set down and turn to next step.


4. Repeat steps 2 and 3 to set the [MONTH] and [DAY].


SETTING DATE & TIME

- Repeat steps 2 and 3 to set the [HOUR] and [MINUTE].
- After the [MINUTE] is set, the LCD will display "donE" first, then display all the settings you have done and then it will turn off.

♥ SELECT THE USER


- When the monitor is off, press and hold **MEM** button to enter user setting mode. The user ID will blink.
- Then press **MEM** button again, select the user ID between user 1 and user 2.
- After selecting the suitable user ID, press **SET** button to confirm. It will display "User ID+donE" and then turn off.


SETTING LANGUAGE

TURN ON OR TURN OFF THE VOICE


1. When the monitor is off, press and hold **MEM** button and **SET** button to enter the mode for turn on or turn off the voice. The trumpet will blink.

2. Press **MEM** button again to select turn on or turn off the voice. It will speak "Voice on" or "Voice off" at the same time.

3. Then press **SET** button to confirm. It will display "donE" to save the setting you have done and then turn off.


FITTING & APPLYING CUFF

PROPER CUFF FIT AND PLACEMENT IS CRITICAL IN OBTAINING ACCURATE BLOOD PRESSURE MEASUREMENTS. PLEASE READ THESE INSTRUCTIONS CAREFULLY AND CONSIDER CONSULTING WITH YOUR PHYSICIAN FOR FURTHER DIRECTION.

Before applying your blood pressure cuff, be sure you have selected the appropriate size cuff.

To determine the proper cuff size, measure the circumference of your left upper arm at the point midway between inside bend of your elbow and your shoulder.

Wide-Range Adult Cuff

fits arm circumference 8.6 – 16.5

If for any reason you are unable to or should not use your left arm, please modify the instructions for cuff application to your right arm. Your physician can help you identify which arm is best for you to take measurements from.

Assembling Your Cuff

Your cuff may need to be 'threaded' through the metal bar to form the cuff cylinder for your arm.

1. Lay the open cuff hook and loop side down, size/instruction wording side up on a flat surface.
2. Keeping the words on the inside of the cuff, thread the end of the cuff through the metal bar, Fig 1. The end of the cuff should be able to fold over onto itself for securing via the hook and loop closure, Fig 2.

Fig. 1


Fig. 2

FITTING & APPLYING CUFF

APPLYING YOUR CUFF

Be seated at a table or desk with your feet flat on the floor. To avoid accidentally pulling on or dropping your monitor, the cuff should not be plugged into the monitor until after the cuff is applied to your arm.

The cuffs are printed with a useful guide for applying the cuff; please consult the following steps for additional direction.

1. Remove any constrictive clothing or jewelry that may interfere with cuff placement.
2. Slide the cuff cylinder onto your left arm and align the artery mark over the brachial artery on the inside of your arm, Fig 3.

The tube should follow the inside length of the arm, Fig 4.

3. The bottom edge of the cuff should be positioned approximately 0.5" above the elbow joint, Fig 5.
4. When the cuff is closed around the arm, the triangle guide imprint should be within the OK range bar marking; if the triangle is outside of the range, a different sized cuff is needed
5. The cuff should fit comfortably, yet snugly around your arm, Fig 6. You should be able to insert one finger easily between your arm and the cuff.


Fig. 3


Fig. 4


Fig. 5


Fig. 6

TAKING YOUR BP READING

HELPFUL TIPS

- Take the measurement in a silent room.
- Rest for 5 minutes before a measurement.
- Wait at least 3 minutes before another measurement. This allows your blood circulation to recover.
- Be relax and do not move and talk during the measurement procedure.
- The cuff should maintain at the same level as the right atrium of the heart.
- Please sit comfortably. Do not cross your legs and keep your feet flat on the ground.
- Keep your back against the backrest of the chair.
- For a meaningful comparison, try to measure under similar conditions. For example, take daily measurements at approximately the same time, on the same arm, or as directed by a physician.

TAKING YOUR BP READING

START THE MEASUREMENT

Be seated at a table or desk with your feet flat on the floor. To avoid accidentally pulling on or dropping your monitor, the cuff should not be plugged into the monitor until after the cuff is applied to your arm.

The cuffs are printed with a useful guide for applying the cuff; please consult the following steps for additional direction.

1. When the monitor is off, press the **START/STOP** to turn on the monitor, and it will finish the whole measurement.


(When the volume is on, it will speak "Starting measurement, please remain still and quiet" and then complete the measurement process.)

LCD Display

Inflating and measuring.

Adjust the zero point.

Display and save the results.
When volume is on, It will broadcast your measuring result.)

2. Press the **START/STOP** to power off, otherwise it will turn off within 1 minute.


MEMORY FUNCTION

RECALL THE RECORDS

1. When the monitor is off, please press **MEM** button to show the recent record. If the records are less than 3 groups, it will display the latest record instead. Take user 1 for example. (When the volume is on, it will read the number for you).
2. Press the **MEM** or **SET** to get the record you want. (When the volume is on, it will read the number for you)

The date and time of the record will be shown alternately.

CAUTION

The most recent record (1) is shown first. Each new measurement is assigned to the first (1) record. All other records are pushed back one digit (e.g., 2 becomes 3, and so on), and the last record (120) is dropped from the list.

MEMORY FUNCTION

DELETE RECORDS

If you did not get the correct measurement, you can delete results by following steps below.

Delete one record

1. Hold pressing **SET** button (except the AVG) for 3 seconds when the monitor is in the memory recall mode, the flash display "dEL one+User ID" will show.
2. Press and hold **START/STOP** button to confirm deleting and the monitor will display "dEL donE+User ID" and then the device will show the latest record. (When the volume is on, it will read the number for you)

Note: To exit out the delete mode press **START/STOP** button shortly before you press and hold **START/STOP** button.

Delete all records

1. Hold pressing **MEM** button for 3 seconds when the monitor is in the memory recall mode, the flash display "dEL ALL+User ID" will show.
2. Press and hold **START/STOP** button to confirm deleting and the monitor will display "dEL donE+USER ID" and then turn off.

Note: To exit out the delete mode press **START/STOP** button shortly before you press and hold "Start/stop" button.

3. If there is no record, the right display will show.

CARE & MAINTENANCE

Proper care and maintenance of your blood pressure monitor will help ensure long life for the product; improper use of the monitor may void the included warranty.

Always disconnect the cuff and tubing from the monitor prior to cleaning and storage.

MONITOR CARE & MAINTENANCE

- Do not use any liquids on the monitor; use a soft, dry cloth to clean the monitor as needed.
- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: -4°F – 131°F (-20°C – 55°C); Humidity: 10%~93%RH (non- condensing).
- This monitor requires 6 hours to cool from the maximum storage temperature between uses until the monitor is ready for its INTENDED USE when the ambient temperature is 68°F(20 °C).
- This monitor requires 6 hours to warm from the minimum storage temperature between uses until the monitor is ready for its INTENDED USE when the ambient temperature is 68°F(20 °C).
- Avoid high temperatures and direct sunlight.
- Always remove the batteries if the unit will not be used for an extended period of time.
- Do not attempt to disassemble your monitor; disassembling your monitor may void the manufacturers warranty.
- Do not attempt to service the monitor while in use or on a patient.
- Do not subject the monitor to strong shocks; take care not to drop the monitor.

CARE & MAINTENANCE

- No component can be maintained by user in the monitor. The circuit diagrams, component part lists, descriptions, calibrations instructions, or other information which will assist the user's appropriately qualified technical personnel to repair those parts of equipment which are designated repairable can be supplied.


CUFF CARE & MAINTENANCE

- The cuff may be spot cleaned with a mild-detergent as needed.
- To disinfect the cuff, wipe the inside (skin contact side) of the cuff using a cloth lightly moistened with 70% - 90% ethyl rubbing alcohol and allow to air dry; it is suggested to disinfect the cuff frequently when used in clinical settings or if the cuff is contaminated in any way.
- Never attempt to iron or dry the cuff; if the cuff becomes wet, allow the cuff to air dry completely before use or storage.
- Do not attempt to remove the tube and bladder from the nylon cuff as this may damage the cuff.
- Take care not to forcefully bend or kink the tubing during storage.
- It is recommended the performance should be checked every 2 years or after a repair. Please contact Customer Care at 866-326-1313.
- It is recommended to clean the cuff every 200 uses.
- The monitor can maintain the safety and performance characteristics for minimum of 10,000 measurements of three years.
- The cuff integrity is maintained after 1,000 open- close cycles of the closure.

HYPERTENSION INDICATOR

WHAT ARE SYSTOLIC PRESSURE AND DIASTOLIC PRESSURE?

When ventricles contract and pump blood out of the heart, the blood pressure reaches its maximum value in the cycle, which is called systolic pressure. When the ventricles relax, the blood pressure reaches its minimum value in the cycle, which is called diastolic pressure.

WHAT IS THE STANDARD BLOOD PRESSURE CLASSIFICATION?

The chart on the right is the standard blood pressure classification published by American Heart Association (AHA).



This chart reflects blood pressure categories defined by the American Hearth Association.			
Blood Pressure Category	Systolic mmHg (upper#)		Diastolic mmHg (lower#)
Normal	Less than 120	and	Less than 80
Prehypertension	120 - 129	and	Less than 80
High Blood Pressure (Hypertension) Stage 1	139 - 139	or	80 - 89
High Blood Pressure (Hypertension) Stage 2	140 or higher	or	90 or higher
Hypertension Crisis (Consult your doctor immediately)	Higher than 180	and/or	Higher than 120

⚠ CAUTION

Please consult a physician if your measuring result falls outside the range. Please note that only a physician can tell whether your blood pressure value has reached a dangerous point.

IRREGULAR PULSE RATE INDICATOR

An irregular pulse rate will be detected if there is an irregular pulse rhythm while measuring systolic and diastolic blood pressure. When measurements were performed, the monitor will record all pulse intervals and calculate the average. If two or more pulse intervals were recorded, and the difference between each interval and the average is larger than $\pm 25\%$ of the average; or if four or more pulse intervals were recorded, and the difference between each interval and the average is larger than $\pm 15\%$ of the average value, the irregular pulse symbol will be displayed along with measurement results.

CAUTION

The appearance of the IHB icon indicates that a pulse irregularity consistent with an irregular pulse rate was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice.

Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

TROUBLESHOOTING

This section includes a list of error messages and frequently asked questions for problems you may encounter with your blood pressure monitor. If the product is not operating as you think it should, check here before arranging for servicing.

Problem	Symptom	Check This	Remedy
No Power	Display does not light up	Batteries are exhausted.	Replace with new batteries.
		Batteries are inserted incorrectly.	Insert the batteries correctly.
		AC Adapter is inserted incorrectly.	Insert the AC adaptor tightly
Low Battery	Display is dim or shows 	Batteries are low.	Replace with new batteries.
Error Code	E 01 shows	The cuff is too tight or too loose.	Refasten the cuff and then measure again.
	E 02 shows	The monitor detected motion, talking or the pulse is too poor while measuring.	Relax for 5 minutes, keep still and then measure again.
	E 03 shows	The measurement process does not detect the pulse signal.	Loosen the clothing on the arm and then measure again.
	E 04	The treatment of the measurement failed.	Relax for 5 minutes and measure again.
	EExx, shows on the display.	A calibration error occurred.	Retake the measurement. If the problem persists, contact the retailer or our customer service department for further assistance. Refer to the warranty for contact information and return instructions.
Warning Message	"Out" shows	Out of measurement range.	Relax for a moment. Refasten the cuff and then measure again. If the problem persists, contact your physician.

DEVICE & SYMBOLS

SAFETY INFORMATION

The signs below might be in the user manual, labeling or other component. They are the requirement of standard and using.

	Refer to instruction manual/booklet to signify that the instruction manual/booklet must be read. NOTE: Symbol is BLUE
	Serial Number
	Direct Current
	General symbol for recovery/recyclable
	Caution Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.
	Type BF Applied Part
	The symbol indicates that the product should not be discarded as unsorted waste but must be sent to separate collection facilities for recovery and recycling.
	Date of Manufacture
	MR Unsafe To identify an item which poses unacceptable risks to the patient, medical staff or other persons within the MR environment.

FCC STATEMENT

NOTE:

POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for U.S.A. only)

This product has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. The product generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If the product does cause harmful interference to radio or television reception, which can be determined by turning the product on and off, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the product and the receiver.
- Connect the product into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

ELECTROMAGNETIC COMPATIBILITY

The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments.

Warning: Don't be near the active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Technical description:

- 1. All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the expected service life.
- 2. Guidance and manufacturer's declaration -electromagnetic emissions and Immunity

Table 1

Guidance and manufacture's declaration - electromagnetic emissions

Emission test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations and flicker emissions IEC 61000-3-3	Comply

ELECTROMAGNETIC COMPATIBILITY

Table 2

Guidance and manufacture's declaration - electromagnetic immunity

Immunity test	IEC60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact; ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact; ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines; ± 1 kV signal input/output; 100 kHz repetition frequency	± 2 kV for power supply lines; Not applica- ble; 100 kHz repetition frequency
Surge IEC61000-4-5	± 0.5 kV, ± 1 kV differential mode; ± 0.5 kV, ± 1 kV, ± 2 kV common mode	± 0.5 kV, ± 1 kV differential mode; Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0% UT; 1 cycle and 70% UT; 25/30 cycles; Single phase: at 0°. 0% UT; 250 / 300 cycle	0% UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0% UT; 1 cycle and 70% UT; 25/30 cycles; Single phase: at 0°. 0% UT; 250 / 300 cycle
Power frequency magnetic field IEC 61000-4-8	30 A/m 50 Hz / 60 Hz	30 A/m 50 Hz / 60 Hz
Conducted RF IEC61000-4-6	3 V; 0,15 MHz – 80 MHz; 6 V in ISM and amateur radio bands; between 0,15 MHz and 80 MHz; 80% AM at 1 kHz	3 V; 0,15 MHz – 80 MHz; 6 V in ISM and amateur radio bands; between 0,15 MHz and 80 MHz; 80% AM at 1 kHz
Radiated RF IEC61000-4-3	10 V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

ELECTROMAGNETIC COMPATIBILITY

Table 3

Guidance and manufacture's declaration - electromagnetic immunity

Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication s equipment)	Test Immunity MHz	Band MHz	Service	Modulation	Modulation (W)	Distance (m)	Immunity Test Level (V/m)
	385	380-390	TETRA 400	Pulse modulation b) 18Hz	1.8	0.3	27
	450	430-470	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0.3	28
	710	704-787	LTE Band 13,17	Pulse modulation b) 217 Hz	0.2	0.3	9
	745						
	780						
	810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18 Hz	2	0.3	28
	870						
	930						
	1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1,3,4,25; UMTS	Pulse modulation b) 217 Hz	2	0.3	28
	1845						
	1970						
	2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	2	0.3	28
	5240	5100-5800	WLAN 802.11a/n	Pulse modulation b) 217 Hz	0.2	0.3	9
	5500						
	5785						

SPECIFICATIONS

Name	Upper Arm Talking Blood Pressure Monitor
Model Number	UAM-901
Display Mode	Digital LCD display V.A. 2.36" x 2.91" (60mm * 74mm)
Measuring Mode	Oscillographic testing mode
Measuring Range	Rated cuff pressure: 0mmHg ~ 299mmHg (0kPa ~ 39.9kPa) Measurement pressure: SYS: 60mmHg ~ 230mmHg (8.0kPa~30.7kPa) DIA: 40mmHg ~ 130mmHg (5.3kPa~17.3kPa) Pulse value: (40 - 199) beat/minute
Accuracy	Pressure: 41°F - 104°F (5°C - 40°C) within±3mmHg (0.4kPa) Pulse value:±5%
Operating Environment	Temperature:+41°F to +104°F (+5°C to +40°C) Relative humidity: 15% to 90%, non-condensing, but not requiring a water vapor partial pressure greater than 50 hPa; Atmospheric pressure: 700 hPa to 1060 hPa
Storage/transportation Environment	Temperature: -4°F to +140°F (-20°C to +60°C) Relative humidity: ≤93%, non-condensing, at a water vapor pressure up to 50 hPa
Arm Measurement	Wide Range: 8.6" - 16.5" (22 ~ 42 cm) (upper arm circumference)
Weight	Approx. 173g (Excluding the batteries and cuff)
External Dimensions	Approx. 3.66" x 5.12" x 1.28" (93mm*130mm*32.5mm)
Accessories	Instruction Manual, Four (4) AAA Batteries
Mode of Operation	Continuous operation
Device Classification	Battery Powered Mode: Internally Powered ME Equipment AC Adaptor Powered Mode: Class II ME Equipment
Degree of Protection	Type BF applied part
Protection against ingress of water	IP21 It means the device could be protected against solid foreign ingress objects of 12,5mm Φ and greater, and against vertically falling water drops.
Software Version	A01

Specifications are subject to change without notice.

5-YEAR LIMITED WARRANTY

Congratulations on your purchase of a Blood Pressure Monitor. Your Blood Pressure Monitor is covered by the following limited warranty commencing upon the date of purchase, and subject to the following terms and conditions:

The Warrantor warrants that its Blood Pressure Monitor will be free from defects in materials and workmanship under normal consumer usage for a period of five years for the original purchaser of the product.

Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the Blood Pressure Monitor for commercial purposes or subjecting the Blood Pressure Monitor to abnormal usage or conditions; or (d) other acts which are not the fault of the Warrantor, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the Blood Pressure Monitor.

If the Blood Pressure Monitor fails to conform to this limited warranty, return the Blood Pressure Monitor postage prepaid to: **Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031**. When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone number. Carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. At our option, the Warrantor will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will be responsible for parts charges, if any, and repair charges not covered by this limited warranty.

EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL THE WARRANTOR BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

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