Blood Pressure Monitor Premium

Model UAM-880 Wide Range Arm Cuff

INSTRUCTION MANUAL

ENGLISH & SPANISH


Please read this instruction manual
completely before operating this unit.


STOP!
**PLEASE ENSURE YOU HAVE ALL OF THE
FOLLOWING COMPONENTS BEFORE USING
YOUR BLOOD PRESSURE MONITOR.**

Monitor


Wide-Range Cuff

Four (4) AA
Batteries,
Preinstalled

Carrying
Case

Instruction Manual
App guide

Toll-Free Customer Care Help Line:

1-866-326-1313

Monday – Friday

8:30 a.m. – 4:30 p.m. CST

www.veridianhealthcare.com

INDEX

Introduction	4
Care & Safety.....	5-9
Blood Pressure Monitor Features	10-11
MyHealth Connected App Guide	12
Battery Installation.....	13
Setting Date & Time.....	14-16
Fitting & Applying Cuff	17-18
Taking Your BP Reading	19-22
Memory Function	23-24
Care & Maintenance	25-26
Hypertension Indicator.....	27
Irregular Pulse Rate Indicator	28
Troubleshooting	29
Device & Label Symbols.....	30
FCC Statement.....	31
Electromagnetic Compatibility.....	32-34
Specifications	35
Warranty	36
Español	37-76
Reference Page.....	74
MyHealth Connect App	76

Zewa®
Manufactured for
Veridian Healthcare
1175 Lakeside Drive
Gurnee, IL 60031

Made in China

#93-1648 12/25
©2025 Veridian Healthcare, LLC

INTRODUCTION

GENERAL DESCRIPTION

Thank you for selecting the Blood Pressure Monitor. The monitor features blood pressure measurement, pulse rate measurement and the result storage. Readings taken by the monitor are equivalent to those obtained by a trained observer using the cuff and stethoscope auscultation method.

This manual contains important safety and care information, and provides step by step instructions for using the product.

Read the manual thoroughly before using the product.

INDICATIONS FOR USE

This Blood Pressure Monitor is digital monitor intended for use in measuring blood pressure and pulse rate with arm circumference ranging from:

Wide-Range Arm Cuff fits arm size: 8.6" – 16.5" (22cm~42cm)

It is intended for adult indoor use only.

The device can be used to detect irregular pulse rate.

MEASUREMENT PRINCIPLE

This product uses the Oscillometric Measuring method to detect blood pressure. Before every measurement, the unit establishes a "zero pressure" equivalent to the air pressure. Then it starts inflating the arm cuff, meanwhile, the unit detects pressure oscillations generated by beat-to-beat pulsatile, which is used to determine the systolic and diastolic pressure, and also pulse rate.

CARE & SAFETY INFORMATION

NOTE: Read all instructions carefully before use. The following basic precautions are needed when using an electrical product.

CAUTION: Failure to read and observe all precautions could result in personal injury or equipment damage.

PRECAUTIONS

- This device is intended for indoor, home use and is not intended for self-use in public areas.
- This device is portable, but it is not intended for use during patient transport.
- This device is not suitable for continuous monitoring during medical emergencies or operations.
- This device is intended for non-invasive measuring and monitoring of arterial blood pressure. It is not intended for use on extremities other than the arm, or for any purpose other than obtaining a blood pressure measurement.
- This device is for adults. Do not use this device on neonates or infants. Do not use it on children and adolescents unless otherwise instructed by a medical professional.
- Consult with your physician before using this monitor if you suffer from the following conditions: common arrhythmias such as premature ventricular beats or atrial fibrillation; peripheral arterial disease; pregnancy; preeclampsia; implantation with electrical devices; undergoing intravascular therapy; arteriovenous shunt or mastectomy. Please note that any of these conditions may affect measurement readings, in addition to patient motion, trembling or shivering.
- Do not use this device for diagnosis or treatment of any health problem or disease. Contact your physician if you have or suspect any medical problem. Do not change your medications without the advice of your physician or health care professional.

CARE & SAFETY INFORMATION

- If you are taking medication, consult your physician to determine the proper time to measure your blood pressure.
- This device may be used only for the intended use described in this manual, the manufacturer shall have no liability for any incidental, consequential, or special damages caused by misuse or abuse.
- Please use the device under the environment which is provided in the user manual. Otherwise, the performance and lifetime of the device will be impacted and reduced.
- The device may require up to 30 minutes to warm up / cool down from the minimum/ maximum storage temperature before it is ready for use.
- The blood pressure monitor and the cuff are suitable for use within the patient environment.
- Do not wash the cuff in a washing machine or dishwasher !
- The device contains sensitive electronic components. To avoid measurement errors, avoid taking blood pressure measurements near a strong electromagnetic field radiated interference signal or electrical fast transient/burst signal.
- Blood Pressure Monitor is intended for use by medical staffs and lay persons, and patient is also an intended user or operator.
- Wireless communication equipment, such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies may cause interference that may affect the accuracy of measurements. A minimum distance of 1 foot (30 cm) should be kept from such devices during a measurement.

CAUTIONS

- Do not attempt to repair the unit yourself if it malfunctions. Only have repairs carried out by authorized service centers.

CARE & SAFETY INFORMATION

- It is recommended that the performance should be checked after repair, maintenance, and every two years of use, by retesting the requirements in limits of the error of the cuff pressure indication and air leakage (testing at least at 50 mmHg and 200 mmHg). Please contact manufacturer or distributor for authorized service personnel.
- Store your device, cuff in a clean and dry place, protect it against extreme moisture, heat, lint, dust and direct sunlight. Never place any heavy objects on it.
- Make sure the rubber tube of the cuff is not squeezed, stretched, or kinked during storage.
- Dispose of accessories, detachable parts, and the device according to the local guidelines.

WARNINGS

- Do not apply the cuff on an arm that has an intravenous drip or a blood transfusion attached.
- Do not kink, fold, stretch, compress, or otherwise deform the tube during measuring, as the cuff pressure might continuously increase, which could prevent blood flow and result injury.
- Taking blood pressure measurements too frequently could disrupt blood circulation and cause injuries.
- Do not apply cuff to areas on patient where skin is delicate or damaged. Check cuff site frequently for irritation.
- Do not place the cuff on the arm of a person whose arteries or veins are undergoing medical treatment, i.e. intra-vascular access or intra-vascular therapy or an arteriovenous (A-V) shunt, which could disrupt blood circulation and cause injuries.
- Do not place the cuff on the arm on the same side of a mastectomy (especially when lymph nodes have been removed). it is recommended to take measurements on the unaffected side.
- Do not wrap the cuff on the same arm to which another monitoring

CARE & SAFETY INFORMATION

device is applied. One or both devices could temporarily stop functioning if you try to use them at the same time.

- Please check (for example, by observation of the limb concerned) that the operation of the device does not result in prolonged impairment of patient blood circulation.
- On the rare occasion of a fault causing the cuff to remain fully inflated during measurement, loosen and remove the cuff immediately. Prolonged high pressure applied to the arm (cuff pressure >300 mmHg or constant pressure >15 mmHg for more than 3 minutes) might lead to bruising and discolored skin.
- Do not use this device with high-frequency (HF) surgical equipment at the same time.
- This device is not used in conjunction with oxygen rich environments, not intended for use with flammable anaesthetics, not intended for use in conjunction with flammable agents.
- Do not touch output of the batteries and the user simultaneously.
- The power cord is considered the disconnect device for isolating this equipment from supply mains. Do not position the equipment so that it is difficult to reach or disconnect.
- Do not use this device if you are allergic to polyester, nylon, or plastic.
- Only use accessories approved by manufacturer. Using unapproved accessories might cause damage to the unit and injure users.
- If you experience discomfort during a measurement, such as pain in the arm or other complaints, press the Power button immediately to release the air from the cuff.
- Do not use the device while under maintenance, or being serviced.

CARE & SAFETY INFORMATION


- The air tube poses a risk of strangulation. Furthermore, the small parts of product and batteries present a choking hazard if swallowed. They should therefore always be kept away from infants/children.
- Sensor degradation or looseness may reduce performance of device or cause other problems.

NOTICE

- You can use this device to take your own measurement, no third-party operator is required.
- At the request of authorized service personnel, circuit diagrams, component part lists, descriptions, and calibration procedures will be made available by the manufacturer or distributor.
- The expected lifetime of the cuff may vary by the frequency of washing, skin condition, and storage state.
- Please report to the manufacturer and the competent authority of the Member State / the FDA in which you are established about any serious incident that has occurred in relation to this device.

BLOOD PRESSURE MONITOR FEATURES


MONITOR COMPONENTS

**PLEASE READ THIS INSTRUCTION MANUAL
COMPLETELY BEFORE OPERATING THIS UNIT.**

BLOOD PRESSURE MONITOR FEATURES

LCD DISPLAY SIGNAL

Symbol	Description	Explanation
SYS.	Systolic blood pressure	High blood pressure
DIA.	Diastolic blood pressure	Low blood pressure
Pul min	Pulse display	Pulse in beats per minute
M D 88/88 PM	Current time	Time (year: month: day: hour: minute)
AVG.3	Average value	The average value of last 3 groups value blood pressure
mmHg	mmHg	Measurement Unit of the blood pressure
Lo +	Low battery indicator	Batteries are running low and need to be replaced
Heart icon with wavy lines	Irregular Pulse Rate Detection	Blood pressure monitor is detecting an irregular pulse rate during measurement
Triangle icon	Blood Pressure level indicator	Indicate the blood pressure level
Heart icon	Pulse Rate	Blood pressure monitor is detecting a pulse rate during measurement
User 1 icon	User 1	Start measurement,save and transmit the measuring results for User 1
User 2 icon	User 2	Start measurement,save and transmit the measuring results for User 2
Hand icon	Motion indicator	Motion may result in an inaccurate measurement
Document icon	Data pending to transmit	Measurement data stored in the equipment
Bluetooth icon	Bluetooth icon	The Bluetooth icon blinks when the Bluetooth is working

MYHEALTH CONNECTED APP

Installing the **MyHealth Connected** app is fast and easy.

Simply download the app, prepare your device, and start your health journey. Consult your healthcare professional before starting a new health regime, changing medication and for aid in interpreting your personal results.

Bluetooth Smart Compatible
Connects to Android and iOS
mobile devices.

1. Download the **MyHealth Connected** app from the App Store (iOS) or Google Play (Android).
2. Launch the app and follow the instructions for account profile set up.
3. Pair your device following the on-screen directions.

Consult the included instruction manual for detailed instructions, tips, cautions and cleaning direction.


THIS BLOOD PRESSURE MONITOR OPERATES WITH OR WITHOUT THE APP. SOME FEATURES, SUCH AS MEMORY STORAGE AND RECALL, CAN ONLY BE ACCESSED THROUGH THE APP.

BATTERY INSTALLATION

INSTALLING AND REPLACING THE BATTERIES

Battery powered mode: 6VDC 4xAA batteries

1. Slide off the battery cover.
2. Install the batteries by matching the correct polarity, as shown.
3. Replace the cover.

Replace the batteries whenever the below happens.

- The  shows
- The display dims
- The display does not light up

CAUTION

- Do not use new and used batteries together.
- Do not use different types of batteries together.
- Do not dispose the batteries in fire. Batteries may explode or leak.
- Remove batteries if the device is not likely to be used for some time.
- Worn batteries are harmful to the environment. Do not dispose with daily garbage.
- Remove the old batteries from the device following your local recycling guidelines.

AC adaptor powered mode:


AC Adapter **NOT INCLUDED**

Type: BLJ06L060100P-U

Input: 100-240V, 50-60Hz, 0.2A max

Output: 6V  1000mA

(Please use the AC adapter which is authorized by the manufacturer!) Please unplug the adaptor to depart from using utility power.

CAUTION


In order to get the best effect and protect your monitor, please use the right batteries and special power adapter which complies with local safety standard.

SETTING DATE & TIME


SETTING DATE, TIME AND MEASUREMENT UNIT

It is important to set the clock before using your blood pressure monitor, so that a time stamp can be assigned to each record that is stored in the memory. (year: 2018—2058 time: 12 H)


1. When the monitor is off, hold pressing  button for 3 seconds to enter the year setting mode.

2. Press the **MEM** button to change the [YEAR].


3. When you get the right year, press  button to confirm and turn to next step.


SETTING DATE & TIME


4. Repeat step 2 and 3 to set the [MONTH] and [DAY].

5. Repeat steps 2 and 3 to set the [HOUR] and [MINUTE].


6. After the [MINUTE] is set, the LCD will display "donE" first, then display all the settings you have done and then turn off.


SETTING DATE & TIME

SELECT THE USER

1. When the monitor is off , press the  button shortly to enter user setting mode.

2. Then press  button again, select the user ID between user 1 and user 2.

3. After selecting the suitable user ID, press **START/STOP** button to start the measurement.


FITTING & APPLYING CUFF

PROPER CUFF FIT AND PLACEMENT IS CRITICAL IN OBTAINING ACCURATE BLOOD PRESSURE MEASUREMENTS. PLEASE READ THESE INSTRUCTIONS CAREFULLY AND CONSIDER CONSULTING WITH YOUR PHYSICIAN FOR FURTHER DIRECTION.

Before applying your blood pressure cuff, be sure you have selected the appropriate size cuff.

To determine the proper cuff size, measure the circumference of your left upper arm at the point midway between inside bend of your elbow and your shoulder.

Wide-Range Adult Cuff fits arm size:
8.6" – 16.5" (22cm~42cm))

If for any reason you are unable to or should not use your left arm, please modify the instructions for cuff application to your right arm. Your physician can help you identify which arm is best for you to take measurements from.

Assembling Your Cuff

Your cuff may need to be 'threaded' through the metal bar to form the cuff cylinder for your arm.

1. Lay the open cuff hook and loop side down, size/instruction wording side up on a flat surface.
2. Keeping the words on the inside of the cuff, thread the end of the cuff through the metal bar, Fig 1. The end of the cuff should be able to fold over onto itself for securing via the hook and loop closure, Fig 2.


Fig. 1


Fig. 2

FITTING & APPLYING CUFF

APPLYING YOUR CUFF

Be seated at a table or desk with your feet flat on the floor. To avoid accidentally pulling on or dropping your monitor, the cuff should not be plugged into the monitor until after the cuff is applied to your arm.

The cuffs are printed with a useful guide for applying the cuff; please consult the following steps for additional direction.

1. Remove any constrictive clothing or jewelry that may interfere with cuff placement.
2. Slide the cuff cylinder onto your left arm and align the artery mark over the brachial artery on the inside of your arm, Fig 3.
The tube should follow the inside length of the arm, Fig 4.
3. The bottom edge of the cuff should be positioned approximately 0.5" above the elbow joint, Fig 5.
4. When the cuff is closed around the arm, the triangle guide imprint should be within the OK range bar marking; if the triangle is outside of the range, a different sized cuff is needed
5. The cuff should fit comfortably, yet snugly around your arm, Fig 6. You should be able to insert one finger easily between your arm and the cuff.


Fig. 3


Fig. 4


Fig. 5


Fig. 6

TAKING YOUR BP READING

HELPFUL TIPS


- Take the measurement in a silent room.
- Rest for 5 minutes before a measurement.
- Wait at least 3 minutes before another measurement. This allows your blood circulation to recover.
- Be relax and do not move and talk during the measurement procedure.
- The cuff should maintain at the same level as the right atrium of the heart.
- Please sit comfortably. Do not cross your legs and keep your feet flat on the ground.
- Keep your back against the backrest of the chair.
- For a meaningful comparison, try to measure under similar conditions. For example, take daily measurements at approximately the same time, on the same arm, or as directed by a physician.

TAKING YOUR BP READING

START THE MEASUREMENT

Be seated at a table or desk with your feet flat on the floor. To avoid accidentally pulling on or dropping your monitor, the cuff should not be plugged into the monitor until after the cuff is applied to your arm.


The cuffs are printed with a useful guide for applying the cuff; please consult the following steps for additional direction.


1. When the monitor is off, press the **START/STOP** button to turn on the monitor, and it will finish the whole measurement, save the desired user. (Take User 1 for example.)

START
STOP

LCD Display


Adjust the zero point.

Inflating and measuring.


Display and save the results.


TAKING YOUR BP READING

2. This device will proceed to data transmission after measurement. The Bluetooth symbol blinks on the LCD indicates data is transmitting.

3. If the data is successfully transmitted, the Bluetooth symbol will light up, and the [document icon] will disappear. If the data transmission fails, the device will turn off after one minute.

4. Press **START/STOP** button to power off, otherwise it will turn off within 1 minute.

TIPS FOR BLOOD PRESSURE MEASURING

- You can press **START/STOP** button at any time to stop measuring during the process of measurement.
- Maximum 60 records are both for **USER 1** and **USER 2**.
- If the measurement result is out of the measurement range (SYS: 60mmHg to 230mmHg; or DIA: 40mmHg to 130mmHg; or Pulse: 40-199 pulse/minute), the LCD will display "out".

TAKING YOUR BP READING

Bluetooth Module No.: LS51802

RF Frequency Range: 2.402 - 2.480 GHz

Output Power Range: -6dBm -4dBm

Supply Voltage: 1.8V - 3.6V

Transmitting Distance: 1-10 meters

CAUTION

- Interference may occur in the vicinity of equipment marked with the following symbol . UAM-880 may interfere with electrical equipment.
- Keep the monitor at least 20 centimeters away from the human body (especially the head) when the data transmission is proceeding after measurement.
- To enable the data transmission function, this product should be paired to Bluetooth end at 2.4 GHz.


HOW TO MITIGATE POSSIBLE INTERFERENCE?

1. The range between the device and BT end should be reasonably close, from 1 meter to 10 meters. Please ensure no obstacles between the device and BT end so as to obtain quality connection and to lower the RF output range.
2. To avoid interference, other electronic devices (particularly those with wireless transmission / Transmitter) should be kept at least 1 meter away from the monitor.



MEMORY FUNCTION

RECALL THE RECORDS

1. When the monitor is off, please press **MEM** button shortly to show the average value of the latest three records. If the records are less than 3 groups, it will display the latest record instead. (Take user 1 for example.)

2. Press the **MEM** button again to rotate the records. The order of the record, date and time will be displayed alternatively.

CAUTION


The most recent record (1) is shown first. Each new measurement is assigned to the first (1) record. All other records are pushed back one digit (e.g., 2 becomes 3, and so on), and the last record (60) is dropped from the list.

MEMORY FUNCTION


DELETE THE RECORDS

You can delete all results of the selected user by following steps below.

1. When device is turned OFF, press the **MEM** button to go into memory mode. Now press the **MEM** button for 3 seconds until the "User ID + dEL ALL" is flashing.


2. Press **MEM** button to confirm deleting all the memories, the LCD will display "User ID + dEL donE" and the monitor will turn off.

3. To exit the delete mode without deleting any records, press **START/STOP** button before pressing **MEM** button to confirm any delete commands.

4. If there is no record, the below display will be shown

CARE & MAINTENANCE

Proper care and maintenance of your blood pressure monitor will help ensure long life for the product; improper use of the monitor may void the included warranty.

Always disconnect the cuff and tubing from the monitor prior to cleaning and storage.

MONITOR CARE & MAINTENANCE

- Do not use any liquids on the monitor; use a soft, dry cloth to clean the monitor as needed.
- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: -4°F – 140°F (-20°C – 60°C); Humidity: ≤93%RH (non- condensing).
- This monitor requires 6 hours to cool from the maximum storage temperature between uses until the monitor is ready for its INTENDED USE when the ambient temperature is 68°F(20 °C).
- This monitor requires 6 hours to warm from the minimum storage temperature between uses until the monitor is ready for its INTENDED USE when the ambient temperature is 68°F(20 °C).
- Avoid high temperatures and direct sunlight.
- Always remove the batteries if the unit will not be used for an extended period of time.
- Do not attempt to disassemble your monitor; disassembling your monitor may void the manufacturers warranty.
- Do not attempt to service the monitor while in use or on a patient.
- Do not subject the monitor to strong shocks; take care not to drop the monitor.
- No component can be maintained by user in the monitor. The circuit diagrams, component part lists, descriptions, calibrations instructions, or other information which will assist the user's appropriately qualified technical personnel to repair those parts of equipment which are designated repairable can be supplied.

CARE & MAINTENANCE


CUFF CARE & MAINTENANCE

- The cuff may be spot cleaned with a mild-detergent as needed.
- To disinfect the cuff, wipe the inside (skin contact side) of the cuff using a cloth lightly moistened with 70% - 90% ethyl rubbing alcohol and allow to air dry; it is suggested to disinfect the cuff frequently when used in clinical settings or if the cuff is contaminated in any way.
- Never attempt to iron or dry the cuff; if the cuff becomes wet, allow the cuff to air dry completely before use or storage.
- Do not attempt to remove the tube and bladder from the nylon cuff as this may damage the cuff.
- Take care not to forcefully bend or kink the tubing during storage.
- It is recommended the performance should be checked every 2 years or after a repair. Please contact Customer Care at 866-326-1313.
- It is recommended to clean the cuff every 200 uses.
- The monitor can maintain the safety and performance characteristics for minimum of 10,000 measurements of three years.
- The cuff integrity is maintained after 1,000 open- close cycles of the closure.

HYPERTENSION INDICATOR

WHAT ARE SYSTOLIC PRESSURE AND DIASTOLIC PRESSURE?

When ventricles contract and pump blood out of the heart, the blood pressure reaches its maximum value in the cycle, which is called systolic pressure. When the ventricles relax, the blood pressure reaches its minimum value in the cycle, which is called diastolic pressure.

WHAT IS THE STANDARD BLOOD PRESSURE CLASSIFICATION?

The chart on the right is the standard blood pressure classification published by American Heart Association (AHA).


This chart reflects blood pressure categories defined by the American Heart Association.			
Blood Pressure Category	Systolic mmHg (upper#)		Diastolic mmHg (lower#)
Normal	Less than 120	and	Less than 80
Prehypertension	120 - 129	and	Less than 80
High Blood Pressure (Hypertension) Stage 1	139 - 139	or	80 - 89
High Blood Pressure (Hypertension) Stage 2	140 or higher	or	90 or higher
Hypertension Crisis (Consult your doctor immediately)	Higher than 180	and/or	Higher than 120

CAUTION

Please consult a physician if your measuring result falls outside the range. Please note that only a physician can tell whether your blood pressure value has reached a dangerous point.

IRREGULAR PULSE RATE INDICATOR

An irregular pulse rate will be detected if there is an irregular pulse rhythm while measuring systolic and diastolic blood pressure. When measurements were performed, the monitor will record all pulse intervals and calculate the average. If two or more pulse intervals were recorded, and the difference between each interval and the average is larger than $\pm 25\%$ of the average; or if four or more pulse intervals were recorded, and the difference between each interval and the average is larger than $\pm 15\%$ of the average value, the irregular pulse symbol will be displayed along with measurement results.

CAUTION

The appearance of the IPR icon indicates that a pulse irregularity consistent with an irregular pulse rate was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

TROUBLESHOOTING

This section includes a list of error messages and frequently asked questions for problems you may encounter with your blood pressure monitor. If the product is not operating as you think it should, check here before arranging for servicing.

Problem	Symptom	Check This	Remedy
No Power	Display does not light up	Batteries are exhausted.	Replace with new batteries.
		Batteries are inserted incorrectly.	Check battery polarity.
		Adapter is inserted incorrectly.	Check adapter and replug.
Low Battery	Display is dim or shows 	The battery is too low.	Replace with new batteries.
Error Code	E 01 shows	The cuff is not wrapped or wrapped incorrectly, or the cuff air plug is loose.	Refasten the cuff and insert air tube plug correctly then measure again.
	E 02 or  shows	Excessive body motion (such as shaking of the arm with the cuff on) or weak Pulse is detected.	Relax for 5 minutes, and then keep still, measure again.
	E 03 shows	Pulse is not detected during measuring.	Loosen the clothing on the arm and then measure again.
	E 04 shows	The measurement failed.	Relax for 5 minutes and measure again.
	EExx, shows on the display.	A calibration error occurred (XX can be some digital symbol, such as 01, 02, etc., if this similar situation appear, all belong to calibration error.)	Retake the measurement. If the problem persists, contact the retailer or our customer service department for further assistance. Refer to the warranty for contact information and return instructions.
	E 12 shows	The connection failed.	Open the app, enable Bluetooth and network, then retry.
Warning Message	"Out" shows	Out of measurement range.	Relax for a moment and then measure again. If the problem persists, contact your physician.

DEVICE & LABEL SYMBOLS

SAFETY INFORMATION

The signs below might be in the user manual, labeling or other component. They are the requirement of standard and using.

	Refer to instruction manual/booklet to signify that the instruction manual/booklet must be read. NOTE: Symbol is BLUE		Type BF Applied Part
	Consult instructions for use or consult electronic instructions for use		Serial Number
	Direct Current		Polarity of d.c. power connector
	Class II Equipment		For indoor use only
	Batch Code		General symbol for recovery/recyclable
	Date of Manufacture		Temperature limit
	Atmospheric pressure limitation		Humidity limitation
	MR Unsafe To identify an item which poses unacceptable risks to the patient, medical staff or other persons within the MR environment.		
	Caution Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.		
	The symbol indicates that the product should not be discarded as unsorted waste but must be sent to separate collection facilities for recovery and recycling.		

FCC STATEMENT

NOTE:

POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for U.S.A. only)

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Caution: The user is cautioned that changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment.

This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

ELECTROMAGNETIC COMPATIBILITY

The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments.

Warning: Don't be near the active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Technical description:

1. All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the expected lifetime.
2. Guidance and manufacturer's declaration-electromagnetic emissions and Immunity.

Table 1

Guidance and manufacture's declaration - electromagnetic emissions

Emission test	Compliance	Emission test	Compliance
RF emissions CISPR 11	Group 1	Harmonic emissions IEC 61000-3-2	Class A
RF emissions CISPR 11	Class B	Voltage fluctuations and flicker emissions IEC 61000-3-3	Comply

ELECTROMAGNETIC COMPATIBILITY

Table 2

Guidance and manufacture's declaration - electromagnetic immunity

Immunity test	IEC60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact; ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact; ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines; ± 1 kV signal input/output; 100 kHz repetition frequency	± 2 kV for power supply lines; Not applicable; 100 kHz repetition frequency
Surge IEC61000-4-5	± 0.5 kV, ± 1 kV differential mode; ± 0.5 kV, ± 1 kV, ± 2 kV common mode	± 0.5 kV, ± 1 kV differential mode; Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0% UT; 1 cycle and 70% UT; 25/30 cycles; Single phase: at 0°. 0% UT; 250 / 300 cycle	0% UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0% UT; 1 cycle and 70% UT; 25/30 cycles; Single phase: at 0°. 0% UT; 250 / 300 cycle
Power frequency magnetic field IEC 61000-4-8	30 A/m 50 Hz / 60 Hz	30 A/m 50 Hz / 60 Hz
Conducted RF IEC61000-4-6	3 V; 0,15 MHz – 80 MHz; 6 V in ISM and amateur radio bands; between 0,15 MHz and 80 MHz; 80% AM at 1 kHz	3 V; 0,15 MHz – 80 MHz; 6 V in ISM and amateur radio bands; between 0,15 MHz and 80 MHz; 80% AM at 1 kHz
Radiated RF IEC61000-4-3	10 V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz

NOTE: U_i is the a.c. mains voltage prior to application of the test level.

ELECTROMAGNETIC COMPATIBILITY

Table 3

Guidance and manufacture's declaration - electromagnetic immunity

Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication s equipment)	Test Immunity MHz	Band MHz	Service	Modulation	Modulation Power (W)	Distance (m)	IEC 60601-1-2 Level (V/m)	Compliance Level (V/m)
	385	380-390	TETRA 400	Pulse modulation b) 18Hz	1.8	0.3	27	27
	450	430-470	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0.3	28	28
	710	704-787	LTE Band 13,17	Pulse modulation b) 217 Hz	0.2	0.3	9	9
	745							
	780							
	810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18 Hz	2	0.3	28	28
	870							
	930							
	1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1,3,4,25; UMTS	Pulse modulation b) 217 Hz	2	0.3	28	28
	1845							
	1970							
	2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	2	0.3	28	28
	5240	5100-5800	WLAN 802.11a/n	Pulse modulation b) 217 Hz	0.2	0.3	9	9
	5500							
	5785							

SPECIFICATIONS

Name	Blood Pressure Monitor Premium
Model Number	UAM-880 — Wide Range Cuff
Display System	LCD V.A. 80 mm x 60 mm
Dimensions	Approx. 102 mm x 143 mm x 73 mm
Weight	Approx. 250g (Excluding the batteries and cuff)
Measuring Mode	Oscillographic testing mode
Mode of Operation	Continuous operation
Measuring Range	Rated cuff pressure: 0 mmHg ~ 299 mmHg; Measurement pressure: SYS: 60 mmHg ~ 230 mmHg DIA: 40 mmHg ~ 130 mmHg; Pulse value: (40-199) beat/minute
Accuracy	Static Pressure: 5°-40°C within ±3mmHg; Pulse value: ±5%
Operating Environment	Temperature: +5°C to +40°C Relative humidity: 15% to 90%, non-condensing, but not requiring a water vapor partial pressure greater than 50 hPa; Atmospheric pressure: 700 hPa to 1060 hPa
Storage Environment	Temperature: -20°C to +60°C; Relative humidity: ≤93%, non-condensing, at a water vapor pressure up to 50 hPa; an atmospheric pressure range of 500 hPa to 1060 hPa
Arm Measurement	Wide Range: 8.6" - 16.5" (22- 42 cm) (upper arm circumference)
Degree of Protection	Type BF applied part
Device Classification	Battery Powered Mode: Internally Powered ME Equipment AC Adaptor Powered Mode: Class II ME Equipment
Protection against ingress of water	IP21 It means the device could be protected against solid foreign objects of 12,5mm Φ and greater, and against vertically falling water drops.
Expected Lifetime	Device: 3 years or 10,000 measurements (may vary based on usage conditions); Cuff: 10,000 uses; Alkaline Battery: Approximately 200-300 times
Types of use/reuse	Multiple patient multiple use
Power Source	6V DC 4 x AA batteries; AC adaptor powered mode: (Adapter Not included) Type: BLJ06L060100P-U; Input: 100-240V, 50-60Hz, 0.2A max Output: 6V <u> </u> 1000mA
Software Version	A01
Accessories	Instruction Manual, App Guide, Four (4) AA Batteries, Case

Specifications are subject to change without notice

LIFETIME LIMITED WARRANTY

Congratulations on your purchase of a Blood Pressure Monitor. Your Blood Pressure Monitor is covered by the following limited warranty commencing upon the date of purchase, and subject to the following terms and conditions:

The Warrantor warrants that its Blood Pressure Monitor will be free from defects in materials and workmanship under normal consumer usage for a period of a lifetime for the original purchaser of the product.

Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the Blood Pressure Monitor for commercial purposes or subjecting the Blood Pressure Monitor to abnormal usage or conditions; or (d) other acts which are not the fault of the Warrantor, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the Blood Pressure Monitor.

If the Blood Pressure Monitor fails to conform to this limited warranty, return the Blood Pressure Monitor postage prepaid to: **Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031**. When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone number. Carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. At our option, the Warrantor will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will be responsible for parts charges, if any, and repair charges not covered by this limited warranty.

EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL THE WARRANTOR BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

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PLEASE FILL-IN INFORMATION FOR FUTURE REFERENCE AND ATTACH YOUR RECEIPT BELOW.

This information is necessary should you need to contact Customer Care in the future.

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measure • track • share

Setting up your health monitor is fast and easy. Simply download the app, prepare your device, and start your health journey. Consult your healthcare professional before starting a new health regime, changing medication and for aid in interpreting your personal results.

Bluetooth Smart Compatible
Connects to Android and iOS mobile devices.

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