R.P.M. PRO™

BLOOD PRESSURE MONITOR

Model 210-2111C

INSTRUCTION MANUAL

ENGLISH & SPANISH


Please read this instruction manual
completely before operating this unit.


STOP!
PLEASE ENSURE YOU HAVE ALL OF THE
FOLLOWING COMPONENTS BEFORE USING
YOUR BLOOD PRESSURE MONITOR


Monitor

Adult Cuff and
XXL Adult Cuff

Four (4) AA
Batteries,
Preinstalled

Instruction Manual
Quick Start Guide

INDEX

Care & Safety Information	4-8
Introduction & Indications for Use	10
Blood Pressure Monitor Features	11-12
Getting Started	13
Battery Installation / AC Power	14
Applying Your Blood Pressure Cuff	15-16
Before You Begin	17
Taking Your Blood Pressure Reading	18-21
Irregular Pulse Rate Detection	22
Memory Recall	23
Care & Maintenance	24-25
Troubleshooting	26
Device & Label Symbols	27
FCC Statement	28
Electromagnetic Compatibility	29-32
Specifications	33
Warranty	34
Instrucciones en Español	35-68

Toll-Free
Customer Care Help Line:
1-866-326-1313
Monday – Friday
8:30 a.m. – 4:30 p.m. CST
www.veridianhealthcare.com

RPM PRO
Manufactured for
Veridian Healthcare
1175 Lakeside Drive
Gurnee, IL 60031

Made in China

#93-1672 12/25
©2025 Veridian Healthcare, LLC

CARE & SAFETY INFORMATION

NOTE: Read all instructions carefully before use. The following basic precautions are needed when using an electrical product.

CAUTION: Failure to read and observe all precautions could result in personal injury or equipment damage.

PRECAUTIONS

- Blood Pressure Monitor is intended to be operated by adults, including medical staffs and lay persons. Adult patients could also be intended users or operators.
- This device is intended for indoor, home use and is not intended for self-use in public areas.
- This device is portable, but it is not intended for use during patient transport.
- This device is not suitable for continuous monitoring during medical emergencies or operations.
- This device is intended for non-invasive measuring and monitoring of arterial blood pressure.
- It is not intended for use on extremities other than the arm, or for any purpose other than obtaining a blood pressure measurement.
- This device is intended for adults. Do not use this device on neonates or infants. Do not use it on children and adolescents unless otherwise instructed by a medical professional.
- Consult with your physician before using this monitor if you suffer from the following conditions: common arrhythmias such as premature ventricular beats or atrial fibrillation; peripheral arterial disease; pregnancy or preeclampsia; implantation with electrical devices; undergoing intravascular therapy; arteriovenous shunt or mastectomy. Please note that any of these conditions may affect measurement readings and patient motion, trembling or shivering.
- If you are taking medication, consult your physician to determine the proper time to measure your blood pressure.
- This device may be used only for the intended purpose described in this manual. The manufacturer shall have no liability for any incidental, consequential, or special damages caused by misuse or abuse.

CARE & SAFETY INFORMATION

- Please use the device under the environment provided in the user manual. Otherwise, the device's performance and lifetime will be impacted and reduced.
- The device may require up to 30 minutes to warm up / cool down from the minimum/ maximum storage temperature before it is ready for use.
- The blood pressure monitor, its adapter, and the cuff are suitable for use within the patient environment.
- Do not wash the cuff in a washing machine or dishwasher.
- The device contains sensitive electronic components. To avoid measurement errors, avoid taking blood pressure measurements near a strong electromagnetic field radiated interference signal or electrical fast transient/burst signal.
- Wireless communication equipment, such as wireless home network devices, mobile phones, cordless telephones, and their base stations, walkie-talkies may cause interference that may affect the accuracy of measurements. A minimum distance of 1 foot (30 cm) should be kept from such devices during a measurement.
- Please choose the appropriate cuff according to your arm circumference and physical health.

CAUTIONS

- Do not attempt to repair the unit yourself if it malfunctions. Only have repairs carried out by authorized service centers.
- It is recommended that the performance should be checked after repair, maintenance, and every two years of use, by retesting the requirements in limits of the error of the cuff pressure indication and air leakage (testing at least at 50 mmHg and 200 mmHg). Please contact manufacturer or distributor for authorized service personnel.
- Store your device, cuff, and adapter in a clean and dry place, protect it against extreme moisture, heat, lint, dust, and direct sunlight. Never place any heavy objects on it.
- Make sure the rubber tube of the cuff is not squeezed, stretched, or kinked during storage.
- Dispose of accessories, detachable parts, and the device according to the local guidelines.

CARE & SAFETY INFORMATION

WARNINGS

- DO NOT self-diagnose the measurement results and start treatment by yourself. The measurement results given by this device is not a diagnosis. ALWAYS consult your doctor for evaluation of the results and treatment.
- NOT adjust medication based on readings from this blood pressure monitor. Take medication as prescribed by your physician. ONLY a physician is qualified to diagnose and treat high blood pressure.
- DO NOT apply the cuff on an arm that has an intravenous drip or a blood transfusion attached.
- DO NOT kink, fold, stretch, compress, or otherwise deform the tube during measuring, as the cuff pressure might continuously increase, which could prevent blood flow and result in injury.
- Taking blood pressure measurements too frequently could disrupt blood circulation and cause injuries.
- DO NOT apply the cuff to areas on patient where the skin is delicate or damaged. Check the cuff site frequently for irritation.
- DO NOT place the cuff on the arm of a person whose arteries or veins are undergoing medical treatment, i.e. intra-vascular access or intra-vascular therapy or an arteriovenous (A-V) shunt,
- which could disrupt blood circulation and cause injuries.
- DO NOT place the cuff on the arm on the same side of a mastectomy (especially when lymph nodes have been removed). it is recommended to take measurements on the unaffected side.
- DO NOT wrap the cuff on the same arm to which another monitoring device is applied. One or both devices could temporarily stop functioning if you try to use them at the same time.
- Warning: Please check (for example, by observation of the limb concerned) that the operation
- of the device does not result in prolonged impairment of patient blood circulation.
- Warning: On the rare occasion of a fault causing the cuff to remain fully inflated during measurement, loosen and remove the cuff immediately. Prolonged high-pressure applied to the arm (cuff pressure >300

CARE & SAFETY INFORMATION

mmHg or constant pressure >15 mmHg for more than 3 minutes) might lead to bruising and discolored skin.

- DO NOT use this device with high-frequency (HF) surgical equipment at the same time.
- This device is not used in conjunction with oxygen rich environments, not intended for use with
- flammable anesthetics, not intended for use in conjunction with flammable agents.
- Excessive cuff tube lengths could cause strangulation if you don't manage them properly.
- DO NOT touch output of the batteries/adaptor and the user simultaneously.
- The power cord is considered the disconnect device for isolating this equipment from supply mains. DO NOT position the equipment so that it is difficult to reach or disconnect.
- DO NOT use this device if you are allergic to polyester, nylon, or plastic.
- Only use accessories approved by manufacturer. Using unapproved accessories might cause damage to the unit and injure users.
- If you experience discomfort during a measurement, such as pain in the arm or other complaints, press the Power button immediately to release the air from the cuff.
- DO NOT use the device while under maintenance, or being serviced.
- Sensor degradation or looseness may reduce performance of device or cause other problems.
- The air tube poses a risk of strangulation. Furthermore, the small parts of product and batteries present a choking hazard if swallowed. They should therefore always be kept away from infants/children.

CARE & SAFETY INFORMATION

NOTICE

- You can use this device to take your own measurements, no third-party operator is required.
- Adapter is specified as a part of ME EQUIPMENT.
- At the request of authorized service personnel, circuit diagrams, component part lists, descriptions, and calibration procedures will be made available by the manufacturer or distributor.
- The expected lifetime of the cuff may vary by the frequency of washing, skin condition, and storage state.
- Please report to the manufacturer and the competent authority of the Member State / the FDA in which you are established about any serious incident that has occurred in relation to this device.

ATTENTION!

This manual is intended for instruction on operating this blood pressure monitor as a stand-alone measurement device. Many features of this device require an App or Assistive Technology to fully access features.

Your Blood Pressure Monitor **MAY HAVE** connectivity to a Smartphone App, Alert Device, Hub or other Assistive Technology.

For instructions on how to connect your Blood Pressure Monitor with other devices, please consult the documentation that was provided with the device.

INTRODUCTION & INDICATIONS FOR USE

It is recommended that you first seek the advice and recommendation of your physician or healthcare professional when using home diagnostic devices, including blood pressure monitors.



This blood pressure monitor is intended for use in measuring blood pressure and pulse rate for users with the arm circumference of 8.6" – 16.5" (22–42cm) or 15.7" – 20.5" (40–52cm), and is intended for adults only, excluding those with the condition of diabetes, pregnancy or preeclampsia.

Automatic blood pressure monitors use the oscillometric method to electronically measure your blood pressure. The monitor detects your blood's movement through the artery in your arm and converts the movements into a reading. The oscillometric method does not require a stethoscope, making the monitor ideal for home use.


Blood pressure readings determined with this device are equivalent to measurements obtained by a trained healthcare professional using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers.

Thank you for purchasing an Automatic Blood Pressure Monitor. With proper care and use, your monitor will provide you with many years of reliable readings.

BLOOD PRESSURE MONITOR FEATURES

BLOOD PRESSURE MONITOR FEATURES

Symbol	Description	Explanation
SYS	Systolic Pressure	Top blood pressure
DIA	Diastolic Pressure	Low blood pressure
PULSE/min	Pulse Display	Pulse beats per minute
mmHg	mmHg	Unit of blood pressure measurement
 + 	Low Battery	Battery replacement needed
88:88 AM	Current Tie	Hour : Minute
	Pulse Rate	Pulse rate is shown
	Irregular Pulse Rate	Appears when an irregular pulse is detected
	Blood Pressure Level Indicator	Indicates blood pressure level
	Network Signal Indication	Indicates signal strength
	Three Consecutive Measurements	Device is using 3-measurement option
	Data Transmit Pending	Transmission failed and data is stored up to 60 readings until the device is connected to cellular internet
MEM 88	Memory Recall	Display is showing memory recall function
 Signal linking...  Sending Data  Measure again  Pls go to a good signal  Contact customer service	Warning Messages	Refer to the trouble shooting section of this manual

GETTING STARTED

BEFORE TAKING A MEASUREMENT

1. Your readings should only be interpreted by your physician or healthcare professional with access to your individual medical history. Regular use of a home blood pressure monitor will allow you to track and record your readings for discussion with your physician.
 2. Conduct your measurement in a quiet place while seated in a relaxed position. Rest for 15 minutes before taking your reading. If an error occurs or you desire to take a second reading, allow 15 minutes between readings for your blood vessels to return to normal.
 3. Avoid smoking, eating, taking medication, alcohol, physical activity or any other stressful activity for 30 minutes prior to taking a reading.
 4. Always remove any jewelry or constrictive clothing that may interfere with the cuff placement.
 5. Keep yourself and the monitor still during measuring; do not talk during the reading.
 6. It is recommended that you take your readings at the same time each day to better monitor any indications in your results.
- Wait at least 30 min required for ME equipment to warm from the minimum storage temperature between uses until it is ready for intended use. At least 30 min required for ME equipment to cool from the maximum storage temperature between uses until it is ready for intended use.
 - Before use, make sure the device functions safely and is in proper working condition.

BATTERY INSTALLATION

THIS UNIT REQUIRES 4 AA BATTERIES, PREINSTALLED.


When using for the first time, remove the pull strip to activate the batteries.

It is necessary to replace the batteries when the "Low Battery" symbol + Replace Batteries' appears on the display or when the display does not turn on after the **START STOP** button is pressed.

The battery compartment is located on the underside of this monitor.

1. Remove battery cover.
2. Insert or replace 4 AA batteries into the battery compartment, ensuring to match the indicated polarity symbols. Always use new batteries.
3. Replace the battery cover.
4. Dispose of batteries according to local disposal and recycling regulations.

It is recommended to remove the batteries if the unit will not be used for an extended period of time.

AC WALL POWER; ADAPTER NOT INCLUDED

AC ADAPTER requirement;

USE AUTHORIZED COMPONENT ONLY (Not Included)

Type: BLJ06L060100P-U

Input: AC100-240V, 0.2A Max, 50/60Hz

Output: 6V  1A


- Adapter should only be used during measurement; unplug from wall when not in use.
- This is not a rechargeable monitor; unit must have AC Adapter or 4x AA batteries for power during use.

FITTING & APPLYING YOUR CUFF

PROPER CUFF FIT AND PLACEMENT IS CRITICAL IN OBTAINING ACCURATE BLOOD PRESSURE MEASUREMENTS. PLEASE READ THESE INSTRUCTIONS CAREFULLY AND CONSIDER CONSULTING WITH YOUR PHYSICIAN FOR FURTHER DIRECTION.

Before applying your blood pressure cuff, be sure you have selected the appropriate size cuff.

To determine the proper cuff size, measure the circumference of your left upper arm at the point midway between inside bend of your elbow and your shoulder.

Adult Cuff fits arm circumference 8.6" – 16.5" (22–42cm)

XXL Adult Cuff fits arm circumference 15.7" – 20.5" (40–52cm)

If for any reason you are unable to or should not use your left arm, please modify the instructions for cuff application to your right arm. Your physician can help you identify which arm is best for you to take measurements from.

Assembling Your Cuff

Your cuff may need to be 'threaded' through the metal bar to form the cuff cylinder for your arm.


Fig. 1

1. Lay the open cuff hook and loop side down, size/instruction wording side up on a flat surface.
2. Keeping the words on the inside of the cuff, thread the end of the cuff through the metal bar, Fig 1. The end of the cuff should be able to fold over onto itself for securing via the hook and loop closure, Fig 2.


Fig. 2

FITTING & APPLYING YOUR CUFF

APPLYING YOUR CUFF

Be seated at a table or desk with your feet flat on the floor. To avoid accidentally pulling on or dropping your monitor, the cuff should not be plugged into the monitor until after the cuff is applied to your arm.

The cuffs are printed with a useful guide for applying the cuff; please consult the following steps for additional direction.

1. Remove any constrictive clothing or jewelry that may interfere with cuff placement.
2. Slide the cuff cylinder onto your left arm and align the artery mark over the brachial artery on the inside of your arm, Fig 3.

The tube should follow the inside length of the arm, Fig 4.

3. The bottom edge of the cuff should be positioned approximately 0.5" above the elbow joint, Fig 5.
4. When the cuff is closed around the arm, the triangle guide imprint should be within the OK range bar marking; if the triangle is outside of the range, a different sized cuff is needed
5. The cuff should fit comfortably, yet snugly around your arm, Fig 6. You should be able to insert one finger easily between your arm and the cuff.


Fig. 3


Fig. 4


Fig. 5




Fig. 6

BEFORE YOU BEGIN

INITIAL SET UP:


Network Connection:

- Every time batteries are installed / replaced, the device will start the network connection process. This will take up to 10-minutes.

Successful Connection:

- If the connection is successful, the LCD will display 'Ok' and the device will turn off after 10-seconds. The monitor is now ready for measurement.

Failed Connection:

- If the device fails to connect to the network within 10-minutes, the display will show FAiL and 'Pls go to a good signal' and then turn off automatically.

USER SETTINGS:

- There is no need to set date/time on this device. Once connected or assistive device, the date and time will automatically sync.

TAKING YOUR BLOOD PRESSURE READING

Please read the preceding portions of this manual prior to taking your first reading.

Reminder, it is important to avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading.

1. Position the monitor on a flat, stable surface with the display in view. Be seated with your feet flat on the floor and legs uncrossed.
2. Insert the cuff tubing connector into the port on the left side of your monitor.
3. Rest your elbow on a solid surface with your palm facing upward. Elevate your arm so that the cuff is at the same level as your heart, Fig 1. Relax your left hand.
4. Press and release the **START/STOP** button to start the measurement process.
5. The unit will show the time and 0 on the display.
6. The cuff will automatically inflate; you will see the numbers increase on the display screen and feel the cuff inflate.
7. When the measurement is complete, the cuff will automatically deflate, then your blood pressure measurement and pulse reading results will show on the display screen at the same time.
8. 'Signal Linking' indicates device is looking for signal.
9. 'Sending Data' indicates the results are transmitting.
10. Ok indicates data was transmitted successfully
11. Press **START/STOP** to turn off monitor.

YOU CAN STOP THE INFLATION OR DEFLATION PROCESS AT ANY TIME BY PRESSING THE START/STOP BUTTON.

TAKING YOUR BLOOD PRESSURE READING

SIMPLE STEPS:


BE SEATED, FEET FLAT ON THE FLOOR, CUFF AT THE SAME LEVEL AS YOUR HEART.

Press **START/STOP**; 0 will show on the display.

The numbers will increase as the cuff inflates.

After measurement, results will appear

- 'Signal Linking' indicates device is looking for signal
- 'Sending Date' indicates the results are transmitting
- Ok indicates data was transmitted successfully

TAKING YOUR BLOOD PRESSURE READING



3X CONSECUTIVE READING OPTION: TRIPLE MEASUREMENT MODE

With the 3x Consecutive Measurement feature, this device will take 3 successive measurements and automatically calculate the average of the 3 readings.

1. Position the monitor on a flat, stable surface with the display in view. Be seated with your feet flat on the floor and legs uncrossed.
2. Insert the cuff tubing connector into the port on the left side of your monitor.
3. Rest your elbow on a solid surface with your palm facing upward. Elevate your arm so that the cuff is at the same level as your heart, Fig 1. Relax your left hand.
4. Press the '**TRIPLE MEASUREMENT**' button to enter the Triple Measurement mode.
5. The LCD will show '1' indicating the first measurement is beginning.
6. After measurement 1 is complete, the cuff will deflate and a timer will show, counting down for 30-seconds.
7. The LCD will show '2' indicating the second measurement is beginning.
8. After measurement 2 is complete, the cuff will deflate and a timer will show, counting down for 30-seconds.
9. The LCD will show '3' indicating the final measurement is beginning.
10. After measurement is complete, the cuff will automatically deflate and the average of all three measurements will appear on the display.


TAKING YOUR BLOOD PRESSURE READING

11. 'Signal Linking' indicates device is looking for signal.
12. 'Sending Data' indicates the results are transmitting.
13. Ok indicates data was transmitted successfully
14. Press **START/STOP** to turn off monitor..

Press **START/STOP**;
0 will show on the
display and 1 to
indicate first reading
is taking place.

The display will count down
between measurements, then
show 2 and 3 as readings
complete.

After measurement, results will appear



- 'Signal Linking' indicates device is looking for signal
- 'Sending Date' indicates the results are transmitting
- Ok indicates data was transmitted successfully

IRREGULAR PULSE RATE DETECTION

An irregular pulse rate will be detected if there is an irregular pulse rhythm while measuring systolic and diastolic blood pressure. When measurements were performed, the monitor will record all pulse intervals and calculate the average. If two or more pulse intervals were recorded, and the difference between each interval and the average is larger than $\pm 25\%$ of the average; or if four or more pulse intervals were recorded, and the difference between each interval and the average is larger than $\pm 15\%$ of the average value, the irregular pulse symbol will be displayed along with measurement results.

CAUTION


This monitor will operate if the icon appears, but the results may be compromised, especially if this symbol appears often. Please consult with your physician or trained healthcare professional for further information regarding an irregular pulse rate and if this symbol appears frequently.

MEMORY RECALL

This monitor stores and recalls up to 60 readings. To access the memory:

1. When the monitor is off, press and release the **MEM** button .
2. The display will show the most recent reading; the reading number (1, 2, 3 etc) and time (alternating) on the lower left of the display.
3. Pres **MEM** to scroll through the records.
4. After 60 records, the oldest will be replaced.
5. Press **START/STOP** to turn off the monitor at any time; it will automatically shut off after 10 seconds.

CARE & MAINTENANCE

Proper care and maintenance of your blood pressure monitor will help ensure long life for the product; improper use of the monitor may void the included warranty.

Monitor Care and Maintenance

- Do not use any liquids on the monitor; use a soft, dry cloth to clean the monitor as needed.
- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: -4°F – 140°F (-20°C – 60°C); Humidity: ≤93% RH, Atmospheric pressure: 500 hPa - 1060 hPa
- Avoid high temperatures and direct sunlight.
- Always remove the batteries if the unit will not be used for an extended period of time.
- Do not attempt to disassemble your monitor; disassembling your monitor may void the manufacturers warranty.
- Do not subject the monitor to strong shocks; take care not to drop the monitor.
- No component can be maintained by user in the monitor. The circuit diagrams, component part lists, descriptions, calibrations instructions, or other information which will assist the user's appropriately qualified technical personnel to repair those parts of equipment which are designated repairable can be supplied.

Cuff Care and Maintenance

- Never attempt to iron or dry the cuff; if the cuff becomes wet, allow the cuff to air dry completely before use or storage.
- Do not attempt to remove the tube and bladder from the nylon cuff as this may damage the cuff.
- Take care not to forcefully bend or kink the tubing during storage.

CARE & MAINTENANCE

- It is recommended the performance should be checked every 2 years or after a repair. Please Contact Customer care at 1-866-326-1313.
- It is recommended to clean the cuff every 200 uses.
- The monitor can maintain the safety and performance characteristics for minimum of 10,000 measurements of three years.
- The cuff integrity is maintained after 1,000 open- close cycles of the closure.

Cleaning

1. Before cleaning, ensure monitor is off and disconnected from any external power sources.
2. Clean the cuff with a soft dampened with mild soapy water until no visible contaminants remain.
3. Rinse the cuff, ensuring no cleaning solution remains.
4. Wipe off with a dry cloth and allow to air dry thoroughly before storage or next use.

Disinfection

1. Before disinfection, ensure monitor is off and disconnected from any external power sources.
2. Disinfect the cuff with a soft dampened with 70% ethyl rubbing alcohol for at least three minutes.
3. Rinse the cuff, ensuring no disinfecting solution remains; wipe with a clean damp towel.
4. Wipe off with a dry cloth and allow to air dry thoroughly before storage or next use.

TROUBLESHOOTING

Issue	Meaning	Review	Solution
No Power	Display does not light up	Batteries are exhausted Batteries are inserted incorrectly Adapter is inserted incorrectly	Replace with new batteries Check battery polarity Check adapter and replug
High Battery	H bAt on display	Adapter voltage is higher than 7.5V	Replace with authorized adapter
Low Battery	Display is dim or  + Replace batteries	The battery is too low	Replace with new batteries
Error Code	E1 + Measure again	The cuff is not wrapped or wrapped incorrectly, or the cuff air plug is loose.	Refasten the cuff and insert air tube plug correctly then measure again.
	E2 + Measure again	Excessive body motion (such as shaking of the arm with the cuff on) or weak pulse is detected during measuring.	Relax for 5 minutes, and then keep still, measure again.
	E3 + Measure again	Pulse is not detected during measuring.	Loosen the clothing on the arm and then measure again.
	E4 + Measure again	The measurement failed.	Relax for 5 minutes and measure again.
	E5 + Measure again	Abnormal communication with server / transmit fail	Attempt to find a new signal; contact Customer Service
	E6 + Pls go to a good signal	Cannot connect to internet	Contact Customer Service
	E8 + Pls go to a good signal	SIM card detection fail (attempted 3 times)	Turn on and off, attempt new measurement
	EExx	Hardware error (xx can be some symbol, such as 01, 02, etc.)	Turn off monitor and measure again. If EExx still appears on the display, please contact customer service.
	Err + UsB	Adapter error (the voltage is too high or too low)	Replace with the authorized adapter.
	Out	Out of measurement range	Relax for a moment. Refasten the cuff and then measure again. If the problem persists, contact your physician.

DEVICE & LABEL SYMBOLS

These symbols may appear on your device, instructions or packaging and may vary by make and model.

Symbol	Meaning	Symbol	Meaning
	Refer to instruction manual/booklet to signify that the instruction manual/booklet must be read. NOTE: Symbol is BLUE		Type BF Applied Part
	Date of manufacture		Serial Number
	Direct Current		Polarity of d.c. power connector
	Class II Equipment		For Indoor Use Only
	Batch Code		Manufacturer
	Temperature Limit		Humidity Limit
	Atmospheric Pressure Limitation		This symbol shows that this device can be expected to comply with part 15 of the FCC Rules
	General symbol for recovery/recyclable		Consult instructions prior to use
	MR Unsafe To identify an item which poses unacceptable risks to the patient, medical staff or other persons within the MR environment.		
	Caution Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.		
	The symbol indicates that the product should not be discarded as unsorted waste but must be sent to separate collection facilities for recovery and recycling.		

FCC STATEMENT

FCC ID: OU9BB-1597-G

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Caution: The user is cautioned that changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC Radiation Exposure Statement:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with a minimum distance of 20cm between the radiator and your body. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

ELECTROMAGNETIC COMPATIBILITY

Guidance and manufacture's declaration - electromagnetic emissions

The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments.

Essential Performance

Accuracy of measuring blood pressure pulse rate

Measurement Range: Systolic pressure: 60-230 mmHg; Diastolic pressure: 40-130 mmHg
Pulse: 40-199 beats/minute

Accuracy Pressure: ± 3 mmHg; Pulse $\pm 5\%$

The Basic Safety of the Blood Pressure Monitor (210-2111C) is as following: Deviation from normal operation that poses an unacceptable risk to the patient or operator.

Warning: Don't be near the active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Technical description:

1. All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the expected lifetime.
2. Guidance and manufacturer's declaration-electromagnetic emissions and Immunity.

Table 1

Guidance and manufacture's declaration - electromagnetic emissions

Emission test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations and flicker IEC 61000-3-3	Comply

ELECTROMAGNETIC COMPATIBILITY

Table 2

Immunity test	IEC60601-1-2 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact; ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact; ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines; ± 1 kV signal input/output; 100 kHz repetition frequency	± 2 kV for power supply lines; Not applicable; 100 kHz repetition frequency
Surge IEC61000-4-5	± 0.5 kV, ± 1 kV differential mode; ± 0.5 kV, ± 1 kV, ± 2 kV common mode	± 0.5 kV, ± 1 kV differential mode; Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0% UT; 1 cycle and 70% UT; 25/30 cycles; Single phase: at 0°. 0% UT; 250 / 300 cycle	0% UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0% UT; 1 cycle and 70% UT; 25/30 cycles; Single phase: at 0°. 0% UT; 250 / 300 cycle
Power frequency magnetic field IEC 61000-4-8	30 A/m 50 Hz / 60 Hz	30 A/m 50 Hz / 60 Hz
Conducted RF IEC61000-4-6	3 V; 0,15 MHz – 80 MHz; 6 V in ISM and amateur radio bands; between 0,15 MHz and 80 MHz; 80% AM at 1 kHz	3 V; 0,15 MHz – 80 MHz; 6 V in ISM and amateur radio bands; between 0,15 MHz and 80 MHz; 80% AM at 1 kHz
Radiated RF IEC61000-4-3	10 V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz

NOTE: U_i is the a.c. mains voltage prior to application of the test level.

ELECTROMAGNETIC COMPATIBILITY

Table 3

Guidance and manufacture's declaration - electromagnetic emissions

Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication s equipment)	Test Immunity MHz	Band MHz	Service	Modulation	Modulation Power (W)	Distance (m)	Immunity Test Level (V/m)	Compliance Level (V/m)
	385	380-390	TETRA 400	Pulse modulation b) 18Hz	1, 8	0, 3	27	27
	450	430-470	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0, 3	28	28
	710	704-787	LTE Band 13,17	Pulse modulation b) 217 Hz	0, 2	0, 3	9	9
	745							
	780							
	810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18 Hz	2	0, 3	28	28
	870							
	930							
	1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1,3,4,25; UMTS	Pulse modulation b) 217 Hz	2	0, 3	28	28
	1845							
	1970							
	2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	2	0, 3	28	28
	5240	5100-5800	WLAN 802.11a/n	Pulse modulation b) 217 Hz	0, 2	0, 3	9	9
	5500							
	5785							

SPECIFICATIONS FOR 4G

SPECIFICATIONS FOR 4G

E-MTC	Supporting Band:	Cat-M1
		B2/4/12/13/25
	Transmitted frequency range	1850—1910MHz, 1710—1755MHz, 699—716MHz, 777—787MHz, 1850—1915MHz
	Transmitted power	20±2dBm
	Throughput	588~1119 kbps
	Latency	≤21ms
	Jitter	≤21ms
	PER	<5%
	Data Integrity Accessibility	Data shall be transmitted correctly and completely Accessibility is high since 4G is broadband

Remark

1. The Quality of Service (QOS) is considered as KPI's here.
2. For the Blood Pressure Monitor, the E-MTC is used to transmit the control data between the EUT and the controller, the manufacture considered 37~1119 kbps is a suitable throughput for the usage.
3. The PER of the wireless system is normally less than 1% for Semi-anechoic chamber, we considered the PER cannot be greater than 10%.
4. Use the wireless module which the EUT used as the test model, instead the EUT.

WARNING

Failure to comply with the following warnings may cause cyber security issues.

DO NOT casually use the blood pressure monitor for others, as it may lead to leakage of personal measurement data.

DO NOT connect the blood pressure at places with poor network or wireless connectivity as it may interrupt data upload and lead to deficient data.

PRODUCT SPECIFICATIONS

Name	Blood Pressure Monitor
Model Number	210-2111C
Display System	LCD 73 mm x 49 mm
Dimensions	Approx. 118 mm x 126 mm x 72 mm
Weight	Approx. 305g (Excluding the batteries and cuff)
Measuring Mode	Oscillographic testing mode
Mode of Operation	Continuous operation
Measuring Range	Rated cuff pressure: 0 mmHg ~ 299 mmHg; Measurement pressure: SYS: 60 mmHg ~ 230 mmHg DIA: 40 mmHg ~ 130 mmHg; Pulse value: (40-199) beat/minute
Accuracy	Static Pressure: 5°C - 40°C within ± 3 mmHg; Pulse value: $\pm 5\%$ Clinical validation: Mean difference within ± 5 mmHg; Standard deviation ≤ 8 mmHg
Operating Environment	Temperature: +5°C to +40°C Relative humidity: 15% to 90%, non-condensing, but not requiring a water vapor partial pressure greater than 50 hPa; Atmospheric pressure: 700 hPa to 1060 hPa
Storage Environment	Temperature: -20°C to +60°C; Relative humidity: $\leq 93\%$, non-condensing, at a water vapor pressure up to 50 hPa; atmospheric pressure range of 500 hPa to 1060 hPa
Arm Measurement	Upper arm circumference: Adult Cuff: 8.6" – 16.5" (22–42cm) XXL Adult Cuff: 15.7" – 20.5" (40–52cm)
Degree of Protection	Type BF applied part
Device Classification	Battery Powered Mode: Internally Powered ME Equipment AC Adaptor Powered Mode: Class II ME Equipment
Protection against ingress of water	IP21 It means the device could be protected against solid foreign objects of 12,5mm Φ and greater, and against vertically falling water drops.
Expected Lifetime	Device: 3 years or 30,000 measurements (may vary based on usage conditions); Cuff: 15,000 uses; Battery: Approximately 200-300 times
Types of use/reuse	Single patient multiple use
Power Source	6V 4x AA batteries
Software Version	A01
Accessories	Instruction Manual, Quick Start Guide; Four (4) AA Batteries

Specifications are subject to change without notice

TWO-YEAR LIMITED WARRANTY

Congratulations on your purchase of a Blood Pressure Monitor. Your Blood Pressure Monitor is covered by the following limited warranty commencing upon the date of purchase, and subject to the following terms and conditions:

The Warrantor warrants that its Blood Pressure Monitor will be free from defects in materials and workmanship under normal consumer usage for a period of two years for the original purchaser of the product.

Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the Blood Pressure Monitor for commercial purposes or subjecting the Blood Pressure Monitor to abnormal usage or conditions; or (d) other acts which are not the fault of the Warrantor, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the Blood Pressure Monitor.

If the Blood Pressure Monitor fails to conform to this limited warranty, return the Blood Pressure Monitor postage prepaid to: **Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031**. When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone number. Carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. At our option, the Warrantor will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will be responsible for parts charges, if any, and repair charges not covered by this limited warranty.

EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL THE WARRANTOR BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

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