SMARTHEART®

Deluxe Pulse Oximeter

Model #11-51DA

INSTRUCTION MANUAL ***ENGLISH AND SPANISH***

**PLEASE READ THESE INSTRUCTIONS
CAREFULLY TO ENSURE PROPER USE AND
SAFE OPERATION OF YOUR PULSE OXIMETER.**

STOP!
PLEASE ENSURE YOU HAVE ALL OF THE
FOLLOWING COMPONENTS BEFORE
USING YOUR PULSE OXIMETER.

Pulse Oximeter, Lanyard &
Storage Case

2 AAA
Batteries

Instruction
Manual

IF YOU ARE MISSING ANY PARTS,
INCLUDING INSERTS OR INSTRUCTION
MANUALS, DO NOT RETURN TO PLACE OF
PURCHASE. CONTACT CUSTOMER CARE AT
866-326-1313.

INDEX

Care & Safety Information	4-7
Introduction	8
Pulse Oximeter Features	9
How Pulse Oximeter Works	9
Display Features.....	10
How to use Pulse Oximeter	11-12
System Settings	13
Display Brightness Setting	13
Battery Installation	14
Lanyard Information	15
Care & Maintenance.....	15
Troubleshooting	16
Device & Label Symbols.....	17
FCC Statement.....	18
Electromagnetic Compatibility	19-23
Specifications.....	24
Warranty	25-26
Instrucciones en Español	27-52

**Toll-Free
Customer Care Help Line:
1-866-326-1313
Monday – Friday
8:30 a.m. – 4:30 p.m. CST**

SmartHeart®
Manufactured for
Veridian Healthcare
5650 Centerpoint Court
Gurnee, IL 60031

Made in China

#93-1737 07/2026
©2026 Veridian Healthcare, LLC

CARE & SAFETY INFORMATION

GENERAL CAUTIONS AND WARNINGS

- Keep unit and lanyard out of reach of children as the included lanyard may present an entanglement or choking hazard to small children. Adult supervision required; never leave children unattended with unit or lanyard.
- To avoid electrical shock keep the unit away from water; never immerse the unit in liquid as this will damage the unit.
- Do not handle the unit with wet hands.
- Do not use the unit if it has been submersed in water or dropped.
- If any abnormality in product function occurs, discontinue use until the unit has been examined; discontinue use and return to the manufacturer for review and replacement under the terms of the one-year warranty.
- Do not shake the unit when in operation.
- Do not use while actively exercising; not suitable for continuous monitoring.
- Do not use during an MRI or CT scan. Do not use while taking a blood pressure measurement

CARE & SAFETY INFORMATION

- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Please follow local ordinances and recycling instructions regarding disposal or recycling of the equipment and batteries.

OPERATING CAUTIONS AND WARNINGS

- Never tie or tape Pulse Oximeter to your finger.
- Explosion hazard: do not use the product in an explosive atmosphere.
- The part expected to come into contact with the human body has been biologically evaluated and the biological safety meets the standard requirements, very few people may have allergic reactions, patients with allergic reactions should stop using it immediately.

CARE & SAFETY INFORMATION

STORAGE CAUTIONS AND WARNINGS

- Never store the Pulse Oximeter the unit direct sunlight, high temperature or humidity.
- Keep the unit out of reach of small children.

CLEANING CAUTIONS AND WARNINGS

- Never immerse the device in water to clean as it may be damaged.
- Do not disassemble, install, or repair the device.
- Do not look directly at the light-emitting devices, which may cause eye damage.
- Follow the included instructions for cleaning the unit after each use; if the unit has not been used for an extended period, the device should be cleaned prior to use.
- The part expected to come into contact with the human body has been biologically evaluated and the biological safety meets the standard requirements, very few people may have allergic reactions, patients with allergic reactions should stop using it immediately.

CARE & SAFETY INFORMATION

CAUTIONS

The equipment is designed to measure the percentage of arterial oxygen saturation of functional hemoglobin. Factors that may degrade pulse oximeter performance or affect the accuracy of the measurement include the following:

- Do not apply the pulse oximeter on the same arm as a blood pressure cuff, arterial catheter or infusion line(s)(IVs)
- Excessive light, such as sunlight or direct home lighting
- Moisture in the equipment
- Finger is beyond recommended size range
- Poor pulse quality
- Venous pulsations
- Anemia or low hemoglobin concentrations
- Cardio green and other intravascular dyes
- Carboxyhemoglobin
- Methemoglobin
- Dysfunctional hemoglobin
- Artificial nails or fingernail polish
- Excessive movement during testing may result in incorrect readings
- Improper wearing position may result in incorrect reading

INTRODUCTION

Please read all instructions carefully prior to first use.


Thank you for purchasing a Pulse Oximeter. With proper care and use, your Pulse Oximeter will provide you with many years of reliable service. The Pulse Oximeter measures user blood oxygen saturation (SpO₂) and pulse rate information for wellness use only in a small, portable, easy-to-use device.

The Pulse Oximeter is intended for use for persons engaging in sporting and aviation activities. Persons who need SpO₂ and pulse rate measurements because of a medical condition should not use the Pulse Oximeter and should consult with their physician.

Some studies have indicated that certain physical traits, including, but not limited to darker skin pigmentation, finger thickness and peripheral perfusion may create variations in readings that make the oxygen level appear higher than it actually is. Users should establish a baseline value when healthy and at rest. It is recommended that users consult a licensed healthcare provider on how to interpret user readings based on physical traits or characteristics.


**PLEASE READ THIS INSTRUCTION MANUAL
COMPLETELY BEFORE OPERATING THIS UNIT.**

PULSE OXIMETER FEATURES


HOW PULSE OXIMETER WORKS

The Pulse Oximeter works by shining two light beams into the small vessels or capillaries of the finger to determine the user's blood oxygen saturation percentage and pulse rate that is then displayed on the Pulse Oximeter's screen. The unit uses the phrase "SpO2" to identify the oxygen saturation level and displays it as a percentage of full capacity (100%). Pulse rate is identified by the phrase "PR bpm" and is displayed as beats per minute.


DISPLAY FEATURES

When in use, the Pulse Oximeter display reflects a 'live' reading of the wearers current SpO2 and PR bpm.

SCREEN DISPLAY ORIENTATION

The Pulse Oximeter display can be adjusted to allow for easy viewing from any angle. While the unit is ON, rotating the finger while the oximeter is on and measuring will automatically change the display orientation, as shown below.

HOW TO USE PULSE OXIMETER

USING YOUR PULSE OXIMETER

The Pulse Oximeter may be used while the user is seated, standing or lying down. The user should not be walking or running during measurements and should take care not to excessively move the arm or hand where the oximeter is attached.

PRIOR TO USE

Keep the sensors clean to help ensure accurate readings. Before and after each use, carefully swab inside the unit with a soft cloth or cotton swab with rubbing alcohol. Do not pour the alcohol directly on or into the unit.


Dirty hands or nail polish (especially dark shades) can also affect reading accuracy. Washing your hands and/or removing any nail polish prior to use helps ensure accurate readings.

Suitable for use on any finger (excluding the thumb). In most instances, the index or middle finger work best.


1. Squeeze the hinged end of the unit to open the measurement side of the unit.
2. Insert finger nail side up into the unit until resistance is felt. The unit should gently close below the first knuckle of the finger, and the sensors

HOW TO USE PULSE OXIMETER

inside the unit should align with the finger nail and pad at the tip of the finger.

3. Press the power button to turn the unit on. The screen will briefly flash 'finger out'.

4. Solid dashes (— —) will appear on the display for a few seconds, then the display will begin to show the oxygen saturation (SpO2) and pulse rate (PR bpm) measurements. The display will vary depending on which display mode orientation is in use. Turning the user's hand will rotate the display for best view.

5. Take care to not move excessively during monitoring. These are active readings and may fluctuate slightly as you observe the display, responding to your heart rate (bps) and other factors.
6. The pulse oximeter should only be used for spot-checking blood oxygen saturation and pulse rate and should not remain on the finger for more than a few minutes at a time.
7. Remove the unit from the finger. The unit will automatically turn off after 10 seconds.
8. Clean the unit as directed before storage.

SYSTEM SETTINGS

This device includes the option to activate an audible beep during pulse measurements.

1. When the device is OFF, press and hold the POWER button to enter the Settings menu.
2. Press and release the POWER button to scroll through the menu.
3. When Beep is highlighted, press and hold POWER to toggle between On and Off.
4. When the desired setting is shown, release the POWER button, then press and release until Exit is highlighted. Press and hold POWER to exit the settings menu.

Settings		
Beep	*	On
Language		EN
Default		
Exit		

DISPLAY BRIGHTNESS SETTING

This device has up to 5 brightness levels.

1. When the device is ON, press and hold the POWER button to enter the brightness menu.
2. Press and release POWER to scroll through the brightness menu. When desired level is shown, release the button; the display will revert to the live read after a few moments.

BATTERY INSTALLATION

BATTERY INSTALLATION AND REPLACEMENT

This unit comes complete with 2 AAA batteries.

It is necessary to replace the batteries when the Low Battery symbol appears on the display or when the display does not turn on after the power button is pressed.

1. Slide off the battery cover.
2. Insert or replace 2 AAA batteries into the battery compartment, ensuring to match the indicated polarity symbols.
3. Replace the battery cover.
4. It is recommended to remove the batteries if the unit will not be used for an extended period of time.


Always use new batteries. Do not mix alkaline, standard or rechargeable batteries. Please dispose of batteries according to local recycling regulations.

LANYARD INFORMATION

QUICK-RELEASE LANYARD


The lanyard features a quick-release clip for easy release while in use.

CHOKING HAZARD-The included lanyard may present an entanglement or choking hazard to small children. Adult supervision required; never leave children unattended with unit or lanyard.

CARE & MAINTENANCE

Before and after each use, carefully swab inside the unit with a soft cloth or cotton swab with rubbing alcohol. Do not pour the alcohol directly on or into the unit.

The unit should be stored in a dry environment, within the ideal temperature range of -4°F–131°F (-20°C–60°C).

TROUBLESHOOTING

SpO2 or PR is not displaying a reading though the power is on

- The user's oxygen level is too low to be measured; consult your physician if the measurement cannot be read
- Reinsert the finger, remain still

SpO2 or PR display readout is unstable and fluctuating highly

- Finger may not be inserted deep enough
- Finger may be trembling, try to remain still

The unit will not turn on

- Ensure the batteries have been properly installed
- Replace the batteries
- If the unit will not turn on after changing batteries, consult Customer Care

Display screen goes blank

- The unit will automatically shut-off after 10 seconds if no finger is inserted
- Check the batteries and replace if needed

DEVICE & LABEL SYMBOLS

These symbols may appear on your device, instructions or packaging and may vary by make and model:

SpO₂ Blood Oxygen Saturation Value	 Serial Number
PRbpm Pulse Rate (bpm)	 Type BP Applied Parts
 Read This Manual	 Environment Protection
 Low Battery	 Manufacturer
 Warning or Caution	 Manufacturer Date
IP22 Protected against dripping water	 Magnetic Resonance Unsafe
 Audio Beep ON	
 Audio Beep OFF	

FCC STATEMENT

1.6 Federal communications commission compliance statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable/mobile exposure condition without restriction.

ELECTROMAGNETIC COMPATIBILITY

This device satisfied the EMC requirements of the international standard IEC 60601-1-2. The requirements are satisfied under the conditions described in the table below. The device is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the instructions for use. Portable and mobile HF communications equipment can affect the device negatively and alter the electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

TABLE 1

**GUIDANCE & MANUFACTURE'S DECLARATION
ELECTROMAGNETIC EMISSION**

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment – guidance
Conducted emission CISPR 11	Not applicable	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Radiated emission CISPR 11	Group 1 Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	

ELECTROMAGNETIC COMPATIBILITY

TABLE 2

GUIDANCE & MANUFACTURE'S DECLARATION ELECTROMAGNETIC IMMUNITY

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance level
Electrostatic discharge(ESD)* IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrostatic transient / burst IEC 61000-4-4	± 2 kV for power supply lines; 100 kHz repetition frequency; ± 1 kV for input/output lines	N/A
Surge IEC 61000-4-5	± 0.5 kV, ± 1 kV differential mode line-line	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT (100 % dip in UT) for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0 % UT (100 % dip in UT) for 1 cycle at 0° 70 % UT (30 % dip in UT) for 25/30 cycles at 0° 0 % UT (100 % dip in UT) for 250/300 cycle at 0°	N/A
Power frequency** (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m, 50/60Hz	30 A/m, 50/60Hz

EMC environment - guidance

***Electrostatic Discharge** -- Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.

****Power frequency** -- magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: UT is the a. c. mains voltage prior to application of the test level.

ELECTROMAGNETIC COMPATIBILITY

TABLE 3

Guidance and manufacture's declaration-electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The user of the device should assure that it is used in such an environment.

Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment guidance
Conducted RF ICE 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms 150 kHz to 80 MHz outside ISM bands	N/A	Portable and mobile RF communication equipment should be used no closer to any part of the Pulse Oximeter, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF ICE 61000-4-3	10 V/m 80 MHz to 2.7 GHz	10 V/m	Recommended separation distance $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80\text{MHz to } 800\text{MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800\text{MHz to } 2.7\text{GHz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

ELECTROMAGNETIC COMPATIBILITY

TABLE 3 -- CONTINUED

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a. The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.
- b. The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.
- c. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.
- d. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

ELECTROMAGNETIC COMPATIBILITY

TABLE 4

RECOMMENDED SEPARATION DISTANCES BETWEEN RF WIRELESS COMMUNICATIONS EQUIPMENT

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between RF wireless communications equipment and the device as recommended below, according to the maximum output power of the communications equipment.

Frequency (MHz)	Maximum Power W	Distance	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
385	1.8	0.3	27	27	RF wireless communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
450	2	0.3	28	28	
710	0.2	0.3	9	9	
745					
780					
810	2	0.3	28	28	Recommended separation distance $E = \frac{6}{d} \sqrt{P}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitter, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
870					
930					
1720	2	0.3	28	28	
1845					
1970					
2450	2	0.3	28	28	
5240	0.2	0.3	9	9	
5500					
5785					

NOTE 1: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

WARNINGS!

This device should not be used in the vicinity or on the top of other electronic equipment such as cell phone, transceiver or radio control products. If you have to do so, the device should be observed to verify normal operation.

• The use of accessories and power cord other than those specified, with the exception of cables sold by the manufacturer of the equipment or system as replacement parts for internal components, may result in increased emissions or decreased immunity of the equipment or system.

PRODUCT SPECIFICATIONS

NAME	Deluxe Pulse Oximeter
MODEL NUMBER	11-51DA
POWER SOURCE	2 "AAA" batteries
SPO2 MEASURING RANGE	70-100%
SPO2 ACCURACY	70-100%, $\pm 2\%$: $\leq 69\%$ no definition
PULSE RATE MEASURING RANGE	30-240 bpm
PULSE RATE ACCURACY	30-100 bpm, ± 2 bpm; 101-240 bpm, $\pm 2\%$
AUTOMATIC SHUT-OFF	After 5 seconds of no indication on the sensors
OPERATION ENVIRONMENT	Temperature 50°F – 104°F (10°C – 40°C); Humidity 15% – 93% Atmospheric Pressure: 700hPa ~1060hPa
STORAGE ENVIRONMENT	Temperature -4°F – 131°F (-20°C – 55°C); Humidity 15% – 93% Atmospheric Pressure: 500hPa ~1060hPa
UNIT DIMENSIONS	60mm x 32mm x 31.4mm (2.36" x 1.26" x 1.24")
UNIT WEIGHT	50g oz. (without batteries)

Specifications are subject to change without notice

TWO-YEAR LIMITED WARRANTY

Congratulations on your purchase of a Pulse Oximeter. Your Pulse Oximeter is covered by the following limited warranty commencing upon the date of purchase for a period of two years, and subject to the following terms and conditions:

The warrantor warrants that its Pulse Oximeter will be free from defects in materials and workmanship under normal consumer usage for one year as the original purchaser owns the product.

Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the Pulse Oximeter for commercial purposes or subjecting the Pulse Oximeter to abnormal usage or conditions; or (d) other acts which are not the fault of the warrantor, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the Pulse Oximeter.

If the Pulse Oximeter fails to conform to this limited warranty, return the Pulse Oximeter to: **Attn: Repair Department, 5650 Centerpoint Court, Gurnee, IL, 60031**. When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone number. Carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. At our option, The warrantor will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will be responsible for parts charges, if any, and repair charges not covered by this limited warranty.

EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE

TWO-YEAR LIMITED WARRANTY

CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL THE WARRANTOR BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

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SmartHeari®
Manufactured for
Veridian Healthcare
5650 Centerpoint Court
Gurnee, IL 60031

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